Avantgarde T, Ti8.9, TiUltra, VR, XXL 2

<input type="checkbox"/> DE	Gebrauchsanweisung	3
<input type="checkbox"/> EN	Instructions for use	55

DE

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5.5.9 Elevating footrest (Fig. 27)	73
5.5.10 Elevating footrest with automatic length adjustment (Fig. 28)	73
5.5.11 Amputation legrest (Fig. 29)	73
5.6 Seat (Fig. 30/31)	73
5.6.1 Seat upholstery with pocket for catheter (Fig. 30).....	73
5.6.2 "The Seat" (Fig. 31).....	73
5.7 Back (Fig. 32–35)	74
5.7.1 Standard back, height-adjustable/adaptable (Fig. 32).....	74
5.7.2 Ergo back (Fig. 33)	74
5.7.3 "The Back" (Fig. 34).....	74
5.7.4 Back stabilizer (Fig. 35).....	74
5.8 Side panels (Fig. 36–39)	75
5.8.1 Plug-on clothing protector side panel (Fig. 36).....	75
5.8.2 Side panels, adjustable, with protection against cold (Fig. 37).....	75
5.8.3 Aluminium clothing protector side panels with neoprene insert (Fig. 38)...	75
5.8.4 Carbon clothing protector side panels (Fig. 39)	75
5.8.5 Desk side panels (Fig. 40).....	76
5.9 Push handles (Fig. 41/42).....	76
5.9.1 Push handles, height-adjustable (Fig. 41)	76
5.9.2 Telescoping push handles (Fig. 42).....	76
5.10 Caster wheel adapter/caster wheel accessories (Fig. 43–46).....	77
5.10.1 Caster wheel adapter for large anterior seat heights (Fig. 43)	77
5.10.2 Caster wheel adapter for low anterior seat heights (Fig. 44).....	77
5.10.3 Caster wheel with quick-release axle (Fig. 45).....	77
5.10.4 Caster wheels (Fig. 46).....	77
5.11 Rear wheels (Fig. 47/48)	78
5.12 Brake/wheel lock systems/wheel lock accessories (Fig. 49/50)	78
5.12.1 Rear wheel with drum brake (not illustrated)	78
5.12.2 Wheel lock, swing-away (Fig. 49).....	78
5.12.3 Wheel lock for one-handed operation (Fig. 50)	79
5.12.4 Wheel lock for user and attendant (not illustrated).....	79
5.12.5 Plug-on wheel lock lever extension (not illustrated)	79
5.13 One-arm drive (Fig. 51)	79
5.14 Wheelbase extension (Fig. 52)	79
5.15 Spoke protector (Fig. 53).....	80
5.16 Shock absorber system (Fig. 54)	80
5.17 Safety belt (Fig. 55)	80
5.18 Tray (Fig. 56)	80
5.19 Further instructions for use	81
6 Adjustment/assembly instructions	82
6.1 Adjusting the wheel lock (Fig. 57)	82

6.2	Changing the rear wheel adapter position in the frame (Fig. 58)	83
6.3	Changing the rear wheel position in the rear wheel adapter (Fig. 59–61)	84
6.3.1	Changing the position in the rear wheel adapter Avantgarde T/Ti 8.9/VR	84
6.3.2	Changing the position in the rear wheel adapter Avantgarde XXL 2	85
6.4	Setting the wheel camber (Fig. 62/63)	85
6.4.1	Setting the rear wheel camber to 0° (Fig. 62 A).....	86
6.4.2	Setting the rear wheel camber to 2° (Fig. 52 B)	86
6.4.3	Setting the rear wheel camber to 4° (Fig. 62 C).....	87
6.5	Adjusting the caster wheel journal angle (Fig. 64/65)	87
6.6	Adjusting the footplate angle (Fig. 66/67).....	88
6.6.1	Adjusting the footplate angle Avantgarde T/VR/XXL 2 (Fig. 66)	88
6.6.2	Adjusting the footplate angle Avantgarde Ti 8.9/Ti Ultra (Fig. 67)	88
6.7	Adjusting the footrest (Fig. 68-69)	88
6.7.1	Adjusting the footrest height Avantgarde T/Ti 8.9/Ti Ultra (Fig. 68).....	88
6.7.2	Adjusting the footrest height Avantgarde VR and XXL 2 (Fig. 69).....	88
6.7.3	Adjusting the height of the footrest for short lower leg lengths (Fig. 70).....	89
6.7.4	Adjusting the height of the elevating footrest with automatic length	89
6.7.5	Adjusting the height of the elevating footrest Avantgarde VR/XXL 2.....	89
6.7.6	Adjusting the angle of the elevating footrest Avantgarde VR/XXL 2	90
6.7.7	Swinging away the elevating Avantgarde FL footrest (Fig. 74).....	90
6.8	Adjusting the brake force – drum brake (Fig. 75)	90
6.9	Tension adjustable seat upholstery (Fig. 76)	91
6.10	Adapting the back upholstery (Fig. 77)	91
6.11	Adjusting the back height (Fig. 78–79).....	92
6.12	Mounting the anti-tipper (Fig. 80–83).....	93
6.13	Attaching the safety belt (Fig. 84/85)	94
7	Maintenance, cleaning, and disinfection	94
7.1	Maintenance	94
7.2	Cleaning and disinfection.....	95
8	Changing tyres (Fig. 86–89)	96
9	Appendices.....	98
9.1	Technical data	98
9.2	Threshold values for wheelchairs transportable by train	102
10	Re-use / disposal	104
10.1	Information on re-use.....	104
10.2	Information on disposal	104

1 General information

INFORMATION

Date of the last update: 2015-12-21

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

INFORMATION

- You can request this document as a PDF file from the Customer Care Center (CCC) at oa@ottobock.com or from the manufacturer's service department (see inside back cover or back page for addresses).
- It is possible to increase the display size of the PDF document.
- For further questions about the instructions for use, please contact the authorised personnel who issued the product to you.

1.1 Preface

With the purchase of the Avantgarde wheelchair for active use, you have selected a quality product for versatile use during everyday life, both indoors and outdoors.

The chapter "Options" presents various additional components of the Avantgarde, which can expand its field of application and improve comfort for the wheelchair user. The chapter "Adjustment/assembly instructions" provides an overview of the ways the wheelchair can be customised to your requirements. Technical changes to the design described in these instructions for use are reserved.

1.2 Indications for use

The Avantgarde wheelchair for active use is designed solely for individual use by persons who are unable to walk or who have a walking impediment, and can be operated either by the patient or by another person.

The Avantgarde wheelchair for active use may only be combined with the options mentioned in these instructions for use and vice versa. Ottobock assumes no liability for combinations with medical products and/or accessories from other manufacturers outside of the modular system.

1.3 Field of application

The versatility of the Avantgarde wheelchair and its modular design make it suitable for patients who have walking impediments or walking disabilities, for example due to:

- Paralysis (paraplegia / tetraplegia)
- Loss of limbs (lower limb amputation)
- Limb defect / deformation
- Joint contractures / defects
- Muscle and nerve disorders

- Diseases such as cardiac or circulatory insufficiency, balance disorders or cachexia as well as for geriatric patients who still have usable residual strength in the upper limbs.

The Avantgarde was specially designed for individuals who are normally able to move in a wheelchair independently.

The following points should also be considered for an individual fitting:

- Body height and weight (maximum load is 125 kg/ 275 lbs; model XXL2: 180 kg/397 lbs)
- Physical and mental condition
- Age of the user
- Home and living conditions
- Environment

1.4 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

1.4.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

1.4.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

1.4.3 Service life

Expected service life: **4 years**.

The design, manufacturing and requirements for the intended use of the product are based on the expected service life. These also include the requirements for maintenance, ensuring effectiveness and the safety of the product.

Using the product beyond the specified expected service life leads to increased residual risk and should only take place subject to the due diligence and deliberations of qualified personnel. If the service life is reached, the user or a responsible attendant should contact the qualified personnel who fitted the product or the manufacturer's servicing department (see inside rear cover or back page for address). Here the user can obtain information about known risks and the current options for refurbishing the product.

1.4.4 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

2 Safety instructions

2.1 Explanation of symbols

 WARNING	Warnings regarding possible risks of severe accident or injury.
 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.
 INFORMATION	Information regarding operation. Information for service personnel.

2.2 General safety instructions

INFORMATION
<p>Please read the instructions for use first! Before using the product, you should become familiar with the handling, function and use of the product. You may compromise your safety if you do not observe these instructions.</p> <p>However, it is impossible to protect against and avoid all unforeseeable circumstances and situations.</p>

 WARNING
<p>Risk of injury as a result of unauthorised use in vehicles for transporting persons with reduced mobility. We recommend that, wherever and whenever possible, users of a wheelchair, mobility base with seating shell or buggy transfer to the seats installed in vehicles for transporting persons with reduced mobility and use the corresponding vehicle restraint systems, as this is the only way to ensure optimal protection of the passengers in case of an accident. The Avantgarde wheelchair for active use is permissible for transporting passengers in wheelchair accessible vehicles when using Ottobock safety components as well as appropriate restraint systems.</p> <p>For more information, please refer to our instructions for use manual "Using your product for transportation in wheelchair accessible vehicles", order no. 646D158.</p>

 WARNING
<p>Serious injuries due to failure to observe the manufacturer's requirements.</p> <ul style="list-style-type: none"> • Using the product beyond the specified expected service life (see chapter "Service life") leads to increased residual risk and should only take place subject to the due diligence and deliberations of qualified personnel. • If the service life is reached, the user or a responsible attendant should contact the qualified personnel who fitted the product or the manufacturer's servicing department (see inside rear cover or back page for address). Here the user can obtain information about known risks and the current options for refurbishing the product.

Danger when overcoming stairs/obstacles

WARNING

Risk of falling in case of lacking assistance. Do not ascend or descend stairs without the assistance of attendants. If devices such as ramps or elevators are available, they should be used. If wheelchair-friendly access is missing, two attendants must carry the wheelchair over the obstacle.

WARNING

Risk of falling as a result of incorrectly set anti-tipper. If only one attendant is available when ascending or descending stairs, an incorrectly set anti-tipper (if mounted at all) can lead to severe falls. Adjust the anti-tipper so that it does not come in the way of the steps during transport. Afterwards, swing the anti-tipper back to its operational position.

WARNING

Risk of falling as a result of lifting the wheelchair incorrectly. Attendants must lift the wheelchair only by parts that are firmly attached (e.g. not at the footrest or rear wheels). If your wheelchair is equipped with height-adjustable push handles, ensure that the clamping levers are always firmly tightened.

WARNING

Risk of tipping when using lifting platforms. Ensure that the anti-tipper (if mounted) is outside the danger-area.

Danger when getting into the wheelchair

CAUTION

Risk of injury as a result of using the footplate for getting into the wheelchair. If possible, the footplates should be flipped up before getting into or out of the wheelchair.

Danger while driving

CAUTION

Risk of accident when driving without experience. This can lead to falls and other dangerous situations. Therefore, familiarize yourself with the new wheelchair by practising on even, straightforward terrain first. Learn how the wheelchair reacts when the centre of gravity shifts; for example on slopes or inclines or when clearing obstacles such as steps and curbs. This should be done only with assistance from another person. Using an anti-tipper is strongly recommended for inexperienced wheelchair users.

CAUTION

Risk of tipping when ascending too steep slopes. When ascending slopes or ramps and when crossing obstacles on upward slopes, always lean your upper body far forward.

CAUTION

Risk of tilting or tipping over as a result of centre of gravity shifts.

- When reaching for objects in front, to the side or behind the wheelchair, do not lean out of the wheelchair too far.
- Prior to activities that require you to bend forward in the wheelchair (e.g. tying your shoes), maximise the stability of the wheelchair. In order to do so, push the wheelchair backwards until the caster wheels turn forward.

⚠ CAUTION

Risk of tipping when driving against obstacles (steps, curbs) without braking. Use your wheelchair properly. Do not "jump" the chair down from higher surfaces.

⚠ CAUTION

Risk of accident as a result of incorrect use of the wheel lock. All brakes acting on the tyres do not serve as service brake but are only designed as parking brake (wheel lock). The wheel locks must not be used as driving brakes for slowing down the wheelchair, as in extreme cases, the abrupt stopping of the wheelchair can lead to falls.

⚠ CAUTION

Risk of accident if the wheelchair starts rolling. Engage the wheel lock to prevent your wheelchair from moving on uneven ground or during transfers (e.g. into a car).

⚠ CAUTION

Risk of accident as a result of driving in the dark. In the dark, wear light clothing or clothing with reflectors in order to improve your visibility. Ensure that the reflectors installed on the sides and rear of the wheelchair are easily visible. We also recommend installing lights.

⚠ CAUTION

Risk of tipping when descending slopes. When descending slopes, do not drive without braking and reduce your speed.

⚠ CAUTION

Reaching into exposed drive parts.

Crushing, pinching due to incorrect handling.

- Do not reach between the rear wheel and wheel lock or rear wheel and side panel when driving the wheelchair.
- Do not reach into the spokes of the rotating rear wheel.

⚠ CAUTION

Risk of falling due to lack of tipping resistance on public transport

- Always observe currently applicable legal requirements when using public transport.
- Always ensure that you are held in place securely when travelling on public transport. To do so, use the wheelchair areas, wheelchair bays and restraint systems provided. Firmly tighten the wheel locks.
- Please note that the anti-tipper, if activated, may be exposed to heavy loads when the public transport vehicle starts to move. In order to avoid damage, the manufacturer recommends that the wheelchair is positioned at a right angle to the direction of travel if no restraint system for passengers with reduced mobility is available.

Danger as a result of incorrect assembly or settings

⚠ CAUTION

Risk of tipping as a result of incorrectly set wheels. Extreme settings (e.g. rear wheels mounted in the foremost position) combined with an unfavourable body posture can cause the wheelchair to tip even on level ground. Avoid extreme settings. For transfemoral amputees it is essential to shift the rear wheels backwards. Use of an anti-tipper is required.

⚠ CAUTION

Risk of tipping as a result of missing or incorrectly mounted anti-tipper. Under no conditions should the anti-tipper assume the function of transport wheels, for example to transport a person in the wheelchair with the rear wheels removed. The anti-tipper must audibly lock in place, before it is able to bear loads. Firm seating must be verified by the user or by an attendant.

⚠ CAUTION

Risk of tipping as a result of changed diameters/mounting positions of the wheels. If the size and position of the caster wheels or the rear wheel size are modified, the caster wheels may flutter when driving at higher speeds. This can lead to blocking of the caster wheels and cause the wheelchair to tip over. For this reason, you should always maintain the factory settings or, if modifications are required, ensure that the wheelchair frame is aligned horizontally (see "Adjustment/assembly instructions").

⚠ CAUTION

Risk of falling as a result of unintentional loosening of the rear wheels. Always make sure that the quick-release axles are correctly set on the rear wheel. The rear wheel must not be removable without the button of the quick-release axle being depressed.

⚠ CAUTION

Risk of pinching, crushing, and falling due to unsecured screw connections.

- After any adjustment/ readjustment work, retighten the mounting screws/nuts firmly. Observe the torque specifications.
- Any time you loosen a screw connection with thread lock, replace it with a new screw connection or secure the old screw connection with medium strength thread locker (e.g. Loctite® 241).

⚠ CAUTION

Risk of falling due to missing self-locking nuts.

Always replace self-locking nuts with new self-locking nuts after disassembly.

Danger as a result of improper use of the tyres**⚠ CAUTION**

Risk of accident as a result of insufficient or excessive tyre pressure. The effectiveness of the wheel lock and the overall driving quality are dependent on adequate air pressure. The wheel locks are only effective with sufficient air pressure and proper adjustment (approx. **5 mm** distance, subject to technical changes).

Before starting to use your wheelchair, check that the tyres are inflated correctly. The required air pressure is printed on the side of the tyre. For rear wheels, it should be **at least 3.5 bar**. With properly inflated rear wheels and even tyre pressure on both sides, your wheelchair is much easier to operate and manoeuvre.

⚠ CAUTION

Risk of accident as a result of bad tyres. Insufficient tread depth of the tyres reduces adhesion. Please note that when you operate your wheelchair in public road traffic, you are subject to road traffic regulations.

Risk of skin damage

CAUTION

Be careful in case of unhealthy skin. Using the wheelchair with unhealthy skin can cause medically indicated complications such as reddening of the skin and pressure sores. Over several hours of use, the skin on the buttocks, back, and rear thighs is subjected to considerable strain. Therefore, always check these areas carefully before using the wheelchair.

No liability will be assumed for health injuries caused by use of this wheelchair with damaged skin.

Danger caused by fire/heat and cold

CAUTION

Risk of burns when near to fire. The seat and back upholstery might catch fire. Keep the wheelchair away from all sources of ignition, especially lit cigarettes.

CAUTION

Be careful in case of extreme temperatures. The wheelchair can heat up significantly when in the sun or in the sauna. In extreme cold, there is a risk of hypothermia.

Risk of hand injuries

CAUTION

Risk of hand injuries. To avoid hand injuries, do not grasp between the rear wheel and the wheel lock when driving your wheelchair.

CAUTION

Risk of pinching between wheelchair components. The wheel lock lever and the side panel or frame part can represent a risk of crushing with various adjustment options.

CAUTION

Risk of friction burns. Slowing down from high speeds or when descending longer slopes tends to heat up the hands and fingers, especially if using aluminium push rings. When using your wheelchair outdoors, wear leather gloves. Gloves provide you with a better grip and protect your fingers from dirt and hot metal.

Warnings regarding wheelchair damage

NOTICE

Damage caused by overload. The maximum load capacity for the Avantgarde wheelchair for active use is **125 kg / 275 lbs** and for the XXL 2 model **180 kg / 397 lbs**.

NOTICE

Damage caused by objects lying upon the device. If devices are transported in the folded state, the weight of objects lying upon the device can cause deformations and result in problems when unfolding the device. Never put heavy objects on a folded rehabilitation device.

NOTICE

Damage caused by wear and tear. If your seat and back upholstery become worn out or damaged, replace it immediately.

NOTICE


Damage caused by incorrect packaging. Only use original packaging for shipment of the product.

2.3 Further notes

INFORMATION

Even in the event of compliance with all applicable guidelines and standards, it is possible that alarm systems (e.g., in department stores) may respond to your product. Should this happen, remove your product from the area where the alarm was triggered.

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3 Delivery/preparing the wheelchair for use

⚠ CAUTION

Risk of pinching. Only grip the components shown when folding or collapsing the wheelchair.

Your wheelchair is usually delivered fully assembled, folded, and with the rear wheels removed. A few simple steps are sufficient to prepare it for use:

1. When folding or unfolding the wheelchair, hold it only at the indicated positions (otherwise there is risk of injury).
2. Loosen the safety strap (Fig. 4). While standing next to the wheelchair, tip it slightly toward you and press on the edge of the seat upholstery closest to you. The wheelchair will unfold up to its entire seat width (Fig. 1).
3. Flip the footplate down (Fig. 2).

INFORMATION When using a single-panel footrest, be sure the support of the footplate is locked into the receiver.

⚠ CAUTION Risk of tipping. The rear wheels have been pre-positioned according to your indications. Before using the wheelchair, this pre-adjusted position must be checked by the wheelchair user for function and stability against tipping with the assistance of skilled staff.

4 Transportation

Flip up the footplate.

Grasp the middle of the seat upholstery in the front and the back and pull up. The wheelchair will fold (Fig. 3).

Use the safety strap to hold the folded wheelchair in position when it is being lifted (Fig. 4).

INFORMATION During reassembly, ensure that the footrest locking device has engaged.

Depending on the wheelchair model, further options are available for making the wheelchair more manageable for car transport etc. With the model »Avantgarde T« with single-panel footrest you can swing away the left part of the footrest with footplate to the outside to reduce the folding with (Fig. 5). The Avantgarde VR and XXL 2 models allow you to swing away and detach the footrests, thus reducing the overall length of the wheelchair.

Transportation is greatly facilitated by the removable quick-release axles on the rear wheels (included in the standard version) (Fig. 6).

Reach into the spokes close to the hub of the rear wheel and press the button on the quick-release axle with your thumb. You can then easily pull off the rear wheel (Fig. 6). Follow the corresponding procedure when remounting.

CAUTION Risk of accident. During assembly, confirm that all bolts engage in the corresponding bore holes and that the quick-release axles of the rear wheels and caster wheels securely lock into their axle housings. The rear wheels and caster wheels must not be removable without the button being depressed! Check the function of the wheel locks.

4.1 Transfer

Transfer into the wheelchair is described below. The Avantgarde is designed for persons who normally move around in the wheelchair independently. Therefore, the independent transfer is described. Should you require an attendant, the steps described here are carried out correspondingly but with support.

Manoeuvre the wheelchairs so that they are positioned at an angle of approx. 45° degrees to each other (Fig. 7). Ensure that the wheel locks do not get caught in each other but remain freely movable. First activate the wheel locks of the wheelchair that is currently in use (Fig. 8).

Then activate the wheel locks of the wheelchair you want to transfer into (Fig. 9). Flip up the footplate (Fig. 10). Move your feet from the footplate to the floor and move ahead in the wheelchair if required (Fig. 11).

Transfer into the second wheelchair. The method illustrated is one possible way to support yourself (Fig. 12). Find the most suitable path for you by practising with an attendant.

Flip down the footplate of the Avantgarde. Position your feet on the footrest of the wheelchair (Fig. 13). Release the wheel locks (Fig. 14). Now you can use the wheelchair.

4.2 Use in vehicles for transporting persons with reduced mobility

4.2.1 Avantgarde T, Ti8.9, TiUltra, VR

⚠ WARNING

Use as a seat for transportation in vehicles for transporting persons with reduced mobility. Risk of injury as a result of unauthorised use in vehicles for transporting persons with reduced mobility. Whenever possible while travelling in vehicles for transporting persons with reduced mobility, passengers should use the seats and vehicle restraint systems already installed in the vehicle. This is the only way to provide optimal protection for the passengers in case of accident. You can use the product as a seat for transporting passengers in wheelchair accessible vehicles when using the relevant Ottobock safety components as well as the appropriate restraint systems.

For further information please refer to our instructions for use „Using your product for transportation in a wheelchair accessible vehicle“, order number 646D158.

4.2.2 Avantgarde XXL2

⚠ WARNING

Unallowable use in vehicles for transporting persons with reduced mobility. Risk of serious injury when using the product as a seat.

- The product has **not** been approved by the manufacturer for use as a seat in vehicles for transporting persons with reduced mobility.
- Always use the seats installed in the wheelchair accessible vehicle with the corresponding vehicle restraint systems.

5 Options

The wheelchair is designed as a modular component system. This means that certain accessories can be adapted to your wheelchair. In the following, we will present you a selection of our variants and accessories which improve the use and comfort of your wheelchair.


We recommend using a seat cushion at all times. The seat upholstery of the Avantgarde wheelchair has a loop texture so that seat cushions can be attached with hook and loop connection. All seat cushions listed on the order form have the required hook texture.

5.1 Anti-tipper (Fig. 15)

Prevents the wheelchair from tipping too far backward. The anti-tipper is strongly recommended for untrained wheelchair users and children, as well as when extreme adjustments have been made to the rear wheel. With slight downward pressure on the anti-tipper, it can be disengaged and swung to the inside (refer to section "Adjustments/assembly instructions").

5.2 Tip-assist (Fig. 16)

Makes it easier for an attendant to tip up the wheelchair, e.g. to cross a step. The tip-assist must be mounted so that the distance to the ground is at least 4 cm.

5.3 Crutch holder with hook and loop strap (Fig. 17)

Allows crutches to be attached to the wheelchair.

5.4 Transport wheels (Fig. 18)

With the rear wheels removed, the wheelchair can still be used as a rolling base; this makes it possible to manoeuvre through narrow passages (e.g. narrow bathroom doorways or airplane aisles).

CAUTION Risk of accident. Note that once the rear wheels are removed, the wheel locks are no longer functional!

5.5 Footrests/footplates + accessories (Fig. 19–29)

5.5.1 Single-panel footrest, narrow (Fig. 19)

Makes the wheelchair more compact (Figure incl. frame pad).

5.5.2 Individual footrests, angle-adjustable (Fig. 20)

Depth = 16 cm.

5.5.3 Footrest for short lower leg lengths (Fig. 21)

For 16 to 31 cm lower leg length.

5.5.4 Footrest bar for Avantgarde Ti 8.9 (Fig. 22)

The footrest is angle-adjustable.

5.5.5 Plastic plate for the footrest of the Avantgarde Ti 8.9 (Fig. 23)

Non-skid surface, colour: black.

5.5.6 Locking device for single-panel footrest (Fig. 24)

Prevents inadvertent loosening of the footplate, especially with spastic users.

5.5.7 Lateral heel blocks for the footplate (Fig. 25)

For 16 cm deep footplates.

5.5.8 Adjustable lateral heel blocks for the footplate (Fig. 26)

Prevent the feet from slipping off the footplate and can be adjusted in their position.

5.5.9 Elevating footrest (Fig. 27)

Allows positioning of the leg at different angles. The adjustment range of the legrest is 15° to 90°.

5.5.10 Elevating footrest with automatic length adjustment (Fig. 28)

A spring mechanism provides for the corresponding length of the swivel arm.

5.5.11 Amputation legrest (Fig. 29)

An alternative to be mounted to an elevating footrest.

5.6 Seat (Fig. 30/31)

5.6.1 Seat upholstery with pocket for catheter (Fig. 30)

Accommodates disposable catheter and offers easy access below the seat.

5.6.2 "The Seat" (Fig. 31)

Variably adjustable seat base to improve positioning.

5.7 Back (Fig. 32–35)

5.7.1 Standard back, height-adjustable/adaptable (Fig. 32)

Height adjustment is possible in increments of 2.5 cm and with setting ranges of 25 to 35 cm, 30 to 40 cm, and 37.5 to 47.5 cm.

5.7.2 Ergo back (Fig. 33)

Continuously height-adjustable option with setting ranges of 25 to 35 cm, 35 to 45 cm, and 40 to 50 cm. Angle-adjustable from -9° to $+15^{\circ}$ in 6° increments.

5.7.3 "The Back" (Fig. 34)

Firm back to support positioning (see 647H491 instructions for use).

5.7.4 Back stabilizer (Fig. 35)

Increases the wheelchair's stability, especially in case of high loads, while the wheelchair remains foldable.

5.8 Side panels (Fig. 36–39)

5.8.1 Plug-on clothing protector side panel (Fig. 36)

The side panel can be adapted to the rear wheel position. In addition, the clothing protector side panel can be removed for transfer, for example.

5.8.2 Side panels, adjustable, with protection against cold (Fig. 37)

5.8.3 Aluminium clothing protector side panels with neoprene insert (Fig. 38)

Many wheelchair users know the problem of clothes coming into contact with the rear wheels and getting dirty. The clothing protector prevents this. It is horizontally and vertically adjustable and can be adapted to the position of the rear wheel.

5.8.4 Carbon clothing protector side panels (Fig. 39)

The side panel has been optimised in terms of weight and has a cold insulation effect.

5.8.5 Desk side panels (Fig. 40)

This side panel can be flipped up and removed, and features a height-adjustable armrest. For getting into and out of the wheelchair, the side panel can be swung to the rear and removed.

5.9 Push handles (Fig. 41/42)

5.9.1 Push handles, height-adjustable (Fig. 41)

Can be adjusted vertically to a more comfortable height for the attendant.

5.9.2 Telescoping push handles (Fig. 42)

Continuously adjustable up to max. 10 cm. These push handles can be adjusted vertically to a more comfortable height for the attendant.

5.10 Caster wheel adapter/caster wheel accessories (Fig. 43–46)

5.10.1 Caster wheel adapter for large anterior seat heights (Fig. 43)

Increases the anterior seat height by 3.5 cm. Also increases the adjustment range of the lower leg lengths. Offers more freedom of movement for the feet due to caster wheels being positioned more exteriorly.

5.10.2 Caster wheel adapter for low anterior seat heights (Fig. 44)

Reduces the anterior seat height by 3.5 cm. Offers more freedom of movement for the feet due to caster wheels being positioned more exteriorly.

5.10.3 Caster wheel with quick-release axle (Fig. 45)

The caster wheel with quick-release axle serves the quick exchange of the caster wheels for different kinds of application (indoors/outdoors).

5.10.4 Caster wheels (Fig. 46)

Different caster wheel models are available.

(45)

(46)

5.11 Rear wheels (Fig. 47/48)

Spox wheels (Fig. 47) and Infinity Ultralight wheels (Fig. 48) can be selected among others.

(47)

(48)

5.12 Brake/wheel lock systems/wheel lock accessories (Fig. 49/50)

5.12.1 Rear wheel with drum brake (not illustrated)

Allows the attendant to comfortably and safely brake the wheelchair. The rear wheels remain removable via the quick-release axle system.

5.12.2 Wheel lock, swing-away (Fig. 49)

Can be swung to the inside of the frame. Makes transfer from the wheelchair easier and, in combination with transport wheels, reduces the overall width of the wheelchair.

5.12.3 Wheel lock for one-handed operation (Fig. 50)

Recommended especially for hemiplegics. It can be activated either on the right or left side and ensures the safe locking of both rear wheels with a pull-wire system.

INFORMATION Please note that the folding wheel lock lever extension must always be attached when releasing or engaging the wheel lock.

5.12.4 Wheel lock for user and attendant (not illustrated)

The wheelchair can also be braked by the attendant. The rear wheels remain removable via the quick-release axle system.

5.12.5 Plug-on wheel lock lever extension (not illustrated)

Facilitates use of the wheel lock.

5.13 One-arm drive (Fig. 51)

Especially for hemiplegics who want to move either by »toddling« or by using the push ring on one side.

5.14 Wheelbase extension (Fig. 52)

This option makes the wheelchair particularly stable against tipping and at the same time increases the turning radius.

INFORMATION For transfemoral amputees the long wheelbase setting is indispensable.


5.15 Spoke protector (Fig. 53)

Protects the fingers from getting caught in the wheel spokes. Also improves the appearance of the wheelchair.

5.16 Shock absorber system (Fig. 54)

The shock absorber system absorbs loading peaks, for example when driving over steps or curbs or on uneven ground.

5.17 Safety belt (Fig. 55)

Prevents the wheelchair user from falling out of the wheelchair.

5.18 Tray (Fig. 56)

Serves as a supporting surface during meals, when working or when playing. The clear material allows visibility of the legs and correction of the sitting posture.

5.19 Further instructions for use

- Attaching loads (e.g. backpacks) can adversely affect stability. Suspending additional loads on the wheelchair is not permitted.
- The recommended overall width for manual wheelchairs in operational condition is 700 mm. This requirement is intended to ensure unimpeded use of, for example, escape routes. Please note that the wheelchair dimensions can exceed the recommended value in versions with very large seat widths (for more information see chapter 9 > “Technical Data”).
- The wheelchairs in this series fully satisfy the minimum technical requirements for wheelchairs transportable by train. However, please note that it is possible that not every individual wheelchair will fulfil the minimum requirements due to the variation in designs (see chapter 9.2 > “Threshold values for wheelchairs transportable by train”).

6 Adjustment/assembly instructions

INFORMATION

The following tools are required for repair and maintenance:

- Allen wrenches, sizes 4 mm, 5 mm, and 6 mm
- Open-end wrenches, sizes 8 mm, 10 mm, 11 mm, 12 mm, 13 mm, 19 mm, and 24 mm
- Set of screwdrivers
- Torque wrench
- Tyre levers
- Tyre pump

Each wheelchair user has their own ideas regarding the wheelchair type and accessories that are correct for him/her, as well as the position in the wheelchair and the handling characteristics. Therefore, the Avantgarde wheelchair for active use offers many options that permit the wheelchair to be adapted according to your requirements. We recommend consulting your dealer or therapist to determine the correct wheelchair settings for you.

Special features of the Avantgarde Ti 8.9/Ti Ultra

The ultralight Avantgarde Ti 8.9 and Ti Ultra wheelchair models for active use have consistently been fabricated with a light-weight construction. For this reason, any adjustment works must be done with utmost care. It is of particular importance that the torque values for the screw connections of high-strength aluminium alloys are observed (see table of torque values for aluminium screw connections).

If you do not have of the proper tools (e.g. torque wrench), please contact your medical supply dealer.

Torque values for aluminium screw connections

Screw diameter in mm	Max. torque in Nm
6	5.5
8	13

Before re-assembly and re-use of the aluminium screw connections, it must be ensured that the threads are not damaged. It is recommended to use new aluminium screws for each (re-) assembly. Please use a liquid thread lock substance for securing the aluminium screws (e.g. Loctite® 241).

6.1 Adjusting the wheel lock (Fig. 57)

To change the position of the rear wheels, first loosen the wheel lock attachment screws and push the wheel lock forward. Then reposition the rear wheels. The wheel lock should then be remounted such that the distance between the tyres and brake block is **max. 5 mm** when the wheel lock is not activated (Fig. 57; subject to technical changes).

⚠ CAUTION The wheel locks are only effective with sufficient air pressure and proper adjustment. The proper air pressure is indicated on the side wall of the tyre, but should at least be **3.5 bar** (350 kPa). To ensure sufficient braking effect, use only original rear wheels with an approved radial excursion of maximally **± 1 mm**.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments! The torque for the attachment screws of the wheel lock is **10 Nm**.

6.2 Changing the rear wheel adapter position in the frame (Fig. 58)

The further back the rear wheel is attached, the longer the wheelbase and the greater the wheelchair stability (refer to section 5.14 "Wheelbase extension").

Attaching the rear wheels in a more forward position (only for advanced wheelchair users!) reduces the load on the caster wheels, which makes the wheelchair more manoeuvrable. The caster wheels can then also be lifted more easily to clear obstacles such as curbs or steps.

To adjust the horizontal position of the rear wheel, first remove the countersunk head screws (Fig. 58, item A). Choose the desired position (Fig. 58, item B) and re-assemble the screws.

Ensure that the left and right rear wheel adapters have the same horizontal position in the frame.


CAUTION Risk of tipping. Please note: If the rear wheels are in a more forward position and the user's body is not appropriately positioned, the user may tip backward – even on level ground!

CAUTION Risk of tipping. For untrained wheelchair users and when there are extreme adjustments to the rear wheel, an anti-tipper is strongly recommended!

CAUTION Risk of tipping. For transfemoral amputees it is essential to shift the rear wheels backwards!

INFORMATION Changing the rear wheel position can also change the angle of the caster axle in relation to the ground. Ensure that this angle is always approx. **90°** wide. The wheel lock must be re-adjusted as well.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The torque for the screws (Fig. 58, item A) is **10 Nm**.

6.3 Changing the rear wheel position in the rear wheel adapter (Fig. 59–61)

Each change of position of the quick-release axle housing in the rear wheel adapter has an influence on the posterior seat height and the centre of gravity! The relation of anterior and posterior seat height represents the individual seat inclination. The higher you attach the rear wheel in the rear wheel adapter, the more the seat pitches down. This has two effects:

1. The wheelchair has a greater tendency to tilt backwards.
2. The user sits deeper in the wheelchair and thus has more stability.

6.3.1 Changing the position in the rear wheel adapter Avantgarde T/Ti 8.9/VR (Fig. 59)


Loosen the nut (Fig. 59, item A) of the quick-release axle housing (Fig. 59, item B) on the inside of the frame, until the tapered washers can be spread apart slightly. Now you can move the quick-release axle housing with nuts and tapered washers into the desired position.

Let the lugs of the tapered washers engage in the new desired orientation and firmly re-tighten the nut (Fig. 59, item A). Make sure that the position of the tapered camber washers is not changed and that both the left and right rear wheel attachment devices are set in the same vertical position.

The distance between the rear wheel and side panel can now be continuously adjusted by adjusting the rear wheel fitting (Fig. 59, item C).



INFORMATION When the rear wheel position is changed, the angle of the caster axle in relation to the ground also changes. Ensure that this angle is always approx. 90° wide. The wheel lock must be re-adjusted as well.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The torque for the hexagon nuts (Fig. 59, item A/C) is **50 Nm**.

6.3.2 Changing the position in the rear wheel adapter Avantgarde XXL 2 (Fig. 60/61)

Loosen and remove the nut (Fig. 60, item A) of the rear wheel fitting (Fig. 60, item B) on the inside of the frame. The possible positions result from defined bore holes in the rear wheel adapter (Fig. 61, item C).


6.4 Setting the wheel camber (Fig. 62/63)

(does not apply to XXL 2 and Ti Ultra, since here only 0° is possible)

The Avantgarde modular system provides a rear wheel adapter which can be used to set the rear wheels at different camber angles (0°, 2°, and 4°). The greater the negative camber angle of the rear wheels, the more manoeuvrable the wheelchair becomes (Fig. 63).

The rear wheel adapter itself is set to a wheel camber of 2° due to its form. Deviations from this setting (0° or 4°) are achieved by rotating the special 2° tapered washers.

6.4.1 Setting the rear wheel camber to 0° (Fig. 62 A)

Position the tapered washer on the outside of the rear wheel adapter with the thicker side pointing upwards into the desired engaging position (you will find -2° engraved on the outside).

Place the second tapered washer (with the thinner side pointing upwards) on the inside of the rear wheel adapter and into the corresponding engaging position (you will find +2° engraved on the inside). The outsides of the two camber washers are now in parallel positions (Fig. 62 A).

Insert the rear wheel fitting and slightly tighten the hexagon nuts (Fig. 59, items A/C) once you have applied the retaining washer from the inside. Check the space between the rear wheel and the side panel. You can adjust the space by rotating the fitting.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The torque for the hexagon nuts (Fig. 59, item A/C) is **50 Nm**.

6.4.2 Setting the rear wheel camber to 2° (Fig. 52 B)

Position the tapered washer on the outside of the rear wheel adapter with the thicker side pointing upwards into the desired engaging position (you will find -2° engraved on the outside).

Place the second tapered washer (with the thinner side pointing upwards) from the outside against the first tapered washer (you will find +2° engraved on the outside). The outsides of the two camber washers are now in parallel positions (Fig. 62 B).

Insert the rear wheel fitting and slightly tighten the hexagon nuts (Fig. 59, items A/C) once you have applied the retaining washer from the inside. Check the space between the rear wheel and the side panel. You can adjust the space by rotating the fitting.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The torque for the hexagon nuts (Fig. 59, item A/C) is **50 Nm**.

6.4.3 Setting the rear wheel camber to 4° (Fig. 62 C)

Position the tapered washer on the outside of the rear wheel adapter with the thinner side pointing upwards into the desired engaging position (you will find +2° engraved on the outside).

Place the second tapered washer (with the thicker side pointing upwards) on the inside of the rear wheel adapter and into the corresponding engaging position (you will find -2° engraved on the inside). The outsides of the two camber washers are now in parallel positions (Fig. 62 C).

Insert the rear wheel fitting and slightly tighten the hexagon nuts (Fig. 59, items A/C) once you have applied the retaining washer from the inside. Check the space between the rear wheel and the side panel. You can adjust the space by rotating the fitting.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The torque for the hexagon nuts (Fig. 59, item A/C) is **50 Nm**.

INFORMATION Please be aware that with a larger negative camber the overall width of the wheelchair increases and the rear wheels are subjected to greater stresses!



6.5 Adjusting the caster wheel journal angle (Fig. 64/65)

After the rear wheels have been mounted in the appropriate position, the caster wheel journal angle must be readjusted. The caster wheel adapter should be as horizontal as possible in relation to the ground to ensure optimal driving behaviour of the wheelchair. The caster wheel adapter allows for continuous angle adjustment.

First, remove the protective caps from the caster wheel adapters (Fig. 64, item A). Next, loosen the two cap screws (Fig. 65, item B) on the inside and then the cap screw on the outside (Fig. 64, item C). Place the level from the tool pouch on the caster wheel adapter (Fig. 65).

Using a large screwdriver, screw the eccentric disk so that the caster wheel adapter moves into horizontal position. When in this position, first tighten the cap screw (Fig. 64, item D) followed by the cap screws (Fig. 64/65, item B).

INFORMATION Make sure to **very firmly** retighten all screws and nuts after making adjustments. The torque of the cap screws (Fig. 64, item B) for the Avantgarde T, VR, and XXL2 is **23 Nm**; for the Avantgarde Ti 8.9 and Ti Ultra with aluminium screws, the torque is **13 Nm**. Tighten the Allen head screw on the eccentric to **8 Nm** (Fig. 65, item C).

6.6 Adjusting the footplate angle (Fig. 66/67)

An angle adjustable footplate is included as standard equipment.

6.6.1 Adjusting the footplate angle Avantgarde T/VR/XXL 2 (Fig. 66)

Loosen the Allen screw at the rear suspension of the footplate with several turns of the wrench (Fig. 66). Adjust the footplate to the desired position and firmly re-tighten the screw.

6.6.2 Adjusting the footplate angle Avantgarde Ti 8.9/Ti Ultra (Fig. 67)

To adjust, slide a 5 mm Allen wrench through the side hole in the footrest bar until the wrench locks into the Allen screw in the footplate. Loosen the screw by several turns, then sharply tap the end of the Allen wrench, which will disengage the clamp connection.

Adjust the footplate to the desired position and firmly re-tighten the screw.

6.7 Adjusting the footrest (Fig. 68-69)

6.7.1 Adjusting the footrest height Avantgarde T/Ti 8.9/Ti Ultra (Fig. 68)

By loosening the screw at the footrest bar (Fig. 68), the footplate can be continuously adapted to your lower leg length and to the thickness of your seat cushion.

INFORMATION Loosen the screw – do not completely remove the screw. There are two threaded holes in the footrest bar. Depending on the desired adjustment, however, it may be necessary to completely remove the screw and use the second hole.


INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The torque for the cap screws is **10 Nm**.

6.7.2 Adjusting the footrest height Avantgarde VR and XXL 2 (Fig. 69)

By loosening the screw (Fig. 69, item A), the footplate can be adapted to your lower leg length and to the thickness of your seat cushion.

Make sure when adjusting the footrest height, that the footrest bar (Fig. 69, item A) is slid into the swivel segment (Fig. 69, item C) at least **40 mm**.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The torque for the screw (Fig. 69, item A) is **8 Nm**.

6.7.3 Adjusting the height of the footrest for short lower leg lengths (Fig. 70)

Loosen the screws (Fig. 70, item A). Move the footplate to the desired position.



INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments. The torque for all screws is **9 Nm**.

6.7.4 Adjusting the height of the elevating footrest with automatic length adjustment (Fig. 71)

First completely remove the screw (Fig. 71, item A) in the footrest bar.

To change the lower leg length, you can now choose one of the upper or lower threaded holes in the footrest bar and screw together the two components. Depending on the adjustment, it may be necessary to remove one or several of the spacer sleeves (Fig. 71, item B).

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments.

6.7.5 Adjusting the height of the elevating footrest Avantgarde VR/XXL 2 (Fig. 72)

By loosening the screw at the footrest bar (Fig. 72), the footplate can be adapted to your lower leg length and to the thickness of your seat cushion.

INFORMATION Loosen the screw – do not completely remove the screw. There are five threaded holes in the footrest bar. Depending on the desired adjustment, however, it may be necessary to completely remove the screw and use another hole.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments.


6.7.6 Adjusting the angle of the elevating footrest Avantgarde VR/XXL 2 (Fig. 73)

The legrest can be adjusted within the range from 15° to 90°.

Pull the footrest upward until it has reached the desired position. The legrest automatically locks in place.

To reposition the legrest, relieve the load on the footrest and press the release lever (Fig. 73, item A) downward to lower the legrest. Bring the legrest to the desired position.

INFORMATION Make sure to always secure the legrest with your hand against uncontrolled falling downward when lowering it.

6.7.7 Swinging away the elevating Avantgarde FL footrest (Fig. 74)

Slightly lift the legrest and swing it to the inside or outside (Fig. 74).

6.8 Adjusting the brake force – drum brake (Fig. 75)


The brake force is adjusted at the adjustment screw (Fig. 75, item A) to achieve optimal braking effect. The brake force is increased by turning the screw counter-clockwise and decreased by turning the screw clockwise.

Loosen the counter nut (Fig. 75, item B) and turn the adjustment screw anti-clockwise until a frictional noise occurs at the wheel. Then tighten the adjustment screw until the frictional noise disappears. The wheel should run freely.

After finishing the adjustment, fasten the adjustment screw by tightening the counter nut (Fig. 75, item B).

INFORMATION Ensure that the brake force of both rear wheels is adjusted equally.

INFORMATION Make sure that the drum brake engages with sufficient force only when the manual brake lever is set to the second ratchet position.

6.9 Tension adjustable seat upholstery (Fig. 76)

If the seat fabric has stretched due to extended use, the upholstery can be re-adjusted to restore optimal pressure distribution of the thighs.

Loosen the hook and loop connection under the seat fabric. Re-tighten the upholstery and secure this position using the hook and loop fastener.

Unfold the wheelchair. The cross brace must be seated fully in the support guides.

6.10 Adapting the back upholstery (Fig. 77)

You can adjust the back upholstery in segments according to your needs.

To do this, remove the back upholstery padding and undo the hook and loop fasteners of the straps (Fig. 77).

Reposition the straps as desired, attach them again with the hook-and-loop fasteners and cover with the padding.

INFORMATION The contours of the adjustable back will result from the user's individual requirements. Please do not make the straps too taut to allow for some play. A well adjusted back provides lasting comfort for the wheelchair user and reduces the risk of secondary damages and pressure zones.

INFORMATION See to it that the user's pelvis is positioned as far posterior in the wheelchair as possible, i.e. between the back tubes. Especially in the lower section (pelvis) the padding must have sufficient sag.


6.11 Adjusting the back height (Fig. 78–79)

With the “Back, height-adjustable/adaptable” option, the height of the back tube can be adjusted. Adjustment is possible within the hole matrix with spaces of 25 mm over a range of 100 mm. The adjustment range is predetermined by the selected back size (175 – 275 mm, 250 – 350 mm, 300 – 400 mm).

1. Preparation: Take off the wheel and secure the wheelchair by jacking it up. Also remove other add-on components (e. g. anti-tipper or desk side panel mountings) that might be installed and connected to the back tubes.
2. Loosen and remove the upper and lower screw on the rear frame on each side (Fig. 78, item A, B).
3. Loosen and remove the screw on the side panel clamp on each side (Fig. 78, item C).
Slide the two back tubes to the desired height (Fig. 79). Adjust both back tubes to the same height.

Please note: Make adjustments only within the range provided for the back tubes (175 – 275 mm, 250 – 350 mm, 300 – 400 mm). Always use 2 screws on each side to fix the back tube.

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments.





6.12 Mounting the anti-tipper (Fig. 80–83)

1. Mount the receiver tube with the two screws (Fig. 80, item 1) on the inner side of the posterior frame part. Tighten the nuts (Fig. 80, item 2) with a torque of **9 Nm**.
2. Insert the expansion spring (Fig. 81, item 4) with the close eyelet into the receiver tube. When the receiver tube's bore hole and the eyelet hole are flush, insert the sleeve (Fig. 81, item 3). Secure the sleeve with the screw (Fig. 82, item 7). Tighten the nut (Fig. 82, item 8) with a torque of **7 Nm**.
3. Hook a cord (Fig. 82, item 10) into the lower eyelet of the expansion spring and push the anti-tipper tube (Fig. 81, item 6) over the cord and expansion spring into the receiver tube.
4. Pull at the cord (Fig. 82, item 10) until the expansion spring extends from the anti-tipper tube (Fig. 81, item 6).
5. Secure the expansion spring using a screwdriver (Fig. 83, item 11) and remove the cord. Put the spring holder (Fig. 83, item 12) onto the lower eyelet of the expansion spring. Remove the screwdriver. The eyelet of the expansion spring will then snap into place.

The anti-tipper can now be adjusted to the position of your rear wheel. After the screw connection has been loosened, the anti-tipper can be horizontally and vertically adjusted. The anti-tipper wheel must at least completely extend from the rear wheel to the back, and the distance between the wheel and ground should not exceed **5 cm** maximum (Fig. 15). Find the appropriate position with the assistance of a helper!

INFORMATION Make sure to firmly retighten all screws and nuts after making adjustments.

6.13 Attaching the safety belt (Fig. 84/85)

Attach the safety belt with the screws provided to the uppermost bore hole of the posterior frame part.

7 Maintenance, cleaning, and disinfection

7.1 Maintenance

Your wheelchair has been provided with the CE marking. The manufacturer herewith guarantees that this medical product as a whole conforms to the requirements of the European Directive for Medical Products 93/42/EEC.

The proper function of the wheelchair, especially of the wheel locks, should be checked before every use. Safety nuts should be used only once. If they have been loosened several times, these nuts must be replaced.

If you notice any defects, please contact your authorised specialist dealer immediately in order to have them corrected. This is particularly important if you notice that the driving characteristics change or the wheelchair becomes less stable. **We also recommend that you have the wheelchair maintained regularly at 12-month intervals by your authorised specialist dealer.** The functions described in the following list must be verified by the user at the specified intervals.

Check	Before each use	Monthly	Quarterly
Function test of the wheel locks/brakes	X		
Sagging of the seat or back upholstery		X	
Check stability of footplate		X	
Visual inspection of wear parts (e.g. tyres, bearings)		X	
Dirt on bearings		X	
Push rings for damages		X	
Air pressure (see indication on tyre)		X	

Check	Before each use	Monthly	Quarterly
Folding mechanism for wear and tear		X	
Spoke tension of the rear wheels			X
Screw connections			X

To ensure smooth operation at all times, the following precautions and cleanings should be carried out. Users may carry these out by themselves.

- Hair and dirt particles generally accumulate between the caster wheel and fork. This can restrict the caster wheels from rotating smoothly. Remove the caster wheels and thoroughly clean the forks and caster wheels using a mild household cleaner.
- The rear wheels feature a quick-release system. To keep this system operational, ensure that no dirt adheres to the quick-release axle or axle housing. The quick-release axle should also be lightly lubricated regularly with resin-free sewing machine oil.
- If your wheelchair gets wet, towel-dry it as soon as possible.
- Do not use the wheelchair in salt water and keep sand or other particles from damaging the wheel bearings.
- Screw connections should be checked frequently, in particular when beginning to use the wheelchair and after any adjustment. If a screw connection comes loose repeatedly, please consult your specialist dealer.

7.2 Cleaning and disinfection

Clean the padding and upholstery with warm water and a washing-up liquid. Remove spots with a sponge or a soft brush. Rinse with clear water and let the components dry.

INFORMATION Do not use any aggressive cleansing agents, solvents, or hard brushes etc.

INFORMATION No wet-cleaning. Components are not machine washable.

INFORMATION For disinfection, water-based agents (e.g. Sagrotan original concentrate) should be used. The manufacturer's instructions for use must be observed.

INFORMATION Before disinfecting, the pads and handles must be cleaned.

INFORMATION Plastic parts, frame parts, and the chassis and wheels can be cleaned with a mild detergent and a damp cloth. Dry thoroughly afterwards.

8 Changing tyres (Fig. 86–89)

Changing flat tyres requires only the necessary tools and minimal technical skills. Users may change tyres themselves if they wish. We recommend carrying an emergency repair kit and tyre pump at all times (if using pneumatic tyres).

Suitable tyre pumps are listed on the order form and are supplied with the product. Alternatively, users can keep a spray can of hardening foam (available at bicycle shops) to fill a tyre in case of punctures.

- In case of a flat tyre, carefully remove the tyre from the rim using the appropriate tools.
- Take care not to damage the rim or the inner tube.
- Repair the inner tube according to the directions in the repair kit, or, replace the damaged inner tube with a new one.
- Before re-mounting the tyre, inspect the inner rim surface and the inside of the tyre for any object that may have caused the puncture.
- Be sure that the rim bands are not damaged. They protect the tube from being damaged by the ends of the spokes.

Mounting

- Place the protective rim band over the valve. Then push the valve through the rim. Screw on the valve nut. The rim band can then be easily applied.
- Ensure that all spoke ends are covered.

Tyres

- Beginning behind the valve, place the lower side of the tyre over the edge of the wheel rim. Pump the tube with a small amount of air until it assumes its round shape, and then insert it into the tyre.
- Check to ensure the inner tube is properly positioned all around and is free of folds; if this is not the case, release some of the air. Now you can easily install the upper side of the tyre starting across from the valve, using both hands and working towards the valve.

Inflating

- Check around the tyre on both sides to ensure that the tube is not clamped between the edge of the tyre and the rim.
- Push the valve back slightly and then pull it out again to ensure that the tyre is properly seated in the valve area.
- Inflate the tyre until it can still be pressed easily with your thumb. If the control line on both sides of the tyre is equally distanced to the rim, the tyre sits centrally. If the tyre does not sit centrally, let out some air and realign the tyre. Then inflate the tyre to maximum operating pressure (see sidewall) and tighten the dust cap.

9 Appendices

9.1 Technical data

Avantgarde T, Ti 8.9, and Ti Ultra – Overall length (in cm) ¹⁾								
	Rear wheel size							
	22"		24"		MTB 24"		26"	
Seat depth	Rear axle position	Front axle position	Rear axle position	Front axle position	Rear axle position	Front axle position	Rear axle position	Front axle position
36 – 42	89	80.5	91.5	83	92.5	84	94	85.5
44 – 46	93	84.5	95.5	87	96.5	88	98	89.5

Avantgarde VR – Overall length (in cm) ¹⁾									
Seat depth	Lower leg length max./min.	22"		24"		MTB 24"		26"	
		Rear axle position	Front axle position	Rear axle position	Front axle position	Rear axle position	Front axle position	Rear axle position	Front axle position
36	55	89	80.5	91.5	83	92.5	84	94	85.5
	32	83	74.5	85.5	77	86.5	78	88	79.5
38	55	91	82.5	93.5	85	94.5	86	96	87.5
	32	85	76.5	87.5	79	88.5	80	90	81.5
40	55	93	84.5	95.5	87	96.5	88	98	89.5
	32	87	78.5	89.5	81	90.5	82	92	83.5
42	55	95	86.5	97.5	89	98.5	90	100	91.5
	32	89	80.5	91.5	83	92.5	84	94	85.5
44	55	97	88.5	99.5	91	100.5	92	102	93.5
	32	91	82.5	93.5	85	94.5	86	96	87.5
46	55	99	90.5	101.5	93	102.5	94	104	95.5
	32	93	84.5	95.5	87	96.5	88	98	89.5
48	55	101	92.5	103.5	95	104.5	96	106	97.5
	32	95	86.5	97.5	89	98.5	90	100	91.5
50	55	103	94.5	105.5	97	106.5	98	108	99.5
	32	97	88.5	99.5	91	100.5	92	102	93.5

Avantgarde T, Ti 8.9, Ti Ultra, and VR – Overall width (in cm) ¹⁾			
Seat width	Standard side panels and clothing protectors	Side panel with height adjustable armrest	Rear wheel with drum brake
32	48.5	51	52
34	50.5	53	54
36	52.5	55	56
38	54.5	57	58
40	56.5	59	60
42	58.5	61	62
44	60.5	63	64
46	62.5	65	66
48	64.5	67	68
50	66.5	69	70

1) Applicable to push ring attachment, narrow, and a 0° camber of the rear wheels (for push ring attachment, wide: +20 mm).

Avantgarde XXL 2 – Overall length (in cm)

Seat depth	Lower leg length max./ min.	24"		26"	
		Rear axle position	Front axle position	Rear axle position	Front axle position
40	55	95.5	85	98	87.5
	32	89.5	79	92	81.5
42	55	97.5	87	100	89.5
	32	91.5	81	94	83.5
44	55	99.5	89	102	91.5
	32	93.5	83	96	85.5
46	55	101.5	91	104	93.5
	32	97.5	85	98	87.5
48	55	103.5	93	106	95.5
	32	97.5	87	100	89.5
50	55	105.5	95	108	97.5
	32	99.5	89	102	91.5
52	55	107.5	97	110	99.5
	32	101.5	91	104	93.5
54	55	109.5	99	112	101.5
	32	103.5	93	106	95.5
56	55	111.5	101	114	103.5
	32	105.5	95	108	97.5

Avantgarde XXL 2 – Overall width (in cm)¹⁾


Seat width	Standard side panels and clothing protectors	Side panel with height adjustable armrest	Rear wheel with drum brake
50	69.5	72	73
52	71.5	74	75
54	73.5	76	77
56	75.5	78	79
58	77.5	80	81
60	79.5	82	83
62	81.5	84	85

Increase of the overall width because of the wheel camber setting of the rear wheels (all models, in cm)

Wheel camber	Rear wheel size		
	22"	24"	26"
0°	0	0	0
2°	4	4	4.5
4°	8	8.5	9

Seat heights — Avantgarde T, Ti 8.9, Ti Ultra, and VR

Anterior seat height — Avantgarde T, Ti 8.9, Ti Ultra, and VR (in cm)¹⁾


	Position	Caster wheel diameter					
		3"	4"	5"	140 mm	6"	7"
Short caster fork	Position 4						
	Position 3		42.4	43.5			
	Position 2	42.3	43.6	44.7	45.7	46.2	
	Position 1	43.6	44.9	46	47	47.5	49
Long caster fork	Position 4			46	47	47.5	49
	Position 3		46.1	47.2	48.2	48.7	50.2
	Position 2	46.1	47.4	48.5	49.5	50	51.5
	Position 1	47.3	48.6	49.8	50.7	51.2	52.7

INFORMATION The anterior seat height and posterior seat height are dependent on the chosen wheel size and mounting position. Please refer to the seat height table!

INFORMATION Seat inclination: 10 cm maximal seat height difference between anterior seat height and posterior seat height (equaling a seat inclination of about 15°).

INFORMATION Without indication of the desired seat inclination, it will be factory-set at approx. 4°.

Posterior seat height — Avantgarde T, Ti 8.9, Ti Ultra, and VR (in cm)¹⁾

Free mounting positions ²⁾	Rear wheel size/tyres			
	22"	24"	24" MTB	26"
1	37	40	41	42
2	38	41	42	43
3	39	42	43	43
4	40	43	44	45
5	41	43.5	45	46
6	42	44	45.5	47
7	43	45	46	48
8	44	46	47	49
9	45	47	48	50
10	45.5	48	49	50.5
11	46	49	50	51
12	47	50	51	52
13	48	51	52	53


INFORMATION The axle can be adjusted horizontally between the front and rear axle position in 6 different positions (A through F). Position A allows a passive centre of gravity, whilst in position F, the most active centre of gravity is reached. Unless otherwise specified, the rear wheel position will have the standard setting of C/7.

1) Measured at 0° seat angle, without seat cushion.

2) Mounting positions can be counted on the posterior brace of the rear wheel attachment device (see sketch).

Seat heights — Avantgarde XXL 2

Anterior seat height — Avantgarde XXL (in cm)¹⁾


	Position	Caster diameter		
		140 mm	7"	8"
Short caster fork	Position 4			
	Position 3			
	Position 2	45		
	Position 1	47	49	
Long caster fork	Position 4	47	49	
	Position 3	48	50	51
	Position 2	49.5	51.5	52.5
	Position 1	50.5	52.5	53.5

INFORMATION The anterior seat height and posterior seat height are dependent on the chosen wheel size and mounting position. Please refer to the seat height table!

INFORMATION Seat inclination: 10 cm maximal seat height difference between anterior seat height and posterior seat height (equaling a seat inclination of about 15°).

INFORMATION Without indication of the desired seat inclination, it will be factory-set at approx. 4°.

Posterior seat height — Avantgarde XXL (in cm)¹⁾

Mounting positions	Posterior seat height	
	24"	26"
1	40	42
2	42.5	44.5
3	44.5	46.5
4	47	49
5	49	51
6	51	53

INFORMATION The posterior seat height is dependent on the chosen wheel size and mounting position in the rear wheel attachment device. Please refer to the seat height table!

INFORMATION The axle can be adjusted horizontally between the front and rear axle position in 7 different positions (A through G). Position A allows a passive centre of gravity, whilst in position G, the most active centre of gravity is reached. Unless otherwise specified, the rear wheel position will have the standard setting of D/4.


1) Measured at 0° seat angle, without seat cushion.

Lower leg length (in cm)²⁾

	Footrest for short lower leg lengths		Short footrest bar (23 cm)		Long footrest bar (32 cm)	
	Min. (with 28 cm seat depth)	Max. (with 28 cm seat depth)	Min. (with 38 cm seat depth)	Max. (with 38 cm seat depth)	Min. (with 42 cm seat depth)	Max. (with 42 cm seat depth)
T (70°)	16	31	33.5	45	43	54
T (80°)	16	31	32	45.5	42	54
Ti (70°)	16	31	33.5	45	43	54
Ti (80°)	16	31	32	45.5	42	54
VR	—	—	32	46	39	55
XXL 2	—	—	32	46	39	55

1) Specified for narrow push ring installation (in case of wide installation: +20 mm)

2) Measured from anterior edge of the seat upholstery to the posterior edge of the footplate (lower leg length minus thickness of the seat cushion used)

9.2 Threshold values for wheelchairs transportable by train

INFORMATION

- The wheelchairs in this series fully satisfy the minimum technical requirements of regulation (EU) No. 1300/2014 for wheelchairs transportable by train. However, not all versions can comply with all threshold values due to different settings.
- With the help of the table that follows, you or the qualified personnel can take measurements and verify whether the concrete wheelchair meets the threshold values.

Feature	Threshold value (according to regulation (EU) No. 1300/2014)
Length [mm]	1200 (plus 50 mm for the feet)
Width [mm]	700 (plus 50 mm on each side for the hands when moving)
Smallest wheels ["]	approx. 3 or greater (according to the regulation, the smallest wheel must be able to accommodate a gap measuring 75 mm horizontal and 50 mm vertical)
Height [mm]	max. 1375; including a 1.84 m large male wheelchair user (95th percentile)
Turning radius [mm]	1500
Maximum weight [kg]	200 (for wheelchair and occupant, including any baggage)

Maximum obstacle height that can be overcome [mm]	50
Ground clearance [mm]	60 (at an upward slope angle of 10°, ground clearance must measure at least 60 mm under the foot rest for going forward at the end of the slope)
Maximum inclination angle on which the wheelchair will remain stable [°]	6 (dynamic stability in all directions) 9 (static stability in all directions, including with wheel lock engaged)

10 Re-use / disposal

10.1 Information on re-use

⚠ CAUTION

Risk due to improper use. Seat padding which come into direct contact with the skin may cause functional or hygienic risks for another patient using one and the same product. They should be replaced in case of re-use.

The wheelchairs for active use of the Avantgarde series are suitable for re-use.

Similar to second-hand machines or cars, products that are being re-used are subject to increased strain. Features and functions must not change in a way that could endanger patients or other persons within the product's life cycle.

In cases of re-use, the corresponding product must first be thoroughly cleaned and disinfected. Subsequently, the product must be examined by an authorised specialist to check the condition, wear, and possible damages. Any worn and damaged components as well as components which do not fit or are not suitable for the user must be replaced.

The service manual includes a service schedule, detailed information, and a list of the required tools.

10.2 Information on disposal

If a product is to be disposed of, all components of the product must be disposed of properly in accordance with the applicable environmental protection regulations in the respective country.

Kundenservice/Customer Service

Europe

Otto Bock HealthCare Deutschland GmbH
Max-Näder-Str. 15 · 37115 Duderstadt · Germany
T +49 5527 848-3455 · F +49 5527 848-1510
healthcare@ottobock.de · www.ottobock.de

Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
F +43 1 5267985
service-admin.vienna@ottobock.com · www.ottobock.at

Otto Bock Adria Sarajevo D.O.O.
Omladinskih radnih brigada 5
71000 Sarajevo · Bosnia-Herzegovina
T +387 33 766200 · F +387 33 766201
obadria@bih.net.ba · www.ottobockadria.com.ba

Otto Bock Bulgaria Ltd.
41 Tzar Boris III Blvd. · 1612 Sofia · Bulgaria
T +359 2 80 57 980 · F +359 2 80 57 982
info@ottobock.bg · www.ottobock.bg

Otto Bock Suisse AG
Luzerner Kantonsspital 10 · 6000 Luzern 16 · Suisse
T +41 41 455 61 71 · F +41 41 455 61 70
suisse@ottobock.com · www.ottobock.ch

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Otto Bock South Africa (Pty) Ltd
Building 3 Thornhill Office Park · 94 Bekker Road
Midrand · Johannesburg · South Africa
T +27 11 312 1255
info-southafrica@ottobock.co.za
www.ottobock.co.za

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Otto Bock HealthCare Canada
5470 Harvester Road
Burlington, Ontario, L7L 5N5, Canada
T +1 289 288-4848 · F +1 289 288-4837
infocanada@ottobock.com · www.ottobock.ca

Otto Bock HealthCare Andina Ltda.
Calle 138 No 53-38 · Bogotá · Colombia
T +57 1 8619988 · F +57 1 8619977
info@ottobock.com.co · www.ottobock.com.co

Otto Bock de Mexico S.A. de C.V.
Prolongación Calle 18 No. 178-A
Col. San Pedro de los Pinos
C.P. 01180 México, D.F. · Mexico
T +52 55 5575 0290 · F +52 55 5575 0234
info@ottobock.com.mx · www.ottobock.com.mx

Otto Bock HealthCare LP
11501 Alterra Parkway Suite 600
Austin, TX 78758 · USA
T +1 800 328 4058 · F +1 800 655 4963
usa.customerservice@ottobockus.com
www.ottobockus.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
Suite 1.01, Century Corporate Centre
62 Norwest Boulevard
Baulkham Hills NSW 2153 · Australia
T +61 2 8818 2800 · F +61 2 8814 4500
healthcare@ottobock.com.au · www.ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co., Ltd.
B12E, Universal Business Park
10 Jiuxianqiao Road, Chao Yang District
Beijing, 100015, P.R. China
T +8610 8598 6880 · F +8610 8598 0040
news-service@ottobock.com.cn
www.ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Unit 1004, 10/F, Greenfield Tower, Concordia Plaza
1 Science Museum Road, Tsim Sha Tsui
Kowloon, Hong Kong · China
T +852 2598 9772 · F +852 2598 7886
info@ottobock.com.hk · www.ottobock.com

Otto Bock HealthCare India Pvt. Ltd.
20th Floor, Express Towers
Nariman Point, Mumbai 400 021 · India
T +91 22 2274 5500 / 5501 / 5502
information@indiaottobock.com · www.ottobock.in

Otto Bock Japan K. K.
Yokogawa Building 8F, 4-4-44 Shibaura
Minato-ku, Tokyo, 108-0023 · Japan
T +81 3 3798-2111 · F +81 3 3798-2112
ottobock@ottobock.co.jp · www.ottobock.co.jp

Otto Bock Korea HealthCare Inc.
4F Agaworld Building · 1357-74, Seocho-dong
Seocho-ku, 137-070 Seoul · Korea
T +82 2 577-3831 · F +82 2 577-3828
info@ottobockkorea.com · www.ottobockkorea.com

Otto Bock South East Asia Co., Ltd.
1741 Phaholyothin Road
Kwaeng Chatuchark · Khet Chatuchark
Bangkok 10900 · Thailand
T +66 2 930 3030 · F +66 2 930 3311
obsea@ottobock.co.th · www.ottobock.co.th

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Ottobock has a certified Quality Management System in accordance with ISO 13485.

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Otto Bock Mobility Solutions GmbH