







ottobock.











Terra Flair

DE Gebrauchsanweisung	6
EN Instructions for use	13
FR Instructions d'utilisation	19
IT Istruzioni per l'uso	26
ES Instrucciones de uso	33
CS Návod k použití	40
TR Kullanma talimatı	47

INFORMATION

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INFORMATION

Last update: 2015-02-20

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Function

The Terra Flair wheelchair cushion has an anatomically shaped rest suspension foam base (see fig. 1, item 1) to provide extra cushioning. This distributes pressure over the largest possible surface area in order to reduce pressure points and provide good postural support.

The foam base is combined with a ROHO inflatable air insert positioned at the rear of the cushion area to provide optimal protection to the most vulnerable, bony prominences (see fig. 1, item 2).

The cushion is protected by a two-way stretch water-impermeable Dartex cover (see fig. 2).

The components of the wheelchair cushion fulfil the following functions:

Component	Figure	Function
Rest suspension foam base for cushioning; anatomically contoured	see fig. 3	Stable leg position due to built-in medial and lateral thigh guides
High-profile ROHO dry air floatation insert*	see fig. 4	Improved immersion and envelopment that ensures optimal pressure redistribution under the most vulnerable bony areas
Dartex fabric cover	see fig. 5	PU-coated fabric with two-way stretch for incontinent users

* The dry air floatation principle enables weight to be distributed evenly without shear forces.

2 Application

2.1 Indications for use

The wheelchair cushion is to be used **exclusively** as a wheelchair accessory for adult users. It provides comfort, postural guidance and pressure redistribution. An appropriately sized wheelchair cushion should be selected based on the user's size and should be used exclusively as an accessory for a correspondingly sized wheelchair.

The wheelchair cushion is particularly suitable for low to moderate positioning and preventing pressure sores in patients with paraplegia, hemiplegia, neurological deficits, and rheumatic nerve and muscle disorders as well as for geriatric patients.

Please note that the most appropriate wheelchair cushion for an individual should always be determined by a suitably qualified or experienced person based on a full risk assessment and used in conjunction with a full programme of care.

2.2 Area of Application

The wheelchair cushion can be used by people weighing up to **150 kg [330 lbs]** with no restrictions. We assume no liability if the maximum user weight is exceeded.

2.3 Contraindications

Seek professional advice before selecting an appropriate wheelchair cushion where there is an existing case or history of pressure sores/decubitus ulcers.

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative Contraindications

If redness of the skin or other problems occur during use, consult a doctor or therapist as soon as possible.

2.4 Further Usage Restrictions

The product is designed for use by **one person** only. Re-use of the product is not permitted due to hygiene considerations.

2.5 Qualification

This product may only be selected and adapted by trained qualified personnel.

3 Safety

3.1 Explanation of Warning Symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General Safety Instructions

 **CAUTION**

Altered centre of gravity resulting from use of the seat cushion

Risk of tipping over or falling due to unstable settings

- ▶ After installing any seating into a wheelchair, always check that the user's sitting position and the wheelchair's centre of gravity are appropriate for the intended use and the user's abilities. Work with the user to find the optimal sitting position.
- ▶ If required, adjust the sitting position or axle position to optimise stability, or alternatively assess the suitability of anti-tip devices.

 **CAUTION**

Lack of testing prior to long-term use

Risk of skin reddening or pressure sore formation due to incorrectly selected seat cushion

- ▶ Any seat cushion should be tried out by the user under medical supervision for several hours prior to long-term use. The skin should be monitored for redness and/or damage during this time.

 **NOTICE**

Use of a wrong air pump

Damage to the product

- ▶ Only use the supplied hand air pump for inflating the ROHO insert.

3.3 Information to be Provided to the Patient/Attendant

 **CAUTION**

Skin damage

Risk of skin reddening during long-term use

- ▶ Do not use the wheelchair cushion until it has been adjusted by a suitably qualified or experienced person.

- ▶ Examine the skin for redness on a regular basis. Skin reddening is a clinical indicator of tissue damage.
- ▶ If redness of the skin or other problems occur during use, stop using the seat cushion immediately. Only start using the seat cushion again after consulting the prescribing doctor or therapist.

⚠ CAUTION

Failure to observe care instructions

Risk of pressure sores due to incorrect skin care

- ▶ All seat cushions were developed by the manufacturer to reduce pressure on the buttocks. A seat cushion alone, however, cannot eliminate sitting pressure or possible pressure sores entirely.
- ▶ When using the seat cushion, always pay attention to correct skin care – especially through a healthy diet, hygiene and regular pressure relief.

⚠ CAUTION

Improper handling of ignition sources

Risk of burns due to seat cushions, padding and upholstery catching on fire

- ▶ The seat cushions, padding and upholstery are not highly flammable, but can catch fire. Therefore the utmost caution must be exercised in the vicinity of an open flame.
- ▶ Keep away from all ignition sources, especially lit cigarettes.

NOTICE

Improper use

Risk of damage to the cushion surface due to user error

- ▶ Do not allow the wheelchair cushion to come into contact with sharp objects. This also applies to animals such as pet cats with sharp claws.
- ▶ Do not expose the wheelchair cushion to excessive heat or open flame.
- ▶ If the wheelchair cushion is expected to come into contact with liquid, such as spilt drinks or episodes of incontinence, always use it in conjunction with a liquid-repellent cover.
- ▶ Only use the official Ottobock incontinence barriers for this product. These barriers have been optimised for maximum pressure redistribution in the wheelchair cushion. Contact your therapist or specialist dealer to obtain a spare Ottobock cover.

4 Delivery

4.1 Scope of Delivery

The product is delivered preinstalled. The scope of delivery includes:

Component	Figure
Wheelchair Seat Cushion	
Extra-cushioning foam base	see fig. 3
High-profile ROHO insert	see fig. 4
Dartex cover	see fig. 5
Accessories	
Hand air pump	see fig. 6
Insert repair kit	see fig. 7
Instructions for use	not illustrated

5 Preparation for Use

INFORMATION

- ▶ The procedures described in this section may not be carried out by the user of the product.
- ▶ These instructions for use must be given to the user along with the seat cushion.

5.1 Positioning

- 1) Place the cushion into the seat with the 'REAR' label facing the back of the wheelchair (see fig. 5).
- 2) The flat surface of the cushion should be facing down and the contoured surface should be facing up.
- 3) The hook-and-loop strap sewn to the cushion is designed to be fastened to the corresponding straps found on most wheelchair seats to allow the cushion to be secured in place.

5.2 Setting Up

- 1) Place the wheelchair cushion on the chair, making sure it is centred and right side up. The valve must be located in the front right corner (see fig. 8).
- 2) Turn the valve anti-clockwise to open (see fig. 9).
- 3) Attach the flexible hose of the pump to the valve (see fig. 10).
- 4) Inflate the cushion until it begins to slightly arch upwards (see fig. 11).
- 5) Close the air valve by twisting in a clockwise direction and then disconnect the pump from the valve (see fig. 12).
- 6) Lower the user down carefully on to the wheelchair cushion. The cushion must be centred below the user.

INFORMATION: Ensure the individual is sitting in the desired sitting position.

- 7) Slide your hand palm down between the cushion surface and the user's buttocks (see fig. 13).

INFORMATION: You should be able to feel the ischial tuberosity (IT).

- 8) Keep your hand under the user's ischial tuberosity and rotate the valve anti-clockwise to slowly let out air (see fig. 9).
- 9) Allow the ROHO insert to sink until your fingertips can only move up and down over a distance of **2 – 2.5 cm** (see fig. 14).
- 10) Close the valve by turning it clockwise (see fig. 12).
- 11) Wiggle your fingers to ensure the correct clearance. If there is **less than 2 cm** clearance between the user and the wheelchair seat, air must be added and the process repeated.

INFORMATION: For all dry air floatation cushions a hand check should be performed daily. The seat cushion needs to be readjusted only when there is a change in air level during a hand check, or after a puncture in the air insert has been repaired.

- 12) Slide the valve into the cover and close the zip.

6 Use

6.1 Positioning

Correct positioning of the seat cushion in the wheelchair must be maintained for daily use:

- The label that says 'REAR' has to face the back (see fig. 5).
- The hook-and-loop strap sewn to the seat cushion serves to secure it in the wheelchair.

6.2 Removing/Putting On the Cover

Removing the cover

The cover can be removed for cleaning of the cushion:

- 1) Unzip the cover along all three sides (see fig. 15).
- 2) Lift up the cover from the seating surface of the cushion (see fig. 16; see fig. 17).

- 3) Pull the tube on the ROHO insert out of the cover (see fig. 18).
- 4) Carefully remove the foam base with the ROHO insert from the cover (see fig. 19).

INFORMATION: See the cleaning instructions for information on washing.

Putting on the cover

- 1) Ensure the ROHO insert is in the correct position.
- 2) Unzip the cover along all three sides (see fig. 15).
- 3) Insert the cushion into the cover with the zipper running along the rear and side edges of the cushion (see fig. 19).
- 4) Push the tube on the ROHO insert through the cover (see fig. 18).
- 5) Pull the cover over the cushion (see fig. 17; see fig. 16).
- 6) Zip up the cover again (see fig. 20).
- 7) Ensure the cover is fitted well and there are no large folds or wrinkles on the top surface.

6.3 Cleaning

NOTICE

Improper cleaning of the wheelchair cushion

Risk of damage to cushion due to user error

- ▶ Never wash the cover at a temperature higher than the recommended as this may cause it to shrink. This could impede the cushion's ability to redistribute pressure or it could damage the fabric.
- ▶ Do not use any bleaching agents for cleaning.
- ▶ Never dry the foam components or inserts in a dryer.

INFORMATION

- ▶ Replace the cover when it shows signs of advanced wear.
- ▶ Should it be necessary to disinfect the cover, add a standard hygiene detergent during washing (e.g. Sagrotan® laundry hygiene detergent).
- ▶ In case of heavy incontinence, an additional incontinence barrier should be used. Please observe the care and cleaning instructions for each incontinence barrier used.

Cover

- 1) Unzip and remove the cover for separate washing.
- 2) Zip up the cover before washing.
- 3) The cover can be washed at 60 °C [140 °F] using a mild, environmentally friendly detergent.
Recommendation: Wash the cover using the delicate cycle at **40 °C [104 °F]** to prevent excessive wear.
- 4) Allow to air dry. When using an automatic clothes dryer, select the delicate cycle at low temperature.

7 Repair

INFORMATION

- ▶ If the ROHO insert is punctured by a sharp object, refrain from using the cushion as far as possible.
- ▶ If it is not possible to repair the ROHO insert yourself, contact your therapist or your specialist dealer.

Removing the ROHO insert

- 1) Unzip and remove the cover as described above (see Page 16).

- 2) Gently ease the ROHO insert and valve away from the foam with your fingertips. Do not pull on the air cells (see fig. 21).

Sealing the hole

- 1) Locate the puncture.
INFORMATION: If the puncture is not visible, inflate the air chamber, close the valve and immerse the insert in water. The puncture should be detectable by the air bubbles escaping from the insert.
- 2) Follow the instructions for repair provided with the repair kit (see fig. 7).
- 3) Re-inflate the ROHO insert with air.
- 4) Immerse the ROHO insert in water to check that the puncture has been sealed.

Inserting the ROHO insert

- 1) Dry the ROHO insert thoroughly.
- 2) Position the ROHO insert so that the valve fits into the valve hole in the foam base (see fig. 22).
- 3) Carefully reinsert the ROHO insert. The edge must be located between the foam and the black cushion base. Do not pull on the air cells (see fig. 21).
- 4) Put the cover back on (see Page 16).
- 5) Set up the cushion for the user again (see Page 16).

8 Legal Information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as class I accessories according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

8.3 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

9 Technical data

General information	
Functional height	70 mm [3 in]
Weight (420 x 420 mm; 16 x 16 in*)	1474 g [3.25 lbs]
Maximum user weight	150 kg [330 lbs]

Foam base	
Material	PUR extra-cushioning foam
Foam density (ISO 845)	100 kg/m ³
Compressive stress of the foam (ISO 3386)	11 kPa

Cover	
Material	Dartex
Maximum washing temperature	60 °C

ROHO insert	
Material	Neoprene rubber

* without incontinence cover

Available sizes – Seat cushions

Seat width in mm [in]*	Seat depth in mm [in]*		
	405 [16]	455 [18]	505 [20]
380 [15]	X		
405 [16]	X	X	X
430 [17]	X	X	X
455 [18]	X	X	X
505 [20]	X	X	X

* Seat width and seat depth of the wheelchair seat

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© Ottobock · 647G877-04-1502
Otto Bock Korea HealthCare Inc.
South East Asia Co., Ltd.
Otto Bock HealthCare GmbH
Otto Bock Mobility Solutions GmbH