



Zenit

EN Instructions for use (user) 3

Table of contents

EN

1	Foreword	5
2	Product description	5
2.1	Function	5
2.2	Product Overview	6
3	Safety	6
3.1	Explanation of warning symbols	6
3.2	Safety Instructions for Use	7
3.3	Further information	10
3.4	Nameplate and warning labels	10
3.4.1	Nameplate.....	10
3.4.2	Warning labels	10
4	Delivery	11
4.1	Scope of delivery	11
4.2	Options	11
4.3	Storage.....	11
4.3.1	Storage during daily use	11
4.3.2	Storage during extended disuse.....	11
5	Preparation for use	11
5.1	Assembly	11
6	Use	12
6.1	Techniques for using a wheelchair	12
6.2	Further instructions for use.....	13
6.3	Getting in and transferring	13
6.4	Legrests	14
6.5	Seat cushion and back upholstery	15
6.6	Side panels	15
6.6.1	Adjusting the armrest height	15
6.7	Push handles.....	15
6.7.1	Adjusting the push handles	15
6.8	Rear wheels	15
6.8.1	Removing/mounting the rear wheels	16
6.9	Caster wheels and caster forks.....	17
6.9.1	Approach in case of stiffness	17
6.10	Wheel locks	17
6.10.1	Using the wheel locks.....	17
6.11	Anti-tipper/tip-assist	18
6.11.1	Activating and deactivating the anti-tipper	19
6.11.2	Using the tip-assist.....	19
6.12	Lap belt	20
6.13	Additional options	20
6.14	Disassembly and transport	20
6.15	Use in vehicles for transporting persons with reduced mobility.....	22
6.16	Care.....	22
6.16.1	Cleaning.....	22
6.16.1.1	Cleaning Seat Cushions	22
6.16.1.2	Cleaning belts	22
6.16.2	Disinfection	23
7	Maintenance and repair	23
7.1	Maintenance.....	23
7.1.1	Maintenance intervals	23
7.1.2	Maintenance tasks	23
7.2	Repair	24
7.2.1	Inner tube, rim tape and tyre replacement	24

8	Disposal	25
8.1	Disposal Information.....	25
9	Legal information.....	25
9.1	Liability	25
9.2	CE Conformity.....	25
9.3	Warranty	26
9.4	Service Life	26
9.5	Trademarks	26
10	Technical data	26
11	Appendices	27
11.1	Threshold values for wheelchairs transportable by train	27

1 Foreword

INFORMATION

Last update: 2016-12-07

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Have yourself instructed by qualified personnel in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

INFORMATION

- ▶ New information regarding product safety and product recalls can be obtained from the Customer Care Center (CCC) at oa@ottobock.com or from the manufacturer's service department (see inside back cover or back page for addresses).
- ▶ You can request this document as a PDF file from the Customer Care Center (CCC) at oa@ottobock.com or from the manufacturer's service department (see inside back cover or back page for addresses). It is possible to increase the display size of the PDF document.
- ▶ For further questions about the instructions for use, please contact the qualified personnel who issued the product to you.

You have purchased a high-quality product which can be put to versatile, daily use at home and outdoors.

In order to exclude injuries of any type, familiarise yourself with the handling, functions and use of the product before using it. These instructions for use provide the necessary information.

Please note the following in particular:

- Attendants should also be instructed in the use of the product by qualified personnel with the help of these instructions for use.
- The product has been adapted to the needs of the user. Further changes may be made only by qualified personnel. We recommend checking the product settings regularly in order to assure an optimum fit over the long term. For growing children and youths in particular, fitting should be performed every six months.
- In case of questions or problems, please consult the qualified personnel that adapted the product or the manufacturer's service (see inside back cover or back page for addresses).
- The product may be combined only with the options listed here. The manufacturer assumes no liability for combinations with third-party medical devices and/or accessories not included in the modular system. Please also observe the information in the section "Liability".
- Service and repairs to the product may only be carried out by qualified personnel. If you have any problems, please contact your specialist dealer. This ensures that any necessary repairs will be made exclusively with Ottobock spare parts.
- Your product may differ from the models shown. In particular, not all the options described in these instructions for use will be installed on your product.
- The manufacturer reserves the right to make technical changes to the model described in these instructions for use.

2 Product description

2.1 Function

The wheelchair is designed solely for individual use by persons who are unable to walk or who have a walking impediment, and can be operated either by the patient or by another person.

The wheelchair can be used on solid ground both indoors and outdoors.

2.2 Product Overview



Zenit

Maximum load capacity: 120 kg

The Zenit is equipped with an open frame made of aluminium, or optionally of carbon. Thanks to the Easy-Go folding system, it can be unfolded with very little force.

1	Back/back upholstery	8	Caster wheel
2	Side panel	9	Rear wheel with push ring
3	Seat/seat pad	10	Rear wheel adapter
4	Scissor wheel lock	11	Folding mechanism
5	Quick-release axle release button	12	Frame
6	Footplate	13	Back tube
7	Caster fork	14	Push handle

3 Safety

3.1 Explanation of warning symbols

WARNING	Warning regarding possible serious risks of accident or injury.
CAUTION	Warning regarding possible risks of accident or injury.
NOTICE	Warning regarding possible technical damage.

3.2 Safety Instructions for Use

Hazards during preparation for use

⚠ WARNING

Independent modification of settings

Serious injuries to the user due to unallowable changes to the product

- ▶ Do not modify the settings established by the qualified personnel.
- ▶ In case of problems with the settings (unsatisfactory seating position, caster wheel wobble, etc.) please contact the qualified personnel who adjusted your product.

⚠ WARNING

Settings outside the safety limits

Tipping over, falling due to adjustment errors because of failure to observe the physical/psychological requirements of the user

- ▶ In certain versions and with certain settings, the wheelchair tends to tip backwards due to its design. This is an intentional feature designed to enable users with the corresponding physical prerequisites to manoeuvre quickly and nimbly.
- ▶ Users with the necessary physical and psychological preconditions can be supplied with such versions/settings. If the physical/psychological preconditions change, the wheelchair may no longer be used with these settings. Inform the responsible qualified personnel immediately in this case.

⚠ WARNING

Improper handling of packaging materials

Risk of suffocation due to neglect of the duty to supervise

- ▶ Packaging materials must be kept out of the reach of children.

Risk of hand injuries

⚠ CAUTION

Pinching on components

Pinching, crushing due to lack of caution in danger areas

- ▶ Do not reach between the rear wheel and wheel lock or rear wheel and side panel when driving the wheelchair.
- ▶ Do not reach into the spokes of the rotating rear wheel.
- ▶ Pay attention not to pinch parts of your body at the wheel lock lever or between the side or frame parts.

⚠ CAUTION

Heat development when braking with push rings

Burns due to insufficient hand protection

- ▶ Wear wheelchair gloves when travelling at high speeds.

Hazards while driving

⚠ CAUTION

Lack of driving experience

Tipping over, falling due to errors in handling the product

- ▶ Practice on level, open ground first.
- ▶ Learn with the support of an assistant how the wheelchair reacts to changes in the centre of gravity.

⚠ CAUTION

Leaning forward in the wheelchair

Tipping over, overturning due to incorrect centre of gravity

- ▶ Do not lean too far out of the wheelchair when trying to pick up objects.
- ▶ Lean well forward with your upper body when climbing gradients, overcoming obstacles on gradients or on ramps. If users cannot lean their upper body forward, accompanying persons must support the wheelchair from behind.

⚠ CAUTION

Risky operation

Falling, tipping over backwards due to approaching obstacles incorrectly

- ▶ Push slowly when crossing obstacles (e.g. steps, curbs) and negotiating uphill or downhill slopes and inclines.
- ▶ Never cross obstacles at an angle. Always approach obstacles head on (at an angle of 90°).
- ▶ Raise the front wheels before crossing obstacles.
- ▶ Avoid collisions with obstacles and dropping off curbs/ledges.
- ▶ Avoid riding cross-country.

⚠ CAUTION

Improper use of the wheel lock

Falling due to abrupt braking, rolling away of the wheelchair, damage to the wheel lock

- ▶ Do not use the wheel lock when travelling.
- ▶ Apply the wheel lock to prevent the wheelchair from moving on uneven ground or during transfers (e.g. into a car).

⚠ CAUTION

Lack of tipping resistance on public transport

Tipping over, falling of the user, damage to the product due to incorrect positioning of the wheelchair

- ▶ Always observe currently applicable legal requirements when using public transport.
- ▶ Always ensure that you are held in place securely when travelling on public transport. To do so, use the wheelchair areas, wheelchair bays and restraint systems provided. Firmly tighten the wheel locks.
- ▶ Please note that the anti-tipper, if activated, may be exposed to heavy loads when the public transport vehicle starts to move. In order to avoid damage, the manufacturer recommends that the wheelchair is positioned at a right angle to the direction of travel if no restraint system for passengers with reduced mobility is available.

⚠ CAUTION

Incorrect handling at level crossings

User may fall or tip over due to driving error

- ▶ Never go over a level crossing diagonally, but always at a right angle to the tracks. If you travel diagonally, the wheelchair casters may get caught in the tracks.
- ▶ Only cross railway systems and railway tracks in the designated areas.
- ▶ Stop before the level crossing and check to the left and right to make sure it is safe to cross.

⚠ CAUTION

Driving in the dark

Risk of collisions with other traffic participants due to lack of lighting

- ▶ Wear bright clothing or clothing with reflectors.
- ▶ Install active lighting on your product.
- ▶ Ensure that the reflectors on the product are clearly visible.

Hazards when overcoming obstacles

WARNING

Overcoming steps and obstacles without assistance

Tipping over, falling of the user due to failure to observe transportation instructions

- ▶ Always have accompanying persons help you negotiate steps and other obstacles.
- ▶ Use available facilities (e.g. access ramps or lifts).
- ▶ If such facilities are not available, have 2 assistants carry you over the obstacle.

WARNING

Improper lifting by attendants

Tipping over, falling of the user due to lifting by removable components

- ▶ Only lift the product by firmly mounted components (e.g. main frame, push handles).
- ▶ Be sure that the clamping levers of the height-adjustable push handles are firmly tightened.

Hazard in case of broken skin

CAUTION

Prolonged contact with broken skin

Reddening of the skin or pressure points due to contamination with germs or incorrectly adapted seat cushion

- ▶ Check before using the product that the skin is not broken in areas subject to prolonged pressure (e.g., buttocks, back and rear side of the thighs).
- ▶ In case of problems, please contact the qualified personnel who adjusted the product.
- ▶ No liability is assumed for injuries caused by the use of the product with broken skin.

Hazards due to fire, heat and cold

CAUTION

Extreme temperatures

Hypothermia or burns through contact with components, failure of components

- ▶ Do not expose the product to any extreme temperatures (e.g. direct sunlight, sauna, extreme cold).

Hazards due to improper use of the product

WARNING

Overloading

Severe injuries if the wheelchair tips over due to overloading, damage to the product

- ▶ Do not exceed the maximum load capacity (see the nameplate and section "Technical data").
- ▶ Please note that certain options and add-on components will reduce the remaining load capacity.

WARNING

Exceeding the service life

Serious injuries due to failure to observe the manufacturer's requirements

- ▶ Using the product beyond the specified expected service life (see Page 26) leads to increased residual risk and should only take place subject to the due diligence and deliberations of qualified personnel.
- ▶ If the service life is reached, the user or a responsible attendant should contact the qualified personnel who fitted the product or the manufacturer's servicing department (see inside rear cover or back page for address). Here the user can obtain information about known risks and the current options for refurbishing the product.

CAUTION

Uncontrolled driving behaviour, unexpected sounds or odours
 Falling, tipping, collision with persons or nearby objects due to defects

- ▶ If any faults, defects or other hazards that can lead to personal injury are detected, the product must be taken out of service immediately. This includes uncontrolled movements as well as sounds that are unexpected or previously not noted or odours that deviate significantly from the state of the product at the time of delivery.
- ▶ Please contact your authorised dealer.

NOTICE

Use under incorrect environmental conditions
 Damage to product due to corrosion or abrasion

- ▶ Do not use the product in salt water.
- ▶ Also avoid – if possible – sand or other dirt particles that may damage the wheel bearings.

3.3 Further information

INFORMATION

Even in the event of compliance with all applicable guidelines and standards, alarm systems (e.g. in department stores) may respond to your product. Should this happen, remove your product from the area where the alarm was triggered.

3.4 Nameplate and warning labels

3.4.1 Nameplate

Label	Meaning	
	A Product; product reference number	
	B CE marking – product safety in accordance with EU directives	
	C Maximum load capacity (see section "Technical data")	
	D Manufacturer information, address, country of origin	
	E Serial number*	
	F Manufacturing date**	
	G Read the instructions for use before using the product.	
	H Only for CLT option: The product has not been approved by the manufacturer for use as a seat in vehicles for transporting persons with reduced mobility	
	I European Article Number/Global Trade Item Number	
	K Product version***	
	The nameplate is located on the crossbrace.	

* YYYY = year of manufacture; WW = week of manufacture; PP = production site; XX...X = sequential production number

** YYYY = year of manufacture; MM = month of manufacture; DD = day of manufacturer

*** NN = country code; VV = country version; C = configuration

3.4.2 Warning labels

Label	Meaning
<p>Achtung! Die Antriebsräder sind Ihren Angaben entsprechend vorpositioniert. Die voreingestellte Position muß vor Inbetriebnahme vom Rollstuhlfahrer mit sichernder Unterstützung durch Fachpersonal auf Kippsicherheit und Funktion überprüft werden.</p> <p>Attention! The rearwheels are positioned according to your specifications. Before putting the wheelchair into operation the pre-adjusted position must be checked for tilting safety and proper functioning by the user supported safely by specialist staff.</p>	<p>(Safety loop on the wheelchair) Caution! The rear wheels have been pre-positioned according to your indications. Before using the wheelchair, this pre-adjusted position must be checked by the wheelchair user for function and stability against tipping with the assistance of qualified personnel.</p>

4 Delivery

4.1 Scope of delivery

The wheelchair is delivered preassembled and with the rear wheels removed packed in a crate.

The scope of delivery includes:

- Preassembled wheelchair
- 2 rear wheels
- Options according to the order
- Instructions for use (user)

4.2 Options

The standard model can be fitted to the user's personal requirements thanks to a large range of options. For use of these options: see Page 12 ff.

4.3 Storage

4.3.1 Storage during daily use

The wheelchair must be stored in a dry place.

4.3.2 Storage during extended disuse

The wheelchair must be stored in a dry place. During longer storage periods, an ambient temperature between **-10 °C and +40 °C** must be maintained.

It is not necessary to disassemble or fold up the wheelchair.

During extended storage, the knee lever wheel lock on wheelchairs with PU tyres (= tubeless tyres) must be released since tyre deformation may otherwise result.

Tyres contain chemical substances that can react with other chemical substances (such as cleaning agents, acids, etc.).

5 Preparation for use

5.1 Assembly

⚠ CAUTION

Exposed pinch points

Crushing, pinching due to incorrect handling

- ▶ When folding the wheelchair out or together, only grip by the specified components.

⚠ CAUTION

Failure to verify readiness for use before putting into operation

Tipping, falling due to incorrect adjustment or installation

- ▶ Before you start using the wheelchair for the first time, check the chosen wheelchair settings with the support of qualified personnel.
- ▶ After every assembly, check for proper mounting of the rear wheels. The quick-release axles must be securely locked in the receiver bushings.
- ▶ Pay particular attention to the stability against tipping, free running of the rear wheels and correct function of the wheel locks.
- ▶ Check the tyre pressure. The correct tyre pressure is printed on the sidewall. Ensure that the pressure is the same in both tyres.

INFORMATION

On the topic of disassembly/transport: see Page 20.

It only takes a few simple steps to prepare the wheelchair (see fig. 2) for use:

- 1) Insert the rear wheels into the quick-release axle mountings (see fig. 3). It must not be possible to remove the quick-release axles after releasing the push-button.
- 2) Unfold the wheelchair by pulling the rear wheel to the side (see fig. 4).
Ensure that the folding mechanism is fully unfolded. Hearing the folding mechanism engage audibly indicates that the folding process is complete.

- 3) With the "Folding backrest" option: Pull the release strap for the folding mechanism of the backrest in the direction of travel of the wheelchair and fold the backrest up (see fig. 5). Ensure that the backrest latches in the top position (see fig. 6).
- 4) If necessary, fold the push handles up (see fig. 7).
- 5) Fit the seat cushion. The seat cushion is secured against sliding by being pressed onto the hook-and-loop closure.



6 Use

6.1 Techniques for using a wheelchair

Information and practical tips on fundamental techniques for using the wheelchair can be found in the brochure "Techniques for using a wheelchair" supplied with the product.

Should the brochure not be included with the product, it can be ordered under reference number 646D60. Further information on obtaining the brochure can be provided by the qualified personnel supplying the wheelchair or by Ottobock Service (addresses: see back inside cover or back cover).

6.2 Further instructions for use

Attaching loads (e.g. backpacks) can adversely affect stability. Therefore, when surmounting obstacles or driving on ramps or inclines, always remove the additional load or activate the anti-tipper.

The load must not exceed 5 kg.

We recommend that you always attach a wheelchair backpack to the push handles with the attachment straps provided.

- The recommended overall width for manual wheelchairs in an operational state is **700 mm**. This specification should ensure unhindered use of escape routes, for example. Please note that the wheelchair dimensions can exceed the recommended value in versions with very large seat widths (for more information see see Page 26 ff.).
- The wheelchairs in this series fully satisfy the minimum technical requirements for wheelchairs transportable by train. However, please note that it is possible that not every individual wheelchair will fulfil the minimum requirements due to the variation in designs (please see see Page 27 for further details).

6.3 Getting in and transferring

⚠ CAUTION

Incorrect handling when getting in

Risk of falling, tipping or rolling away due to incorrect handling

- ▶ Activate the wheel lock each time before you get into or out of or transfer to or from the wheelchair.
- ▶ Always get into the wheelchair from the side, where possible.
- ▶ Never step on the footrests/footplates when getting in and out.
- ▶ Do not support yourself on the wheel lock when getting in/out.

⚠ CAUTION

Incorrect caster wheel position when leaning forward in the wheelchair

Tipping over, falling due to incorrect caster wheel positioning

- ▶ Prior to activities that require you to bend forward in the wheelchair (e.g. tying your shoes), maximise the stability of the wheelchair.
- ▶ In order to do so, push the wheelchair backwards until the caster wheels turn forward.

Getting into the wheelchair in combination with transfer from one wheelchair to another is described below.

The wheelchair is designed for persons who normally move around independently in a wheelchair. Hence independently getting in/transferring is described here. If an assistant is required, the steps described here should be performed analogously with support:

- 1) Manoeuvre the wheelchairs such that they are positioned at an angle of approx. 45° to each other (see fig. 8). Ensure that the wheel locks do not become entangled.
- 2) Engage the wheel locks of the wheelchair currently in use.
- 3) On the other wheelchair, engage the wheel locks and fold up the footrest (footplate) (see fig. 9).
- 4) Move your feet from the footplate onto the floor and move forward in the wheelchair (see fig. 10).
- 5) Make the transfer to the second wheelchair (see fig. 11). The method illustrated is one possible way to support yourself.

INFORMATION Find the most suitable method for yourself by practicing with the help of an assistant.

- 6) Finally fold the footplate down, place your feet on the footrest and release the wheel locks.
- The wheelchair can now be used.



6.4 Legrests

⚠ CAUTION

Improper use of the legrests

Tipping over, falling due to errors in handling the product

- ▶ Fold up the footplates before getting in.

The legrests support the user's legs.

The height of the legrest has been adjusted by qualified personnel to the length of the user's lower legs.

The angle of the footrest has been set by the qualified personnel so that it allows the ankles to rest in a comfortable position.

The depth of the footplate is **140 mm**.



6.5 Seat cushion and back upholstery

⚠ WARNING

Seat cushion and back padding may ignite

Burns due to user error

- ▶ The seat and backrest upholstery as well as padding fulfil the requirements for flame resistance according to the norms ISO 8191-2 or DIN EN 1021-2. However, they may still ignite if fire is handled improperly or negligently.
- ▶ Keep away from all ignition sources, especially lit cigarettes.

⚠ CAUTION

Wear of the seat and back upholstery

Loss of functionality due to prohibited continued use

- ▶ Have the seat and back upholstery replaced immediately in case of damage.

The product is equipped with seat and back upholstery. Hook-and-loop fasteners on the seat upholstery are used to attach the seat cushion. The seat cushion ensures pressure relief during use of the wheelchair. It was chosen by qualified personnel according to the needs of the user.

6.6 Side panels

⚠ CAUTION

Pinching at the side panels

Pinching, crushing due to lack of caution in danger areas

- ▶ Pay attention not to pinch parts of your body between the side panel or frame parts.

The side panels protect the user and his/her clothing from getting dirty.

If armrests have been installed on the wheelchair, they offer the user additional support for the forearms.

6.6.1 Adjusting the armrest height

The armrest height can be changed without the use of tools.

- 1) Hold the push-button on the inside of the side panel down.
- 2) Put the armrest in the desired position.
- 3) Release the push-button. The armrest locks into place automatically.

6.7 Push handles

The push handles make pushing the wheelchair easier for the attendant.

Some of the push handle designs allow the height to be adjusted to suit the needs of the attendant.

6.7.1 Adjusting the push handles

If the "height-adjustable push handle" option is used, the push handles can be adjusted in order to make pushing easier for the attendant.

- 1) Release the clamping lever.
 - 2) Adjust the height of the push handle.
 - 3) Lock the clamping lever securely.
- Both push handles must be adjusted to the same height.

6.8 Rear wheels

⚠ WARNING

Incorrect fitting of the rear wheels

Tipping, falling over of the user due to wheels coming off

- ▶ After each assembly, verify the proper fit of the rear wheels. The quick-release axles must be firmly locked in the wheel attachment.

⚠ CAUTION**Reaching into exposed drive parts**

Crushing, pinching due to incorrect handling

- ▶ Do not reach between the rear wheel and wheel lock or rear wheel and side panel when driving the wheelchair.
- ▶ Do not reach into the spokes of the rotating rear wheel.

⚠ CAUTION**Defective tyres**

Accidents/falling due to poor traction, reduced braking force or lack of manoeuvrability

- ▶ Maintain sufficient tyre pressure. Note the information in the section "Technical data" or on the tyre sidewall.
- ▶ Ensure that the pressure is the same in both tyres.
- ▶ Ensure that the tyres have sufficient tread depth.

⚠ CAUTION**Heat development when braking with push rings**

Burns due to insufficient hand protection

- ▶ Wear wheelchair gloves when travelling at high speeds.

The wheelchair is moved, steered, braked and stopped with the help of the push rings on the rear wheels.

To make transport easier, rear wheels with quick-release axle can be taken off the wheelchair.

6.8.1 Removing/mounting the rear wheels**Rear wheels with quick-release axle**


- 1) Release the wheel lock.
- 2) Grip the area between the spokes near the hub with your fingers.
- 3) Use your thumb to press the push-button on the quick-release axle.
- 4) Pull off/push on the rear wheel. The quick-release axles must not be able to be removed after the push-button has been released.

**Rear wheels with lock-release assist**

The lock-release assist allows the rear wheels to be installed even with heavily restricted movement of the arms and hands.

- 1) Open the lock-release assist according to the grip setting (see fig. 15= push down; see fig. 16 = pull up). The lock is open when the bar is horizontal (90° position).
- 2) Pull off / push on the rear wheel.
- 3) Close the lock-release assist (see fig. 15, see fig. 16 = pull up or push down).

INFORMATION: The wheels must fit on to the wheel attachments without there being any play after closure.



6.9 Caster wheels and caster forks

⚠ CAUTION

Incorrect caster wheel position when leaning forward in the wheelchair

Tipping over, falling due to incorrect caster wheel positioning

- ▶ Prior to activities that require you to bend forward in the wheelchair (e.g. tying your shoes), maximise the stability of the wheelchair.
- ▶ In order to do so, push the wheelchair backwards until the caster wheels turn forward.

The combination of caster wheels and caster forks ensures the ability to hold a straight line and navigate bends securely.

The caster wheels and caster forks have been chosen by the qualified personnel according to the needs of the users.

6.9.1 Approach in case of stiffness

In case of stiffness, the caster wheel axles should be cleaned and oiled.

Oiling the caster axle

- 1) Remove dirt (such as hairs) from the caster axle between the caster wheel and caster fork.
- 2) Lubricate the caster axle between the caster wheel and caster fork with a few drops of thin, resin-free oil (sewing machine oil).

6.10 Wheel locks

The wheel locks secure the parked wheelchair against rolling away.

6.10.1 Using the wheel locks

⚠ CAUTION

Improper use of the wheel lock

Falling due to abrupt braking, rolling away of the wheelchair, crushing the hands

- ▶ Do not use the knee lever wheel lock as a driving brake.
- ▶ Always engage the wheel lock on both sides.
- ▶ When parking on uneven ground or during transfers (e.g. into a car), secure the wheelchair by activating the knee lever wheel lock.
- ▶ Do not reach between the rear wheel and the knee lever wheel lock when driving the wheelchair.
- ▶ Ensure the knee lever wheel lock is properly adjusted (approx. **5 mm** gap to the tyres, technical changes reserved) and the tyre pressure is sufficient.


Activating/deactivating the knee lever wheel lock

- 1) Push the handle of the knee lever wheel lock forward (see fig. 17).
 - The wheel lock bolt secures the wheel.
- 2) Pull the handle of the knee lever wheel lock up.
 - The wheel lock bolt releases the wheel.



Activating/deactivating the scissor wheel lock (Outfront scissor wheel lock)

- 1) Reach under the seat and pull the handle of the scissor wheel lock rearwards from the side (see fig. 18).
→ The wheel lock bolt secures the wheel.
- 2) Pull the handle of the scissor wheel lock forwards to unfold (see fig. 19).
→ The wheel lock lever releases the wheel.



6.11 Anti-tipper/tip-assist

⚠ WARNING

Anti-tipper not activated

Tipping over, falling of the user due to improper operation of a safety device

- ▶ Ensure that an installed anti-tipper is activated before travelling over obstacles and on inclines.
- ▶ The anti-tipper has to engage audibly prior to use. Proper engagement must be checked by the user or an accompanying person.
- ▶ For transfemoral amputees, the use of an anti-tipper is highly recommended.

⚠ WARNING

Incorrectly adjusted anti-tipper

Tipping over, falling of the user due to errors in handling the product

- ▶ If only 1 accompanying person helps in overcoming steps, the accompanying person must first adjust the anti-tipper so that it cannot collide with the steps during transport.
- ▶ The accompanying person must adjust the anti-tipper again correctly after overcoming the steps.

⚠ WARNING

Incorrectly adjusted anti-tipper

Risk of falling as a result of incorrectly adjusted anti-tipper.

- ▶ The anti-tipper may only be adjusted by qualified personnel.

The anti-tipper prevents the wheelchair from tipping backwards when overcoming obstacles and going uphill. It is set for a maximum ground clearance of **50 mm** and so that the anti-tipper rollers, as a minimum, project fully beyond the largest diameter of the rear wheels.

The tip-assist makes overcoming obstacles easier for the attendant.

6.11.1 Activating and deactivating the anti-tipper

A mounted anti-tipper must always be active.

When overcoming obstacles (such as steps and curbs that are not lowered) **upwards**, the anti-tipper rollers contact the ground.



When overcoming obstacles (such as steps and curbs that are not lowered) **downwards**, the anti-tipper has to be deactivated by the user or attendant to avoid damaging the anti-tipper.

Activation

- 1) The anti-tipper is pressed down with the hand or foot (see fig. 20, item 1; see fig. 21, item 1).
- 2) Swing the anti-tipper back and allow it to engage (see fig. 20, item 2; see fig. 21, item 2).

Deactivation

- 1) From above, press on the anti-tipper with the hand or foot until it disengages (see fig. 22, item 1; see fig. 23, item 1).
- 2) Rotate the anti-tipper forward by **180°** and release it (see fig. 22, Pos. 1; see fig. 23, Pos. 1).



6.11.2 Using the tip-assist

- 1) At an obstacle, place one foot on the tip-assist and push down.
- 2) The wheelchair can then be easily tipped by pressing down on the push handles at the same time.

6.12 Lap belt

⚠ CAUTION

Incorrectly adjusted lap belt

Incorrect user posture, illness, falling of the user due to installation/adjustment errors

- ▶ Do not modify the settings established by the qualified personnel. In case of problems with the adjustments (unsatisfactory sitting position) please contact the qualified personnel who fitted the product.
- ▶ The lap belt has to fit closely but not too tightly so the user is not injured. It should be possible to slide two fingers comfortably between the belt and thigh.
- ▶ Have the settings of the belt system checked regularly. Adjustments may be required due to the growth of the user or because of changes in the course of the disease or different clothing.

The lap belt prevents the user from slipping and supports positioning.

It is installed on the product by qualified personnel if needed and adapted to the requirements of the user.


Information about subsequent acquisition and mounting is provided by the qualified personnel that handed over the product.

Using the lap belt

- 1) Open the belt using the metal snap / buckle.
- 2) Place the user in an upright, 90° seated position (if physiologically possible). Ensure that the back is up against the backrest padding (if physiologically possible).
- 3) Fasten the belt by means of the metal snap or buckle.
- 4) The lap belt should be at an angle of about **60°** to the seat surface. The belt strap should run across the thighs in front of the pelvic bones (see fig. 25).

Possible errors

- The lap belt is positioned above the pelvis of the user in the area of the soft part of the stomach.
- The user does not sit upright in the seat.
- If the lap belt is too loose, the user can shift/slide out to the front.
- During the installation/adjustment, the safety belt is led over parts of the seating system (e.g. over armrests or seat pads). This causes the lap belt to lose its retaining function.



6.13 Additional options

The product may be equipped with additional options.

The options are firmly mounted to the product by qualified personnel or the manufacturer and are pre-adjusted by qualified personnel at delivery.

6.14 Disassembly and transport

NOTICE

Deformation when folded

Damage to the product, problems unfolding due to unallowable loads

- ▶ Never place heavy objects on the folded product.

The wheelchair must be prepared for transport in a passenger vehicle.

- 1) Apply the wheel lock.
- 2) Remove the seat cushion.
- 3) With the "Folding backrest" option: Pull the release strap for the folding mechanism of the backrest upwards (see fig. 26) and fold the backrest up (see fig. 27).
- 4) Pull the release strap for the folding mechanism (see fig. 28).
- 5) Fold the wheelchair (see fig. 29).
- 6) Fasten the safety strap (see fig. 30).
- 7) Insert four fingers between the spokes of the rear wheel and press the knob of the quick-release axle with your thumb. Pull off the rear wheel outwards (see fig. 3).

INFORMATION: For the "lock-release assist" option, see the section "Rear wheels".



6.15 Use in vehicles for transporting persons with reduced mobility

⚠ WARNING

Improper use in vehicles for transporting persons with reduced mobility

Risk of serious injury when using the product as a seat

- ▶ The product has **not** been approved by the manufacturer for use as a seat in vehicles for transporting persons with reduced mobility. Please refer to the prohibition symbol on the nameplate (see below).
- ▶ In vehicles for transporting persons with reduced mobility, only use the seats installed in the vehicle with the corresponding personal restraint systems.
- ▶ Information on the current status of our measures is available from your specialist dealer.

6.16 Care

6.16.1 Cleaning

- 1) Clean the padding and upholstery with warm water and a mild detergent.
- 2) Remove any spots with a sponge or a soft brush.
- 3) Rinse with clear water and let the components dry.

Important information on cleaning

- Do not use any aggressive cleaners, solvents or hard brushes etc.
- Clean plastic parts, frame parts and the chassis and wheels with a mild cleaner and a damp cloth. Dry thoroughly afterwards.

6.16.1.1 Cleaning Seat Cushions

Standard seat cushion

- ▶ Clean the seat cushion according to the attached care label.

Seat cushion with foam core

- 1) Open the zipper and remove the foam core (see fig. 31).
- 2) Clean the cover according to the attached care label.
- 3) Wash all foam parts by hand at **40°C** using a mild, environmentally friendly detergent. Allow to air dry.
- 4) Reinsert the foam core with proper alignment. Close the zipper.



6.16.1.2 Cleaning belts

Cleaning a belt system with metal buckle

INFORMATION

Observe the washing recommendations on the product and the information in the corresponding instructions for use provided for the product.

- Belts with metal snaps may **not be washed in washing machines** as the penetration of water could cause corrosion and subsequent malfunctions.
- Clean the belt straps by gently dabbing them with warm soapy water (with some disinfectant) or carefully wiped with a dry, clean, absorbent cloth.

Cleaning a belt system with plastic buckle

- Depending on the model, straps with plastic buckles can be washed in the washing machine between **40 °C** and **60 °C**.
- Recommendation: Use a laundry bag or net and mild detergent.

INFORMATION

Alternatively, the belt straps can be cleaned by gently dabbing them with warm soapy water (with some disinfectant) or carefully wiped with a dry, clean, absorbent cloth.

Additional cleaning instructions

- Allow the belts to air dry. Ensure that the belts and pads are completely dry before installation.
- Do not expose the belts to direct heat (e.g. sunshine, stove or radiator).
- Do not iron or bleach the belts.

6.16.2 Disinfection

- Thoroughly clean the pads and handles before disinfecting.
- Wipe all parts of the wheelchair with a disinfectant.

Important information about disinfecting

- If the product is used by more than one person, using a conventional disinfectant is required.
- Only use colourless water-based disinfectants. Observe the instructions for use provided by the manufacturer.

7 Maintenance and repair

7.1 Maintenance

- The function of the product should be checked before each use.
- The product should not be used if defects are noted. This applies in particular in case of instability of the product or altered driving characteristics as well as problems with the user's seating position or the stability of the seat. Inform the qualified personnel promptly for the rectification of defects.
- This also applies if loose, worn, bent or damaged components, cracks in the frame or broken frame components are identified.
- Some maintenance tasks can be carried out to a certain extent by the user at home (see the sections "Maintenance intervals" and "Maintenance tasks").
- The manufacturer also recommends regular maintenance **every 12 months** by authorised, qualified personnel.
- Failure to maintain the product can lead to serious or life-threatening injury to the user of the product.
- Service and repairs may only be carried out by authorised, qualified personnel or the manufacturer. This will ensure that only Ottobock spare parts are used for repairs.

7.1.1 Maintenance intervals

The functions described below must be checked by the user or an attendant at the specified intervals:

Inspection task	Before each use	Monthly	Quarterly
Functional test of the wheel locks	X		
Sagging of the seat or back upholstery		X	
Stability of the legrests		X	
Visual inspection of wear parts (e.g. tyres, bearings)		X	
Soiling of bearings		X	
Damage to the push ring		X	
Air pressure (information see tyre sidewall)		X	
Wear and tear of the folding device		X	
Check of the spoke tension on the rear wheels			X
Check of all bolted connections			X

7.1.2 Maintenance tasks

To ensure smooth operation at all times, users or attendants with some technical skills can maintain some parts of the product:

- Bolted connections must be periodically checked for tightness, especially during the initial period of use or after adjustments have been made to the wheelchair. If a screw connection loosens repeatedly, contact the specialist dealer promptly.
- Hair and dirt particles generally accumulate between the caster wheel and fork. This can restrict the caster wheels from rotating smoothly. For this purpose, lubricate the caster axle between the caster wheel and caster fork with a few drops of thin, resin-free oil (sewing machine oil).
CAUTION: Do not remove the caster wheel yourself. If the caster wheel continues to not rotate smoothly, please contact the qualified personnel.
- The rear wheels are equipped with a quick-release axle system as standard equipment. To keep this system operational, ensure that no dirt adheres to the quick-release axle or receiver bushing. Periodically lubricate the quick-release axle very lightly with thin, resin-free oil (sewing machine oil).
- If the wheelchair becomes wet, it should be rubbed dry again.

7.2 Repair

7.2.1 Inner tube, rim tape and tyre replacement

⚠ CAUTION

Improper tyre replacement

Injuries to the user due to incorrect installation, product damage

- ▶ No person is permitted to sit in the wheelchair during tyre replacement
- ▶ Before removing a wheel, support the product so it cannot tip over.
- ▶ Always replace the tyres in pairs. Two differently worn tyres affect the directional stability of the wheelchair.

ℹ INFORMATION

When driving outdoors, always carry a repair kit and tyre pump (when using pneumatic tyres) in case of emergency.

Suitable tyre pumps are listed on the order form and are supplied with the product. An alternative is tyre foam, which fills your tyre and then hardens (available from bicycle shops, etc.).

Repairing flat tyres requires only the necessary tools and users may change tyres themselves if they wish:

Removal and preparing for installation

- 1) Carefully remove the tyre from the rim using appropriate tools.
INFORMATION: Take care not to damage the rim or the inner tube.
- 2) Unscrew the valve nut from the valve and remove the tube.
- 3) Repair the tube according to the directions in the repair kit or replace it with a new tube.
- 4) Before fitting the tyre again, inspect the rim bed and tyre inner wall for foreign objects. This could have caused the puncture.
- 5) Before installing the tube, check that the rim band is in proper condition. The rim band protects the tube from being damaged by the ends of the spokes.



Replacing the rim band (only when necessary)

- 1) If the rim band needs to be replaced, remove it from the rim.
- 2) Install the new rim band on the inside of the rim, making sure the valve opening is in the right position.

- 3) Glue the rim band in place if this is intended. Ensure that all spoke ends are covered.

Installing the tube and tyre

- 1) Behind the valve, push one side of the tyre over the edge of the rim.
- 2) Slightly inflate the tube until it starts to assume its round shape.
- 3) Unscrew the valve nut from the tube and push the valve through the valve opening in the rim.
- 4) Insert the tube into the tyre.
- 5) Mount the other side of the tyre on the rim, starting from the position across the valve. Ensure that the tube is not pinched between the tyre and rim during this process.



Inflating the tube

- 1) Ensure that the valve is positioned perpendicularly for proper positioning of the tube and tyre in the region of the valve.
- 2) Firmly screw on the valve nut.
- 3) Inflate the tube so that the tyre can still be pressed in easily with your thumb.
INFORMATION: If the circumferential lines on the two sides of the tyre are both at an even distance from the rim, the tyre is centred. If not, let some air out and realign the tyre.
- 4) Inflate the tube to the maximum pressure specified by the tyre manufacturer (see information printed on the tyre sidewall).
- 5) Firmly screw the valve cap onto the valve.

8 Disposal

8.1 Disposal Information

Return the product to the specialist dealer for disposal.

All components of the product must be disposed of properly in accordance with the respective national environmental regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

9.3 Warranty

Further information on the warranty terms and conditions can be obtained from the qualified personnel that has fitted this product or the manufacturer's service (see inside back cover for addresses).

9.4 Service Life

Expected service life: **4 years.**

The design, manufacturing and requirements for the intended use of the product are based on the expected service life. These also include the requirements for maintenance, ensuring effectiveness and the safety of the product.

Using the product beyond the specified expected service life leads to increased residual risk and should only take place subject to the due diligence and deliberations of qualified personnel.

If the service life is reached, the user or a responsible attendant should contact the qualified personnel who fitted the product or the manufacturer's servicing department (see inside rear cover or back page for address). Here the user can obtain information about known risks and the current options for refurbishing the product.

9.5 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

10 Technical data

INFORMATION

- ▶ Much of the technical data below is given in mm. Please note that product settings – unless otherwise specified – cannot be adjusted in the mm range but only in increments of approx. **0.5 cm** or **1 cm**.
- ▶ Note that the values achieved during adjustment may deviate from the values specified below. The deviation can be **±10 mm** and **±2°**.

Zenit – general information

	Zenit
Maximum load capacity [kg]	120 ¹⁾
Weight [kg]	from approx. 8.9 ²⁾
Transport weights [kg] (for seat width of 440 mm)	Frame: 6.5 ²⁾ Rear wheel 24": 1.3 ²⁾
Seat width [mm]	320 – 460
Anterior seat height [mm]	380 – 540
Rear seat height [mm]	370 – 500
Seat depth [mm]	360 – 500
Overall length [mm]	640 – 1040
Overall width [mm]	495 – 710
Overall height	620 – 1000
Folding size	330
Back height	250 – 500
Lower leg length	280 – 510
Front frame angle	75°/85°
Min. tyre pressure [bar]	7 ³⁾

¹⁾ Depending on selected option: 100 kg / 120 kg

²⁾ The specified weights vary according to the selected options and model.

³⁾ Varies according to tyre option; see the print on the tyre wall

11 Appendices

11.1 Threshold values for wheelchairs transportable by train

INFORMATION

- ▶ The products in this series fully satisfy the minimum technical requirements of regulation (EU) No. 1300/2014 regarding train accessibility for people with disabilities. However, not all versions can comply with all threshold values due to different settings.
- ▶ With the help of the table that follows, you or the qualified personnel can take measurements and verify whether the specific product in question meets the threshold values.

Feature	Threshold value (according to Directive (EU) No. 1300/2014)
Length [mm]	1200 (plus 50 mm for the feet)
Width [mm]	700 (plus 50 mm on each side for the hands when moving)
Smallest wheels ["]	approx. 3 or greater (according to the directive, the smallest wheel must be able to overcome a gap measuring 75 mm horizontally and 50 mm vertically)
Height [mm]	max. 1375; including a 1.84 m large male wheelchair user (95th percentile)
Turning radius [mm]	1500
Maximum weight [kg]	200 (for wheelchair and user, including any baggage)
Maximum obstacle height that can be overcome [mm]	50
Ground clearance [mm]	60 (at an upward slope angle of 10°, ground clearance must measure at least 60 mm under the footrest for going forward at the end of the slope)
Maximum inclination angle on which the wheelchair will remain stable [°]	6 (dynamic stability in all directions) 9 (static stability in all directions, also when wheel lock engaged)



A series of horizontal lines spanning the width of the page, providing a template for writing. The lines are evenly spaced and extend from the left margin to the right margin.



A series of horizontal lines spanning the width of the page, providing a template for writing.

Ihr Fachhändler | Your specialist dealer



Otto Bock Mobility Solutions GmbH
Lindenstraße 13 · 07426 Königsee-Rottenbach/Germany
www.ottobock.com

