

Controller
Boîtier de commande
Bedieningskast

000690820.02

VR2

PG Drives Technology

User manual
Notice d'utilisation
Gebruikershandleiding

English

© 2016 Sunrise Medical

All rights reserved.

The information provided herein may not be reproduced and/or published in any form, by print, photoprint, microfilm or any other means whatsoever (electronically or mechanically) without the prior written authorization of Sunrise Medical.

The information provided is based on general data concerning the construction known at the time of the publication of this manual. Sunrise Medical executes a policy of continuous improvement and reserves the right to changes and modifications.

The information provided is valid for the product in its standard version. Sunrise Medical can therefore not be held liable for any damage resulting from specifications of the product deviating from the standard configuration. Illustrations contained in the manual may deviate from the configuration of your product.

The information made available has been prepared with all possible diligence, but Sunrise Medical cannot be held liable for any errors contained in the information or the consequences thereof.

Sunrise Medical accepts no liability for loss resulting from work executed by third parties.

Names, trade names, trademarks etc. used by Sunrise Medical may not, as per the legislation concerning the protection of trade names, be considered as being available.

Contents

1. Introduction	5
1.1 This user manual	5
1.2 For further information	5
1.3 Symbols used in this manual	5
2. The controller	6
2.1 VR2 controller	6
3. Driving the wheelchair with the controller	8
3.1 Switching the controller on or off	8
3.2 Driving the wheelchair	8
3.3 Speed	8
4. Operating the electrical adjustment options	9
4.1 Selecting the desired adjustment option	9
5. Lights	10
6. Troubleshooting	10
7. Locking the controller	12
8. Technical specifications	13
9. Technical diagrams	14
9.1 Technical diagram	14
9.2 Technical diagram battery charger	15

1. Introduction

1.1 This user manual

This user manual will help you to use and maintain the controller of your power wheelchair safely. This user manual is a supplement to Sunrise Medical's general wheelchair user manual. When necessary this user manual refers to other manuals as shown below:

-  Wheelchair: Refers to the general wheelchair user manual.
-  Battery charger: Refers to the user manual for the battery charger.

Read this user manual and the other user manuals referred to carefully before using the product. If one of the user manuals was not included with your wheelchair, please contact your dealer immediately. In addition to this user manual, there is also a service manual for qualified specialists.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu or alternatively is available on request in large text.

1.2 For further information

Please contact your local, authorised Sunrise Medical dealer if you have any questions regarding the use, maintenance or safety of your wheelchair. In case there is no authorised dealer in your area or you have any questions, contact Sunrise Medical either in writing or by telephone.

For information about product safety notices and product recalls, go to www.SunriseMedical.eu

1.3 Symbols used in this manual

Note!

Pointing out possible problems to the user.

Caution!

Advice for the user to prevent damage to the product.

Warning!

Warnings for the user to prevent personal injury.

Not following these instructions may result in physical injury, damage to the product or damage to the environment!

2. The controller

2.1 VR2 controller

A controller will usually have three basic functions:

- Driving and steering a wheelchair
- Operating electrical seat adjustments
- Charging the wheelchairs batteries

There are many different control systems for wheelchairs on the market. If the controller on your wheelchair does not resemble the one in figure 1, 2 or 3 contact your dealer.

1. VR2 drive only
2. VR2 Lights
3. VR2 Seats&Lights

VR2 by PG Driving Technologies is a collective name for the entire control system of the wheelchair. The wheelchair is operated by means of a controller, which includes the following components:


Figure 1

Figure 2

Figure 3


Figure 4 From the left to the right ; VR2 drive only, VR2 Lights, VR2 Seats & Lights

Part	Function
A. Joystick	In the 'drive' mode: driving and steering In the 'adjustment options' mode: <ul style="list-style-type: none"> • Left/right to select the adjustment options • Front/back to select the adjustment mode
B. On/off switch	Switching the controller on or off
C. Horn	Warning signal with sound
D.1 Battery indicator	Displays the power level of the battery
D.2 Maximum speed	Displays the maximum speed limit as set by the user
D.3 Adjustment option	Displays the selected adjustment options
E. "Mode" button	Changing between the 'driving' and the 'adjustment options' mode
F. 1 Speed regulator	Reduce driving speed (slower)
F. 2 Speed regulator	Increase driving speed (faster)
G. Charge connector	Input for the battery charger
H. Charge plug for the battery charger	Connector for the battery charger
I. Lights button	Switching the lights on or off
J. Hazard lights	Warning signal with lights
K.1 Direction indicator left	Switches the left direction indicator on or off
K.2 Direction indicator right	Switches the right direction indicator on or off

3. Driving the wheelchair with the controller

3.1 Switching the controller on or off

To be able to drive or operate the electronic adjustment options of the wheelchair, the controller must be switched on. Press the on/off button (B in figure 4).

3.2 Driving the wheelchair

Driving an electric wheelchair is done by operating a joystick. Move the joystick forwards and the wheelchair will also move forwards. Steer left and right and the wheelchair will turn.

3.3 Speed

The maximum speed can be controlled by the speed regulator on the controller (F1 en F2 in figure 4). The maximum speed is displayed by the LED lights above the speed regulators (D2 figure 4). Speed can be controlled with the joystick while driving. If the joystick is moved a little, the wheelchair will move more slowly.

4. Operating the electrical adjustment options

Not every wheelchair has been equipped with electronic adjustment options. We make a distinction between four different adjustment options to the seating system:

- Tilt adjustments
- High/low adjustments
- Backrest adjustments
- Legrest adjustments

To use electrical adjustment options, the VR2 Seats&Lights controller is required (maximum of two options).


Figure 5 Electrical adjustment options

4.1 Selecting the desired adjustment option

1. Switch on the controller
2. Press on the 'mode' button (E in figure 4) to select the 'adjustment options' mode. The controller will now be in the 'adjustment options' mode. It is also possible to follow the status.

The joystick is used to select and operate the electronic adjustment options.

3. Move the joystick to the left or the right to select the desired adjustment option. The selected adjustment option will be visible on the display screen.
4. Moving the joystick forwards and/or backwards will activate the selected electronic adjustment option (see table 6). Move the joystick forwards or backwards until the desired adjustment option has been attained.

	Move joystick backwards	Move joystick forwards
Tilt adjustment	The entire chair will tilt backwards	The entire chair will tilt forwards
High/low adjustments	The entire chair will be raised	The entire chair will be lowered
Backrest adjustments	The backrest will tilt backwards	The backrest will tilt forwards
Legrest adjustments	The legrest angle will increase, the footplate will be raised	The legrest angle will decrease, the footplate will be lowered

Table 1 Electrical adjustments

5. To return to the 'drive' mode: Press on the 'mode' button to select the 'drive' mode.

Note: If you adjust the seat height by using the high/low option or use the electrical tilt adjustment 0 - 45°, the speed will be reduced due to safety reasons.

5. Lights

Not every wheelchair has been equipped with lights. You can activate the following functions only by using VR2 Lights and/or VR2 Seat&Lights:

- Lights (I in figure 4)
- Hazard Lights (J in figure 4)
- Direction indicators (K1 & K2 in figure 4)

6. Troubleshooting (Fig. 6)


If the wheelchair will not function while the batteries are fully charged, check the following points before consulting your dealer:

- Switch the controller off and then switch it on again. Check to see if the malfunction has been solved.
- Check if the free wheel switch was switched to Drive.
- Check if the joystick was in the 0 position when the controller was switched on. In other words, the joystick must not be moved when the controller is being switched on or off.

Malfunctions list

An extended malfunction list can be found in the appendix of the service manual (for qualified specialists only). The service manual can be found on www.SunriseMedical.eu

Fig. 6 On-board Diagnostics:

	An excessive voltage has been applied to the control system. This is usually caused by a poor battery connection. Check the battery connections
	The parking brakes have a bad connection. Check the parking brake and the motor connections. Make sure that the control system connections are secure
	A control system fault is indicated. Make sure the control system connections are secure.
	A joystick fault is indicated. Make sure that the joystick is in the centre position before switching on the control system
	The wheelchair is being prevented from driving by an external signal. One possibility is that the battery charger is plugged in.
	The right hand motor has a short circuit to a battery connection. Contact your Sunrise Medical Authorised Dealer
	The right hand motor has a bad connection. Check the connections to the right hand motor.
	The left hand motor has a short circuit to a battery connection. Contact your Sunrise Medical Authorised Dealer.
	The left hand motor has a bad connection. Check the connections to the left hand motor.
	The batteries need charging, or there is a bad connection to the battery. Check the connections to the battery. If the connections are good, try charging the batteries
	S=Speed indicator LED's A communication fault is indicated. Make sure that the joystick cable is securely connected and not damaged
	A=Actuator LED's An actuator trip is indicated. If more than one actuator is fitted, check which actuator is not working. Check the actuator wiring.

7. Locking the controller

To lock the controller;

- While the control system is switched on, depress and hold the on/off button.
- After 1 second the control system will beep. Now release the on/off button
- Deflect the joystick forwards until the control system beeps.
- Deflect the joystick in reverse until the control system beeps.
- Release the joystick, there will be a long beep.
- The wheelchair is now locked.

To unlock the controller;

- If the control system has switched off, press the on/off button.
- Deflect the joystick forwards until the control system beeps.
- Deflect the joystick in reverse until the control system beeps.
- Release the joystick, there will be a long beep.
- The wheelchair is now unlocked.

8. Technical specifications

Supply Voltage:	24 Vdc
Operating Voltage:	16Vdc to 35Vdc
Peak Voltage:	35Vdc
Reverse Battery Voltage:	40Vdc
PWM Frequency:	20kHz \pm 0.5%
Brake Voltage:	12Vdc or 24Vdc
Brake Current:	100mA min. 1A max.
Charger Connector:	Use only Neutrik NC3MX
Batt. Charging Current:	12Arms max.
Actuator Current:	12A max.
Maximum Drive Current:	VR2 50 50A VR2 60 60A VR2 70 70A VR2 90 90A
Moisture Resistance:	IPx4
Operating Temperature:	-25°C to +50°C
Storage Temperature:	-40°C to +65°C

EMC tested on sample wheelchair:

Susceptibility:	Tested at 30V/m to EN12184 (1999) and ANSI/ RESNA requirements
Emissions:	To EN55022 Class B
ESD:	IEC801 part 2

9. Technical diagrams

9.1 Technical diagram

The technical diagrams can also be found on the cover of the specific electronic component.


Figure 7


Figure 8


Figure 9

9.2 Technical diagram battery charger

The controller's standard configuration includes a '3-pin connection'. Ensure that the battery charger is properly connected so that the 'negative pole' and the 'inhibit' are connected, enabling the system to prevent the wheelchair from moving when the battery is being charged.


Figure 10

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