ottobock.

**M1 Intro · M2 Effect · M3 Hemi · M4 XXL ·
M5 Comfort · M6 Junior · Asia**

DE Gebrauchsanweisung	3
EN Instructions for Use	43

DE

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EN**Instructions for Use for the Start Series**

Contents	Page
1 General Information	45
1.1 Preface.....	45
1.2 Intended Use	45
1.3 Field of Application.....	45
1.4 Liability	46
1.5 CE Conformity	46
1.6 Trademarks.....	46
1.7 Technical Service	46
2 Safety Instructions	47
2.1 Explanation of Symbols.....	47
2.2 General Safety Instructions	47
2.3 Warning Symbols and Type Plates.....	51
2.4 Further information	51
3 Delivery and Preparation for Use	52
3.1 Storage.....	52
3.2 Preparation for Use.....	52
4 Transportation	53
4.1 Transfer.....	54
4.2 Use in a wheelchair accessible vehicle	56
5 Options	57
5.1 Anti-tipper.....	57
5.2 Swinging Anti-tipper.....	58
5.3 Transport Wheels	60
5.4 Crutch Holder with Hook and Loop Strap.....	60
5.5 Plug-on Side Panels.....	61
5.6 Desk Side Panels (standard/stepped/height-adjustable)	61
5.7 Plug-on Wheel Lock Lever Extension	61
5.8 Rear Wheels with Drum Brake.....	62
5.9 Wheel Lock for One-handed Operation	62
5.10 Tray	62
5.11 Lap Belt.....	63
5.12 Elevating Footrest.....	63
5.13 Push Handles, Height-adjustable.....	63
5.14 Wheelbase, Long	63

5.15	<i>Back Upholstery, Adaptable</i>	64
5.16	<i>Folding Back</i>	64
5.17	<i>Angle-adjustable Back</i>	64
5.18	<i>Spoke Protector</i>	64
5.19	<i>Amputation Legrest</i>	64
5.20	<i>“The Back”</i>	65
5.21	<i>“The Seat”</i>	65
6	Adjustment / Assembly Instructions	65
6.1	<i>Adjusting the Seat Height and Seat Inclination</i>	65
6.2	<i>Changing the Wheelbase</i>	66
6.3	<i>Setting the Wheel Camber – Only Start M5 Comfort</i>	67
6.4	<i>Caster Trailing Adjustment</i>	69
6.5	<i>Adjusting the Wheel Lock</i>	69
6.6	<i>Adjusting the Brake Force – Drum Brake</i>	70
6.7	<i>Adjusting the Footplate Angle</i>	70
6.8	<i>Adjusting the Lower Leg Length</i>	71
6.9	<i>Adjusting the Angle of the Elevating Footrest</i>	72
6.10	<i>Adjusting the Height of the Armrest</i>	72
6.11	<i>Adjusting the Seat Depth</i>	73
6.12	<i>Adjusting the Back Height</i>	73
6.13	<i>Seat Upholstery</i>	74
6.14	<i>Adapting the Back Upholstery</i>	74
6.15	<i>Changing the Back Angle</i>	74
6.16	<i>Push Handles, Height-adjustable</i>	75
7	Maintenance, Cleaning and Disinfection	76
7.1	<i>Maintenance</i>	76
7.2	<i>Cleaning and Disinfection</i>	77
8	Changing Tyres	77
9	Technical Data	79
10	Re-use / Disposal	82
10.1	<i>Information on Re-Use</i>	82
10.2	<i>Information on Disposal</i>	82

1 General Information

INFORMATION

Date of the last update: 2014-03-20

- Please read this document carefully.
- Follow the safety instructions.

1.1 Preface

With the purchase of a lightweight wheelchair of the Start series, you have selected a quality product for versatile use during everyday life, both indoors and outdoors. Before you start using your wheelchair, please make sure you read the chapter “Delivery / Preparing the Wheelchair for Use” and the safety instructions.

The chapter “Accessories” presents various additional components of the Start wheelchair, which can expand its field of application and improve comfort for the wheelchair user. The chapter “Adjustment/Assembly Instructions” provides an overview of the ways the wheelchair can be customised to your requirements. Technical changes to the design described in these instructions for use are reserved.

1.2 Intended Use

The Start lightweight wheelchair is designed solely for individual use by persons who are unable to walk or who have a walking impediment, and can be operated either by the patient or by another person.

The Start lightweight wheelchair may only be combined with the options mentioned in these Instructions for Use and vice versa. Ottobock assumes no liability for combinations with medical devices and/or accessories from other manufacturers outside of the modular system.

1.3 Field of Application

The versatility of the wheelchair and its modular design make it suitable for patients who have walking impediments or walking disabilities, for example due to:

- Paralysis (paraplegia / tetraplegia)
- Loss of limb (leg amputation)
- Limb defect / deformation
- Joint contractures / defects
- Muscle and nerve disorders
- Diseases such as cardiac or circulatory insufficiency, balance disorders or cachexia as well as for geriatric patients who still have usable residual strength in the upper limbs.

The Start lightweight wheelchair was specially designed for individuals who are normally able to move in a wheelchair independently.

The following points should also be considered for an individual fitting

- Body height and weight (maximum load capacity: 125 kg / M4 XXL: 160 kg / M6 Junior: 90 kg)
- Physical and mental condition
- Age of the user
- Home and living conditions
- Environment

1.4 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

The manufacturer is not responsible for damages caused by components and spare parts not approved by the manufacturer. Repairs must be carried out exclusively by authorised dealers or by the manufacturer.

1.5 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This product has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

1.6 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

1.7 Technical Service

Service and repairs on the Ottobock product may only be carried out by an authorised dealer. Should any problems arise, please contact the supplier of your product. Any necessary repairs will be made exclusively with authentic Ottobock spare parts there. Your product requires regular maintenance / regular service (see 647G345=GB Service Instructions).

Your authorised Ottobock dealer:

2 Safety Instructions

2.1 Explanation of Symbols

 WARNING	Warnings regarding possible risks of severe accident or injury.
 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.
 INFORMATION	Information regarding operation. Information for service personnel.

2.2 General Safety Instructions

 INFORMATION
<p>Please read the instructions for use first! Before using the product, familiarise yourself with the handling, functioning, and use of the product. You may compromise your safety if you do not observe these instructions.</p> <p>However, it is impossible to protect against and avoid all unforeseeable circumstances and situations.</p>

Danger when overcoming stairs/obstacles

 WARNING
<p>Risk of falling in case of lacking assistance. Do not ascend or descend stairs without the assistance of attendants. If devices such as ramps or elevators are available, they should be used. If wheelchair-friendly access is missing, two attendants must carry the wheelchair over the obstacle.</p>
 WARNING
<p>Risk of falling as a result of incorrectly set anti-tipper. If only one attendant is available when ascending or descending stairs / undropped kerbs, an incorrectly set anti-tipper can lead to severe falls. If an anti-tipper is mounted, adjust the anti-tipper to the disengaged position, so that it does not interfere with the stairs / undropped kerbs during ascending or descending. Afterwards, bring the anti-tipper back to its operational position.</p>
 WARNING
<p>Risk of falling as a result of lifting the wheelchair incorrectly. Attendants must lift the wheelchair only by parts that are firmly attached (e.g. not at the footrest or rear wheels). If your wheelchair is equipped with height-adjustable push handles, ensure that the clamp levers are always firmly tightened.</p>
 WARNING
<p>Risk of tipping when using lifting platforms. Ensure that the anti-tipper (if mounted) is outside the danger-area.</p>

Danger when getting into the wheelchair

 CAUTION
<p>Risk of injury due to incorrect boarding.</p> <ul style="list-style-type: none"> • Activate the wheel lock each time before you get into, out of, or transfer to or from the wheelchair. • If possible, the footplates should be flipped up before getting into or out of the wheelchair. • Always get into the wheelchair from the side, where possible.

Danger while driving**⚠ CAUTION**

Risk of accident when driving without experience. This can lead to falls and other dangerous situations. Therefore, familiarise yourself with the new wheelchair by practising on even, straightforward terrain first. Learn how the wheelchair reacts when the centre of gravity shifts; for example on slopes or inclines or when clearing obstacles such as steps and curbs. This should be done only with assistance from another person. Using an anti-tipper is strongly recommended for inexperienced wheelchair users.

⚠ CAUTION

Risk of tipping when ascending too steep slopes. When ascending slopes or ramps and when crossing obstacles on upward slopes, always lean your upper body far forward.

⚠ CAUTION

Risk of tilting or tipping over as a result of centre of gravity shifts.

- When reaching for objects in front, to the side or behind the wheelchair, do not lean out of the wheelchair too far.
- Prior to activities that require you to bend forward in the wheelchair (e.g. tying your shoes), maximise the stability of the wheelchair. In order to do so, push the wheelchair backwards until the steering casters turn forward.

⚠ CAUTION

Risk of tipping when driving against obstacles (steps, curbs) without braking. Use your wheelchair properly. Do not “jump” the wheelchair down from higher surfaces.

⚠ CAUTION

Risk of accident as a result of incorrect use of the wheel lock. All brakes acting on the tyres do not serve as service brake but are only designed as parking brake (wheel lock). The wheel locks must not be used as driving brakes for slowing down the wheelchair, as in extreme cases, the abrupt stopping of the wheelchair can lead to falls.

⚠ CAUTION

Risk of accident if the wheelchair starts rolling. Engage the wheel lock to prevent your wheelchair from moving on uneven ground or during transfers (e.g. into a car).

⚠ CAUTION

Risk of accident as a result of driving in the dark. In the dark, wear light clothing or clothing with reflectors in order to improve your visibility. Ensure that the reflectors installed on the sides and rear of the wheelchair are easily visible. We also recommend installing active illumination.

⚠ CAUTION

Risk of tipping when descending slopes. When descending slopes, do not drive without braking and reduce your speed.

⚠ CAUTION

Risk of falling due to lack of tipping resistance on public transport

- Always observe currently applicable legal requirements when using public transport.
- Always ensure that you are held in place securely when travelling on public transport. To do so, use the wheelchair areas, wheelchair bays and restraint systems provided. Firmly tighten the wheel locks.

- Please note that the anti-tipper, if activated, may be exposed to heavy loads when the public transport vehicle starts to move. In order to avoid damage, the manufacturer recommends that the wheelchair is positioned at a right angle to the direction of travel if no restraint system for passengers with reduced mobility is available.

Danger as a result of incorrect assembly or settings

⚠ CAUTION

Risk of tipping as a result of incorrectly set wheels. Extreme settings (e.g. rear wheels mounted in the foremost position) combined with an unfavourable body posture can cause the wheelchair to tip even on level ground. For transfemoral amputees it is essential to shift the rear wheels backwards. Use of an anti-tipper is required.

⚠ CAUTION

Risk of tipping as a result of missing or incorrectly mounted anti-tipper. Under no conditions should the anti-tipper assume the function of transport wheels, for example to transport a person in the wheelchair with the rear wheels removed. The anti-tipper must audibly lock in place, before it is able to bear loads. Firm seating must be verified by the user or by an attendant.

⚠ CAUTION

Risk of tipping as a result of changed diameters/mounting positions of the wheels. If the size and position of the casters or the rear wheel size are modified, the casters may flutter when driving at higher speeds. This can lead to blocking of the casters and cause the wheelchair to tip over. For this reason, you should always maintain the factory settings or, if modifications are required, ensure that the wheelchair frame is aligned horizontally (see "Adjustment / Assembly Instructions").

⚠ CAUTION

Risk of falling as a result of unintentional loosening of the rear wheels. Always make sure that the quick-release axles are correctly set on the rear wheel. It must not be possible to remove the rear wheel unless the button on the quick-release axle is pushed.

Danger as a result of improper use of the tyres

⚠ CAUTION

Risk of accident as a result of insufficient or excessive tyre pressure. The effectiveness of the wheel lock and the overall driving quality are dependent on adequate air pressure. The wheel locks are only effective with sufficient air pressure and proper adjustment (approx. **4 mm** distance, subject to technical changes).

Before starting to use your wheelchair, check that the tyres are inflated correctly. The required air pressure is printed on the side of the tyre. For rear wheels, it should be **at least 3.5 bar**. With properly inflated rear wheels and even tyre pressure on both sides, your wheelchair is much easier to operate and manoeuvre.

⚠ CAUTION

Risk of accidents due to bad tyres. Insufficient tread depth of the tyres reduces adhesion. Please note that when you operate your wheelchair in public road traffic, you are subject to road traffic regulations.

Risk of skin damage

⚠ CAUTION

Possible complications in case of unhealthy skin. Using the wheelchair with unhealthy skin can cause medically indicated complications such as reddening of the skin and pressure sores. Over several hours of use, the skin on the buttocks, back, and rear thighs is subjected to considerable strain. Therefore, always check these areas carefully before using the wheelchair. No liability will be assumed for health injuries caused by use of this wheelchair with damaged skin.

Danger caused by fire / heat and cold

⚠ CAUTION

Risk of burns in the proximity of fire. The seat and back upholstery might catch fire. Keep the wheelchair away from all sources of ignition, especially lit cigarettes.

⚠ CAUTION

Be careful in case of extreme temperatures. The wheelchair can heat up significantly when in the sun or in the sauna. In extreme cold, there is a risk of hypothermia.

Risk of hand injuries

⚠ CAUTION

Risk of hand injuries. To avoid hand injuries, do not grasp between the rear wheel and the wheel lock when driving your wheelchair.

⚠ CAUTION

Risk of pinching between wheelchair components. The wheel lock lever and the side panel or frame part can represent a risk of crushing with various adjustment options.

⚠ CAUTION

Risk of friction burns. Slowing down from high speeds or when descending longer slopes tends to heat up the hands and fingers, especially if using aluminium push rings. When using your wheelchair outdoors, wear leather gloves. Gloves provide you with a better grip and protect your fingers from dirt and hot metal.

Warnings regarding wheelchair damage

NOTICE

Risk of damage caused by overload. The maximum load for modules **M1, M2, M3** and **M5** is **125 kg**. The maximum load for module **M4** is **160 kg** and for module **M6** **90 kg**.

NOTICE

Damage caused by objects lying upon the wheelchair. If appliances are transported in the folded state, the weight of objects lying upon the appliance can cause deformations and result in problems when unfolding the appliance. Never put heavy objects on a folded rehab device.

NOTICE

Damage caused by wear and tear. If the upholstery gets worn or damaged, replace it immediately. With new seat upholstery, you will automatically receive new lamellar plugs for the seat tube which you should use to replace the previously installed old lamellar plugs.

NOTICE

Risk of damage caused by incorrect packaging. Only use original packaging for shipment of the product.

2.3 Warning Symbols and Type Plates

Label	Meaning
	<ul style="list-style-type: none"> A Type designation B Article number of the manufacturer C Maximum load capacity (see section “Technical Data”) D Manufacturer information / address E Serial number F International article number G Read the instructions for use prior to using the product. H CE symbol – product safety according to EU guidelines
	<p>Attention! Read the Instructions for Use prior to using the product.</p>

2.4 Further information

INFORMATION

Even in the event of compliance with all applicable guidelines and standards, alarm systems (e.g. in department stores) may respond to your product. Should this happen, remove your product from the area where the alarm was triggered.

3 Delivery and Preparation for Use

3.1 Storage

The wheelchair must be stored in a dry place.

The ambient temperature must be maintained between -10 °C and +40 °C during transportation and storage.

During extended storage, the knee lever wheel lock on wheelchairs with PU tyres must be released since tyre deformation may otherwise result.

INFORMATION

Tyres contain chemical substances that can react with other chemical substances (such as cleaning agents, acids, etc.).

3.2 Preparation for Use

⚠ CAUTION

Risk of pinching. Only grip the components shown when folding or collapsing the wheelchair.

Your wheelchair is usually delivered fully assembled, folded, and with the rear wheels removed. A few simple steps are sufficient to prepare it for use:

1. Attach the wheels. To attach the rear wheels reach into the spokes with four fingers and press the button of the quick-release axle with your thumb (Fig. 1). Then attach the rear wheels to the fittings and the axles will automatically engage.
2. While standing next to the wheelchair, tip it slightly toward you. With flat hands, press the crossbrace apart. Make sure that the crossbrace has locked in the support brackets (Fig. 2).
3. Pull the back upholstery flap up and attach it to the wheelchair's seat upholstery (hook and loop connection) (Fig. 3). Insert the seat cushion.
4. If applicable, insert the footrests into the receivers (Fig. 4). Finally, fold down the footplates.

⚠ CAUTION

Accident risk. When remounting, be sure the quick-release axle securely locks into the quick-release axle housing! The rear wheel must not be removable when the button is not being pressed!

⚠ CAUTION

Risk of tipping. The rear wheels have been pre-positioned according to your indications. Before using the wheelchair, this pre-adjusted position must be checked by the wheelchair user for function and stability against tipping with the assistance of skilled staff.

4 Transportation

Flip up the footplate and remove the seat cushion.

Grab the seat upholstery at the front and the back in the middle and pull it upwards: the wheelchair will fold (Fig. 5).

Once the wheelchair is folded, secure the strap of the folding fixture with the snap fastener. This will prevent the wheelchair from unfolding unexpectedly (Fig. 6).

For more manageable transportation, e.g. in cars, the footrests can be swung outwards or be entirely removed in just a few easy steps.

To swing away or remove the footrest, pull the handle to the rear (Fig. 7). The footrest with footplate can then be swung away to the inside or outside by 90° (Fig. 8), or be swung away and removed by pulling upward.

INFORMATION

During reassembly, ensure that the footrest locking device has engaged.

Transportation is greatly facilitated by the removable quick-release axles on the rear wheels (included in the standard version).

To remove the rear wheels, press the quick-release axle button, as described in section 3. The rear wheels can then be removed. (Fig. 1).

⚠ CAUTION

Accident risk. During assembly, confirm that all bolts engage in the corresponding bore holes and that the quick-release axles of the rear wheels and casters securely lock into their axle housings. The rear wheels and casters must not be removable without the button being depressed! Check the function of the wheel locks.

4.1 Transfer

Transfer into the wheelchair is described below. The lightweight wheelchair is designed for persons who normally move around in the wheelchair independently. Therefore, the independent transfer is described. Should you require an attendant, the steps described here are carried out correspondingly but with support.

Manoeuvre the wheelchairs so that they are positioned at an angle of approx. 45 ° degrees to each other (Fig. 9). Ensure that the wheel locks do not get caught in each other but remain freely movable. First activate the wheel locks of the wheelchair that is currently in use (Fig. 10).

Then activate the wheel locks of the wheelchair you want to transfer into (Fig. 11). Move your feet from the footrest to the floor and move ahead in the wheelchair if required (Fig. 12).

Flip up the footplates of the wheelchair (Fig. 13). Transfer into the second wheelchair. The method illustrated is one possible way to support yourself (Fig. 14). Find the most suitable path for you by practising with an attendant.

Fold the footplates of the wheelchair down again. Position your feet on the footplates of the wheelchair (Fig. 15). Release the wheel locks (Fig. 16). Now you can use the wheelchair.

4.2 Use in a wheelchair accessible vehicle

⚠ WARNING

Use as a seat for transportation in a wheelchair accessible vehicle

Risk of injury as a result of unauthorised use in a wheelchair accessible vehicle. Whenever possible while travelling in a wheelchair accessible vehicle, passengers should use the seats and vehicle restraint systems already installed in the vehicle. This is the only way to provide optimal protection for the passengers in case of accident. The models M1, M2, M3, M5 and M6 can be used as passenger seats in wheelchair accessible vehicles when used in conjunction with Ottobock safety components and the appropriate restraint systems.

For further information please refer to our instructions for use "Using Your Wheelchair/Seating Shell Mobility Base or Buggy as a Seat for Transportation in a Wheelchair Accessible Vehicle", order number 646D158.

⚠ WARNING

Risk of accident and injury due to improper use as seat during transportation in a wheelchair accessible vehicle.

Use in a wheelchair accessible vehicle is not permissible if the following options are used:

- 30° back angle adjustment
- folding back.

A list of further options not permitted in a wheelchair accessible vehicle can be found in our instructions for use "Using Your Wheelchair/Seating Shell Mobility Base or Buggy as a Seat for Transportation in a Wheelchair Accessible Vehicle", order number 646D158.

5 Options

The wheelchair is designed as a modular component system. This means that certain accessories can be adapted to your wheelchair. We recommend using a seat cushion at all times. The seat upholstery of the wheelchair of the Start series has a loop texture so that seat cushions can be attached with hook and loop connection. All seat cushions listed on the order form have the required hook texture.

Furthermore, we will present you a selection of our variants and accessories which improve the use and comfort of your wheelchair.

5.1 Anti-tipper (Fig. 17)

WARNING

Risk of falling as a result of incorrectly set anti-tipper. The anti-tipper wheel must at least completely extend from the rear wheel to the back. The space between the floor and lower edge of the anti-tipper wheels must not exceed 5 cm.

After making adjustments, the anti-tipper must snap in place in the new position with an audible click.

CAUTION

Risk of tipping. Never drive over obstacles (steps and kerbs) at a slant angle. Always approach obstacles straightly (at an angle of 90°).

The anti-tipper allows the wheelchair user to safely overcome thresholds and dropped kerbs without the help of an attendant. It prevents the wheelchair from tipping too far backward.

The anti-tipper can be adjusted in length without the need for tools, depending upon the seat height and seat angle (Fig. 17, item A).

Mounting of the anti-tipper is strongly recommended for untrained wheelchair users as well as when extreme adjustments have been made to the rear wheel.

5.1.1 Using the Anti-tipper

a) Negotiating thresholds / dropped kerbs

When ascending an obstacle, the anti-tipper wheels hit the ground and prevent the wheelchair from tipping over to the rear.

b) Negotiating steps / undropped kerbs

WARNING

Risk of falling. Please note that an attendant is required to assist with negotiating steps and undropped kerbs.

Before tilting the wheelchair in order to get over undropped kerbs or steps, the attendant must disengage the anti-tipper by rotating it to an upwards position. This prevents the anti-tipper from getting damaged.

To disengage the anti-tipper, press the spring button of the tripod spring and turn up the anti-tipper by 180° (Fig. 17, item B) until it locks in place.

Now the wheelchair can be moved over the obstacle forwards/backwards.

INFORMATION

Please note that it is necessary to re-engage the anti-tipper after getting over the obstacle. To do so, press the spring button of the tripod spring and turn down the anti-tipper by approx. 180° until it locks in place (Fig. 17, item B).

5.2 Swinging Anti-tipper (Fig. 18–22)**⚠ WARNING**

Risk of falling as a result of incorrectly set swinging anti-tipper. The space between the floor and anti-tipper ends must not exceed 5 cm. This basic setting may only be made by qualified personnel.

⚠ CAUTION

Risk of tipping. Never drive over obstacles (steps and kerbs) at a slant angle. Always approach obstacles straightly (at an angle of 90°).



The swinging anti-tipper allows the wheelchair user to overcome thresholds and dropped kerbs without the help of an attendant. It prevents the wheelchair from tipping and after overcoming the obstacle swings back to the rear.

If necessary, adjust the length of the swinging anti-tipper by loosening the marked screw connections (Fig. 18, items A/B). To change the height setting of the anti-tipper, loosen the screw connection (Fig. 18, item B) and use another bore hole.

⚠ WARNING

Risk of falling as a result of incorrectly set anti-tipper. These basic settings may only be made by qualified personnel.

When making these basic settings, make sure that the outer arm of the swinging anti-tipper protrudes beyond the largest diameter of the tyre (Fig. 18, item C).

5.2.1 Using the Swinging Anti-tipper


a) Negotiating thresholds / dropped kerbs

To negotiate thresholds and dropped kerbs, the arms of the anti-tipper should hang freely towards the ground (Fig. 18).

When ascending an obstacle, the arms of the anti-tipper hit the ground and prevent the wheelchair from tipping over to the rear (Fig. 19).

When descending an obstacle, the front arm of the anti-tipper must protrude beyond the edge of the obstacle to enable the anti-tipper to swing out of the way.

When continuing to move the wheelchair (forwards/backwards), the anti-tipper rotates around its own axis (Fig. 20, item A). Once the obstacle has been overcome, the anti-tipper swings back into the start position.

b) Negotiating steps / kerbs

⚠ WARNING

Risk of falling. Please note that an attendant is required to assist with negotiating steps and undropped kerbs.


Before tilting the wheelchair in order to get over undropped kerbs or steps, the attendant must disengage the anti-tipper by rotating it to an upwards position (Fig 21). This prevents the anti-tipper from getting damaged.

To do so, loosen the thumb screw (Fig. 22, item A) from the swinging zone (Fig. 22, item B) and rotate the anti-tipper approximately 180 degrees. Then screw the thumb screw into the central locking zone (Fig. 22, item C).

Now the wheelchair can be moved over the obstacle forwards/backwards.

INFORMATION

Please note that it is necessary to swing the anti-tipper down after clearing the obstacle to reactivate it! To do so, loosen the thumb screw on the anti-tipper (Fig. 22, item A) and rotate the anti-tipper approximately 180°. Then firmly retighten the thumb screw into the swinging zone (Fig. 22, item B).

5.3 Transport Wheels (Fig. 23)

⚠ CAUTION

Accident risk. Please note that once the rear wheels are removed, the wheel locks are no longer functional!

With the rear wheels removed, the wheelchair can still be used as a rolling base; this makes it possible to manoeuvre through narrow passages (e.g. narrow bathroom doorways or airplane aisles).


The transport wheels have to be mounted to the rear wheel adapters. Remove the attachment screws of the rear wheel adapter one after the other, replace them with the longer M6x80 oval head screws provided and tighten the screw connections with the hexagon nuts provided (Fig. 23, item A). The correct torque for the screw connection is **10 Nm**.

Now attach the transport wheels and spacer sleeves provided to the protruding screw ends and secure the transport wheels with the cap nuts provided (Fig. 23, item B).

Finally screw the plastic brackets to the central unit (Fig. 23, item C) using the cap nuts/spacer washers provided.

5.4 Crutch Holder with Hook and Loop Strap (Fig. 24)

Allows crutches to be attached to the wheelchair.




5.5 Plug-on Side Panels (Fig. 25)

These side panels can be removed for getting into or out of the wheelchair. To remove the side panels, pull them out of the side panel attachment device (Fig. 25, item A). Moreover, the armrest can be adjusted in height (see Section 6.10).

5.6 Desk Side Panels (standard/stepped/height-adjustable; Fig. 26)

These side panels can be swung away to the rear or removed for transfer. To flip up the side panel, manipulate the release lever behind the side panel tube (Fig. 26 A/B/C, item A). Moreover, the armrest can be adjusted in height (see Section 6.10).

5.7 Plug-on Wheel Lock Lever Extension (Fig. 27)

The extension makes operation of the wheel lock easier for users with limited hand function.

5.8 Rear Wheels with Drum Brake (Fig. 28)

The drum brakes allow the attendant to comfortably and safely brake the wheelchair. The rear wheels remain removable via the quick-release axle system.

5.9 Wheel Lock for One-handed Operation (Fig. 29)

This wheel lock is recommended especially for hemiplegics. It can be activated either on the right or left side and ensures the safe locking of both rear wheels with a pull-wire system.

INFORMATION

Please note that the wheel lock lever extension must always be attached when releasing or engaging the wheel lock.

5.10 Tray (Fig. 30)

Serves as a supporting surface during meals, when working or when playing. The clear material allows visibility of the legs and correction of the sitting posture.

5.11 Lap Belt (Fig. 31)

The lap belt prevents the user from slipping and supports positioning.

5.12 Elevating Footrest (Fig. 32)

Allows positioning of the leg at different angles. Upon manipulation of the release lever which is attached to the swivel segment the footrest angle can be adjusted in small increments.

5.13 Push Handles, Height-adjustable (Fig. 33)


Can be adjusted vertically to a more comfortable height for the attendant.

5.14 Wheelbase, Long (Fig. 34)

This option makes the wheelchair particularly stable against tipping; also increases the turning radius.

INFORMATION

For transfemoral amputees the long wheelbase setting is indispensable.

5.15 Back Upholstery, Adaptable (Fig. 35)

This back upholstery can be adjusted in segments.

5.16 Folding Back (Fig. 36)

The folding back reduces the folding size and facilitates transportation in a car.

5.17 Angle-adjustable Back (Fig. 37)

The design allows to recline the backrest up to an angle of 30 degrees. The pull cable is used to unlock the adjustable back.

5.18 Spoke Protector (not illustrated)

Protects the fingers from getting caught in the wheel spokes. Also improves the appearance of the wheelchair.

5.19 Amputation Legrest (Fig. 38)

Allows positioning of residual limbs at different angles. Upon manipulation of the release lever which is attached to the swivel segment the amputation legrest angle can be adjusted in small increments.

5.20 "The Back" (Fig. 39)

Firm back to support positioning (see 647H491 Instructions for Use).

5.21 "The Seat" (Fig. 40)

Seat base with variable setting possibilities to improve positioning (see 647H450 Instructions for Use).

6 Adjustment / Assembly Instructions

INFORMATION

The following tools are required for repair and maintenance:

- Allen wrenches, sizes 4 mm and 5 mm
- Torx socket wrenches, sizes x 30 (for models from 2007)
- Open-end wrenches, sizes 19 and 24 mm
- Socket wrenches, sizes 10 and 19 mm
- Phillips head screwdriver
- Tyre levers
- Torque wrench

Each wheelchair user has their own ideas regarding the wheelchair type and accessories that are correct for him / her, as well as the position in the wheelchair and the handling characteristics. Therefore, the Start wheelchair offers many options that permit the wheelchair to be adapted according to your requirements. We recommend consulting your dealer or therapist to determine the correct wheelchair settings for you.

6.1 Adjusting the Seat Height and Seat Inclination (Fig. 41/42)

The higher you attach the rear wheel to the frame, the more the seat inclines downward. This has two effects:

- First, the wheelchair has a greater tendency to tilt backwards.
- Secondly, the user sits deeper in the wheelchair and thus has more stability.

A careful and synchronized height adjustment of both the rear wheels and the casters will allow to adapt the seat height to the user's individual requirements.

Start M1 Intro, M2 Effect, M3 Hemi and M6 Junior: These versions offer the possibility to change the vertical position of the rear wheel by shifting the rear wheel adapter along the frame. Loosen the two screw connections (Fig. 41, item A) and remove the screws. Move the adapter to the desired position and firmly retighten the screws.

For fine adjustment, loosen the nut of the fitting on the inside of the rear wheel adapter (Fig. 41, item B) and move the fitting along the slot hole. The setting must always be equal on both sides.



Start M4 XXL and M5 Comfort: Loosen the 4 screw connections which connect the rear wheel adapter to the central unit (Fig. 42). The rear wheel adapter can then be shifted up or down. The setting must always be equal on both sides.

INFORMATION

When the rear wheel position is changed, the angle of the caster axle in relation to the ground also changes. Ensure that this angle is always approx. 90° wide (see Section 6.4). The wheel lock must be re-adjusted as well (see Section 6.5).

INFORMATION

Make sure to firmly retighten all screws and nuts after making adjustments. The proper torque for the screw connections of the rear wheel adapter is **10 Nm** for screw connection A, **50 Nm** for screw connection B of the fitting, and **8 Nm** for screw connection C (see Fig. 41/42).

6.2 Changing the Wheelbase (Fig. 43/44)

Shifting the rear wheels backwards lengthens the wheelbase and thus provides for greater stability of the wheelchair. Attaching the rear wheels in a more forward position (only for advanced wheelchair users!) reduces the load on the casters, which makes the wheelchair more manoeuvrable. The caster wheels can then also be lifted more easily to clear obstacles such as curbs or steps.

Start M1 Intro, M2 Effect, M3 Hemi and M6 Junior: Loosen the two screws of the rear wheel adapter, turn it round and mount the rear wheel adapter in the desired front/back position (Fig. 43/44). The setting must always be equal on both sides.

Start M4 XXL and M5 Comfort: Loosen the nut of the fitting on the inside of the rear wheel adapter and mount the fitting in the desired more forward/backward position (Fig. 44). The setting must always be equal on both sides.

The distance between rear wheel and side panel can be continuously adjusted by adjusting the counter nut on the outer side of the rear wheel adapter.

When changing the adjustments originally made in the factory, two other assemblies of the wheelchair must be readjusted, refer to Sections 6.4 and 6.5.

⚠ CAUTION

Risk of tipping. Please note: If the rear wheels are in a more forward position and the user's body is not appropriately positioned, the user may tip backward – even on level ground!

⚠ CAUTION



Risk of tipping. For untrained wheelchair users and when there are extreme adjustments to the rear wheel, an anti-tipper is strongly recommended.

⚠ CAUTION

Risk of tipping. For transfemoral amputees it is essential to shift the rear wheels backwards!

INFORMATION

Make sure to firmly retighten all screws and nuts after making adjustments. The proper torque for the screw connection A of the rear wheel adapter is **10 Nm** and for the screw connection B of the fitting **50 Nm** (see Fig. 39/40).


6.3 Setting the Wheel Camber – Only Start M5 Comfort (Fig. 45–48)

The rear wheel adapters of the Start M5 Comfort allow for different camber angle settings of the rear wheels (0° and 2.5°). The greater the negative camber angle of the rear wheels, the more manoeuvrable the wheelchair becomes (Fig. 45).

The rear wheel adapter itself has a wheel camber of 0° or 2.5°, depending on the arrangement of the wedge-shaped camber washers. The wheel camber can be changed by turning round the mounted camber washers. Proceed as follows:

First disassemble the nuts on the rear wheel fitting.

Setting the rear wheel camber to 0°: If applicable, remove the camber washer from the inside of the rear wheel adapter and place it against the washer on the outside of the adapter such that the thinner side of the washer is next to the thicker side of the other washer (Fig. 46).


Setting the rear wheel camber to 2.5°: Position the camber washer with the thinner side pointing upwards into the desired engaging position (Fig. 47) on the outside of the rear wheel adapter. Place the second camber washer with the thicker side pointing upwards into the corresponding engaging position on the inside of the rear wheel adapter.

Further proceeding: Once you have applied the retaining washer from the inside (Fig. 48, arrow), slightly tighten the nuts on the rear wheel fitting. Check the space between the rear wheel and the side panel. You can adjust the space by rotating the fitting.

To finish, tighten the nuts of the rear wheel fitting.

INFORMATION

Observe the values indicated in section 6.1.

6.4 Caster Trailing Adjustment (Fig. 49)

After the rear wheels have been mounted in the appropriate position, the caster wheel journal angle must be readjusted. The threaded axle should be as vertical as possible in relation to the ground to ensure optimal driving behaviour of the wheelchair. The caster adapter allows for continuous angle adjustment.

Loosen the two screws between caster adapter and frame (Fig. 49 A/B). Now you can move the caster adapter along the frame tube. The level (included in delivery) will help you position the threaded axle as vertically to the ground as possible. Finally, retighten all screws.

INFORMATION

The correct torque for the screws is **8 Nm**.

The correct torque for the screws in case of mounted foot space enlargement (only M3) is **10 Nm**.

6.5 Adjusting the Wheel Lock (Fig. 50/51)

To change the position of the rear wheels, first loosen the attachment screws of the clamp adapters of the wheel lock and push the wheel lock forward (Fig. 50). Then reposition the rear wheels. The wheel lock should then be remounted such that the distance between the tyres and brake block is max. **4 mm** when the wheel lock is not activated (Fig. 51; subject to technical changes).

⚠ CAUTION

Accident risk. The wheel locks are only effective with sufficient air pressure and proper adjustment. The proper air pressure is indicated on the side wall of the tyre, but should at least be **3.5 bar**. To ensure sufficient braking effect, use only original rear wheels with an approved radial excursion of maximally $\pm 1 \text{ mm}$.

INFORMATION

Make sure to firmly retighten all screws and nuts after making adjustments!

INFORMATION

Alternatingly tighten the screws of the clamp adapter with a torque of **16 Nm**.

6.6 Adjusting the Brake Force – Drum Brake (Fig. 52)

The brake force is adjusted at the adjustment screw to achieve optimal braking effect. The brake force is increased by turning the adjustment screw (Fig. 52, item A) counter-clockwise and decreased by turning the screw clockwise.

Turn the adjustment screw counter-clockwise until frictional noise occurs at the turning wheel. Then tighten the adjustment screw until the frictional noise disappears. The wheel should run freely.

After finishing the adjustment, fasten the adjustment screw by tightening the counter nut (Fig. 52, item B).

INFORMATION

Ensure that the brake force of both rear wheels is adjusted equally.

INFORMATION


Make sure that the drum brake engages with sufficient force only when the manual brake lever is set to the second ratchet position.

6.7 Adjusting the Footplate Angle (Fig. 53, item A)

Prior to adjusting the footplate angle, loosen the screw connections as indicated in the figure. Then set the footplate to the desired angle.

INFORMATION

Make sure to re-tighten the loosened screw connections. The proper torque is **6 Nm**.

6.8 Adjusting the Lower Leg Length (Fig. 54/55)

By loosening the set screw, the footplate can be adapted to your lower leg length and to the thickness of your seat cushion (standard footrest = Fig. 54; elevating footrest = Fig. 55).

Make sure during adjustment that the footrest bar is slid into the swivel segment at least **40 mm**. A mark on the footrest bar also indicates the minimum distance the bar must be inserted during installation.

INFORMATION

After having made adjustments, retighten the set screw to a torque of **6 Nm**.

6.9 Adjusting the Angle of the Elevating Footrest (Fig. 56)

To adjust the footrest angle, turn the release lever until it hits the stop.

While propping up the footrest bar, set the desired angle. Then carefully let the release lever turn back. At the next free position, the footrest will snap into place.

6.10 Adjusting the Height of the Armrest (Fig. 57–60)

The armrest can be adjusted in height by bringing it into different positions:

Desk side panel, height-adjustable: For height adjustment of the side panel, press the lock button. Simultaneously, pull the armrest into the desired position (Fig. 57).

Desk side panel (standard/stepped): For height adjustment of the armrest, change the position of the spacers (Fig. 58/59). There are three possible height settings (both spacers are positioned above/below the side panel tube or one spacer is positioned above the tube and the other spacer below the tube). To bring the armrest into the highest position, remove the two oval head screws completely from the tube bend. Then place the two spacers above the tube bend between armrest and side panel tube. Re-attach the oval head screws. Make sure to firmly retighten all screws and nuts after making adjustments.

Plug-on side panel: For height adjustment of the side panel you have to adjust the height of the clothing protector. Press the counter-sunk lock button using a ball point pen / key / screwdriver, etc. (Fig. 60). Move the clothing protector into the desired position.

6.11 Adjusting the Seat Depth (Fig. 61)

Undo the hook and loop connections and remove the back padding/seat cushion.

Loosen all screw connections that connect the central unit to the frame tubes and remove the screws (Fig. 61). The central unit can then be shifted along the frame tubes. Depending on the design, 3 different seat depth settings are possible in increments of 3 cm.

INFORMATION

Make sure to firmly retighten all screws and nuts after making adjustments. The proper torque is **8 Nm**.

Fasten the hook and loop attachment again between the back padding and the seat upholstery.



6.12 Adjusting the Back Height (Fig. 62)

The back height is adjustable by **7.5 cm**. The back tube which is integrated into the rear frame allows for the adjustment of the back height in increments of 2.5 cm.

To do so, remove all screws from the vertical attachment device of the central unit (Fig. 62). Set the desired back height.

INFORMATION

Make sure to firmly retighten all screws and nuts after making adjustments. The proper torque for screw connection A is **8 Nm** and for screw connection B (rear wheel adapter) **10 Nm** (see Fig. 62). When using the maximal back height setting, an additional screw connection with plastic spacers (included in the 481D53=ST170 Conversion Kit for back tube, standard) has to be mounted in the uppermost bore hole of the central unit (Fig. 62, item C).

6.13 Seat Upholstery (not illustrated)

The tension of the seat upholstery on your wheelchair can be readjusted to compensate for sagging.

6.14 Adapting the Back Upholstery (Fig. 63)

If you have selected the option “adaptable back upholstery”, you can adjust the back upholstery in segments according to your needs.

To do this, remove the back upholstery padding and undo the hook-and-loop fasteners of the straps (Fig. 63).

Reposition the straps as desired, attach them again with the hook-and-loop fasteners and cover with the padding.

6.15 Changing the Back Angle (Fig. 64)

⚠ CAUTION

Risk of tipping due to missing anti-tipper. When using the 30° back angle adjustment and in case of short wheelbase, 2 anti-tippers (one on each side) must be mounted and in a functional position; in case of long wheelbase, at least one anti-tipper must be mounted and in a functional position.

For therapeutic reasons, individual adjustment of the preset angle between the seat and backrest may be required. The angle can be adjusted between 90° and 120°.

Pull the release cable of the back release until the angle adjustment is released. Bring the backrest into the desired position and let go of the release cable. Make sure that the lock bolts safely engage on both sides.

6.16 Push Handles, Height-adjustable (Fig. 65)

To adjust the suitable height of the push handles, use the adjustable clamp lever to loosen the clamp screw.

INFORMATION

Be sure to snugly retighten the clamp levers after all adjustments!

7 Maintenance, Cleaning and Disinfection

7.1 Maintenance

Your product has been provided with the CE marking. The manufacturer herewith guarantees that this medical device as a whole conforms to the requirements of the European Directive for Medical Devices 93/42/EEC.

The proper function of the wheelchair, especially of the wheel locks and brakes, should be checked before every use. Safety nuts should be used only once. If they have been loosened several times, these nuts must be replaced.

If you notice any defects, please contact your authorised specialist dealer immediately in order to have them corrected. This is particularly important if you notice that the driving characteristics change or the wheelchair becomes less stable. **We also recommend that you have the wheelchair maintained regularly at 12-month intervals by your authorised specialist dealer.** The functions described in the following list must be verified by the user at the specified intervals.

Check	Before each use	Monthly	Quarterly
Function test of the wheel locks/brakes	X		
Sagging of the seat and back upholstery		X	
Setting of the caster journal bearing		X	
Check stability of footplate		X	
Visual inspection of wear parts (e.g. tyres, bearings)		X	
Dirt on bearings		X	
Push rings for damages		X	
Air pressure (see indication on tyre)		X	
Folding mechanism for wear and tear		X	
Spoke tension of the rear wheels			X
Screw connections			X

To ensure smooth operation at all times, the following precautions and cleanings should be carried out. Users may carry these out by themselves:

- Hair and dirt particles generally accumulate between the caster wheel and fork. This can restrict the caster wheels from rotating smoothly. Remove the casters and thoroughly clean the forks and casters using a mild household cleaner.
- The rear wheels feature a quick-release system. To keep this system operational, ensure that no dirt adheres to the quick-release axle or axle housing. The quick-release axle should also be lightly lubricated from time to time with resin-free sewing machine oil.
- If your wheelchair gets wet, towel dry it as soon as possible.
- Do not use the wheelchair in salt water and keep sand or other particles from damaging the wheel bearings.

- Screw connections should be checked frequently, in particular when beginning to use the wheelchair and after any adjustment. If a screw connection comes loose repeatedly, please consult your specialist dealer.

7.2 Cleaning and Disinfection

Clean the padding and upholstery with warm water and a washing-up liquid. Remove spots with a sponge or a soft brush. Rinse with clear water and let the components dry.

INFORMATION

Do not use any aggressive cleansing agents, solvents, or hard brushes etc.

INFORMATION

No wet-cleaning. Components are not machine washable. For disinfection, water-based agents (e.g., Sagrotan original concentrate) should be used. The manufacturer's instructions for use must be observed.

INFORMATION

Before disinfecting, the pads and handles must be cleaned. Plastic parts, frame parts as well as chassis and wheels can be cleaned with a mild detergent. Dry thoroughly afterwards.

8 Changing Tyres (Fig. 66–69)

Changing flat tyres requires only the necessary tools and minimal technical skills. Users may change tyres themselves if they wish. We recommend carrying an emergency repair kit and tyre pump at all times (if using pneumatic tyres).

Suitable tyre pumps are listed on the order form and are supplied with the product. Alternatively, users can keep a spray can of hardening foam (available at bicycle shops) to fill a tyre in case of punctures.

- In case of a flat tyre, carefully remove the tyre from the rim using the appropriate tools.
- Take care not to damage the rim or the inner tube.
- Repair the inner tube according to the directions in the repair kit, or, replace the damaged inner tube with a new one.
- Before re-mounting the tyre, inspect the inner rim surface and the inside of the tyre for any object that may have caused the puncture.
- Be sure that the rim bands are not damaged. They protect the tube from being damaged by the ends of the spokes.

Mounting

- Place the protective rim band over the valve. Then push the valve through the rim. Screw on the valve nut. The rim band can then be easily applied.
- Ensure that all spoke ends are covered.

Tyres

- Beginning behind the valve, place the lower side of the tyre over the edge of the wheel rim. Pump the tube with a small amount of air until it assumes its round shape, and then insert it into the tyre.
- Check to ensure the inner tube is properly positioned all around and is free of folds; if this is not the case, release some of the air. Now you can easily install the upper side of the tyre starting across from the valve, using both hands and working towards the valve.

Inflating

- Check around the tyre on both sides to ensure that the tube is not clamped between the edge of the tyre and the rim.
- Push the valve back slightly and then pull it out again to ensure that the tyre is properly seated in the valve area.
- Inflate the tyre until it can still be pressed easily with your thumb. If the control line on both sides of the tyre is equally distanced to the rim, the tyre sits centrally. If the tyre does not sit centrally, let out some air and realign the tyre. Then inflate the tyre to maximum operating pressure (see sidewall) and tighten the dust cap.

9 Technical Data

Overall length (in cm)

Module	Rear wheel size					
	20x1 3/8"		22x1 3/8"		24x1 3/8"	
	minimal ¹⁾	maximal ²⁾	minimal ¹⁾	maximal ²⁾	minimal ¹⁾	maximal ²⁾
M1	–	–	–	–	94.5	116.0
M2, M3	–	–	92.0	113.5	94.5	116.0
M4	–	–	87.0	106.0	89.5	109.0
M5	–	–	88.0	107.0	90.5	110.0
M6	80.0	101.5	82.5	104.0	85.0	106.5
M2 Asia/ Am Asia	83.5	105.0	86.0	107.5	88.5	110.0

¹⁾ measured with min. seat depth, smallest standard lower leg length (30 cm), short/active wheelbase

²⁾ measured with max. seat depth, largest standard lower leg length (47 cm), long/passive wheelbase

Overall width (in cm)¹⁾

Seat width	Standard rear wheel max.	Rear wheel with drum brake max.
28.0	49.5	52.0
30.5	51.5	54.5
33.0	54.5	57.0
35.5	57.0	59.5
38.0	59.5	62.0
40.5	62.0	64.5
43.0	64.5	67.0
45.5	67.0	69.5
48.0	69.5	72.0
50.5	72.0	74.5
53.0	74.5	77.0
55,5	77.0	79.5
58,0	79.5	82.0

¹⁾ Applicable to push ring attachment, narrow (for push ring attachment, wide: +2 cm) and a 0° camber of the rear wheels. When using height-adjustable side panels: all overall widths +2 cm.

Start M1, M2, M4, M5, M6, Asia

**Anterior seat height (in cm)
without foot space enlargement**

Caster fork, short		
Position		Setting range ¹⁾
Caster diameter	5"	40–44
	5.5"	42–44
	6"	43–45
	7"	45–47
Caster fork, long		
Position		Setting range
Caster diameter	5"	44–47
	5.5"	44–48
	6"	45–49
	7"	47–50
	8"	48–51

Start M3

**Anterior seat height (in cm)
with foot space enlargement**

Caster fork, short		
Position		Setting range ¹⁾
Caster diameter	5"	38–40
	5.5"	39–41
	6"	40–42
	7"	42–43
Caster fork, long		
Position		Setting range
Caster diameter	5"	41–44
	5.5"	42–45
	6"	42–45
	7"	44–47
	8"	45–48

1) Specified at 0° seat angle, without seat cushion

- Anterior seat height - posterior seat height ≤ 3.5 cm
- Casters in dependency of the anterior seat height
- Caster fork in dependency of the anterior seat height and casters
- Rear wheel size in dependency of the posterior seat height

Start M1, M2, M3, M6, Asia

Posterior seat height (in cm)

Posterior seat height	Setting range	
Rear wheel size	20"	37–43
	22"	38–45
	24"	41–49


Start M4, M5

Posterior seat height (in cm)


Free mounting positions	Rear wheel size	
	22"	24"
1	38	41
2	41	44
3	44	47
4	46	50

Rear wheel adapter

M1, M2, M3, M6

M4, M5

10 Re-use / Disposal

10.1 Information on Re-Use

⚠ CAUTION

Risk due to improper use. Seat paddings that come directly into contact with the skin can cause functional and hygienic risks if the product is used by another person. They must be replaced for re-use.

The product is suitable for re-use.

Similar to second-hand machines or cars, products that are being re-used are subject to increased strain. Features and functions must not change in a way that could endanger patients or other persons within the product's life cycle.

Based on market observations and the current state of technology, the manufacturer has calculated that the product can be used for a period of 4 years, provided that it is used properly and that the service and maintenance instructions are observed. Periods during which the product is stored at the dealer or with the third party payer are not included in this period. It should be clearly pointed out, however, that the product is reliable far beyond this defined period of time, provided that it is cared for and maintained appropriately.

In cases of re-use, the corresponding product must first be thoroughly cleaned and disinfected. Subsequently, the product must be examined by an authorised specialist to check the condition, wear and possible damages. Any worn and damaged components as well as components which do not fit or are not suitable for the user must be replaced.

The service manual includes a service schedule, detailed information and a list of the required tools.

10.2 Information on Disposal

If a product is to be disposed of, all components of the product must be disposed of properly in accordance with the applicable environmental protection regulations in the respective country.

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Otto Bock South Africa (Pty) Ltd
Building 3 Thornhill Office Park · 94 Bekker Road
Midrand · Johannesburg · South Africa
T +27 11 312 1255
info-southafrica@ottobock.co.za
www.ottobock.co.za

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Otto Bock HealthCare Canada
5470 Harvester Road
Burlington, Ontario, L7L 5N5, Canada
T +1 289 288-4848 · F +1 289 288-4837
infocanada@ottobock.com · www.ottobock.ca

Otto Bock HealthCare Andina Ltda.
Calle 138 No 53-38 · Bogotá · Colombia
T +57 1 8619988 · F +57 1 8619977
info@ottobock.com.co · www.ottobock.com.co

Otto Bock de Mexico S.A. de C.V.
Prolongación Calle 18 No. 178-A
Col. San Pedro de los Pinos
C.P. 01180 México, D.F. · Mexico
T +52 55 5575 0290 · F +52 55 5575 0234
info@ottobock.com.mx · www.ottobock.com.mx

Otto Bock HealthCare
Two Carlson Parkway North, Suite 100
Minneapolis, MN 55447 · USA
T +1 763 553 9464 · F +1 763 519 6153
usa.customerservice@ottobockus.com
www.ottobockus.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
Suite 1.01, Century Corporate Centre
62 Norwest Boulevard
Baulkham Hills NSW 2153 · Australia
T +61 2 8818 2800 · F +61 2 8814 4500
healthcare@ottobock.com.au · www.ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co., Ltd.
B12E, Universal Business Park
10 Jiuxianqiao Road, Chao Yang District
Beijing, 100015, P.R. China
T +8610 8598 6880 · F +8610 8598 0040
news-service@ottobock.com.cn
www.ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Unit 1004, 10/F, Greenfield Tower, Concordia Plaza
1 Science Museum Road, Tsim Sha Tsui
Kowloon, Hong Kong · China
T +852 2598 9772 · F +852 2598 7886
info@ottobock.com.hk · www.ottobock.com

Otto Bock HealthCare India
Behind FairLawn Housing Society
St. Gregorios Lane, Sion Trombay Road
Chembur, Mumbai. 400071 · India
T +91 22 2520 1268 · F +91 22 2520 1267
information@indiaottobock.com · www.ottobock.in

Otto Bock Japan K. K.
Yokogawa Building 8F, 4-4-44 Shibaura
Minato-ku, Tokyo, 108-0023 · Japan
T +81 3 3798-2111 · F +81 3 3798-2112
ottobock@ottobock.co.jp · www.ottobock.co.jp

Otto Bock Korea HealthCare Inc.
4F Agaworld Building · 1357-74, Seocho-dong
Seocho-ku, 137-070 Seoul · Korea
T +82 2 577-3831 · F +82 2 577-3828
info@ottobockkorea.com · www.ottobockkorea.com

Otto Bock South East Asia Co., Ltd.
1741 Phaholyothin Road
Kwaeng Chatuchark · Khet Chatuchark
Bangkok 10900 · Thailand
T +66 2 930 3030 · F +66 2 930 3311
obsea@ottobock.co.th · www.ottobock.co.th

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Ottobock has a certified Quality Management System in accordance with ISO 13485.

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Otto Bock HealthCare GmbH
Otto Bock Mobility Solutions GmbH