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We at SUNRISE MEDICAL have been awarded the ISO-9001 certificate, which affirms the quality of our products at every stage, from R & D to production. This products meet the requirements in accordance with EC guidelines. Options or accessories shown are available at extra cost.

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If you have any queries about the use, maintenance or safety of your seating system, please contact your local approved Sunrise Medical service agent. If you do not know of an approved dealer in your area or have any other questions please write or telephone:

**Sunrise Medical**

**Thorns Road**

**Brierley Hill**

**West Midlands**

**DY5 2LD**

**England**

**Phone: 0845 605 66 88**

**Fax: 0845 605 66 89**

**www.SunriseMedical.co.uk**

## Contents

<b>1.0 User information .....</b>	<b>4</b>
1.1 This user manual .....	4
1.2 For further information .....	4
1.3 Symbols used in this manual.....	4
<b>2.0 Safety .....</b>	<b>5</b>
2.1 Seating related symbols and labels used on the product..	5
2.2 Safety: Temperature .....	5
2.3 Safety: Moving parts.....	6
2.4 Safety: Choking hazard .....	6
2.5 Safety: Trays .....	6
<b>3.0 General description .....</b>	<b>7</b>
<b>4.0 Setting up the seating system.....</b>	<b>8</b>
4.1 Pre-set seat width.....	8
4.2 Pre-set armrests (standard SedeoErgo) .....	9
4.3 Preset Leg Lower Leg Length (standard Sedeo Ergo) .....	10
4.4 Preset Leg Rests and foot plates (Sedeo Ergo UP) .....	11
4.4.1 Pre-set lower leg length (Fig 4.10): .....	11
4.4.2 Depth adjustment of the knee fixation (Fig 4.12- 4.13):	11
4.4.3 Width adjustment of the knee fixation (Fig 4.12- 4.13):	11
4.4.4 Height adjustment of the knee fixation (Fig 4.12- 4.13):	11
4.5 Preset seat depth, (Fig. 4.14) .....	12
4.5.1 Sedeo Ergo (Standard), (Fig. 4.14) .....	12
4.5.2 Sedeo Ergo (UP).....	12
4.6 Swing-away controller bracket .....	13
<b>5.0 Seating positioning accessories .....</b>	<b>14</b>
5.1 Lap strap / seating positioning belt.....	14
5.2 Stand Up Functions (Q700-UP M).....	16
5.2.1 Chest strap Adjustments: .....	16
<b>6.0 Using the seating system.....</b>	<b>18</b>
6.1 Getting in and out of the wheelchair (Transferring) .....	18
6.2 (Powered) seating options .....	20
<b>7.0 Maintenance and cleaning .....</b>	<b>21</b>
7.1 General Maintenance.....	21
7.2 Cleaning your seating system:.....	21
7.3 Hygiene measures when being re-used:.....	21
7.4 General cleaning instructions: .....	22
7.5 Cleaning Accessories: .....	22
7.6 Storage:.....	22
7.7 Disposal: .....	22
<b>8.0 Technical specifications: .....</b>	<b>23</b>
8.1 Applicable norms / standards .....	23
8.2 Technical specifications: SEDEO Ergo seating system .....	25
8.3 Technical specifications: SEDEO Ergo UP seating system ..	25
<b>9.0 Warranty .....</b>	<b>26</b>

Dealer signature and stamp

## 1.0 User information

Congratulations on choosing a Sunrise Medical wheelchair. Sunrise Medical's high-quality mobility products are designed to enhance independence and make your everyday life easier.

### 1.1 This user manual

This user manual will help you to use and maintain your wheelchair safely. This entire user manual for this wheelchair consists of three booklets:

- The general wheelchair user manual
- The user manual for the seating system (this booklet)
- The user manual for the controller.

When necessary this general user manual refers to other manuals as shown below:

 General manual:	Refers to the general user manual for the wheelchair.
 Controller:	Refers to the user manual for the controller.
 J3 Back:	Refers to the user manual for the J3 Back Rest.
 Battery charger:	Refers to the user manual for the battery charger.
 Service manual:	Refers to the service manual for qualified personnel.

**Do not use your wheelchair until this entire manual (and all other booklets) has been read and understood!**

If one of the user manuals was not included with your wheelchair, please contact your dealer immediately. In addition to this user manual, there is also a service manual for qualified specialists.

### 1.2 For further information

Please contact your local, authorized Sunrise Medical dealer if you have any questions regarding the use, maintenance or safety of your wheelchair. In case there is no authorized dealer in your area or you have any questions, contact Sunrise Medical either in writing or by telephone.

For information about product safety notices and product recalls, go to [www.SunriseMedical.eu](http://www.SunriseMedical.eu)

### 1.3 Symbols used in this manual

 **DANGER**

**Potential risk of injury or serious injury or death**

 **WARNING!**

**Potential risk of injury.**


 **CAUTION!**

**Potential damage to equipment.**

**NOTE:**

**General user advice.**

Not following these instructions may result in physical injury, damage to the product or damage to the environment!

If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.co.uk](http://www.SunriseMedical.co.uk)

or alternatively is available on request in large text.

Fig. 1.1

## 2.0 Safety

Follow the instructions carefully next to these warning symbols! Not paying careful attention to these instructions could result in physical injury or damage to the wheelchair or the environment. Wherever possible, safety information is provided in the relevant chapter.

### 2.1 Seating related symbols and labels used on the product

The signs, symbols and instructions affixed to the wheelchair comprise part of the safety facilities. They must never be covered or removed. They must remain present and clearly legible throughout the entire lifespan of the wheelchair.

Replace or repair all illegible or damaged signs, symbols and instructions immediately. Please contact your dealer for assistance.

Fig 2.1: Warning – Danger of finger entrapment.

Fig 2.2: Warning – Ensure that the seat interface quick release is screwed in tight before using your wheelchair, especially on a slope.

Fig 2.3: Warning – Do not drive your wheelchair on a slope with the backrest reclined and/or the seating in standing mode.

### 2.2 Safety: Temperature

#### WARNING!

- Avoid physical contact with the wheelchair's motors at all times. Motors are continuously in motion during use and can reach high temperatures. After use, the motors will cool down slowly. Physical contact could cause burns. Allow the motors after using at least 30 min. to cool down.
- If you do not use the wheelchair, ensure that it is not exposed to direct sunlight for lengthy periods of time. Certain parts of the wheelchair, such as the seat, the back and the armrests can become hot if they have been exposed to full sunlight for too long. This may cause burns or allergic reactions to the skin.


Fig. 2.1

Fig. 2.2

Fig. 2.3

### 2.3 Safety: Moving parts

**DANGER!**

A wheelchair has moving and rotating parts. Contact with moving parts may result in serious physical injury or damage to the wheelchair. Contact with the moving parts of the wheelchair should be avoided.

- Wheels (turning and castor)
- Powered tilt in space adjustment
- Powered seat elevator option
- Powered recline adjustment
- Swing-away bracket of the controller

### 2.4 Safety: Choking hazard

**DANGER!**

This mobility aid uses small parts which under certain circumstances may present a choking hazard to young children.

### 2.5 Safety: Trays

#### Safety: Vent tray

**DANGER!**

Using a vent tray will affect the stability and overall weight of your wheelchair.

- Make sure that the batteries and the ventilator are securely fastened to the tray before use.
- Familiarise yourself with the increased size of the wheelchair before driving to prevent potential collisions.
- When activating the recline backrest make sure that the area around the vent tray is clear.
- Be aware that the vent tray assembly can cause a finger trap hazard when the backrest is reclined.
- The vent tray is not designed to hold anything other than the vent unit and its batteries, it is not designed to hold an oxygen bottle etc. Use of these items may cause your wheelchair to become unstable.
- The vent tray must be installed and maintained by a Sunrise Medical authorised dealer.

#### Safety: Swing away tray:

**DANGER!**

- The maximum weight allowed for the tray is 2.5kg.
- Do not overload the tray, this could cause the tray to break or could cause the chair to become unstable.
- Do not leave lit cigarettes or other heat sources on the tray as this could cause the tray to deform and mark.
- Ensure that all extremities and clothing are clear when positioning the tray for use.

### 3.0 General description

This wheelchair is equipped with Sedeo Ergo seating system. The SedeoErgo seating system provides good support for the head, torso, arms and legs. Optimal body support ensures maximum independence for the wheelchair user.

The wheelchair must be adjusted by a qualified specialist. Incorrect settings may have negative consequences for the seating posture.

Fig. 3.1

**Main components:**

1. Headrest
2. Joystick
3. Armrest
4. Side cushion
5. Seat cushion
6. Leg rest
7. Back upholstery

## 4.0 Setting up the seating system

The average wheelchair user does not exist. This is why Sunrise Medical wheelchairs can be adjusted according to specific needs of the user.

In this, we make a distinction between wheelchair settings and wheelchair seating adjustment options. Settings are carried out once, the first time the wheelchair is going to be used. They must be executed by qualified specialist, unless explicitly stated otherwise. Seating adjustment options can be carried out by the user. These adjustment options do not require the use of tools.

### 4.1 Pre-set seat width

#### Pre-set the armrest width:

- Loosen the two screws (5.0mm Allen key) on the traversing tube as shown in (Fig. 4.1 - A).
- Move the armrest receiver brackets to the desired position.
- Tighten the screws firmly prior to use.

#### ⚠ CAUTION!

- Please check that the new arm rest position does not interfere with any seat positioning.

#### Pre-set the seat plate width

Lift the seat cushion off the Velcro® strips to access the seat plates.

- Use a 5.0 mm Allen Key, (Hex Wrench).
- Adjust the seat plate width by loosening the four screws, circled in Fig. (4.2).
- Slide the plates outwards or inwards to achieve the required seat width.
- The seat width is indicated by the scale that runs alongside each slot (Fig.4.2).
- Tighten the four screws securely.

#### Pre-set the back plate width

Lift the back cushion off the Velcro® strips to access the back plate (Fig.4.3).

- Use a 5.0 mm Allen Key, (Hex Wrench).
- Loosen all of the 12 screws located within the slots.

#### Define the back plate width

- Slide the side plates in or out to the required width.
- Tighten the 8 screws securely (white marked in Fig. 4.3)

#### Positioning/centering the back plate width

- Slide the total to the left and the right to align with the marker pin.
- Use the measuring scale embedded on the back plates and the marker pin for reference.
- Tighten all remaining 4 screws securely (black marked in Fig. 4.3).
- Fit a backrest cushion onto the back plate.

#### ⚠ CAUTION!

After finishing seat width adjustments, make sure that the armrest can be flipped up without colliding with the backrest. Not doing so can lead to damaging the product severely.

📖 J3 Back, (Fig.4.3.1):

NOTE: The Jay 3 Backrest can be used in conjunction with the Sedeo Ergo seating system. Please refer to the J3 Back Owner's Manual for fitting and adjustment.


Fig. 4.1


Fig. 4.2

Fig. 4.3

Fig. 4.3.1

## 4.2 Pre-set armrests (standard SedeoErgo)

### Pre-set armrest height

The following section outlines armrest height adjustment. In the case of either fine or gross adjustment, it is a two-part process:

PART 1: Adjust the armrest height using the rear pivot.

PART 2: Adjust arm pad angle using the front pivot.

#### PART 1 (Fig. 4.4)

Up to 145 mm of armrest height adjustment can be achieved using the 5.0 mm fine height adjustment screw to fine tune the armrest's height.

- Use the 5.0 mm Allen key that is stored on the inside of the left hand arm rest.
- Turn the screw, (A), at the upper rear pivot and your armrest will move up or down depending on the rotation direction, (1) + (2).

#### PART 2 (Fig. 4.5)

- When your armrest is at the desired height, use the same tool to turn the screw, (B), below the lower front pivot of the armrest.
- This process will set the arm pad angle to the required setting, (2) + (3).

### Pre-set the arm pad fore and aft: (Fig. 4.6)

- Loosen the 2 bolts at the bottom of the arm pad rail.
- Slide the arm pad along the channel nuts.
- At this point some inward and outward angle adjustment is also possible.
- Move the front bolt from front hole to the slot hole, (A).
- Tighten the bolts securely when the angle is set.
- When the correct position is reached tighten the 2 bolts securely with the 5.0 mm Allen Key, (Hex Wrench).
- Place the tool in its original space so it can be used for further adjustments in the future.

### WARNING!

All screws should be checked periodically for tightness. A loose arm pad could slip down when pressed.


Fig. 4.4


Fig. 4.5

Fig. 4.6

### 4.3 Preset Leg Lower Leg Length (standard Sedeo Ergo)

#### **⚠ WARNING!**

- Be aware of your environment to make sure you do not injure your legs when Leg Rests are extended.
- Always ensure that the Leg Rests or foot plates do not come into contact with the casters before driving the wheelchair.
- Leg Rests are not to be used for lifting or carrying the wheelchair under any circumstances.
- As with all moving parts be careful not to trap your fingers, (Fig.4.7).

#### **Pre-set foot plate height and angle**

The power centre mount footrest allows angle adjustment from 90° to 180° level to ground and a 185 mm powered length extension. The left and right footplate can be adjusted by an angle of +/- 5.0° and has a mechanical height range adjustment of 200 mm to 460 mm measured from the footplate to the seat pan.

#### **Pre-set foot plate height:**

- Undo the 2 bolts with a 5.0 mm Allen key, (A-Fig.4.8).
- Slide the footplate receiver up or down to the height required.
- Tighten the 2 bolts securely.
- Repeat for the other side.

#### **Pre-set foot plate angle:**

- Flip up the footplate so that you can reach the angle setting screw, (B-Fig.4.9).
- Use a 4.0 mm Allen key to adjust the bolt in or out depending on the angle required.

#### **⚠ WARNING!**

- After adjustment and with the seat in its normal drive position: Ensure there is a gap of, (>25 mm.), between the floor and the footplate at all times when maneuvering the Foot plate up or down.
- Perform this test with the user sitting in the chair loading the leg rest.

Using the chair with the leg rest set too low might result in personal injury or damage to the wheelchair and its surroundings.


Fig. 4.7

Fig. 4.8


Fig. 4.9

## 4.4 Preset Leg Rests and foot plates (Sedeo Ergo UP)

### 4.4.1 Pre-set lower leg length (Fig 4.10):

- Loosen bolts and take them out (2x bottom ones) at both sides of the foot plates. (C).
- Put footplates in desired position.
- Tighten bolts (2x bottom ones). (C).

Pre-set footplates individually (Fig 4.11):

They height of the footplates can be individually set by:

- Loosen 4 bolts (A)
- Turn screw B (fig 4.10) to adjust footplate height.
- Tighten 4 bolts (A)

### 4.4.2 Depth adjustment of the knee fixation (Fig 4.12 - 4.13):

The depth of the knee-pads can be adjusted by turning the knob (D).

### 4.4.3 Width adjustment of the knee fixation (Fig 4.12 - 4.13):

- Loosen bolts (E).
- Put knee-pads in desired position.
- Tighten bolts (E).

### 4.4.4 Height adjustment of the knee fixation (Fig 4.12 - 4.13):

2 ways:

Adjust position of knee fixation receiver (2 positions).

- Loosen bolts (4x). (F).
- Put knee fixation receiver in desired position.
- Tighten bolts (4x). (F).

Flip around the knee-pads (2 positions). These are mounted off-centered → 2 positions.

- Loosen bolts (E).
- Flip the knee pads around.
- Tighten bolts (E).


Fig. 4.10


Fig. 4.11


Fig. 4.12


Fig. 4.13

**4.5 Preset seat depth, (Fig. 4.14)**


**⚠ CAUTION!**

When a bigger seat depth is required ensure that the cables are lengthened before the mechanical adjustment takes place.

Be aware. As described below, there is a different procedure for adjusting the seat depth in case of Sedeo Ergo with an standing function (Sedeo Ergo UP).

**⚠ DANGER!**

We strongly recommend that you contact your Sunrise Medical approved service agent to adjust the seat depth, because adjusting the seat depth will affect the centre of gravity.


**4.5.1 Sedeo Ergo (Standard), (Fig. 4.14)**

- Use a 5.0 mm Allen Key, (Hex Wrench).
- Loosen the 2 screws, (A), on the bottom of the outer seat plate, but do not take the screws out as these prevent the shift pieces from pulling out.
- The back assembly can now be shifted backwards or forwards to adjust the seat depth.
- Tighten the screws.

**4.5.2 Sedeo Ergo (UP)**

- Use a 5.0 mm Allen Key, (Hex Wrench).
- Loosen the 2 screws on both sides (A), on the bottom of the outer seat plate, and take the bolts out (Fig 4.14).
- The back assembly can now be shifted backwards or forwards to adjust the seat depth.
- Tighten the screws.
- In addition adjust the armrest linkage simultaneously corresponding the marker holes (B in Fig 4.15).

**⚠ DANGER!**

Adjusting seat depth on a Q700 UP requires additional steps.

- Besides the seat depth adjustment as described in 4.5.2, the armrest linkage (fig 4.15) needs to be adjusted simultaneously.
- Use the corresponding markings on the seat depth in parallel to the armrest linkage.
- Not doing so will lead to mechanical clashes and results in damaging the product.

## 4.6 Swing-away controller bracket

The joystick is mounted with a mechanism that allows the joystick to be locked in a forward or retracted position, and can accommodate a mid-line mount so the joystick can be located in the center of the wheelchair.

### Swing-away, (Fig. 4.16).

- To Swing-away the joystick, push outward on the inside of joystick.
- Push the joystick away from the front of the armrest until it locks into the outside swing-away position.
- To return to forward position, push the back of the joystick, forward until it locks into place.
- To Swing-away the joystick to the inside, push inward on the outside of joystick.
- Push the joystick away from the front of the armrest until it locks into the inside swing-away position.
- To return to forward position, push the back of the joystick, forward until it locks into place.


Fig. 4.16

### Pre-setting depth of the controller, (Fig. 4.17)

1. Loosen (A).
2. Slide (B) into desired position.
3. Tighten (A).


Fig. 4.17

There are 2 ways to pre-set the height of the controller.

### Pre-setting height of the controller (1), (Fig. 4.18).

1. Loosen (A).
2. Slide (C) into desired position.
3. Tighten (A).


Fig. 4.18

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
Fig. 4.19

## 5.0 Seating positioning accessories

### 5.1 Lap strap / seating positioning belt



#### DANGER!

- This product is only to be used to position a single person in a wheelchair.
- Lap straps are not suitable for transit purposes, approved occupant restraint systems must be used.
- Failure to heed these warnings may result in severe injury or death.
- Ensure that the carer or attendant is trained in the correct operation of the belts.
- Untrained carers or attendants may lead to delays, if an emergency situation occurs.

#### Mounting a lap strap / seating positioning belt

Feed the bracket ends of the straps under the backrest brace bar as shown, using the lower fixing hole. Ensure that the adjuster buckles can be accessed and the strap is not twisted, (A-Fig. 5.2).

- Put a washer on the bolt.
- Pass the bolt through the lap strap fixing bracket.
- Put the next washer on to the bolt, (Fig. 5.3).
- Screw the bolt into the lowest hole in the frame. (B-Fig.5.2)
- Tighten using a 5.0mm hex wrench to a torque value of 10Nm, (C-Fig.5.2).
- Adjust the lap strap to suit, leaving no more than a hand's width gap for comfort and safety.
- The hand clearance should be with the lap strap under normal tension and not allow large gaps or loops, (Fig. 4.7).
- Generally, the lap strap should be fixed so that the straps sit at an angle of approximately 45° (Fig. 4.8), and when correctly adjusted should not allow user to slip down in the seat.
- Place the strap loosely across the seat with the opening end of the buckle facing to the right for a left handed person and to the left for a right-handed person. Pass the other ends of the strap through the gap between the backrest posts and the backrest. (Fig. 5.4)

## Positioning an person with a lap strap

### DANGER!

- Always make sure that the lap strap is correctly secured and adjusted prior to use.
- Too loose a strap could cause the user to slip down and cause serious injury.
- Check lap strap and securing components at regular intervals for any signs of fray or damage. Replace if necessary.
- When servicing, check for correct operation of the release buckle and for any signs of wear on the material or plastic brackets.
- Regular Checks/Activities:
- As with all positioning components, there is a need for corrective adjustments as the person changes their seating position over time.
- Check the belts regularly for correct fitting, to ensure the safety and comfort of the user.

### NOTE:

- Standard Sling - Lap strap movement is restricted by upholstery.
- Contoured back - Use universal bottom bracket as per recline back method.

Fig. 5.5


Fig. 5.6


Fig. 5.7

## 5.2 Stand Up Functions (Q700-UP M)

### WARNING!

- The chest strap may only be assembled and adapted by an authorized dealer.
- Before operating the stand-up function the chest strap and the positioning belt, supplied as standard has to be buckled in accordance with these user instructions.
- Always also use the knee fixation for the stand-up procedure.

### DANGER!

Due to the risk of choking it is extremely hazardous to use the belts without using the knee fixation.

### 5.2.1 Chest strap Adjustments:

#### Setting the height of the chest strap, (Fig.5.8).

**NOTE:** A 5.0 mm Hexagon wrench, (Allen Key), is required.

#### Positioning:

- Seat the user into the chair in the correct seating position.
- Loosen the Allen screw to be able to position the chest strap at the desired height.
- Tighten the Allen screws to a torque value of 9 Nm (+/- 10%).

### DANGER!

- To prevent choking hazards the belt may not be adjusted too high or too low relative to the users body.
- Make sure the rider does not slide in the wheelchair seat. If this occurs, the rider may suffer chest compression or suffocate due to pressure from the chest strap.
- Never use the chest strap as a patient restraint.
- Never use the chest strap on a user who is comatose or agitated.

### CAUTION!

Ensure the belt position does not interfere with any other device attached to the wheelchair.

Fig.5.8

**Adaptation/Adjustment:****Seat the user into the wheelchair.**

- Feed the buckle on the chest strap through the loop, (Fig.5.9).
- Pull the chest strap so that it fully supports the users upper body firmly and comfortably, (Fig.5.10).
- Ensure that chest strap is not too tight.
- Close the buckle until a loud click is heard, (Fig.5.11).
- The D-Ring facilitates easy adjustment for the user, (Fig.5.12).
- Ensure the chest strap fully supports the user and that there is enough room to fit a hands width between the user and the support. (Fig 5.13).

**Safety:**

In addition to the chest strap latches, the chest strap incorporates two further safety features that prevent the unintended opening of the belt.

Each of the chest strap latches is secured with an overlaying safety Velcro strap. In addition to these, on each side of the chest strap there is a securing strap which is guided through the D-Ring adding another degree of safety.

 **WARNING!**

- Never use the stand-up function without the knee fixation, chest strap and positioning belt attached and safely locked.
- Make sure the rider does not slide in the wheelchair seat. If this occurs, the rider may suffer chest compression or suffocate due to pressure from the supports.
- Never use positioning supports as a patient restraint.
- Never use positioning supports on a user who is comatose or agitated.

Fig.5.9

Fig.5.10

Fig.5.11

Fig.5.12

Fig.5.13

## 6.0 Using the seating system

### 6.1 Getting in and out of the wheelchair (Transferring)

Sunrise Medical recommend that you consult your healthcare professional for assistance in developing your personal front or side transfer technique to best suit your needs and avoid any personal injury.

#### DANGER!

##### To prevent a fall:

- Ensure controller is switched off during transfers to avoid unintentional movement. If you fail to do so you may touch the joystick and cause your system to move when you do not expect it.
- Ensure that the free wheel lever has been set to 'drive',  General manual.
- Rotate the front casters so that they point forward.
- Work with your health care professional to learn all safe methods for transfers.
- Make sure motor brakes are engaged. This keeps the system from moving when you transfer.
- It is dangerous to transfer on your own. It requires good balance and agility. Be aware that there is a point during every transfer when the seating system seat is not below you.
- Learn how to position your body and how to support yourself during a transfer.
- Have someone help you until you are sure you can do a safe transfer on your own.
- Do not stand on the footrest. Putting your full weight on the footrest may lead to wheelchair damage and to hazardous situations, including the wheelchair tipping up.
- Make sure your feet do not "hang up" or get caught in the space between the footrests.

##### Preparing for a forward transfer (Centre Mount Leg Rest):

#### DANGER!

Move your wheelchair as close as possible to the seat you intend to transfer into.

On the centre mount leg rest, flip both foot plates up out of the way, (Fig.6.1 & Fig.6.2 - A).


Fig. 6.1


Fig. 6.2

## Preparing for side transfer:

**DANGER!**

- Move your wheelchair as close as you can to the seat you are transferring to. If possible, use a transfer board.
- Transfer as far back onto the seat surface as you can. This will reduce the risk that you will miss the seat or fall.

### Flip up armrest

- On the side you wish to transfer from; flip the armrest all the way up until it goes into its mechanical stop. This frees your space for side transfer. (Fig. 6.4).

### Remove side cushion

- To remove the side cushion, pull the plastic lever attached to the aluminum seat side rails backwards, (Fig.6 .5 - A).
- Swing the back part of the cushion outwards to an angle of approximately 10°.
- Swing the back part of the cushion outwards to an angle of approximately 10° (Fig 6.6).
- Lift the cushion out from the front pin.

Fig. 6.4

Fig. 6.5

Fig. 6.6

## 6.2 (Powered) seating options

The wheelchair can be equipped by the following powered options to adjust your seating position. These options can be operated by the user without the use of tools.

### Powered adjustable seating options (Fig. 6.7 - Fig. 6.9)

How to operate these powered seating options can be found in  CONTROLLER.

1. Powered tilt function
2. Powered seat elevator
3. Powered backrest recline
4. Powered elevating leg rests

#### DANGER!

- Use of the powered adjustable seating options may influence the stability of the wheelchair.
- Only use the powered adjustable seating options when the wheelchair is on a flat surface! Do Not attempt to operate any power seating option while positioned on an incline.
- Please stay clear of any powered actuator while the component is in motion. Sunrise power mechanisms move over large ranges of travel. Users should be aware of their surroundings while components are in motion.
- While driving the wheelchair in the reclined, elevated, or tilted position use caution.
- Do Not attempt to operate any power seating option while the wheelchair is in motion.
- Do not attempt to operate the lift or tilt power seating options around children.
- Through programming it is possible to reverse the direction of all powered seating functions, ensure you know which direction your seat is going to move before operating.
- Please be aware that the lift/tilt modules present a trap hazard. Make sure that during operation the lift/tilt modules are free from all clothing, hands, feet and other extremities, to prevent injury.
- Do not drive on ramps or slopes with the seating in the elevated position. Before attempting to climb or descend a slope, return the wheelchair to a standard seated drive position, with the backrest upright.

This power system is designed to automatically reduce the maximum driving speed to a creep mode when certain power seating limits are reached. Depending on the configuration of the chair there might be multiple creep speeds.


Fig. 6.7


Fig. 6.8 Powered adjustable seating options


Fig. 6.9 Powered adjustable seating options

## 7.0 Maintenance and cleaning

### 7.1 General Maintenance

- All fasteners should be checked monthly for wear, such as loose bolts or broken components.
- Loose fasteners should be re-tightened according to the installation instructions.
- Check all straps monthly for fraying, ripped seams or other indications of excessive wear damage. Discontinue use if damage is found.
- Chest straps should be replaced at the first indication of damage and/or excessive wear.
- If a broken or loose component is found, discontinue use immediately and contact your authorised Sunrise Medical supplier for replacement.
- Check all Velcro fastening straps for correct adhesion when pressed together.
- Ensure that any contamination, such as fluff, hair, etc is removed from the Velcro straps. Such contamination may affect adhesion.

### 7.2 Cleaning your seating system:

The wheelchair should be wiped over once per week with a slightly damp, not wet, cloth and any fluff or dust that has accumulated around the actuators should be blown or dusted away.

#### CAUTION!

Make sure that you dry all parts of your wheelchair if it becomes wet or damp after cleaning or if it is used in a wet or damp atmosphere.

#### DANGER!

It is important that should the wheelchair be used by more than one person it is cleaned thoroughly to ensure there is no cross infection.

### 7.2.1 Cleaning and inspection of your seating

#### WARNING!

Tears, dents, wearing or slackening of upholstery particularly near to metal could result in poor posture or lower levels of comfort and pressure relief.

#### **Cleaning instructions for seating:**

You can wash all parts of the covers with a gentle-wash detergent at 40°C. You can spin-dry the covers, but do not dry them in a dryer. You can remove all parts of the covers independently of each other and wash them separately. Take out the foam inlays prior to washing and close the Velcro fasteners!

Where upholstery cannot be removed, clean regularly to prevent build up or soiling. Clean with a damp soapy cloth however disinfectants may be used in dilution as specified by their manufacturer. Ensure surfaces are rinsed well with clean water and dried thoroughly.

### 7.3 Hygiene measures when being re-used:

Prior to the wheelchair being re-used, it must be carefully prepared. All surfaces which come into contact with the user must be treated with a disinfection spray.

To do this, you must use a suitable disinfectant for rapid alcohol-based disinfection of medical products and medical devices.

Please take into account the manufacturer's instructions for the disinfectant you are using.

#### CAUTION!

- Do not use solvents, bleaches, abrasives, synthetic detergents, wax polishes or aerosols.
- Disinfectants may be used in dilution as specified by their manufacturer.
- Ensure surfaces are rinsed with clean water and dried thoroughly.

## 7.4 General cleaning instructions:

### CAUTION!

- You should remove stains on your seat as soon as possible.
- After a long period of use, you should clean the cover of your seat with a standard commercially available dry foam cleaning agent.
- You should always clean the whole surface and not just individual spots in order to avoid unattractive edges. The longer you wait the more difficult it will be to remove them.
- Avoid powerful rubbing with aqueous solutions. This will roughen up cover fabric.
- Dirt stains (eg. beer, blood, cola, red wine etc) are best removed with the dry foam cleaning agents or with a mild action detergent.
- Please follow the respective instructions for use when treating the covers.
- The covers should be allowed to dry for at least 48 hours after cleaning with foam cleaning agents.
- Grease stains (eg. ballpoint pen, lipstick, chewing gum etc). Use a proprietary stain remover. Carefully rub the soiled location on the surface only with a clean cloth impregnated with the stain remover.
- Use only a very small quantity of stain remover.
- Test a small inconspicuous area first.
- The upholstery material must not be saturated, since the stain remover might then attack it and destroy it.
- Do not use a hose or a pressure washer to clean your chair.

### WARNING!

- Always read the label on any commercial or domestic cleaning substances.
- Always follow the instructions carefully.

## 7.5 Cleaning Accessories:

- All parts/ accessories such as swing away trays should be cleaned with a damp cloth.
- All lateral supports, headrest, armrests, side- guards, calf pads, lap belts and knee pads should be cleaned with a damp cloth.

### Cleaning controls:

Should the control of your wheelchair become soiled or dirty, it can be wiped with a damp cloth with a dilute disinfectant until clean.

### Cleaning controls:

Should the control of your wheelchair become soiled or dirty, it can be wiped with a damp cloth with a dilute disinfectant until clean.

### DANGER!

- This is important should the wheelchair be used by more than one person to ensure there is no cross infection.
- Ensure that wafer boards, joysticks (all variants), head arrays and switches (all variants) are cleaned with a mild disinfectant and a cleaning cloth to avoid any cross infection possibilities. Following removal from chair regularly wash the sip and puff mouth piece and tube to maintain cleanliness and functionality.
- Ensure control is switched off before cleaning.

## 7.6 Storage:

To keep the upholstery and other parts of the wheelchair in good condition, store in a dry environment as close to room temperature as possible.


### CAUTION!

- Never store your wheelchair in direct sunlight or in a damp/ outdoor environment.
- Direct sunlight might bleach upholstery / plastic parts and cause metal parts to become hot.

## 7.7 Disposal:

The symbols below mean that in accordance with local laws and regulations your product should be disposed of separately from household waste. When this product reaches the end of its life, take it to the local collection point designated by local authorities. The separate collection and recycling of your product at the time of disposal will help conserve natural resources and ensure that it is recycled in a manner that protects the environment.

Ensure you are the legal owner of the product prior to arranging for the product disposal in accordance with the above recommendations and national requirements.

## 8.0 Technical specifications:

### 8.1 Applicable norms / standards

This product complies with the regulations and guidelines for medical aids and carries a CE symbol.

The Sedeo Ergo meets the requirements and standards below. These are checked by independent institutions.

Standard	Definition / description	Test Dummy Weight
EU guideline 93/42 EEC	Applicable as mentioned in Appendix 1	
“NEN-EN 12182 (2014) Class B”	“The requirements from EN12182 (2014) as mentioned in NEN-EN 12184 (2014): Assistive products for persons with disability - General requirements and test methods”	“SEDEO Ergo”
	Q 700 M	160 kg
“NEN-EN 12184 (2014) Class B”	Electrically powered wheelchairs, scooters and their chargers - Requirements and test methods	SEDEO Ergo
	Q 700 M	160 kg
ISO 7176-8	Requirements and test methods for impact, static and fatigue strengths July 1998	N/A
ISO 7176-9	Climate tests for electric wheelchairs	N/A
ISO 7176-14	Requirements and test methods for control systems for electric wheelchairs 1997	N/A
ISO 7176-16	Requirements for resistance to ignition of upholstered parts May 1997	N/A
ISO 7176-19 (2008)	“Wheelchairs - Part 19: Wheeled mobility devices for use as seats in motor vehicles: The Q-series meet the crash test requirements as described in ISO 7176-19”	SEDEO Ergo
	Q 700 M	104 kg

The Sedeo Ergo UP meets the requirements and standards below. These are checked by independent institutions.

Standard	Definition / description	Test Dummy Weight
EU guideline 93/42 EEC	Applicable as mentioned in Appendix 1	
“NEN-EN 12182 (2014) Class B”	“The requirements from EN12182 (2014) as mentioned in NEN-EN 12184 (2014): Assistive products for persons with disability - General requirements and test methods”	“SEDEO Ergo UP”
	Q 700 Up-M	120 kg
“NEN-EN 12184 (2014) Class B”	Electrically powered wheelchairs, scooters and their chargers - Requirements and test methods	SEDEO Ergo UP
	Q 700 Up-M	120 kg
ISO 7176-8	Requirements and test methods for impact, static and fatigue strengths July 1998	N/A
ISO 7176-9	Climate tests for electric wheelchairs	N/A
ISO 7176-14	Requirements and test methods for control systems for electric wheelchairs 1997	N/A
ISO 7176-16	Requirements for resistance to ignition of upholstered parts May 1997	N/A
ISO 7176-19 (2008)	“Wheelchairs - Part 19: Wheeled mobility devices for use as seats in motor vehicles: The Q-series meet the crash test requirements as described in ISO 7176-19”	SEDEO Ergo UP
	Q 700-Up M	75 kg

## 8.2 Technical specifications: SEDEO Ergo seating system

	Metric system values		Imperial system values & alternative	
	Minimum	Maximum	Minimum	Maximum
Seat angle adjustment				
Powered tilt with COG adjustment optional	0° - 30°		0° - 30°	
Powered tilt with COG adjustment optional	-5° - 25°		-5° - 25°	
Powered tilt with COG adjustment optional	0° - 50°		0° - 50°	
Powered tilt with COG adjustment optional	-5° - 45°		-5° - 45°	
Effective seat depth	400 mm	560 mm	15.7"	22.0"
Effective seat width	400 mm	560 mm	15.7"	22.0"
Seat to floor height: Powered tilt & seat elevator (300 mm)				
pos. 1 60 Ah Batteries 0° seat angle	430 mm	730 mm	16.9"	28.7"
pos. 2 60 or 80 Ah batteries 0° seat angle	445 mm	745 mm	17.5"	29.3"
pos. 3 60 or 80 Ah batteries 0° seat angle	460 mm	760 mm	18.1"	29.9"
pos. 4 60 or 80 Ah batteries 0° seat angle	470 mm	770 mm	18.5"	30.3"
pos. 5 60 or 80 Ah Batteries 0° seat angle	480 mm	780 mm	18.9"	30.8"
Armrest to seat distance	240 mm	340 mm	9.4"	13.4"
Backrest angle	90° - 175°		90° - 175°	
Powered recline with anti shear				
Backrest height	560 mm	700 mm	22.0"	27.6"
Legrest angle (powered centermount)	90° - 180°		90° - 180°	
Powered centermount				
Lower leg length	340 mm	520 mm	13.4"	20.5"
Powered centermount				

## 8.3 Technical specifications: SEDEO Ergo UP seating system

	Metric system values		Imperial system values & alternative	
	Minimum	Maximum	Minimum	Maximum
Seat angle adjustment	0°	22°	0°	22°
Powered tilt with COG adjustment				
Effective seat depth	400 mm	540 mm	15.7"	21.0"
Effective seat width	400 mm	500 mm	15.7"	19.7"
Seat to floor height: Powered tilt & seat elevator (300 mm)				
pos. 1 - 60 or 80 Ah batteries	450 mm	750 mm	17.7"	29.5"
pos. 2 - 60 or 80 Ah batteries	465 mm	765 mm	18.3"	30.1"
pos. 3 - 60 or 80 Ah batteries	480 mm	780 mm	18.9"	30.7"
Backrest angle	85°	175°	85°	175°
Powered recline with anti shear				
Backrest height	560 mm	610 mm	22.0"	24.0"
Legrest angle (powered centermount)	90°	180°	90°	180°
Powered centermount				
Lower leg length	320 mm	420 mm	12.6"	16.5"
Powered centermount				

## 9.0 Warranty

### THIS GUARANTEE DOES NOT AFFECT YOUR LEGAL RIGHTS IN ANY WAY.

Sunrise Medical\* provides a guarantee, as set out in the warranty conditions, for products to its customers covering the following.

#### Warranty conditions:

1. Should a part or parts of the product require repair or replacement as a result of a manufacturing and/or material fault within 24 months, then the affected part or parts will be repaired or replaced free of charge. The warranty will only cover manufacturing defects
2. To enforce the warranty, please contact Sunrise Medical Customer Service with the exact details of the nature of the difficulty. Should you be using the product outside the area covered by the Sunrise Medical customer service agent, repairs or replacement will be carried out by another agency as designated by the manufacturer. The product must be repaired by a Sunrise Medical designated Customer Service agent, (dealer).
3. For parts, which have been repaired or exchanged within the scope of this warranty, we provide a warranty in accordance with these warranty conditions for the remaining warranty period for the product in accordance with point 1.
4. For original spare parts which have been fitted at the customer's expense, these will have a 12 months guarantee, (following the fitting), in accordance with these warranty conditions.
5. Claims from this warranty shall not arise, if a repair or replacement of a product or a part is required for the following reasons:
  - a. Normal wear and tear, which include but is not limited to the following parts where fitted; batteries, armrest pads, upholstery, tyres, brakes shoes, ferrules, etc.
  - b. Any overloading of the product, please check the EC label for maximum user weight.
  - c. The product or part has not been maintained or serviced in accordance with the manufacturer's recommendations as shown in the user instructions and/or the service instructions.
  - d. Accessories have been used which are not specified as original accessories.
  - e. The product or part having been damaged by neglect, accident or improper use.
  - f. Changes/modifications have been made to the product or parts, which deviate from the manufacturer's specifications.
  - g. Repairs have been carried out, before our Customer Service has been informed of the circumstances.
6. This guarantee is subject to the law of the country in which the product was purchased from Sunrise Medical"

\* Means the Sunrise Medical facility from which the product was purchased.

#### Additional Notes For Australia Only:

- i. For goods provided by Sunrise medical Pty Ltd in Australia, our goods come with a guarantee by Sunrise Medical that cannot be excluded under Australian Consumer Law.
- ii. You are entitled to a replacement or refund for a major failure and for compensation for any foreseeable loss or damage.
- iii. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.
- iv. The benefits to you given by this warranty are in addition to your other rights and remedies under a law in relation to the goods to which the warranty relates.

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