



SMART E-MOBILITY

Folding Powered Wheelchair

FreedomChair Model A06



**A lightweight and compact model for both indoor and outdoor users up to 17 stones.
Integrated anti-tip feature and dual battery capacity for a range up to 18 miles.**

Features

Seat Width: 17 inches

User Weight Limit: 17 stone / 110kg

Drive Wheels: 10"

Motors: 200W

Batteries: 1 or 2 can be fitted with up to 9 miles/battery

Joystick controlled, FreedomChairs are versatile powered chairs for both outdoor use and, being highly manoeuvrable, excellent indoors too. They fold incredibly easily and quickly and fit into most car boots and take up very little space indoors. They are great travel chairs for days out, cruises etc and a battery test certificate is available for airline travel. With ISO, CE & FDA approval, the chairs are superbly made and finished. FreedomChairs can be fitted with multiple batteries, doubling or even tripling the range up to 27 miles!

Specifications

Net weight: 22.8kg

Maximum weight capacity: 245lbs/110kg

Seat dimensions: 17" x 17"/43cm x 43cm

Maximum speed: Up to 4 mph/6.5 kph

Length: Open: 37.8"/96cm - Folded: 23.2"/58.9cm

Width: Open: 23.2"/58.9cm - Folded: 27.2"/69cm

Height: Open: 36.8"/93.5cm - Folded: 13.4"/34cm

Suspension: Spring absorber

Driving Range:

With 1 battery: up to 9 miles/15km

With 2 batteries: up to 18 miles/30km

Battery capacity: 24V 10AH per battery

Drive wheels: 10" PU tyre

Front wheels: 8" PU tyre

Motors: 2 x 200W brushless motors with electromagnetic braking

Gradient: 8 degrees
