

A200

The A200 is the ideal wheelchair for indoors. Given its small overall width and compact design, it fits easily through narrow doors and corridors and is easy to manoeuvre. The modular design allows for easy transport between the place of work and home, for example.

The A200 is quick to disassemble and fits in almost all cars. Different adjustment options, such as back angle adjustments or seat tilt, allow the A200 to be customised for individual needs.

What is special about me?

- Ideal for indoor use
- Compact and manoeuvrable
- Low overall width of 570mm
- Easy to transport thanks to its modular design and low weight of individual components



A200

Standard Model


- 6km/h
- 28 Ah (C5) batteries
- Driving range approx. 15km
- VR2 control unit
- Control panel holder, standard
- Armrests with clothing protector
- Standard seat, small or large
- Seat cushion
- 12.5" drive wheels, pneumatic tyres
- Frame colour, silver metallic
- 8" Puncture proof casters

Options

- Puncture proof tyres
- Mechanical back angle adjustment
- Swing-away control panel
- Contour seat
- Curb climbing assist
- Assembly kit for headrests

Technical Data

Speed	6km/h
Battery Capacity	28Ah (C5) AGM Batteries
Range under test Conditions	Approx 15km
Control unit	VR2
Max user weight	100kg
Max. obstacle negotiation height	15mm
Climbing ability	12%
Length incl. footplates	1000mm
Width	570mm
Turning radius	870mm
Seat height (without cushion) At 3° seat tilt	430/480mm
Seat width	340 – 480mm
Seat depth	340 – 500mm
Back angle	-9° to +30°
Seat inclination	0° / 3° / 6°
Back height	450 / 500 / 550mm
Armrest height	225 / 350mm
Lower leg length	250 – 470mm (280 – 540mm with contour seat)
Weight when empty	66kg
Charging time	12hr



A200 disassembled



A200 in a car



Control system



Curb climber