



Adventus



Designed for active users

WHAT DO I GET WHEN I ORDER ADVENTUS?

- One Adventus
- User manual
- Two-year warranty

Checklist when ordering

- Check that the wheelchair has a reinforced footrest in $\varnothing 16$, $\varnothing 18$ or $\varnothing 20$ mm tube diameter
- Check that the user has at least 40 mm of space between their feet (on the foot rest)
- Check that the c/c distance between foot-rest tubes is no less than 60 and no greater than 135 mm



Please feel free to contact us on +46 (0)345 408 80 or at info@decon.se



Adventus is an accessory that makes it much easier for you to use your wheelchair over uneven terrain, e.g. in the garden, on the beach, on footpaths. It will also be much more comfortable on normal surfaces, as Adventus gives a smoother ride.

You can easily mount Adventus yourself using a quick connector that clamps around the footrest.

If you want to remove it temporarily, it is possible to hang Adventus on a holder (accessory) that can be fitted on the back of the wheelchair.



TECHNICAL DATA

Length	approx. 450 mm front foot rest
Height	approx. 350 mm from floor
Wheel dimension	10"
Weight	2.45 kg

ITEM NUMBER

Adventus

ADV3100

Black (standard)

ACCESSORIES



Holder

ADV3050

A holder for Adventus is available for fitting on the back of the wheelchair.

Bracket set for folded frame

ADV5000

A new bracket set for folded frame will be available in winter 2018/19

TYRES



	ETRTO	Dimension	Item number
10"	54-152	10 x 2	DM1013

OTHER COLOURS



Adventus can be delivered in other colours. Contact Decon for more info.



 **decon**

Södra Ekeryd 119 314 93 Hyltebruk

Telephone: +46 (0)345 40880 Fax: +46 (0)345 40895 E-mail: info@decon.se www.decon.se

Decon® is not responsible for any changes made after printing.

© All rights are reserved for Decon®. With reservation for printing and typographical errors, discontinued products and misleading images.