

PRIDE
MOBILITY PRODUCTS LTD.


LIVE YOUR BEST®

Travel Scooter Range 2020





"I really like the Zero Turn Technology! It handles like a 3-wheeler with the stability of a 4-wheeler. I really recommend it!"



PRIDE[®]

MOBILITY PRODUCTS LTD.

LIVE YOUR BEST[®]

Pride Mobility Products is a world leader in the design, development and manufacture of mobility products - powerchairs, scooters and lift chairs - for people with disabilities and mobility impairments. With flagship product ranges such as the Jazzy Range, the British manufactured Hudson Range and our innovative Zero Turn technology – Pride Mobility are committed to helping improving people's lives through quality, innovative products.

QUEST



18stone
Weight Capacity

Up To
4mph
Maximum Speed


820mm
Turning Radius

WEIGHT CAPACITY
18 stone (115 kg)
MAX SPEED
4 mph (6.4 kmh)
TURNING RADIUS
820 mm (32.2")
OVERALL LENGTH
930 mm (36.6")
OVERALL WIDTH
485 mm (19.09")
FRONT/REAR TYRES
Front: 8" Rear: 7"
TOTAL WEIGHT WITHOUT BATTERIES
3 stone 6 lbs (23 kg)
BATTERIES REQUIRED
(1) 11.5 AH (Li-ion battery)



Introducing our lightest four wheeled travel scooter yet, at only 4 stone 3 lbs the Aluma Lite is the perfect travel scooter.


WEIGHT CAPACITY
18 stone 13 lbs (120 kg)
MAX SPEED
4 mph (6.4 kmh)
TURNING RADIUS
1385 mm (54.5")
OVERALL LENGTH
1022 mm (40.25")
OVERALL WIDTH
490 mm (19.25")
FRONT/REAR TYRES
Front: 8" Rear: 8"
TOTAL WEIGHT WITHOUT BATTERIES
4 stone 5 lbs (28.8 kg)
BATTERIES REQUIRED
(1) 12 AH (Li-ion battery)



18stone 13lbs
Weight Capacity

Up To
4mph
Maximum Speed

1385mm
Turning Radius

Available on
 Motability



- Red
- Blue
- Black



21 stone 6 lbs
Weight Capacity

Up To
4mph
Maximum Speed

1385mm
Turning Radius

WEIGHT CAPACITY
21 stone 6 lbs (136 kg)
MAX SPEED
4 mph (6.4 kmh)
TURNING RADIUS
1385 mm (54.4")
OVERALL LENGTH
1022 mm (40.2")
OVERALL WIDTH
490 mm (19.2")
FRONT/REAR TYRES
Front: 8" Rear: 8"
TOTAL WEIGHT WITHOUT BATTERIES
5 stone 8 lbs (36.9 kg)
BATTERIES REQUIRED
(2) 12 AH (deep cycle)

GO GO ELITE TRAVELLER®

Available on
 Motability

WEIGHT CAPACITY
21 stone 6 lbs (136 kg)
MAX SPEED
4 mph (6.4 kmh)
TURNING RADIUS
1120 mm (44.125")
OVERALL LENGTH
1003 mm (39.5")
OVERALL WIDTH
495 mm (19.5")
FRONT/REAR TYRES
Front: 7" Rear: 8"
TOTAL WEIGHT WITHOUT BATTERIES
5 stone 6 lbs (36 kg)
BATTERIES REQUIRED
(2) 12 AH (deep cycle)




Blue

Red

21 stone 6 lbs
Weight Capacity

Up To
4 mph
Maximum Speed

1120 mm
Turning Radius

Available on
 Motability



Compact half back,
fold flat seat

Wraparound delta tiller

Disassembles into
5 pieces

High level charging point

LED headlight

Blue

Red

Optional 17 AH
batteries

Ergonomically
designed tiller
adjustment

Seat post
suspension

Off-board battery charging

19stone 9lbs

Weight Capacity

Up To

4mph

Maximum Speed

1120mm

Turning Radius

WEIGHT CAPACITY

19 stone 9 lbs (125 kg)

MAX SPEED

4 mph (6.4 kmh)

TURNING RADIUS

1120 mm (44.125")

OVERALL LENGTH

1003 mm (39.5")

OVERALL WIDTH

495 mm (19.5")

FRONT/REAR TYRES

Front: 7"
Rear: 8"

TOTAL WEIGHT
WITHOUT BATTERIES

5 stone 6 lbs (35.8 kg)

BATTERIES REQUIRED

(2) 12 AH (deep cycle)



Available on
 Motability

WEIGHT CAPACITY
21 stone 6 lbs (136 kg)
MAX SPEED
4 mph (6.4 kmh)
TURNING RADIUS
1162 mm (44.75")
OVERALL LENGTH
1041 mm (41")
OVERALL WIDTH
521 mm (20.5")
FRONT/REAR TYRES
Front: 8" Rear: 8"
TOTAL WEIGHT WITHOUT BATTERIES
6 stone 2 lbs (39.6 kg)
BATTERIES REQUIRED
(2) 12 AH (deep cycle)



- Blue
- Yellow
- Orange
- Black
- Pink
- Red

21 stone 6 lbs
 Weight Capacity

Up To
4 mph
 Maximum Speed

1162mm
 Turning Radius

Available on


GO GO
ELITE TRAVELLER® PLUS



23stone 3lbs
Weight Capacity

Up To
4mph
Maximum Speed

1162mm
Turning Radius

WEIGHT CAPACITY

23 stone 3 lbs (147 kg)

MAX SPEED

4 mph (6.4 kmh)

TURNING RADIUS

3 wheel: 876 mm
4 wheel: 1162 mm

OVERALL LENGTH

3 wheel: 1016 mm
4 wheel: 1026 mm

OVERALL WIDTH

540 mm (21.25")

FRONT/REAR TYRES

Front: 8"
Rear: 9"

TOTAL WEIGHT
WITHOUT BATTERIES

3 wheel: 39.25 kg
4 wheel: 39.5 kg

BATTERIES REQUIRED

(2) 17 AH (deep cycle)

GO GO

ELITE TRAVELLER *SPORT*

Available on
 Motability

WEIGHT CAPACITY
23 stone 3 lbs (147 kg)
MAX SPEED
4 mph (6.4 kmh)
TURNING RADIUS
1283 mm (50.5")
OVERALL LENGTH
1092 mm (43")
OVERALL WIDTH
540 mm (21.25")
FRONT/REAR TYRES
Front: 9" Rear: 9"
TOTAL WEIGHT WITHOUT BATTERIES
6 stone 9 lbs (44 kg)
BATTERIES REQUIRED
(2) 18 AH (deep cycle)



- Blue
- Red


23stone 3lbs
 Weight Capacity

Up To
4mph
 Maximum Speed

1283mm
 Turning Radius



Our patent-pending, intelligent-turning *iTurn Technology™* ensures effortless navigation in small spaces and around corners.

Available on  Motability



- Blue
- Pewter
- White
- Opal
- Green
- Black



21 stone 6 lbs
Weight Capacity

Up To
4mph
Maximum Speed

972mm
Turning Radius

WEIGHT CAPACITY
21 stone 6 lbs (136 kg)
MAX SPEED
4 mph (6.4 kmh)
TURNING RADIUS
972 mm (38.25")
OVERALL LENGTH
1162 mm (45.75")
OVERALL WIDTH
690 mm (27")
FRONT/REAR TYRES
Front: 7" Rear: 9"
TOTAL WEIGHT WITHOUT BATTERIES
8 stone 5 lbs (54 kg)
BATTERIES REQUIRED
(2) 22 AH (deep cycle)

REVO 2.0

Available on


WEIGHT CAPACITY
26 stone 7 lbs (170.09 kg)
MAX SPEED
4 mph (6.4 kmh)
TURNING RADIUS
1320 mm (52")
OVERALL LENGTH
1193 mm (47")
OVERALL WIDTH
533.4 mm (21")
FRONT/REAR TYRES
Front: 9" Rear: 9"
TOTAL WEIGHT WITHOUT BATTERIES
9 stone 4 lbs (60 kg)
BATTERIES REQUIRED
(2) 35 AH (deep cycle)

Fold-flat seat with adjustable armrests

Convenient tiller charging point

LED headlight

Effortless disassembly

Storage compartment

Pneumatic tyres with solid tyres available

True Blue

Grey Street

Monoshock rear suspension

Comfort-Trac Suspension

Alloy style wheels

26stone 7lbs
 Weight Capacity

Up To
4mph
 Maximum Speed

1320mm
 Turning Radius



ACCESSORIES



XLR USB Charger & Phone Mounts

Charge your phone on the go!
This can be used while your powerchair or scooter is in motion.
Compatible with most mobile phones.



Baskets & Bags

Our handy and versatile range of additional large baskets and storage bags make transporting your essentials hassle free.



Rear Oxygen Tank Holder

Available with most models, featuring 2 thumb screws for safe transportation.



Rear View Mirror

Optional rear view mirror for enhanced visibility and rider safety.



Crutch & Walker Holders

Quick and easy access for your cane, crutch or walking stick with our selection of holders.



Foldable Travel Bags

Light and durable travel bags for the Quest, iGo, iGo Plus and iGo Fold.



Cup Holders

Satisfy your thirst on the move with our handy drinks holders, compatible with most of our range of products.



Replacement Keys

We offer a replacement key service.

*Please contact your local Pride dealer for our full range of accessories and further information



PRIDE[®]

MOBILITY PRODUCTS LTD.

LIVE YOUR BEST[®]

Also from Pride:



Powerchair Range



Large Scooter Range



Lift Chair Range

Pride Mobility LTD
 32 Wedgwood Road
 Bicester, Oxon
 OX26 4UL, England

Tel: 01869 324600
 Email: sales@pride-mobility.co.uk
 Web: www.pride-mobility.co.uk



Pride Mobility LTD reserve the right to alter, at any time, without prior notice, the design, specification, packaging or price without the incurring any obligation.