

Avantgarde DV & DS

The Avantgarde DV impresses with its light weight, but this does not mean it has lost any of its proven advantages. The Avantgarde DV's removable and swing-away leg supports lets you drive up close to objects. The leg supports can be folded in or out for this purpose. Removing the leg supports also makes self-propulsion possible.

The Avantgarde DS stands out for its extremely small folding size. At the same time, it offers active wheelchair users' numerous options. With its compact design and integrated leg supports, the Avantgarde DS has the lowest weight from 8.7kg in the Avantgarde series.

Avantgarde DV and DS Standard Model

- Aluminium tube frame with double cross brace
- Back tubes, different height settings possible
- Driving range approx. 15km
- 14 frame colour options
- Foldable, with central double cross brace
- Choice of frame offsets
- Aluminium or Carbon side panels
- Ultralight rims with quick release axle
- Integrated footrest holders
- Calf strap



Avantgarde DV & DS

Options

- Choice of side panels
- Choice of push handles
- Scissor wheel lock
- Shock absorber
- Various push rings
- Variety of back upholsteries
- Choice of caster wheels and forks
- Various rear wheel rims



Fold down push handles



Scissor Wheel locks



Single panel foot rest



Swing away footplate

Technical Data

Seat width	320 – 520mm
Seat depth	360 - 540mm
Front seat height	380 – 540mm
Rear seat height	360 – 510mm
Back height	250 – 550mm
Lower leg length	160 – 510mm
Overall width	490 – 725mm
Overall length	720 – 1080mm
Folding size	From 27mm
Max user weight	140kg
Weight of chair	From approx. 8.7kg
Front frame angle	70° / 75° / 80°