

Avantgarde XXL

The Avantgarde XXL gives outstanding stability combined with easy operation, with a timeless design. The strengths of the Avantgarde XXL are the robust base structure which offers reliability in any situation while the use of high quality materials keeps the weight of the chair as low as possible at approx. 17kg.

This wheelchair is designed for people with the need for a heavy duty chair. The quality and functionality of the Avantgarde XXL makes this chair the most stable folding wheelchair that can hold a maximum weight limit of 180kg and seat width of 62cm.

Standard Model

- Aluminium tube frame with double cross-brace
- Choice of frame offsets
- Single-panel aluminium leg support
- Back with different height setting options, adaptable upholstery
- Aluminium side panels
- Swing-away leg support holders
- Modular design
- Reinforced hollow rims
- Knee lever wheel lock
- Calf band
- Choice of 17 frame colours



Avantgarde XXL

Options

- Back stabiliser bar
- Choice of push rings
- Choice of cushions
- Anti-tipper
- Tip-assist
- Selection of footplates



Attendant break



Backrest stabilisation bar



Individual swing away
footrests



Wheel lock

Technical Data

Approx. weight of chair	17kg
Max user weight	180kg
Back height	250mm – 500mm
Seat width	460mm – 620mm
Seat depth	400mm – 560mm