

The **Bentley** offers up to 20° of tilt from the front of the seat. This allows the user to keep their **feet on the ground** while the chair is fully tilted. The rigid frame provides **superior stability** and maneuverability.



### Exclusive features

- ✓ Foot and hand propel in tilt position with no knee-rise
- ✓ Short wheelbase for maneuverability
- ✓ Simple center of gravity adjustments

### Available options *See order form*

- ✓ Transit tie downs WC-19 & ISO-7176-19
- ✓ 30° manual dynamic recline
- ✓ Power tilt



**Bentley LT** 20° knee pivot tilt chair with a lightweight aluminum frame

**Bentley** 20° knee pivot tilt chair with a rugged steel frame

The lightweight design of the **Bentley LT** combines the benefits of tilt with a high degree of mobility independence.

**YOUR PICK**

The rugged steel frame of the standard **Bentley** makes it a great choice for heavier users or those with high tone/agitation.

## Bentley LT



DYNAMIC TILT  
Up to 20°



WEIGHT CAPACITY  
250 lbs



SEAT WIDTH  
14" to 22"



SEAT DEPTH  
16" to 22"



SEAT HEIGHT  
13" to 20"



BACK SEAT ANGLE  
90° to 122°



OVERALL WIDTH  
+9.5" to 10.5"



CHAIR WEIGHT  
33 lbs

visit [pdgmobility.com/blt](https://pdgmobility.com/blt)  
for more information



## Bentley

visit [pdgmobility.com/bentley](https://pdgmobility.com/bentley)  
for more information



DYNAMIC TILT  
Up to 20°



WEIGHT CAPACITY  
Up to 450 lbs



SEAT WIDTH  
14" to 32"



SEAT DEPTH  
16" to 22"



SEAT HEIGHT  
13" to 20"



BACK SEAT ANGLE  
90° to 122°



OVERALL WIDTH  
+9.5" to 10.5"



CHAIR WEIGHT  
51 lbs

