

# Stayflex™ MD Series Chest Support

Offers firm trunk control with extra freedom in the shoulders. Dual-zone construction features materials that optimally balance stretch and resistance, promoting contact of the thoracic and lumbar spine with the back support while allowing safe, controlled movement. A great solution for users with lordosis and others needing strong sternal support.

Patented, dual-zone stretch construction prevents harness from shifting up toward the neckline, minimizing risk of injury.



MD (multi-directional) Pull Straps can be configured as front-pull or rear-pull at the time of the fitting to meet user needs.



Swivel Buckles on the lower straps and Side-Release Buckles on the top straps facilitate harness removal.



Optional zipper makes removal even easier. Available on both Standard and Narrow Width styles.

## Part Numbers



Size	Standard Width with Zipper	Standard Width without Zipper	Narrow Width with Zipper	Narrow Width without Zipper
XS	-	SH351XS	-	-
S	SH350S	SH351S	-	-
M	SH350M	SH351M	SH360M	SH361M
L	SH350L	SH351L	SH360L	SH361L
XL	SH350XL	SH351XL	SH360XL	SH361XL

Add suffix -B5 or -B6 to end of part number for selected end-fitting; see options at right.

## End-Fitting Options

Includes four end-fittings. Choose between high-strength plastic Flat-Mounts or Cinch-Mounts and add suffix -B5 or -B6 to the end of the part number.

Plastic Flat-Mount™



- B5

Plastic Cinch-Mount™



- B6

**CAUTION:** Whenever an anterior trunk support is used, a properly-adjusted pelvic support should also be worn to prevent sliding.

206.405.4555 or 800.547.5716 | sales@bodypoint.com | [bodypoint.com](http://bodypoint.com)  
View our product video library: [YouTube.com/BodypointInc](https://www.youtube.com/BodypointInc)

© 2021 Bodypoint. All Rights Reserved. BMM327 2021.3

