

# Boogie Family



Mobile and stationary  
standing frame for children and adolescent  
Thanks to the standing device at eye level



The Mobile



The Stationary



The Tilttable

# Boogie Fix

The stationary standing aid



max. user height  
**150 cm**



weight capacity  
**60 kg**



dimensions  
**67 x 54 cm,**  
**73 x 58 cm,**  
**80 x 64 cm**

# Boogie Swing

The mobile standing aid



max. user height  
**150 cm**



weight capacity  
**60 kg**



dimensions  
**67 x 54 cm,**  
**73 x 58 cm,**  
**80 x 64 cm**



centre column tilt  
**Max. 20° or. 30°**



## Overview

- Stationary standing device
- for children and adolescents from 80 cm up to 150 cm user height
- The Boogie concept supports the therapeutic concept of (early) mobility in an upright pose
- All support functions incl. the foot plate can be set in precise steps or even steplessly in 3 dimensions

## Overview

- tiltable standing device
- for children and adolescents from 80 cm up to 150 cm user height
- The centre column can be adjusted steplessly with a gas pressure spring up to 30° forward
- The Boogie Swing supports the therapeutic concept of (early) mobility in an upright pose

# Boogie Drive



## The tiltable mobile standing aid

German aid index N°: 28.29.01.3010



max. user height  
**150 cm**



camber  
**0°, 3°, 6°, 9°, 12°**



frame width  
**36/ 40/ 44cm**  
(+ 4 cm)



centre column tilt  
**Max. 15°**



wheel size  
**24", 26", 28",  
30", 32", 36"**



weight capacity  
**60 kg**



## Overview

- growable standing device
- for children and adolescents from 80 cm up to 150 cm user height
- the centre column can be adjusted in angle steplessly with a gas pressure spring or telescope from 0° to +15°
- the ideal shoulder-driving wheel-position as well as the camber can be precisely set to the physical possibilities of the user

Unique people, unique solutions.

# One name - three models

## The family at a glance

	<i>Boogie Fix</i>	<i>Boogie Swing</i>	<i>Boogie Drive</i>
frame width	54 / 58 / 64	54 / 58 / 64	36 / 40 / 44 (+4 cm)
weight capacity	60 kg	60 kg	60 kg
user height			
size 1	80 - 110 cm	80 - 110 cm	80 - 110 cm
size 2	100 - 130 cm	100 - 130 cm	100 - 130 cm
size 3	120 - 150 cm	120 - 150 cm	120 - 150 cm
tiltable		✓	✓
independently mobile			✓

## Equipment examples

All equipment options are on the order form of the corresponding standing aid




pluggable posterior pad




foot pedal to trigger the gas pressure spring



Adjustable knee truss pads



undivided foot plate




therapy table



tilting adjustment of the centre column



divided foot plate with heel cup



chest truss pads



Posterior truss pad with crank

And much more ...