



Bariatric Rollator Support

Strongly engineered and functionally designed exclusively for bariatric users.

With a weight capacity of 200 kg our Rollator Support is one of the strongest walkers.

It is both functional and has a light and discrete design. It fits through doorways and folds easily for transport.

It includes a large seat, extra strong tubes and 8" massive tires that provide extra security, stability and comfort.

Bariatric Rollator Support has a large support surface. When the user sits down heavily, the back arch prevents that the weight is placed too far – beyond the support surface. In this way, the risk of accidents is prevented significantly.

The critical point for rollators for bariatric users are the brakes. Our Rollator Support has "feather touch" hand brakes that provide superior control and ensure the walker is easy to maneuver and to hold the weight.

Rollator Support comes with a basket and the handles are fitted with reflectors. Rollator Support is available in two sizes.

- Extremely strong and stable lightweight
- Easy to fold and maneuver
- "Feather touch" hand brakes
- Fits through doorways
- The underarm support helps to walk in a more upright position and minimize stress on wrist

MAX 200 kg



Bariatric Rollator Support Specifications

Bariatric Rollator Support	
Maximum weight capacity is	200 kg
Total length	85 cm
Total width	71 cm
Width between handles	61 cm
Height adjustable	100-130 cm (2.5 cm interval)
Seat width	51 cm
Seat depth	35.5 cm
Seat opening	58 cm
Seat height	50 cm
Weight	17 kg
Weight capacity in basket	5 kg
Resting chair	Yes
Foldable	Yes
Massive wheels	Yes
Wheels	8" massive
Frame	Powder-lacquered aluminum
Seat/backrest	Polyurethane
ISO - classification	ISO 12 06 06 - 11

Ideal use for...

- Bariatric users.

Model	Varenummer
XXL-Rehab Rollator Support including basket	0173-050-000