



Bariatric Stand Tall Walker

Strongly engineered and functionally designed exclusively for bariatric users.

Bariatric Stand Tall Walker is a multi-functional aid. It is used for walking and mobilising from seated to standing. Stand Tall Walker is safe to use, gives independence and has an IV option.

Bariatric Stand Tall Walker has two built-in handgrips on the side which allow the user to stand up in small steps from the bed or chair being in full control of the situation. This minimises the risk of falling and therefore also minimises the risk of injury to caregiver or patient himself. Once upright the user can walk around.

The underarm-padding is adjustable in width and angle giving optimal support to the user. This function is supported in full length of the bar.

In addition, it has many handgrip-combinations and easy adjustment, everything to make the Stand Tall Walker easy to maneuver.

The brakes are mounted onto the bar, useful in stops of course, but also in speed control and turning. The brakes can be locked, i.e. when getting up or sitting down.

Using 32 mm tubing ensures an almost forever lasting frame, but also gives a much more simple look and it can be cleaned in few seconds.

IV holder and holder for urine-bag are optional.

Our Bariatric Stand Tall Walker is mainly designed for use in nursing homes and in hospitals.



MAX 325 kg



Product video
Scan the code and see
the product in use

- Adjustable underarm-padding
- Height adjustable

Stand Tall Walker

Specifications

Bariatric Stand Tall Walker	
Maximum weight capacity is	325 kg
Total length	76 cm
Total width	80-92 cm
Internal width between armrests	59-69 cm
Adjustable height (3 cm interval)	100-130 cm
Weight	20 kg
Massive wheels	Yes
Frame	Powder-lacquered chromium molybdenum steel
Hand-grips	PU foam
Underarm-padding	PU foam
Front wheel	100 mm
Rear wheel	125 mm
ISO - classification	ISO 12 06 12 - 06

Model	Width cm	Length cm	Order No
XXL-Rehab Stand Tall	80-92	76	0180-076-000

Accessories			
IV holder for XXL-Rehab Stand Tall	100	27	0180-999-008

Ideal use for.

- Bariatric users in hospitals or nursing homes.

