

Powerful Performance

Exceptional Style



You're active and always on the go. You need a scooter that is up to the task! Appeal to your sporty side with the Colt Executive. Enjoy premium suspension, bright lighting and other sophisticated features, all in one sporty package.



Specifications

Weight Capacity	28 st 8 lbs (400 kg)
Maximum Speed ¹	Up to 8 mph
Ground Clearance ³	1333 mm (5.25")
Turning Radius ³	1973 mm (76.25")
Overall Length ³	1428 mm (56.25")
Overall Width ³	686 mm (27")
Max. Obstacle Climb	120 mm (4.72")
Drive Tyres	330 mm front & rear - pneumatic
Suspension	Front & rear suspension
Total Weight Without Batteries ^{3, 4}	108 kg (238.3 lbs)
Heaviest Piece When Disassembled ³	Main Frame: 90.3 kg (199 lbs)
Drivetrain	24 volt DC motor
Dual Breaking System	Regenerative & electromechanical
Battery Requirements ³	2 x 100 Ah deep cycle
Battery Weight ⁴	27.4 kg (60.5 lbs)
Battery Charger	8-amp off-board
Colour	Titanium Grey, Black & Yellow

Features

- Comfortable high back seat
- Integrated USB charging point for devices
- Wrap around LCD delta tiller
- Hydraulic sealed breaking system
- High level tiller charger
- Front & rear LED lighting package
- Convenient storage in tiller

Also available in:

- Yellow
- Black
- Titanium Grey



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

² Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.