



COMFORTABLE OUTDOOR PERFORMANCE
WITH LUXURIOUS FEATURES



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

Pride®
Mobility Products Ltd.
— LIVE YOUR BEST® —



Tel: 01869 324600

Email: sales@pride-mobility.co.uk

Web: www.pride-mobility.co.uk

(Version 1: April 2017)

Features and Specification

Introduction

The Colt Pursuit features a powerful drivetrain and large 33 cm (13") pneumatic wheels to handle rugged outdoor terrain. Luxury touches like a wraparound delta tiller and a new LCD console combine to give the Colt Pursuit a sleek, sporty style and a maximum level of comfort.



Available on



Available Colours

The Colt Pursuit comes in two optional colours:



Red



Blue



Technical Specifications

Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
Stone	8 Mph (12.9 km/h)	9.2 cm (3.6 inch's)	167.5 cm (66 inch's)	135 cm (53 inch's)	63.5 cm (25 inch's)	Front: 33 cm Rear: 33 cm	15.3 Stone 97.5 k.g.	43 k.g. Front Section	Hi-Back Seat 21.77 k.g	x2 50 Ah Battery

Warranty

1. Life time limited warranty on frame
2. Two year limited warranty on drive train, motor and brake
3. One-year limited warranty Electronic Components, battery and Charger.

Notes

Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. Battery weight may vary based on manufacturer.

Tel: 01869 324600

Email: sales@pride-mobility.co.uk

Web: www.pride-mobility.co.uk

Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST[®] —