

E-Drive Blackline, Electric drive unit, specifications:

- E-Drive, onboard programmable or programmable with Bluetooth Windows laptop/tablet:
- Speed 0-6km/h, braking distance, turning speed, acceleration, engine power
- 3 programs pre-programmed, Soft, Standard, Sport
- Program and joystick individually programmable
- Motors maximum 2 x 120W
- Wheel self-shutdown time 10, 60 minutes or not shutdown
- Weight per wheel 24" only 6.4 kg
- Max. user weight 125 kg, optional 150 kg "Heavy Duty"
- Normal charge about 2-3 hours, Refresh charge 7-8 hours
- Maximum camber 3 degrees
- 3-in-1: drive yourself with pushrims, drive yourself with joystick or drive with attendant control
- Standard USB-Anschluss im Joystick zum Laden von z.B. Handy, 5V, 500 mA maximal



! It is mandatory to mount an adjustable anti-tilt kit when mounting the e-drive!

Company name:
Name of the customer/Number:

Institute: TÜV Süd, Hanover
Certification: EN12184:2014



Mounting for:

Manufacturer:
Model:
Year of construction:

CE Conforms to the EU directive for:
MDR 2017/745 class 1



Crash test according to ISO 7176-19

E-Drive-Plus

- MED E-Drive Blackline wheels with quick-release axle, aluminum push rims, nickel-metal hydride battery, charger, battery bag, joystick holder can be folded away **£ 5295.00**
- MED3* Bracket* for wheelchair ! Fill in the "Mount for" field!
- HD Upcharge for Heavy Duty, up to 150 kg user weight **£ 400.00**

Anti-tipper

- RL-Kipps Wheelchair will have its own anti-tippers left and right supplied by the wheelchair manufacturer
- XOF- Anti-tippers, bracket on the wheel (16",20",22",24") **£ 120.00**
- TB-2000* Anti-tippers with jacking function*, not for all models available, see www.decon.se/en/adapters/ **£ 352.00**
- TBB1

Wheel size selection

- 316** Wheel size 16" Flex insert, no push rims, (No HD version) **£ 200.00**
- 320 Wheel size 20" **£ 200.00**
- 322 Wheel size 22" **£ 200.00**
- 324 Wheel size 24" **£ 300.00**
- 324-M Wheel size 24" with off-road mountain bike tires **£ 300.00**
- 325 Wheel size 25" **£ 200.00**
- 326 Wheel size 26" **£ 200.00**

Choice of push rims/spoke guards/tyres

- Aluminum push rims
- Upcharge stainless steel push rims **£ 64.00**
- Upcharge poly push rims compact **£ 165.00**
- Upcharge push rim foam compact **£ 95.00**
- Upcharge Max Grip push rims **£ 415.00**
- A.Cover Push rim cover for std aluminum and stainless steel **£ 85.00**
- SG Spoke protection, transparent **£ 105.00**
- Tire standard X x 1"
- MED7692- Tires Marathon plus (22", 24", 25", 26") x 1" **£ 95.00**
- DM2412 Tire tread coarse 24 x 1 3/8"
- MS2 Flex insert (20" x 1" – 26 x 1") (standard on 16") **£ 95.00**
- MEM7691- PU tires (20", 22", 24", 26") **£ 60.00**

Selection control

- Without standard joystick (selection of attendant control required)
- Mount the joystick on the right
- Mount the joystick on the left
- MED1200 Upcharge for dual attendant control with bracket **£ 870.00**
- Ø Push handle (01) 19, 20, 22mm (02) 25,26mm (03) 28, 30, 32mm
- MED7660 Upcharge for parallel swivel arm joystick **£ 120.00**
- Detachable joystick bracket on the wheelchair frame mounted, not available for all models
-
- Holder for joystick as attendant control

Programming

- MEM7611 Software with dongle cable and Bluetooth (Dealer Only) **£ 250.00**
- Joystick Soft-mode
- Joystick Standard-mode
- Joystick Sport-/Aktiv-mode

Battery (bag) & charger

- MEDL Upcharge lithium-ion battery with charger **£ 995.00**
- Battery bag NIMH/Li-Ion mounting on the back or seat
- Mount special bag under the seat (not with Li-ion)
- MED12V Converter for charging in boat, caravan or Camper, Output 230V AC 50Hz 150W **£ 230.00**

Notes:

Offer - Order - Print - Save:

*Not for all available models see www.decon.se/en/adapters/ Please contact us if the bracket you need is not listed. ** User weight max. 100 kg.

