

Eco Buggy

The compact Eco Buggy is just right for the quick, easy transportation of children aged 3 to 8 years. The Eco Buggy is especially agile and manoeuvrable.

What is special about me?

- Higher load capacity compared to conventional buggies
- Umbrella folding for fast transportation
- Lightweight and very compact when folded
- Seat depth adjustable
- Strong recline mechanism
- Swing away legrests that remain attached
- Multiple harness options
- Caster locks for improved road steering



Eco Buggy

Standard Model

- 1 size
- Colour, blue
- Lightweight folding buggy with seat and back upholstery
- Back angle adjustment
- Standard lap belt
- Height-adjustable, divisible, swing-away footrest
- Wheel locks
- Dual casters with caster lock
- Folding storage bag

Technical Data

Age	3 – 8 years
Max user weight	50kg
Seat depth	30cm
Seat width	35cm
Backrest height	55cm
Space between seat edge and footrest	20 – 34cm (max length requires caster wheel adjustment)
Folding size (LxWxH)	115 x 42 x 29cm
Wheel Ø front / rear	17 – 21cm
Total weight	9kg
Back angle	90° to 110°

Options

- 5-point belt
- Padding
- Sun canopy
- Buggy board



Individual footplate



Sun canopy



Umbrella folding