

Eco Buggy

The compact Eco Buggy is just right for quick, easy transportation of children aged 3 to 8 years. The Eco Buggy is especially agile and manoeuvrable.

What is special about me?

- Higher load capacity compared to conventional buggies, and with a footrest
- Minimum equipment for fast transportation
- Lightweight and very compact with folded

Standard Model

- 1 size
- Colour: blue
- Lightweight folding buggy with seat and back upholstery
- Back angle adjustment
- Standard lap belt
- Height-adjustable, divisible, swing-away footrest
- Wheel locks
- Dual casters with caster lock
- Folding storage bag

Options

- 5-point belt
- Padding
- Anchor brackets for transportation in a wheelchair accessible vehicle
- Sun canopy
- Buggy board

Quality for life



Technical Data

Age	3 – 8 years
Max. user weight	50kg
Seat depth	30cm
Seat width	35cm
Backrest height	55cm
Space between seat edge and footrest	20-34cm (max length requires caster wheel adjustment)
Folding size (LxWxH)	115x42x29cm
Wheel Ø front/rear	17 – 21cm
Total weight	9kg
Back angle	90° to 110°