

EGO



Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm
Seat height	Front from 40 to 58 cm Rear from 38 to 45 cm (adjustable)
Footplate distance	From 35 to 54 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	Aluminium 85° - 90° - 100° Carbon 85° - 90° - 95°
Total length	88 cm approx. (std conf.)
Total width	Seat width + 17 cm (with 0° camber)
Encumbrance	28 cm approx. (with 0° camber)
Backrest angle	From 78° to 94° to the seat (adjustable)
Camber	0° - 2°
Weight without rear wheels	Starting from 7,2 kg approx.
Wheelchair weight	Starting from 9,8 kg approx.
Maximum load	125 kg