




Date	UK (English)	
User	Company stamp	
End user weight (kg)		Folding frame wheelchair
Language user manual	Responsible	£ 4064,00
		Retail price - VAT not included

Wheelchair measures						
Seat width (LS)	☐ 33	☐ 36	☐ 39	☐ 42	☐ 45	☐ 48
The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1 cm wider.						
Seat depth (PS)	☐ 35	☐ 37,5	☐ 40	☐ 42,5	☐ 45	☐ 47,5
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.						
Footplate distance (DP)	 (cm) ± 0.5 cm - adjustment steps 1.5cm					
(Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm.)						
Setting	☒ Standard ☐ Prudential ☐ Active					
Backrest Height (TS)	Std.Value 33,5 (min 24.5-max 47) adjustment steps 1.5cm					
camber	☒ 0° ☐ 2°					
Front height STD 49						
Rear height STD 42						

Front frame	
☒ Aluminium standard frame Aluminium frame  ☒  ☐ (max 5" caster)  ☐ (max 4" caster)	Colour aluminium frame ☒ Bright silver ☐ Polished black ☐ Opaque black ☐ Bright blue ☐ Bright anthracite ☐ Bright orange ☐ Bright fuchsia ☐ Bright red ☐ Bright lime ☐ Bright blue china ☐ Pearl white

Backrest ☒ Adjustable velcro straps and padded transpiring upholstery

Handle tube ☒ Tubes with short handle ☐ SPORT tubes Without push handles (Lightest configuration) ☐ Folding handles	 02958 252,00
---	--

Casters ☒ 4" Soft ☐ 3" Sport (Lightest configuration) ☐ 5" Soft ☐ 6" Soft (Available only with 100° angle) ☐ 6" Pneumatic (Available only with 100° angle) ☐ 3" Active Soft (wide) (Aluminium alloy rim) ☐ 4" Active Sport (narrow) (Aluminium alloy rim) ☐ 4" Active Soft (wide) (Aluminium alloy rim) ☐ 5" Active Soft (wide) (Aluminium alloy rim) ☐ 5" Pneumatic	<table border="0"> <tr> <td>06824</td> <td>126,00</td> </tr> <tr> <td>01398</td> <td>126,00</td> </tr> <tr> <td>03943</td> <td>126,00</td> </tr> <tr> <td>07518</td> <td>126,00</td> </tr> <tr> <td>10844</td> <td>63,00</td> </tr> </table>	06824	126,00	01398	126,00	03943	126,00	07518	126,00	10844	63,00
06824	126,00										
01398	126,00										
03943	126,00										
07518	126,00										
10844	63,00										

Support front fork ☒ Standard support front fork ☐ Front fork support V Design 2.0 Exclude the installation of : Quick release front forks	 
---	--

Date	Order	User	Address
-------------------	--------------------	-------------------	----------------------

Special Options

☐ Rear wheels plate extension	01950	177,00
☐ lockingband/strip to close/fold the wheelchair	06833	23,00
☐ Drawstring to unlock the back with a single action	06834	126,00

Colour set: Rear wheels and forks - To be matched only with: Rear wheels Ultra-light Progeo - Standard! (22 - 24)
 ☆ Black 22" 24" 25" 26"

☐ Bright silver 22" 24"	00833	126,00
☐ Bright anthracite 24"	00833	126,00
☐ Bright orange 22" 24"	00833	126,00
☐ Bright fuchsia 22" 24"	00833	126,00
☐ Bright blue 22" 24"	00833	126,00
☐ Bright red 22" 24"	00833	126,00
☐ Bright lime 22" 24"	00833	126,00
☐ Pearl white 22" 24"	00833	126,00

Rear wheels
 ☆ 24" Ultra-light Progeo - Standard!
☐ 22" Ultra-light Progeo - Standard!
☐ 25" Ultra-light Progeo - Standard!
☐ 26" Progeo - Standard!

(colour: black)
(colour: black)

☐ 24" ULTRA-LIGHT PROGEO KEY SPORT	04383	876,00
☐ 24" Ultra-light "Spinergy LX"	01685	876,00
☐ 24" Ultra-light "Spinergy LXL" ☐ Black hub and white spokes ☐ Hub and spokes black	07846	876,00
☐ 24" Ultra-light "Spinergy Flex rim" Integrated aluminium/rubber handrim	02270	1.374,00
 ☐ 25" Ultra-light "Spinergy LX"	 02101	 876,00
☐ 25" Ultra-light "Spinergy LXL" ☐ Black hub and white spokes	07847	876,00
☐ 26" Ultra-light "Spinergy LX"	02102	876,00
 ☐ 24" Off-road (aluminium handrim black)	 01392	 224,00

Rear wheel tyre
 ☆ High pressure tyre (22"-24"-25"-26") (Lightest configuration)
☐ High pressure profile tyre (24")

☐ Marathon Plus (24"-25") with anti-puncture protection	01669	63,00
☐ Anti-puncture - Colour Black (24")	11083	126,00
☐ Anti-puncture - Colour Grey (22")	04181	98,00

Handrim	Wheel distance	☆ 3 cm	□ 4 cm		
☆ Aluminium	(22"-24"-25"-26")	(Lightest configuration)			

☐ Aluminium black SOFT-TOUCH (24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170	126,00
☐ Chromed Steel (24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684	148,00
☐ Titanium (24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016	275,00
☐ Ultra Grepp (24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017	315,00
☐ Max Grepp (24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00834	349,00
☐ Ergo Para (24")		00632	349,00
☐ TOP FRICTION (24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023	399,00
Oval countour with a rubber insert			
 ☐ NOVA H (24"-25")	(24" handrim, on 25" wheel is installed with spacers)	 11025	 263,00
Ergonomic black profile			
 ☐ NOVA D (24"-25")	(24" handrim, on 25" wheel is installed with spacers)	 11026	 286,00
Ergonomic profile with rubber coating (ideal for quadriplegics)			
 ☐ SURGE LT (25")	Oval countour with a rubber insert	 06427	 360,00
☐ Quad (24")		00693	200,00
☐ Silicon cover for standard aluminium handrim	Available for measures (22" - 24" - 26")	04050	160,00

Wheel locks
 ☆ Bent lever push to lock

☐ Straight lever push to lock
☐ Straight lever pull to lock

☐ Sport brake




 (Lightest configuration)

06428
 126,00

Date Order User Address

Automatic Closure footplate

☆ Footplate COMPACT, automatic Closure

(Lightest configuration)



11529


no list

Minimum footplate distance 35cm.

RAKE AND DEPTH ADJUSTABLE

☐ Internal☐ External**Footplate - See reference tables N°3**

Warning!! All flip-up footrests, need a stretch calf band and standard fork-support

Weight increase approx. 500gr.☐ Plastic footrest. rake adjustable

Minimum footplate distance 35cm.

☆ 2/3 Internal ☐ 2/3 External

11525


no list

☐ One piece footplate - Flip-up!. rake and depth adjustable

Minimum footplate distance 35cm.

☐ Internal

☆ 2/3 Internal

☐ External

11527

no list

☐ One piece aluminium footplate - Flip-up!. rake and depth adjustable

Minimum footplate distance 35cm.

☐ Internal

☆ 2/3 Internal

☐ 2/3 External

11528

no list

☐ One piece Carbon fibre footplate, rake adjustable, swing away, flip-up.

Minimum footplate distance 35cm.

☆ Internal


☐ External

03291

275,00

Mudguard**Mudguard - FIXED**

☆ Carbon fibre mudguard

☐ Carbon fibre Straight sideguard☐ Extra charge for : REMOVABLE VERSION

08921

126,00

Armrests☐ Tubular armrests

(Swing away, removable and height adjustable armrest)




00694

252,00

Removable armrest with ABS mudguard☐ Desk armrests PRO GEO

(Removable and height adjustable)



09956

224,00

☐ Removable armrest with carbon-fiber straight sideguard

09079

224,00

Accessories☐ Pad for front frame protectors

07058

52,00

☐ Backrest-seat connection cloth

10241

17,00

☐ Removable storage bag (with quick coupling)

09811

126,00

Date Order User Address

Optional 1

☐ Transit wheels			00836	148,00
☐ Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance) ☐ Right ☐ Left		00023 00024	126,00 126,00
☐ Tipping aid	☐ Right ☐ Left		00025 01046	74,00 74,00
☐ Crutches support	☐ Right ☐ Left		00838 00837	74,00 74,00

Optional 2

☐ Tetraclip	(Quick release axle for quadriplegic)		03947	177,00
☐ Safety pelvic belt			03972	63,00
☐ Harness			04103	126,00
☐ Quick release front forks	Available with: Standard support front fork		00030	148,00
☐ Push handles for companion	(Removable and height adjustable)		01048	189,00
☐ Lumbar support	(to insert in the backrest upholstery)		00013	40,00
☐ Swing away lateral support	Not available with Posture backrest PHYSIO ☐ Right ☐ Left		02097 02098	252,00 252,00
☐ Spoke guards Black	(22"-24")		01043	148,00
☐ Spoke guards Transparent	(22"-24")		00835	148,00
☐ Spoke guards Progeo design	(22"-24")		00019	177,00
☐ Spoke guards Dragon	(22"-24")		02106	177,00

Additional Kit

☐ 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)	03378	498,00
☐ 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)	09357	876,00

Date Order User Address




PHYSIO - Posture backrest
(Std. Backrest upholstery not provided)**Removable with quick release**
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

- ☐ **AIR** Carbon shell
not available for seat width 48

BACKREST HEIGHT		28		34		40	
Shell's shape		Lateral retaining wing A=high - R=regular					
7,5 cm		A	R	A	R	A	R
							
		☐	☐	☐	☐	☐	☐

1.047,00

- ☐ **POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape	Lateral retaining wing B=regular		
			
	☐	☐	☐

1.013,00

Configuration DEEP with side support PAD (height adjustable)

- ☐ Aluminum lateral with padded, breathable upholstery.



Small (h 13cm)

☐ Right
10272

29,00

☐ Left
10272

29,00


Large (h 20cm)

☐ Right
10273

40,00

☐ Left
10273

40,00

SWING AWAY - Lateral supportFixing and adjustment hardware complete with protection sleeve
Aluminum lateral with padded, breathable upholstery.

- ☐ Large (18x12)

05647 183,00

- ☐ Small (14x9)

05643 183,00

- ☐ Lateral support extension

Quantity (number pcs).....

10271 69,00

- ☐ Lateral support extension

Quantity (number pcs).....

- ☐ Small (14x9)

05642 183,00

- ☐ Large (18x12)

05646 183,00

- ☐ Headrest adjustable - Double ball mounting system

Available only with backrest Physio
Not available for backrest Physio AIR carbon shell

- ☐ Large (34x14)

05648 252,00

- ☐ Small (28x14)

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts. AL7 1AN United Kingdom - Phone: +44 1484 722888

Date	Order	User	Address
------------	-------------	------------	---------------

Note

Wheelchair : Progeo EGO - Rev. 2023.1	Prices - VAT excluded to be added to the price list	Send this page only if any optional have been chosen	
		☐ standard	(supplied if no requests are indicated)
		☐ on request	(without extra charge)
		☐ optional	(cost to be added at the basic price)

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

REFERENCE TABLES

1		Front height			
With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster					
FORK		Casters			
		3"	4"	5"	6"
SMALL		max 46	max 47	max 48	
MEDIUM		max 50	max 51	max 52	max 54
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.					

2	Rear height			
Rear wheels				
22"	24"	25"	26"	
min 35	min 38	min 39	min 40	
max 43	max 45	max 46	max 48	

3	Footplate option	Footplate distance
	Footplate COMPACT, automatic Closure	min 36
	Plastic footrest	min 36
	One piece footplate - Flip-up!	min 36
	One piece aluminium footplate - Flip-up!	min 36
	One piece Carbon fibre footplate	min 36
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm.)		
NOTE! Lower footplate distances are possible only after our valuation		
Calculated with the best configuration and with the footplates positioned on inside.		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all colors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)