

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888

Date	UK (English)	
User	Company stamp	Folding frame wheelchair
End user weight (kg)		£ 4186,00
Language user manual	Responsible	Retail price - VAT not included

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm						
Seat width (LS)	<input type="checkbox"/> 33	<input type="checkbox"/> 36	<input type="checkbox"/> 39	<input type="checkbox"/> 42	<input type="checkbox"/> 45	<input type="checkbox"/> 48
The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1 cm wider.						
Seat depth (PS)	<input type="checkbox"/> 35	<input type="checkbox"/> 37,5	<input type="checkbox"/> 40	<input type="checkbox"/> 42,5	<input type="checkbox"/> 45	<input type="checkbox"/> 47,5
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.						
Footplate distance (DP) (cm) ±0.5cm - adjustment steps 1.5cm (Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)					
Setting	<input checked="" type="checkbox"/> Standard <input type="checkbox"/> Prudential <input type="checkbox"/> Active					
Backrest Height (TS)	Std. Value 33,5(min 26-max 47) adjustment steps 1.5cm					
camber	<input checked="" type="checkbox"/> 0° <input type="checkbox"/> 2°					
Front height STD 49					
Rear height STD 42					

Front frame		
☆ Aluminium standard frame		
 <input checked="" type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>
☆	(max 5" caster)	(max 4" caster)

Frame colour			
Frame colour (without extra charge)			
<input checked="" type="checkbox"/> Opaque black	<input type="checkbox"/> Pearl white	<input type="checkbox"/> Matt orange	
<input type="checkbox"/> Tekno anthracite	<input type="checkbox"/> Matt Midnight blue	<input type="checkbox"/> Matt grey	<input type="checkbox"/> Matt red
Frame colour on request			
<input type="checkbox"/> Bright silver	<input type="checkbox"/> Bright fuchsia	<input type="checkbox"/> Bright blue china	12579 105,00
<input type="checkbox"/> Bright blue	<input type="checkbox"/> Bright lime	<input type="checkbox"/> Bright red	

Backrest
☆ Adjustable velcro straps and padded transpiring upholstery

Handle tube
☆ Tubes with short handle
<input type="checkbox"/> SPORT tubes Without push handles <i>Lightest configuration</i>
<input type="checkbox"/> Folding handles Not available with: Push handles for companion
 02958 260,00

Casters			
☆ 4" Soft			
<input type="checkbox"/> 3" Sport	<i>Lightest configuration</i>		
<input type="checkbox"/> 5" Soft			
<input type="checkbox"/> 6" Soft	(Available only with 100° angle)		
<input type="checkbox"/> 6" Pneumatic	(Available only with 100° angle)		
<input type="checkbox"/> 3" Active Soft (wide)	(Aluminium alloy rim)	06824	130,00
<input type="checkbox"/> 4" Active Sport (narrow)	(Aluminium alloy rim)	01398	130,00
<input type="checkbox"/> 4" Active Soft (wide)	(Aluminium alloy rim)	03943	130,00
<input type="checkbox"/> 5" Active Soft (wide)	(Aluminium alloy rim)	07518	130,00
<input type="checkbox"/> 5" Pneumatic		10844	65,00

Support front fork
☆ Standard support front fork
<input type="checkbox"/> Front fork support V Design 2.0 Exclude the installation of : Quick release front forks

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

Date	Order	User	Address
------------	-------------	------------	---------------

Special Options		
<input type="radio"/> Rear wheels plate extension	01950	182,00
<input type="radio"/> lockingband/strip to close/fold the wheelchair	06833	24,00
<input type="radio"/> Drawstring to unlock the back with a singleaction	06834	130,00


Colour set: includes rear wheel hub and rim and front forks		
Only standard Progeo 22"/24" wheels have coloured hubs and rims		
☆ Black	22" 24" 25" 26"	
<input type="radio"/> Bright silver	22" 24"	00833 130,00
<input type="radio"/> Bright lime	22" 24"	00833 130,00
<input type="radio"/> Bright blue	22" 24"	00833 130,00
<input type="radio"/> Bright red	22" 24"	00833 130,00
<input type="radio"/> Bright fuchsia	22" 24"	00833 130,00
<input type="radio"/> Pearl white	22" 24"	00833 130,00
<input type="radio"/> Tekno anthracite	22" 24"	00833 130,00
<input type="radio"/> Matt Midnight blue	22" 24"	00833 130,00
<input type="radio"/> Matt grey	22" 24"	00833 130,00
<input type="radio"/> Matt red	22" 24"	00833 130,00

Rear wheels		
☆ 24" Ultra-light Progeo - Standard!		
<input type="checkbox"/> 22" Ultra-light Progeo - Standard!		
<input type="checkbox"/> 25" Ultra-light Progeo - Standard!	(colour: black)	
<input type="checkbox"/> 26" Progeo - Standard!	(colour: black)	
<input type="radio"/> 24" ULTRA-LIGHT PROGEO KEY SPORT	<i>Lightest configuration</i>	04383 902,00
<input type="radio"/> 24" Ultra-light "Spinergy LX"		01685 902,00
<input type="radio"/> 24" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white spokes <input type="checkbox"/> Hub and spokes black	07846 902,00
<input type="radio"/> 24" Ultra-light "Spinergy Flex rim"	Integrated aluminium/rubber handrim	02270 1.415,00
<input type="radio"/> 25" Ultra-light "Spinergy LX"		02101 902,00
<input type="radio"/> 25" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white spokes	07847 902,00
<input type="radio"/> 26" Ultra-light "Spinergy LX"		02102 902,00
<input type="radio"/> 24" Off-road (aluminium handrim black)		01392 231,00

Rear wheel tyre		
☆ High pressure tyre	(22"-24"-25"-26")	<i>Lightest configuration</i>
<input type="checkbox"/> High pressure profile tyre	(24")	
<input type="radio"/> Marathon Plus	(24"-25")	with anti-puncture protection
<input type="radio"/> Anti-puncture - Colour Black	(24")	01669 65,00
Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"		
<input type="radio"/> Anti-puncture - Colour Grey	(22")	04181 101,00

Handrim	Wheel distance	☆ 3 cm	□ 4 cm		
☆ Aluminium	(22"-24"-25"-26")		<i>Lightest configuration</i>		
<input type="radio"/> Aluminium black SOFT-TOUCH	(24"-25")		(24" handrim, on 25" wheel is installed with spacers)	09170	130,00
<input type="radio"/> Chromed Steel	(24"-25")		(24" handrim, on 25" wheel is installed with spacers)	01684	152,00
<input type="radio"/> Titanium	(24"-25")		(24" handrim, on 25" wheel is installed with spacers)	00016	283,00
<input type="radio"/> Ultra Grepp	(24"-25")		(24" handrim, on 25" wheel is installed with spacers)	00017	324,00
<input type="radio"/> Max Grepp	(24"-25")		(24" handrim, on 25" wheel is installed with spacers)	00834	359,00
<input type="radio"/> TOP FRICTION	(24"-25")		(24" handrim, on 25" wheel is installed with spacers)	11023	411,00
Oval countour with a rubber insert					
<input type="radio"/> NOVA H	(24"-25")		(24" handrim, on 25" wheel is installed with spacers)	11025	271,00
Ergonomic black profile					
<input type="radio"/> NOVA D	(24"-25")		(24" handrim, on 25" wheel is installed with spacers)	11026	295,00
Ergonomic profile with rubber coating (ideal for quadriplegics)					
<input type="radio"/> SURGE LT	(25")		Oval countour with a rubber insert	06427	371,00
<input type="radio"/> Quad	(24")			00693	206,00
<input type="radio"/> Silicon cover for standard aluminium handrim			Available for measures (22"- 24" - 26")	04050	165,00

Wheel locks		
☆ Bent lever push to lock		
<input type="checkbox"/> Straight lever push to lock		
<input type="checkbox"/> Straight lever pull to lock		
<input type="radio"/> Sport brake	<i>Lightest configuration</i>	06428 130,00
<input type="radio"/> Brake with extension lever push to lock		01949 65,00
<input type="radio"/> Brake with extension lever pull to lock		03372 65,00



Date Order User Address

Automatic Closure footplate

☆ Footplate COMPACT, automatic Closure

Lightest configuration



11529

Minimum footplate distance 35cm.

RAKE AND DEPTH ADJUSTABLE


Internal

External

Footplate - See reference tables N°3

Warning!! All flip-up footrests, need a stretch calf band and standard fork-support

Weight increase approx. 500gr.



Plastic footrest. rake adjustable

Minimum footplate distance 35cm.

☆ 2/3 Internal 2/3 External




11525

One piece footplate - Flip-up!. rake and depth adjustable

Minimum footplate distance 35cm.

Internal ☆ 2/3 Internal External



11527

One piece aluminium footplate - Flip-up!. rake and depth adjustable

Minimum footplate distance 35cm.

Internal ☆ 2/3 Internal 2/3 External



11528

One piece Carbon fibre footplate, rake adjustable, swing away, flip-up.

Minimum footplate distance 35cm.

☆ Internal External



03291

283,00

Mudguard

Mudguard - FIXED

☆ Carbon fibre mudguard

Carbon fibre Straight sideguard



Extra charge for : REMOVABLE VERSION

08921

130,00

Armrests

Tubular armrests

(Swing away, removable and height adjustable armrest)



00694

260,00

Removable armrest with ABS mudguard

Desk armrests PROGEO

(Removable and height adjustable)



09956

231,00

Removable armrest with carbon-fiber straight sideguard



09079

231,00

Accessories

Pad for front frame protectors

07058

54,00

Backrest-seat connection cloth

10241

18,00

Removable storage bag (with quick coupling)





09811



130,00

Date	Order	User	Address
------------	-------------	------------	---------------

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

Optional 1				
<input type="radio"/> Transit wheels			00836	152,00
<input type="radio"/> Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance)		00023	130,00
	<input type="radio"/> Right		00024	130,00
	<input type="radio"/> Left			
<input type="radio"/> Tipping aid			00025	76,00
	<input type="radio"/> Right		01046	76,00
	<input type="radio"/> Left			
<input type="radio"/> Crutches support			00838	76,00
	<input type="radio"/> Right		00837	76,00
	<input type="radio"/> Left			

Optional 2				
<input type="radio"/> Tetraclip	(Quick release axle for quadriplegic)		03947	182,00
<input type="radio"/> Safety pelvic belt			03972	65,00
<input type="radio"/> Harness			04103	130,00
<input type="radio"/> Quick release front forks	Available with: Standard support front fork		00030	152,00
<input type="radio"/> Push handles for companion	(Removable and height adjustable)		01048	195,00
<input type="radio"/> Lumbar support	(to insert in the backrest upholstery)		00013	41,00
<input type="radio"/> Swing away lateral support	Not available with Posture backrest PHYSIO		02097	260,00
	<input type="radio"/> Right		02098	260,00
	<input type="radio"/> Left			
<input type="radio"/> Spoke guards Black	(22"-24")		01043	152,00
<input type="radio"/> Spoke guards Transparent	(22"-24")		00835	152,00
<input type="radio"/> Spoke guards Progeo design	(22"-24")		00019	182,00
<input type="radio"/> Spoke guards Dragon	(22"-24")		02106	182,00


Additional Kit				
<input type="radio"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)		03378	513,00
<input type="radio"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)		09357	902,00

Date Order User Address

PHYSIO - Posture backrest (Std. Backrest upholstery not provided)


Removable with quick release
Suitable for tubes with diameter 18+26mm
Depth adjustment (4cm)
Angle adjustment ($\pm 18^\circ$)

- AIR** Carbon shell
not available for seat width 48

BACKREST HEIGHT	28	34	40			
Shell's shape 	Lateral retaining wing A=high - R=regular					
	A	R	A	R	A	R
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

1.078,00

- POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52			
Shell's shape 	Lateral retaining wing R=regular					
	R	R	R			
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

1.043,00

Configuration DEEP with side support PAD (height adjustable)

- Aluminum lateral with padded, breathable upholstery.



Small (h 13cm)



Right
10272

30,00



Left
10272

30,00

Large (h 20cm)



Right
10273

41,00




Left
10273

41,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Aluminum lateral with padded, breathable upholstery.



- Large (18x12)

05647 188,00

- Small (14x9)

05643 188,00

- Lateral support extension

Quantity (number pcs).....

- Lateral support extension

Quantity (number pcs)..... 10271 71,00

- Small (14x9)

05642 188,00

- Large (18x12)

05646 188,00

- Headrest adjustable - Double ball mounting system

Available only with backrest Physio Posture



- Large (34x14)

05648 295,00

- Small (28x14)

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888



customer_support_uk@permobil.com - https://www.permobil.com/en-gb

Date	Order	User	Address
------------	-------------	------------	---------------

Note

Wheelchair : Progeo EGO - Rev. 2024.1	Prices - VAT excluded <small>to be added to the price list</small>	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2" style="text-align: center; padding: 2px;">Send this page only if any optional have been chosen</td> </tr> <tr> <td style="padding: 2px;"><input type="checkbox"/> standard</td> <td style="padding: 2px;"><small>(supplied if no requests are indicated)</small></td> </tr> <tr> <td style="padding: 2px;"><input type="checkbox"/> on request</td> <td style="padding: 2px;"><small>(without extra charge)</small></td> </tr> <tr> <td style="padding: 2px;"><input type="checkbox"/> optional</td> <td style="padding: 2px;"><small>(cost to be added at the basic price)</small></td> </tr> </table>	Send this page only if any optional have been chosen		<input type="checkbox"/> standard	<small>(supplied if no requests are indicated)</small>	<input type="checkbox"/> on request	<small>(without extra charge)</small>	<input type="checkbox"/> optional	<small>(cost to be added at the basic price)</small>
Send this page only if any optional have been chosen										
<input type="checkbox"/> standard	<small>(supplied if no requests are indicated)</small>									
<input type="checkbox"/> on request	<small>(without extra charge)</small>									
<input type="checkbox"/> optional	<small>(cost to be added at the basic price)</small>									

REFERENCE TABLES

1		Front height			
With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster					
FORK		Casters			
		3"	4"	5"	6"
SMALL 		max 46	max 47	max 48	
MEDIUM 		max 50	max 51	max 52	max 54
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.					

2		Rear height			
Rear wheels					
		22"	24"	25"	26"
min		min	min	min	min
35		38	39	40	40
max		max	max	max	max
43		45	46	48	48

3		Footplate option	Footplate distance
		Footplate COMPACT, automatic Closure	min 36
		Plastic footrest	min 36
		One piece footplate - Flip-up!	min 36
		One piece aluminium footplate - Flip-up!	min 36
		One piece Carbon fibre footplate	min 36
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)			
NOTE! Lower footplate distances are possible only after our valuation			
Calculated with the best configuration and with the footplates positioned on inside.			

4		Rear wheels
Model	Features	
22" Ultra-light Progeo - Standard!	(not available in anthracite color)	
24" Ultra-light Progeo - Standard!	(available for all colors)	
24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)	
24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)	
24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)	
24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)	
24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)	
25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)	
25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)	
26" Progeo - Standard!	(colour: black)	
26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)	