



# Pull Device

**E | F35**



## ORDER FORM

Billing Address Or Customer-ID		Shipping Address	
C.-ID:		C.-ID:	
Name:		Name:	
Street:		Street:	
Town/PC:		Town/PC:	
E-Mail:		E-Mail:	
Phone:		Phone:	

Identical To Billing Address

Order Information		User Weight	
<input type="checkbox"/> Order	<input type="checkbox"/> Calculation	<input type="checkbox"/> Built-4-Me®	
Purchase Order No.:		<b>Max. user weight: 110 kg</b>	
Art. No.: 477000X	Reimbursement Code: xxx	<b>User Weight:</b>	<b>kg</b>
Order Date: Day	Month	<b>Base Price    £2,500</b>	
	Year		

Technical Data Overview			
<b>Max. user weight</b>	110 kg	<b>Max. speed forward</b>	6 kph; optional 10 kph or 15 kph
<b>Front frame width</b>	260 – 480 mm Frame width measured from center frame tube to center frame tube at the mounting position. Use only the measured dimension and not the seat width	<b>Wheel size</b>	6,5" solid tyre
<b>Total weight with subframe</b>	from 9 kg According to ISO 7176-4.	<b>Battery type</b>	Li-Ion, 36 V, 4 Ah, intelligent battery management system
<b>Total weight with battery</b>	7,5 kg	<b>Brake type</b>	Dual-brake-system: mechanical drum-brake and separate proportional electrical motor brake
<b>Charger</b>	36 V, 2 A	<b>Operating Voltage</b>	36 Volt
<b>Range of action</b>	15 km max	<b>Motor power</b>	6,5" 250 watt (Nominal)
<b>Colour</b>	black	<b>Climbing power</b>	up to 11% (6°)

Base		£
<input checked="" type="checkbox"/>	<b>Pull-device F35</b> 250 W (Nominal) motor, 4 Ah Li-ion battery, folding handles with lock in folded position, forward- and reverse gear, assistive cruise control, dual-brake-system (drum-brake, e-brake), mudguard, front-light	2,500
Ordering/Bundling		£
<input type="checkbox"/>	EMD090213 Without wheelchair order	0
<input type="checkbox"/>	EMD090217 With wheelchair order, not mounted to the wheelchair Wheelchair order on separate order form. Separate shipment of wheelchair and F35.	0
<input type="checkbox"/>	EMD090218 With wheelchair order, mounted to the wheelchair Wheelchair-order on separate order form. Shipment of readily mounted wheelchair/F35-unit	0
Docking		£
<input type="checkbox"/>	EMD010027 Docking for fixed-front wheelchairs	0
<input type="checkbox"/>	EMD010028 Docking for swing-away wheelchairs	0
Docking-Clamps		
Always check installation space for the docking clamps before ordering. Depending on configuration and option, a particular configuration of a chair might not offer enough space.		
<input type="checkbox"/>	EMD010030 Docking-clamps for oval tubes (Quickie Xenon <sup>2</sup> FF, Helium, Nitrum)	0
<input type="checkbox"/>	EMD010031 Docking-clamps for Quickie Xenon <sup>2</sup> SA	0
<input type="checkbox"/>	EMD010029 Docking-clamps for round tube frames (adaptable to 19.5 mm, 23 mm, 25 mm, 28.6 mm, 30mm tubes & QS5 X) Always check combination agreements for non-Sunrise wheelchairs.	0
<input type="checkbox"/>	EMD010062 Docking-clamps for round tube-frames 32 mm	0
Adjustability (Please order the smaller version when a distance is measured between mentioned values)		£
Frame width, from 260mm to 480mm, measured from center frame tube to center frame tube at the mounting position. Use only the measured dimension and not the seat width to check compatibility.		
<input type="checkbox"/>	FW26 Frame width 260 mm	0
<input type="checkbox"/>	FW27 Frame width 270 mm	0
<input type="checkbox"/>	FW28 Frame width 280 mm	0
<input type="checkbox"/>	FW29 Frame width 290 mm	0
<input type="checkbox"/>	FW30 Frame width 300 mm	0
<input type="checkbox"/>	FW31 Frame width 310 mm	0
<input type="checkbox"/>	FW32 Frame width 320 mm	0
<input type="checkbox"/>	FW33 Frame width 330 mm	0
<input type="checkbox"/>	FW34 Frame width 340 mm	0
<input type="checkbox"/>	FW35 Frame width 350 mm	0
<input type="checkbox"/>	FW36 Frame width 360 mm	0
<input type="checkbox"/>	FW37 Frame width 370 mm	0
<input type="checkbox"/>	FW38 Frame width 380 mm	0
<input type="checkbox"/>	FW39 Frame width 390 mm	0
<input type="checkbox"/>	FW40 Frame width 400 mm	0
<input type="checkbox"/>	FW41 Frame width 410 mm	0
<input type="checkbox"/>	FW42 Frame width 420 mm	0
<input type="checkbox"/>	FW43 Frame width 430 mm	0
<input type="checkbox"/>	FW44 Frame width 440 mm	0
<input type="checkbox"/>	FW45 Frame width 450 mm	0
<input type="checkbox"/>	FW46 Frame width 460 mm	0
<input type="checkbox"/>	FW47 Frame width 470 mm	0
<input type="checkbox"/>	FW48 Frame width 480 mm	0

## Speed Setting

£

<input type="checkbox"/>	EMD000025	6 kph	0
<input type="checkbox"/>	EMD000027	10 kph Check your local laws and regulations for max. permitted speed in your country. Speed-settings > 6 km/h can only be ordered via Fax with a signed waiver. Max. speed-settings can only be altered by Sunrise Medical.	0
<input type="checkbox"/>	EMD000028	15 kph Check your local laws and regulations for max. permitted speed in your country. Speed-settings > 6 km/h can only be ordered via Fax with a signed waiver. Max. speed-settings can only be altered by Sunrise Medical.	0

## Display

£

<input checked="" type="checkbox"/>		Display: Battery charging level, current speed, acceleration mode, speed level, odometer, estimated remaining range	0
-------------------------------------	--	---	---

## Accessories

£

<input checked="" type="checkbox"/>		LED front light	0
<input checked="" type="checkbox"/>		Fender	0
<input checked="" type="checkbox"/>	EMD010033	Padding for sub-frame (enclosed)	0
<input checked="" type="checkbox"/>	EMD090206	Depth stops including safety locks for docking tubes	0
<input type="checkbox"/>	EMD090206	Smart Phone Holder	35
<input type="checkbox"/>	EMD090208	Car charging convertor (12V to 230V)	60
<input type="checkbox"/>	EMD110155	Additional charger	130
<input type="checkbox"/>	EMD090165	Theft-protection, cable-lock	20

## Docking-clamps



// EMD010029 Docking-clamps for round tubes frames (adaptable to 19.5 mm, 23 mm, 25 mm, 28.6 mm & 30 mm tubes & QS5 X)

## Accessories



// EMD090206 Smart phone holder



// EMD090208 Car charging convertor (12V to 230V)



// EMD090207 Depth stops including safety locks for docking tubes

## Letter of Acknowledgement

We herewith advise, that using the Empulse F35 with max. speed-settings higher than 6 kph is only permitted on public roads in certain countries and otherwise may only be used on private property.

In the absence of approval, in accordance with the road traffic laws, the Empulse F35 with speed-settings > 6 kph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the Empulse F35 speed settings higher than 6 kph (10/15 kph). When you're travelling with your F35, always check the specific legislation in the respective country.

When placing an order for a max. speed-setting > 6 kph, please confirm with your signature below, that you have understood the above information and that you only intend the usage of your F35 on private property.

**I herewith confirm, that I am aware of the information above, that I will follow strictly the country specific legislation.**

**Date**

**Customer Signature**

Day

Month

Year

## Comments

F35

£

<b>Base Price</b>	<input type="text"/>
Total Of Additional Equipment	<input type="text"/>
Packaging And Transport	<input type="text"/>
Exclusive VAT	<input type="text"/>
VAT	<input type="text"/>
<b>Total Price</b>	<input type="text"/>

Date

Signature

Day

Month

Year

Send E-Mail

