

## WHEELDRIVE

IMPORTANT NOTE FOR CUSTOMERS: MUST SELECT CHOICE OF MOUNTING.

IMPORTANT NOTE FOR SUNRISE CUSTOMER SERVICES:  
WHEN PROCESSING ORDER, SET TO PARTIAL SHIP-NO.



Power assistance *and*  
continuous drive with  
unique dual rim concept



Picture shown may deviate from the actual product.

## ORDER FORM

### Invoice address or customer number:

Customer number:  
Name:  
Street:  
Town: Postal code:  
Country:  
E-mail:  
Tel:

### Delivery address: Deliver to Invoice address

Customer number:  
Name:  
Street:  
Town: Postal code:  
Country:  
E-mail:  
Tel:

### Order information

Order  Calculation  Built-4-Me  
Order date:  
Your reference:

Base price **£5,800** (complete set, ready to go)

## TECHNICAL DATA OVERVIEW

### WheelDrive // WHEELCHAIR POWER ASSIST

Max. load capacity	130 kg
Total width compared to manual wheels (set)	43 mm
Weight of batteries (set)	3.6 kg
Total weight including battery (one wheel)	11.5 kg
Total weight excluding battery (one wheel)	9.7 kg
Wheel diameter	24"
Max. camber	3°

Max. speed forward (assist rim)	10 km/h
Max. speed forward (drive rim)	6 km/h
Max. range (assist rim)	20 km
Max. range (drive rim)	12 km
Max. kerb climbing, driving backwards	50 mm
Charging time (full charge)	Approx. 2 hours

Intended use and environment: For people with limited mobility; self-propelling; indoor and outdoor use. Wheelchair type class B

ORDER FORM WHEELDRIVE

WHEELDRIVE POWER ASSIST WHEELS (including anti-tip, battery and charger)		£
<input checked="" type="checkbox"/> 1019498	WheelDrive with 24" Marathon Plus Evolution tyre, pneumatic, puncture proof, right	2900
<input checked="" type="checkbox"/> 1019499	WheelDrive with 24" Marathon Plus Evolution tyre, pneumatic, puncture proof, left	2900
<b>WheelDrive mounting (Choose applicable mounting set)</b>		
<input type="checkbox"/> 9012381	Mounting set, quick-release axle diameter 12 mm (25-65mm)	180
<input type="checkbox"/> 9012380	Mounting set, quick-release axle diameter 12.7 mm (25-65mm) - Fits all QUICKIE Wheelchairs	
<input type="checkbox"/> 9011416	Mounting set, quick-release axle diameter 12 mm, short	180
<input type="checkbox"/> 9011417	Mounting set, quick-release axle diameter 12.7 mm, short	180

ACCESSORIES		£
<input type="checkbox"/> 1019491	Spoke guards 24", set	350
<input type="checkbox"/> 1018562	Push rim cover 24", silicone, set (for assist rim)	90
<input type="checkbox"/> 1020824	Wheeldrive, battery 3.8 Ah 1pcs	750
<input type="checkbox"/> 1020824	Wheeldrive, battery 3.8 Ah 1pcs	750

**COMMENTS**

WHEELDRIVE	£
<b>Base price</b> (WheelDrive power assist wheels complete with mounting set)	
Option price	
Price excluding VAT	
VAT	
<b>Total price including VAT</b>	

