



From the first concept to final system, every aspect of the Leggero Enzo was designed to maximize both aesthetics and performance.

ENZO™

Shown with Stealth Products TRU-Comfort seating system

- 45° Posterior Tilt ◆
- Optional 5° Anterior Tilt ◆
- Folding Frame with Excellent Rigidity ◆
- Rigid Frame Option ◆
- Numerous Wheel and Caster Choices ◆
- Hand or Foot Tilt Options ◆
- Push to Lock, Hub Lock, or Drum Brakes ◆
- Adjustable Depth Clamping Seat Rails ◆
- Adjustable Height Clamping Back Canes ◆

Facebook.com/
quantumrehabUK



Web: quantumrehab.co.uk



Smooth Tilt Operation



Adjustable Rails



Adjustable Axle Position



Folds Compactly

01869 324600 | sales@quantumrehab.co.uk

Prices and specifications are subject to change without notice. ©2019 Leggero LLC

ENZO™



	Folding & Rigid
Capacity	185 lbs / 84 kg
Seat Depth	11-20" / 28-51 cm
Seat Width	12-18" / 30-46 cm
Seat to Floor	13.5-20" / 34-51 cm
Tilt Angle	-5° to 45°
Hip Angle	90° to 110°

The Enzo's cast aluminum frame and precise machining provide unparalleled rigidity in both folding and rigid configurations.



Impressively Stable Folding Frame



Easy to Operate Foot Tilt Option

For more information visit www.quantumrehab.co.uk or call 01869 324600