



RETAIL

ESSENCE SPP CUSHION

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Email _____

Phone _____ Fax _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

PO No _____

Order Type _____ Quote _____ Order _____

DETAILS

This cushion offers skin protection and positioning with a Coolcore, 4-way stretch, breathable cover, a urethane liner, and an ABS flat rigidiser.

Low profile Tru-Comfort 2, lateral and pre-ischial shelf. Contoured, HD molded foam base, coccyx and sacral relief.



Weight capacities
 Up to 20" - 158.75 kg
 22" and Above - 204.11 kg
 Overall weight: 1.5 kg - 2.2 kg
 Overall height 3.5"
 15x15" up to 22x24"
 Odd sizes in 15"
 Custom sizes available

CUSHION

			ADD WIDTH	ADD DEPTH	
Essence SPP, No Wedge. 15"-20" W.		ST-ESSSPP	<input type="text"/>	<input type="text"/>	£210
Essence SPP, No Wedge. 22" W.		ST-ESSSPP	<input type="text"/>	<input type="text"/>	£265
Essence SPP, Small Wedge. 15"-20" W.		ST-ESSSPP	<input type="text"/>	<input type="text"/>	£265
Essence SPP, Small Wedge. 22" W.		ST-ESSSPP	<input type="text"/>	<input type="text"/>	£295
Essence SPP, Med Wedge. 15"-20" W.		ST-ESSSPP	<input type="text"/>	<input type="text"/>	£265
Essence SPP, Medium Wedge. 22" W.		ST-ESSSPP	<input type="text"/>	<input type="text"/>	£295
Essence SPP, Large Wedge. 15"-20" W.		ST-ESSSPP	<input type="text"/>	<input type="text"/>	£265
Essence SPP, Large Wedge. 22" W.		ST-ESSSPP	<input type="text"/>	<input type="text"/>	£295
Cushion Cover	Coolcore Cover.	TC2SPP-NBLK	<input type="text"/>	<input type="text"/>	STD
	Mesh Cover.	ESSSPP-MBLK	<input type="text"/>	<input type="text"/>	FOC
Fluid Proof Liner.	Polyurethane.	PUL1709080	<input type="text"/>	<input type="text"/>	£48
Extra Covers And Liners.	Coolcore Cover.	ESSSPP-NBLK	<input type="text"/>	<input type="text"/>	£70
	Mesh Cover.	ESSSPP-MBLK	<input type="text"/>	<input type="text"/>	£44
	Polyurethane Liner.	PUL1709080	<input type="text"/>	<input type="text"/>	£48

		Width										
		10	12	14	15	16	17	18	19	20	21	22
Depth	10											
	12											
	14											
	15											
	16											
	17											
	18											
	19											
	20											
	21											
	22											
24												

Select your size, then complete the Part Number by adding the Width first then Depth.

ACM CONTOURED RIGIDISER SUITABLE FOR SLING SEATS

15"W x 16"D. <i>CUS142467-1616.</i>	£56	17"W x 18"D. <i>CUS142466-1818.</i>	£56	19"W x 18"D. <i>CUS142465-2018.</i>	£56
15"W x 17"D. <i>CUS142467-1616.</i>	£56	17"W x 19"D. <i>CUS142466-1818.</i>	£56	19"W x 19"D. <i>CUS142465-2018.</i>	£56
15"W x 18"D. <i>CUS142467.</i>	£56	17"W x 20"D. <i>CUS142466.</i>	£56	19"W x 20"D. <i>CUS142466.</i>	£56
15"W x 19"D. <i>CUS142467.</i>	£56	17"W x 21"D. <i>CUS142466.</i>	£56	19"W x 21"D. <i>CUS142465.</i>	£56
15"W x 20"D. <i>CUS142467.</i>	£56	17"W x 22"D. <i>CUS142466.</i>	£56	19"W x 22"D. <i>CUS142465.</i>	£56
16"W x 16"D. <i>CUS142467.</i>	£56	18"W x 16"D. <i>CUS142466-1816.</i>	£56	20"W x 18"D. <i>CUS142465-2018.</i>	£56
16"W x 17"D. <i>CUS142467-1616.</i>	£56	18"W x 17"D. <i>CUS142466-1816.</i>	£56	20"W x 19"D. <i>CUS142465-2018.</i>	£56
16"W x 18"D. <i>CUS142467.</i>	£56	18"W x 18"D. <i>CUS142466-1818.</i>	£56	20"W x 20"D. <i>CUS142466.</i>	£56
16"W x 19"D. <i>CUS142467.</i>	£56	18"W x 19"D. <i>CUS142466-1818.</i>	£56	20"W x 21"D. <i>CUS142465.</i>	£56
16"W x 20"D. <i>CUS142467.</i>	£56	18"W x 20"D. <i>CUS142466.</i>	£56	20"W x 22"D. <i>CUS142465.</i>	£56
17"W x 16"D. <i>CUS142467-1616.</i>	£56	18"W x 21"D. <i>CUS142466.</i>	£56		
17"W x 17"D. <i>CUS142467-1616.</i>	£56	18"W x 22"D. <i>CUS142466.</i>	£56		

For use on manual wheelchairs only.

NOTES & SPECIAL REQUESTS

Special Order Essence SPP Cushion.

£450

Please detail your special order cushion request.

All special order cushion requests must be approved with the Quantum Custom Department.

Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply.

Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available for sketch details.


FRONT VIEW



RIGHT SIDE VIEW



TOP VIEW



PRICE

TOTAL

2



Stealth Products

32 Wedgwood Road, Bicester,
Oxon, OX26 4UL

Phone: 01869 324600

Fax: 01869 323070

Email: sales@quantumrehab.co.uk

Web: www.quantumrehab.co.uk

ESSENCE SPP - 02/11/23 V1



V1: 14/06/2022