

Date	UK (English)	Folding frame wheelchair
User	Company stamp	
End user weight (kg)		£ 2359,00
Language user manual	Responsible	Retail price - VAT not included

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm	
Seat width (LS)	<input type="checkbox"/> 33 <input type="checkbox"/> 36 <input type="checkbox"/> 39 <input type="checkbox"/> 42 <input type="checkbox"/> 45 <input type="checkbox"/> 48
Seat depth (PS)	<input type="checkbox"/> 35 <input type="checkbox"/> 37,5 <input type="checkbox"/> 40 <input type="checkbox"/> 42,5 <input type="checkbox"/> 45 <input type="checkbox"/> 47,5
Footplate distance (DP) (cm) ±0.5cm - adjustment steps 1.5cm (Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)
Setting	☆ Standard <input type="checkbox"/> Prudential <input type="checkbox"/> Active <input type="checkbox"/> Extreme
Backrest Height (TS)	Std. Value 36(min 30-max 54) adjustment steps 1.5cm
Crossbrace (black)	☆ Single <input type="checkbox"/> Double (Not seat width 33-36)
Front height (STD 49) - Rear height (STD 44)	

Rear frame	
Aluminium frame	<div style="display: flex; align-items: flex-start;"> <div> <p>☆ Aluminium standard frame - Same colour of front frame</p> <p>The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one.</p> <p>camber ☆ 0° <input type="checkbox"/> 1,5° <input type="checkbox"/> 3°</p> </div> </div>
Aluminium frame	<div style="display: flex; align-items: flex-start;"> <div> <p><input type="radio"/> Aluminium frame DYNAMIC 2.0 - (COLOUR: BLACK) 09955 324,00</p> <p>The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.</p> <p>With rake adjustable backrest (±9°) Exclude the installation of : Brake handles for companion (Drum Brake)</p> <p>camber ☆ 0° <input type="checkbox"/> 1,5° <input type="checkbox"/> 2° <input type="checkbox"/> 4°</p> </div> </div>
Titanium frame TMAX	<div style="display: flex; align-items: flex-start;"> <div> <p>Front frame + Titanium frame T MAX 06858 708,00</p> <p>The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one.</p> <p>Brushed Titanium frame - not varnished Weight reduction approx. 450gr</p> <p>camber ☆ 0° <input type="checkbox"/> 1,5° <input type="checkbox"/> 3°</p> </div> </div>

Front frame	
<input type="checkbox"/> Short frame - Ideal with Seat depth (PS) 35 - 37,5 <input type="checkbox"/> Medium frame - Ideal with Seat depth (PS) 40 - 42,5 <input type="checkbox"/> Long frame - Ideal with Seat depth (PS) 45 - 47,5	Frame Angle ☆ 110° <input type="checkbox"/> 100°

Frame colour	
<p>Frame colour (without extra charge)</p> <p>☆ Opaque black <input type="checkbox"/> Tekno anthracite <input type="checkbox"/> Matt grey <input type="checkbox"/> Pearl white <input type="checkbox"/> Matt Midnight blue <input type="checkbox"/> Matt red</p> <p>Frame colour on request</p> <p><input type="radio"/> Bright silver <input type="radio"/> Bright fuchsia <input type="radio"/> Bright blue china <input type="radio"/> Bright blue <input type="radio"/> Bright lime <input type="radio"/> Bright red</p>	<p>12579 105,00</p>

Backrest
☆ Adjustable velcro straps and padded transpiring upholstery

Handle tube
☆ Tubes with standard handle <input type="radio"/> Folding handles 02958 260,00
Not available with: Push handles for companion, Brake handles for companion (Drum Brake)

Date	Order	User	Address
------	-------	------	---------

Colour set: Includes rear wheel hub and rim, rear wheel plates and front forks			
Only standard Progeo 22"/24" wheels have coloured hubs and rims			
☆ Black	22" 24" 25" 26"		
<input type="radio"/> Bright silver	22" 24"	00833	130,00
<input type="radio"/> Bright lime	22" 24"	00833	130,00
<input type="radio"/> Bright blue	22" 24"	00833	130,00
<input type="radio"/> Bright red	22" 24"	00833	130,00
<input type="radio"/> Bright fuchsia	22" 24"	00833	130,00
<input type="radio"/> Pearl white	22" 24"	00833	130,00
<input type="radio"/> Tekno anthracite	22" 24"	00833	130,00
<input type="radio"/> Matt Midnight blue	22" 24"	00833	130,00
<input type="radio"/> Matt grey	22" 24"	00833	130,00
<input type="radio"/> Matt red	22" 24"	00833	130,00







Casters			
☆ 5" Soft			
<input type="checkbox"/> 4" Soft			
<input type="checkbox"/> 3" Sport	<i>Lightest configuration</i>		
<input type="checkbox"/> 6" Soft	(Available only with 110° angle)		
<input type="checkbox"/> 6" Pneumatic	(Available only with 110° angle)		
<input type="checkbox"/> 7" Soft	(Available only with 110° angle)		
<input type="radio"/> 3" Active Soft (wide)	(Aluminium alloy rim)	06824	130,00
<input type="radio"/> 4" Active Sport (narrow)	(Aluminium alloy rim)	01398	130,00
<input type="radio"/> 4" Active Soft (wide)	(Aluminium alloy rim)	03943	130,00
<input type="radio"/> 5" Active Soft (wide)	(Aluminium alloy rim)	07518	130,00
<input type="radio"/> 5" Pneumatic		10844	65,00







Rear wheels			
☆ 24" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 22" Ultra-light Progeo - Standard!			
<input type="checkbox"/> 25" Ultra-light Progeo - Standard!	(colour: black)		
<input type="checkbox"/> 26" Progeo - Standard!	(colour: black)		
<input type="radio"/> 24" ULTRA-LIGHT PROGEO KEY SPORT	<i>Lightest configuration</i>	04383	902,00
<input type="radio"/> 24" Ultra-light "Spinergy LX"		01685	902,00
<input type="radio"/> 24" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white spokes <input type="checkbox"/> Hub and spokes black	07846	902,00
<input type="radio"/> 24" Ultra-light "Spinergy Flex rim"	Integrated aluminium/rubber handrim	02270	1.415,00
<input type="radio"/> 25" Ultra-light "Spinergy LX"		02101	902,00
<input type="radio"/> 25" Ultra-light "Spinergy LXL"	<input type="checkbox"/> Black hub and white spokes	07847	902,00
<input type="radio"/> 26" Ultra-light "Spinergy LX"		02102	902,00
<input type="radio"/> 24" Off-road (aluminium handrim black)		01392	231,00
<input type="radio"/> 22" with brake handles for companion (drum brake) (aluminium handrim - Wheel distance 3cm)	Available only with camber 0° Not available with Frame Dynamic 2.0°	03385	371,00
<input type="radio"/> 24" with brake handles for companion (drum brake) (aluminium handrim - Wheel distance 3cm)	Available only with camber 0° Not available with Frame Dynamic 2.0°	03384	371,00
<input type="radio"/> One arm drive - (24" rear wheels Aluminium handrim High pressure Tyre) - camber 0°	<input type="radio"/> Right <input type="radio"/> Left	02103	413,00

Rear wheel tyre			
☆ High pressure tyre	(22"-24"-25"-26")	<i>Lightest configuration</i>	
<input type="checkbox"/> High pressure profile tyre	(24")		
<input type="radio"/> Marathon Plus	(24"-25")	with anti-puncture protection	01669 65,00
<input type="radio"/> Anti-puncture - Colour Black	(24")		11083 130,00
Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"			
<input type="radio"/> Anti-puncture - Colour Grey	(22")		04181 101,00

Date	Order	User	Address
------	-------	------	---------






Handrim	Wheel distance ☆ 3 cm	□ 4 cm		
☆ Aluminium	(22"-24"-25"-26")	<i>Lightest configuration</i>		
<input type="radio"/> Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170	130,00
<input type="radio"/> Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684	152,00
<input type="radio"/> Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016	283,00
<input type="radio"/> Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017	324,00
<input type="radio"/> Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00834	359,00
<input type="radio"/> TOP FRICTION Oval countour with a rubber insert	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023	411,00
<input type="radio"/> NOVA H Ergonomic black profile	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025	271,00
<input type="radio"/> NOVA D Ergonomic profile with rubber coating (ideal for quadriplegics)	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026	295,00
<input type="radio"/> SURGE LT	(25")	Oval countour with a rubber insert	06427	371,00
<input type="radio"/> Quad	(24")		00693	206,00
<input type="radio"/> Silicon cover for standard aluminium handrim		Available for measures (22"- 24" - 26")	04050	165,00

Wheel locks			
<input checked="" type="radio"/> Bent lever push to lock			
<input type="checkbox"/> Straight lever push to lock			
<input type="checkbox"/> Straight lever pull to lock			
<input type="radio"/> Sport brake			06428 130,00
<input type="radio"/> Brake with extension lever push to lock			01949 65,00
<input type="radio"/> Brake with extension lever pull to lock			03372 65,00





Footplate - See reference tables N°3			
<input checked="" type="radio"/> Plastic footrest, rake adjustable (Minimum footplate distance 33cm.)	☆ 2/3 Internal □ 2/3 External		11525
<input type="checkbox"/> Automatic Closure footplate , rake and depth adjustable (Minimum footplate distance 30cm.)	<input type="checkbox"/> Internal ☆ 2/3 Internal □ External		11526
<input type="radio"/> Minimum footplate distance 25cm. Available only with Multi-function and rotative system -	03945 195,00		
<input type="checkbox"/> One piece footplate - Flip-up! rake and depth adjustable (Minimum footplate distance 30cm.)	<input type="checkbox"/> Internal ☆ 2/3 Internal □ External		11527
<input type="radio"/> Minimum footplate distance 25cm. Available only with Multi-function and rotative system -	03945 195,00		
<input type="checkbox"/> One piece aluminium footplate - Flip-up! rake and depth adjustable (Minimum footplate distance 30cm.)	<input type="checkbox"/> Internal ☆ 2/3 Internal □ 2/3 External □ External		11528
<input type="radio"/> One piece Carbon fibre footplate rake adjustable, swing away, flip-up. (Minimum footplate distance 30cm.)	☆ Internal □ External		03291 283,00
<input type="radio"/> Footplate COMPACT, automatic Closure, rake and depth adjustable Minimum footplate distance 35cm.	☆ Internal □ External		09957 260,00




Mudguard			
Mudguard - FIXED			
<input checked="" type="radio"/> ABS mudguard			
<input type="checkbox"/> ABS Straight sideguard			
<input type="radio"/> Carbon fibre mudguard			03946 324,00
<input type="radio"/> Carbon fibre Straight sideguard			08918 324,00
Mudguard - REMOVABLE			
<input type="radio"/> Removable ABS mudguard			08921 130,00
<input type="radio"/> Removable ABS fibre straight sideguard			08921 130,00
<input type="radio"/> Removable carbon fibre mudguard			08919 448,00
<input type="radio"/> Removable Carbon fibre straight sideguard			08920 448,00

Date	Order	User	Address
------	-------	------	---------

Armrests			
<i>Tubular armrests with standard support combined with anti-tipper wheels exclude the possibility of installing the transit wheels.</i>			
<input type="radio"/> Tubular armrests	(Swing away, removable and height adjustable armrest) <input type="radio"/> Standard Support (Not with frame Dynamic) <input type="radio"/> Steel support with safety locking system		04828 182,00 00694 260,00
Removable armrest with ABS mudguard			
<input type="radio"/> Desk armrests PROGEO	(Removable and height adjustable) Not available with Frame Dynamic 2.0 + camber 4°		09956 231,00
<input type="radio"/> Flip-up armrests PROGEO	(Removable and height adjustable) Not available with Frame Dynamic 2.0		12277 224,00
<input type="radio"/> Removable armrest with ABS straight sideguard			08881 231,00
<input type="radio"/> Removable armrest with carbon-fiber straight sideguard			08882 467,00

Accessories			
<input type="radio"/> Pad for front frame protectors			07058 54,00
<input type="radio"/> Backrest-seat connection cloth			10241 18,00
<input type="radio"/> Removable storage bag (with quick coupling)			09811 130,00

Optional 1			
<input type="radio"/> Transit wheels			00836 152,00
<input type="radio"/> Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance) <input type="radio"/> Right <input type="radio"/> Left		00023 130,00 00024 130,00
<input type="radio"/> Tipping aid	<input type="radio"/> Right <input type="radio"/> Left		00025 76,00 01046 76,00
<input type="radio"/> Crutches support	<input type="radio"/> Right <input type="radio"/> Left		00838 76,00 00837 76,00

Optional 2			
<input type="radio"/> Tetraclip	(Quick release axle for quadriplegic)		03947 182,00
<input type="radio"/> Rear wheels plate extension	(Not with frame Dynamic 2.0)		01950 182,00
<input type="radio"/> Handybike adaptor			00506 76,00
<input type="radio"/> Push handles for companion	(Removable and height adjustable)		01048 195,00
<input type="radio"/> Harness			04103 130,00
<input type="radio"/> Safety pelvic belt			03972 65,00
<input type="radio"/> Tray table (removable)			00508 283,00
<input type="radio"/> Lumbar support	(to insert in the backrest upholstery)		00013 41,00
<input type="radio"/> Super Head rest	(Folding with the wheelchair)		04049 467,00
<input type="radio"/> Swing away lateral support	<input type="radio"/> Right <input type="radio"/> Left		02097 260,00 02098 260,00
<input type="radio"/> Not available with Posture backrest PHYSIO			
<input type="radio"/> Spoke guards Black	(22"-24")		01043 152,00
<input type="radio"/> Spoke guards Transparent	(22"-24")		00835 152,00
<input type="radio"/> Spoke guards Progeo design	(22"-24")		00019 182,00
<input type="radio"/> Spoke guards Dragon	(22"-24")		02106 182,00


Additional Kit			
<input type="radio"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)		03378 513,00
<input type="radio"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)		09357 902,00

Date	Order	User	Address
------------	-------------	------------	---------------


PHYSIO - Posture backrest - Std. Backrest upholstery not provided

Removable with quick release
 Suitable for tubes with diameter 18÷26mm
 Depth adjustment (4cm)
 Angle adjustment (±18°)

- AIR** Carbon shell
not available for seat width 48

BACKREST HEIGHT	28	34	40	1.078,00																	
Shell's shape	Lateral retaining wing A=high - R=regular																				
	<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:50%; text-align: center;">A</td> <td style="width:50%; text-align: center;">R</td> </tr> <tr> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	A	R			<input type="checkbox"/>	<input type="checkbox"/>	<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:50%; text-align: center;">A</td> <td style="width:50%; text-align: center;">R</td> </tr> <tr> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	A	R			<input type="checkbox"/>	<input type="checkbox"/>	<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:50%; text-align: center;">A</td> <td style="width:50%; text-align: center;">R</td> </tr> <tr> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	A	R			<input type="checkbox"/>	<input type="checkbox"/>
A	R																				
<input type="checkbox"/>	<input type="checkbox"/>																				
A	R																				
<input type="checkbox"/>	<input type="checkbox"/>																				
A	R																				
<input type="checkbox"/>	<input type="checkbox"/>																				



- POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52	1.043,00								
Shell's shape	Lateral retaining wing R=regular											
	<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:100%; text-align: center;">R</td> </tr> <tr> <td style="text-align: center;"></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	R		<input type="checkbox"/>	<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:100%; text-align: center;">R</td> </tr> <tr> <td style="text-align: center;"></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	R		<input type="checkbox"/>	<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:100%; text-align: center;">R</td> </tr> <tr> <td style="text-align: center;"></td> </tr> <tr> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	R		<input type="checkbox"/>
R												
<input type="checkbox"/>												
R												
<input type="checkbox"/>												
R												
<input type="checkbox"/>												

Configuration DEEP with side support PAD (height adjustable)


- Alluminum lateral with padded, breathable upholstery.



<p>Small (h 13cm)</p> 	<p>Large (h 20cm)</p> 	<p><input type="checkbox"/> Right 10272 30,00</p> <p><input type="checkbox"/> Left 10272 30,00</p>	<p><input type="checkbox"/> Right 10273 41,00</p> <p><input type="checkbox"/> Left 10273 41,00</p>
---	---	--	--

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
 Alluminum lateral with padded, breathable upholstery.

	<ul style="list-style-type: none"> <input type="radio"/> Large (18x12) 05647 188,00 <input type="radio"/> Small (14x9) 05643 188,00 <input type="radio"/> Lateral support extension Quantity (number pcs)..... 10271 71,00 <input type="radio"/> Lateral support extension Quantity (number pcs)..... <input type="radio"/> Small (14x9) 05642 188,00 <input type="radio"/> Large (18x12) 05646 188,00
--	---

- Headrest adjustable - Double ball mounting system Large (34x14) **05648 295,00**
- Available only with backrest Physio Posture Small (28x14)





Date	Order	User	Address
------	-------	------	---------

Note

Wheelchair : Progeo EXELLE - Rev. 2024.2	Prices - VAT excluded to be added to the price list	<p style="text-align: center; font-size: small;">Send this page only if any optional have been chosen</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%; text-align: center;">☆</td> <td style="width: 70%;">standard</td> <td style="width: 20%; font-size: x-small;">(supplied if no requests are indicated)</td> </tr> <tr> <td style="text-align: center;">□</td> <td>on request</td> <td style="font-size: x-small;">(without extra charge)</td> </tr> <tr> <td style="text-align: center;">○</td> <td>optional</td> <td style="font-size: x-small;">(cost to be added at the basic price)</td> </tr> </table>	☆	standard	(supplied if no requests are indicated)	□	on request	(without extra charge)	○	optional	(cost to be added at the basic price)
☆	standard	(supplied if no requests are indicated)									
□	on request	(without extra charge)									
○	optional	(cost to be added at the basic price)									

REFERENCE TABLES

1	Front height	
With closed frame (95°) max 5" caster		
Casters	Front height	
	min.	max
3"	41	52
4"	42	53
5"	44	54
6"	48	55
7"	50	56
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.		

2	Rear height				
	Rear frame	Rear wheels			
		22"	24"	25"	26"
STANDARD	Straight plate 	min. 37 max 45	min. 39 max 47	min. 40 max 48	min. 41 max 49
	Shaped plate 	min. 36 max 49	min. 38 max 51	min. 39 max 52	min. 40 max 53
	Rear wheels plate extension with Tubular armrests 	min. 37 min. 39 max 46	min. 39 min. 41 max 48	min. 40 min. 43 max 50	min. 41 min. 44 max 51
DYNAMIC 2.0 		min. 36 max 49	min. 38 max 51	min. 39 max 52	min. 40 max 53

3	Footplate option	Footplate distance	
		Frame Angle 110°	Frame Angle 100°
	Plastic footrest	min. 33	min. 31
	Automatic Closure footplate	min. 30	min. 28
	Automatic Closure footplate , with multifunction rotative system	min. 25	min. 23
	One piece footplate - Flip-up!	min. 30	min. 28
	One piece footplate - Flip-up! , with multifunction rotative system	min. 25	min. 23
	One piece aluminium footplate - Flip-up!	min. 30	min. 28
	One piece Carbon fibre footplate	min. 30	min. 28
Calculated with the best configuration and with the footplates positioned on inside.			
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)			

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all coulors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)