

## EXELLE VARIO



<b>Seat width</b>	33 - 36 - 39 - 42 - 45 - 48 cm
<b>Seat depth</b>	35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm
<b>Seat height</b>	Front from 39 to 58 cm (adjustable) Rear from 38 to 51 cm (adjustable)
<b>Footplate distance</b>	From 32 to 51 cm (adjustable)
<b>Backrest height</b>	From 30 to 54 cm (adjustable)
<b>Front frame angle</b>	95° - 110° - 120°
<b>Total length</b>	90 cm approx. (std conf.) 70 cm approx. (without footplates) (std conf.)
<b>Total width</b>	Seat width + 17 cm (with 0° camber)
<b>Encumbrance</b>	26 cm approx. (with 0° camber)
<b>Backrest angle</b>	90° to the seat (with standard rear frame) adjustable from 81° to 99° to the seat (with Dynamic rear frame)
<b>Camber</b>	0° - 1,5° - 3° (with standard rear frame) 0° - 1,5° - 2° - 4° (with Dynamic rear frame)
<b>Weight without rear wheels</b>	Starting from 8,0 kg approx.
<b>Wheelchair weight</b>	Starting from 10,6 kg approx.
<b>Maximum load</b>	125 kg