




BUFFALO

BARIATRIC ROLLATOR

BUFFALO ROLLATOR



The Mobilex Buffalo has a specially designed frame made from aluminium, with a reinforced cross brace making it one of the strongest rollators on the market today. Although very strong, it is also lightweight at just 7.5kgs* and offers a high degree of stability and ease of movement.

The push handles are height-adjustable with easy to read numerical markings and a simple thumb turn lock. It features a textile seat, padded back strap, walking stick holder and a unique locking system that stops it reopening once folded.

The large 200mm (8") smooth-running wheels have soft tyres for enhanced comfort, and the rear wheels are equipped with brakes for optimum speed and movement control. These are also lockable for parking. The mesh bag has a load capacity of 5kgs and can be removed and used as a convenient shopping bag.

Specifications	MB55	MB62
Overall Width	68cm	68cm
Overall Length	66cm	66cm
Product Weight	7.9kgs	8.2kgs
Width Folded	22.5cm	22.5cm
Seat Height	55cm	62cm
Seat Width	56cm	56cm
Height of Push Handles	66cm – 86cm	74cm – 102cm
Max. User Weight	200kgs	200kgs
Wheels	200mm x 35mm	200mm x 35mm
Max. Bag Load	5kgs	5kgs

*Weight of Buffalo MB55

OPTIONAL ACCESSORIES



WATERPROOF BAG



PADDED SEAT



SOFT WHEELS



ONE HANDED BRAKE



SPEED REDUCTION WHEELS

Sold and maintained by: