



LEOPARD

ALUMINIUM ROLLATOR

LEOPARD ROLLATOR



FRAME COLOUR



Red



Blue



Green



Anthracite

55cm only available in anthracite
62cm available in all frame colours

The Mobilex Leopard has a frame made of aluminium, making it lightweight but also very strong. Although light, it offers a high degree of stability coupled with extreme ease of movement.

The push handles are height-adjustable with an easy to read numerical marking and a simple thumb turn lock. It features a textile seat and padded back strap, and a locking system that stops it reopening once folded.

The large 8" smooth-running wheels are equipped with 2 brakes for optimum speed control that are also lockable for parking. The detachable bag has a load capacity of 5kgs and can also be used as a convenient shopping bag.

The Rollator is tested and approved for a user weight of up to 150kgs.

Specifications	MLP55	MLP62
Overall Width	60cm	60cm
Overall Length	66cm	66cm
Product Weight	7.5kgs	7.8kgs
Width Folded	22.5cm	22.5cm
Seat Height	55cm	62cm
Seat Width	49cm	49cm
Height of Push Handles	66cm – 86cm	74cm – 102cm
Max. User Weight	150kgs	150kgs
Wheels	200mm x 35mm	200mm x 35mm
Max. Bag Load	5kgs	5kgs

OPTIONAL ACCESSORIES



WATERPROOF BAG



LARGE BACKSTRAP



PADDED SEAT



SOFT WHEELS



ONE HANDED BRAKE



NON-SLIP TRAY

Sold and maintained by: