



RETAIL

GLACIAL SP CUSHION

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Email _____

Phone _____ Fax _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____


PO No _____

Order Type _____ Quote _____ Order _____

DETAILS

This cushion offers skin protection and a Coolcore 4-way stretch, breathable cover.

New SIM G base, fluid resistant overlay, top and sides.
TPU heat-sealed, low profile lightweight cooling gel topper



Weight capacities
79.37 kg: 10"-14" width
158.75 kg: 15"-21" width,
181.43 kg: 22" width+
Overall weight: 1.5kg
Overall height 2.5"
10x10" up to 24x24"
Odd sizes in 15" up to 21"
Custom sizing not available

CUSHION

		ADD WIDTH		ADD DEPTH		Width												
						10	12	14	15	16	17	18	19	20	21	22	24	
Glacial SP. 10"-20" W.	ST-GLSP	<input type="checkbox"/>	<input type="checkbox"/>		£190	10												
Glacial SP. 22"-24" W	ST-GLSP	<input type="checkbox"/>	<input type="checkbox"/>		£210	12												
Cushion Cover.	Coolcore.	GLSP-NBLK	<input type="checkbox"/>	<input type="checkbox"/>	STD	14												
	Mesh.	GLSP-MBLK	<input type="checkbox"/>	<input type="checkbox"/>	FOC	15												
	Stealth Tek.	GLSP-TBLK	<input type="checkbox"/>	<input type="checkbox"/>	FOC	16												
Fluid Proof Liner.	Urethane.	URL1709070	<input type="checkbox"/>	<input type="checkbox"/>	£40	17												
	Polyurethane.	PUL1709030	<input type="checkbox"/>	<input type="checkbox"/>	£48	18												
Extra Covers And Liners.	Coolcore Cover.	GLSP-NBLK	<input type="checkbox"/>	<input type="checkbox"/>	£70	19												
	Mesh Cover.	GLSP-MBLK	<input type="checkbox"/>	<input type="checkbox"/>	£44	20												
	Stealth Tek Cover.	GLSP-TBLK	<input type="checkbox"/>	<input type="checkbox"/>	£50	21												
	Urethane Liner.	URL1709070	<input type="checkbox"/>	<input type="checkbox"/>	£40	22												
	Polyurethane Liner.	PUL1709030	<input type="checkbox"/>	<input type="checkbox"/>	£48	23												
						24												

Select your size, then complete the Part Number by adding the Width first then Depth.

ADD SOLID BASE

ABS FLAT RIGIDISER

10"W x 10"D. <i>DWR1010E188A.</i>	£34	16"W x 20"D. <i>DWR1010E191K.</i>	£34	20"W x 20"D. <i>DWR1010E193K.</i>	£34
10"W x 12"D. <i>DWR1010E188C.</i>	£34	17"W x 17"D. <i>DWR1010E234H.</i>	£34	20"W x 22"D. <i>DWR1010E193M.</i>	£34
10"W x 14"D. <i>DWR1010E188E.</i>	£34	17"W x 18"D. <i>DWR1010E234H.</i>	£34	22"W x 16"D. <i>DWR1010E193G.</i>	£34
12"W x 10"D. <i>DWR1010E189A.</i>	£34	17"W x 19"D. <i>DWR1010E234J.</i>	£34	22"W x 18"D. <i>DWR1010E194I.</i>	£34
12"W x 12"D. <i>DWR1010E189C.</i>	£34	17"W x 20"D. <i>DWR1010E234H.</i>	£34	22"W x 20"D. <i>DWR1010E194K.</i>	£34
12"W x 14"D. <i>DWR1010E189E.</i>	£34	18"W x 16"D. <i>DWR1010E192G.</i>	£34	22"W x 22"D. <i>DWR1010E194M.</i>	£34
14"W x 10"D. <i>MEC1804152-1410.</i>	£34	18"W x 18"D. <i>DWR1010E192I.</i>	£34	22"W x 24"D. <i>DWR1010E193O.</i>	£34
14"W x 12"D. <i>DWR1010E190C.</i>	£34	18"W x 20"D. <i>DWR1010E192K.</i>	£34	24"W x 18"D. <i>DWR1010E236I.</i>	£34
14"W x 14"D. <i>DWR1010E190E.</i>	£34	18"W x 22"D. <i>DWR1010E192M.</i>	£34	24"W x 19"D. <i>DWR1010E236I.</i>	£34
14"W x 16"D. <i>DWR1010E190G.</i>	£34	19"W x 17"D. <i>MEC1804152-1917.</i>	£34	24"W x 20"D. <i>DWR1010E236K.</i>	£34
15"W x 15"D. <i>MEC1804152-1515.</i>	£34	19"W x 18"D. <i>MEC1804152-1917.</i>	£34	24"W x 21"D. <i>DWR1010E236K.</i>	£34
15"W x 16"D. <i>MEC1804152-1515.</i>	£34	19"W x 19"D. <i>MEC1804152-1919.</i>	£34	24"W x 22"D. <i>DWR1010E236M.</i>	£34
15"W x 17"D. <i>DWR1010E232H.</i>	£34	19"W x 20"D. <i>MEC1804152-1919.</i>	£34	24"W x 23"D. <i>DWR1010E236M.</i>	£34
16"W x 16"D. <i>DWR1010E191G.</i>	£34	20"W x 16"D. <i>DWR1010E193G.</i>	£34	24"W x 24"D. <i>DWR1010E236O.</i>	£34
16"W x 18"D. <i>DWR1010E191I.</i>	£34	20"W x 18"D. <i>DWR1010E193I.</i>	£34		

ACM CONTOURED RIGIDISER SUITABLE FOR SLING SEATS

16"W x 16"D. <i>CUS142467-1616.</i>	£56	17"W x 20"D. <i>CUS142466.</i>	£56	19"W x 19"D. <i>CUS142465-2018.</i>	£56
16"W x 18"D. <i>CUS142467.</i>	£56	18"W x 16"D. <i>CUS142466-1816.</i>	£56	19"W x 20"D. <i>CUS142465.</i>	£56
16"W x 20"D. <i>CUS142467.</i>	£56	18"W x 18"D. <i>CUS142466-1818.</i>	£56	20"W x 18"D. <i>CUS142465-2018.</i>	£56
17"W x 17"D. <i>CUS142466-1816.</i>	£56	18"W x 20"D. <i>CUS142466.</i>	£56	20"W x 20"D. <i>CUS142465.</i>	£56
17"W x 18"D. <i>CUS142466-1818.</i>	£56	18"W x 22"D. <i>CUS142466.</i>	£56	20"W x 22"D. <i>CUS142465.</i>	£56
17"W x 19"D. <i>CUS142466-1818.</i>	£56	19"W x 18"D. <i>CUS142465-2018.</i>	£56		

For use on manual wheelchairs only.

NOTES & SPECIAL REQUESTS

PRICE

TOTAL

2



Stealth Products

32 Wedgwood Road, Bicester,
Oxon, OX26 4UL

Phone: 01869 324600

Fax: 01869 323070

Email: sales@quantumrehab.co.uk

Web: www.quantumrehab.co.uk

GLACIAL SP - 15/01/24 V1



V1: 14/06/2022