

Enhance Your Independence!



GO CHAIR



The Go Chair® is #1 for portability! With feather-touch disassembly, you can take your Go Chair to all your favorite places and enjoy greater independence.

PRIDE[®]
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST*

www.pride-mobility.co.uk

Specifications


Weight Capacity	21 st 6 lbs (136 kg)
Maximum Speed ¹	Up to 3.7 mph
Ground Clearance ³	42 mm (1.6")
Turning Radius ³	6451 mm (25.4")
Overall Length ³	8382 mm (33") with front rigging
Overall Width ³	5740.4 mm (22.6")
Max. Obstacle Climb	30 mm (1.18")
Drive Tyres	228.6 mm (9") - solid
Suspension	Limited
Total Weight Without Batteries ^{3, 5}	30.39 kg (66.9 lbs)
Heaviest Piece When Disassembled ³	16.32 kg (35.9 lbs)
Drivetrain	Rear-wheel drive, two motor
Dual Breaking System	"Intelligent Braking" electronic regenerative, disc park brake
Battery Requirements ³	2 x 18 Ah batteries 2 x 22 Ah batteries
Battery Weight ⁴	14.06 kg (31 lbs)
Battery Charger	18 Ah: 2-amp, off-board 22 Ah: 3.5-amp, off-board
Colour	Ruby red, amber orange, rose quartz, sapphire blue, citrine yellow, pearl white

Features

- Disassembles into five pieces
- Large comfortable foot plate
- Height & width adjustable armrests
- Dual swing-away under seat storage
- Five point swivel seat with latch release

Also available in:

- Ruby Red
- Amber Orange
- Rose Quartz
- Sapphire Blue
- Citrine Yellow
- Pearl White



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

² Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.