

The Quickest Way To Travel



GO GO
ELITE TRAVELLER® PLUS



It's perfect for traveling, PLUS so much more! The Go-Go Elite Traveller® Plus offers easy disassembly and includes an ergonomic tiller for one-handed operation. Today's fast-paced world is no match for this stylish and portable scooter.

PRIDE[®]
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST*

www.pride-mobility.co.uk

Specifications


	4-Wheel	3-Wheel
Weight Capacity	23 st 3 lbs (147 kg)	23 st 3 lbs (147 kg)
Maximum Speed ¹	Up to 4 mph	Up to 4 mph
Ground Clearance ³	38 mm (1.5")	38 mm (1.5")
Turning Radius ³	1162 mm (45.75")	876 mm (34.5")
Overall Length ³	1023 mm (40.25")	1016 mm (40")
Overall Width ³	540 mm (21.25")	540 mm (21.25")
Max. Obstacle Climb	65 mm (2.5")	65 mm (2.5")
Drive Tyres	Front: 203 mm (8") - solid Rear: 229 mm (9") - solid	Front: 203 mm (8") - solid Rear: 229 mm (9") - solid
Suspension	Limited	Limited
Total Weight Without Batteries ^{3, 4}	39.5 kg (87 lbs)	39.2 kg (86.5 lbs)
Heaviest Piece When Disassembled ³	15.9 kg (35 lbs)	14.7 kg (32.5 lbs)
Drivetrain	24 volt DC motor	24 volt DC motor
Dual Breaking System	Regenerative & electromechanical	Regenerative & electromechanical
Battery Requirements ³	2 x 18 Ah deep cycle	2 x 18 Ah deep cycle
Battery Weight ⁴	18 Ah - 13.2 kg (29 lbs)	18 Ah - 13.2 kg (29 lbs)
Battery Charger	2 amp off-board	2 amp off-board
Colour	Red & blue	Red & blue

Features

- Feather-touch disassembly with one hand
- Front and rear LED lighting
- Wrap around delta tiller
- Seat post suspension
- Available in both 3 & 4 wheel
- Exclusive & stylish low profile tyres
- Convenient front basket
- Disassembles into 5 lightweight pieces

Also available in:

- Red
- Blue



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

² Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.