

Sleek, Sporty Undeniably Stylish



GO GO
ELITE TRAVELLER *SPORT*



It's time to make your mark on this world! The Go-Go® Elite Traveller Sport offers simple disassembly and an ergonomic delta tiller, perfect for one-handed operation. Enjoy the versatility of this sporty scooter, as well as bright LED lighting.

PRIDE[®]
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST*

www.pride-mobility.co.uk

Specifications

| | |
|--|--|
| Weight Capacity | 23 st 3 lbs (147 kg) |
| Maximum Speed ¹ | Up to 4 mph |
| Ground Clearance ³ | 38 mm (1.5") at motor |
| Turning Radius ³ | 1283 mm (50.5") |
| Overall Length ³ | 1092 mm (43") |
| Overall Width ³ | 540 mm (21.25") |
| Max. Obstacle Climb | 50mm (2") |
| Drive Tyres | 229 mm front & rear - solid |
| Suspension | - |
| Total Weight Without Batteries ^{3, 4} | 44 kg (100 lbs) |
| Heaviest Piece When Disassembled ³ | 20 kg (44 lbs) |
| Drivetrain | 24 volt DC motor |
| Dual Breaking System | Regenerative & electromechanical |
| Battery Requirements ³ | 2 x 18 Ah deep cycle or 2 x 20 Ah deep cycle |
| Battery Weight ⁴ | 18 Ah - 12.5 kg (27.4 lbs) 20 Ah - 20 kg (44 lbs) |
| Battery Charger | Off-board |
| Colour | Red & blue |

Features

- Feather-touch disassembly with one hand
- Front and rear LED lighting
- Wrap around delta tiller
- Extra 2" of space for added leg room
- Includes 2 sets of interchangeable shrouds
- Exclusive & stylish low profile tyres
- Convenient front basket

Available in:

- Red
- Blue



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

² Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.