



Helio C2 ORDER FORM - UK

P.O. number:	Date:	Ship to:
Company Name:	Phone:	Name:
Purchasing contact:		Address:
Clinician/Sales rep:		City:
Email:		Country:
Marked for:	* Required	Phone:
		VAT not included

MA083 English Manual (STD)

MODEL								GBP			
1	WC05	Ultralight Folding Carbon Fiber Helio C2 Wheelchair (120 kg - 265 lb weight capacity)						£3 250	<input type="checkbox"/>		
2	WC09	HD Heavy Duty package (159 kg - 350 lb weight capacity)						£395	<input type="checkbox"/>		
3	WC10	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber						£35	<input type="checkbox"/>		
SEAT WIDTH ³ - cm (inches)											
4	SW	Seat width:	35,6 (14") <input type="checkbox"/>	38,1 (15") <input type="checkbox"/>	40,6 (16") <input type="checkbox"/>	43,2 (17") <input type="checkbox"/>	45,7 (18") ¹ <input type="checkbox"/>	48,3 (19") ¹ <input type="checkbox"/>	NC	<input type="checkbox"/>	
5	SW	Seat width:	50,8 (20") ¹ <input type="checkbox"/>			53,3 (21") ² <input type="checkbox"/>		55,9 (22") ² <input type="checkbox"/>		£350	<input type="checkbox"/>
SEAT DEPTH											
6	SD	Seat depth:	35,6 (14") <input type="checkbox"/>	38,1 (15") <input type="checkbox"/>	40,6 (16") <input type="checkbox"/>	43,2 (17") <input type="checkbox"/>	45,7 (18") <input type="checkbox"/>	48,3 (19") <input type="checkbox"/>	NC	<input type="checkbox"/>	
7	SD	Seat depth:	50,8 (20") <input type="checkbox"/>						£350	<input type="checkbox"/>	

1 - Available with HD option

2 - Only available with HD option

[3. Width calculator](#)

NEWTON ULTRALIGHT PACKAGE										
8	PK01	Includes: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels						Save £???	£???	<input type="checkbox"/>

C2 FRONT Seat-to-floor Height - Please select your choice

Forks / Position		7,62cm (3") caster			10,16cm (4") caster			12,7cm (5") caster			15,24cm (6") caster			17,78cm (7") caster			20,32cm (8") caster			
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			
		STD	+2,54 cm 1"	+5,08cm 2"	STD	+2,54 cm 1"	+5,08cm 2"	STD	+2,54 cm 1"	+5,08cm 2"	STD	+2,54 cm 1"	+5,08cm 2"	STD	+2,54 cm 1"	+5,08cm 2"	STD	+2,54 cm 1"	+5,08cm 2"	
4" 10,2 cm	P1	33cm / 13"	35,6cm / 14"	38,1cm / 15"	34,3cm / 13,5"	36,8cm / 14,5"	39,4cm / 15,5"													
	P2	34,3cm / 13,5"	36,8cm / 14,5"	39,4cm / 15,5"	35,6cm / 14"	38,1cm / 15"	40,6cm / 16"	36,8cm / 14,5"	39,4cm / 15,5"	41,9cm / 16,5"										
5" 12,7 cm	P1																			
	P2				35,6cm / 14"	38,1cm / 15"	40,6cm / 16"	36,8cm / 14,5"	39,4cm / 15,5"	41,9cm / 16,5"										
	P3	35,6cm / 14"	38,1cm / 15"	40,6cm / 16"	36,8cm / 14,5"	39,4cm / 15,5"	41,9cm / 16,5"	38,1cm / 15"	40,6cm / 16"	43,2cm / 17"	39,4cm / 15,5"	41,9cm / 16,5"	44,5cm / 17,5"							
7" 17,8 cm	P4	36,8cm / 14,5"	39,4cm / 15,5"	41,9cm / 16,5"	38,1cm / 15"	40,6cm / 16"	43,2cm / 17"	39,4cm / 15,5"	41,9cm / 16,5"	44,5cm / 17,5"	40,6cm / 16"	43,2cm / 17"	45,7cm / 18"	41,9cm / 16,5"	44,5cm / 17,5"	47cm / 18,5"				
	P1										41,9cm / 16,5"	44,5cm / 17,5"	47cm / 18,5"	43,2cm / 17"	45,7cm / 18"	48,3cm / 19"				
	P2								41,9cm / 16,5"	44,5cm / 17,5"	47cm / 18,5"	43,2cm / 17"	45,7cm / 18"	48,3cm / 19"	44,5cm / 17,5"	47cm / 18,5"	49,5cm / 19,5"			
	P3								43,2cm / 17"	45,7cm / 18"	48,3cm / 19"	44,5cm / 17,5"	47cm / 18,5"	49,5cm / 19,5"	45,7cm / 18"	48,3cm / 19"	50,8cm / 20"	47cm / 18,5"	49,5cm / 19,5"	52,1cm / 20,5"
	P4	41,9cm / 16,5"	44,5cm / 17,5"	47cm / 18,5"	43,2cm / 17"	45,7cm / 18"	48,3cm / 19"	44,5cm / 17,5"	47cm / 18,5"	49,5cm / 19,5"	45,7cm / 18"	48,3cm / 19"	50,8cm / 20"	47cm / 18,5"	49,5cm / 19,5"	52,1cm / 20,5"	48,3cm / 19"	50,8cm / 20"	53,3cm / 21"	
P5	43,2cm / 17"	45,7cm / 18"	48,3cm / 19"	44,5cm / 17,5"	47cm / 18,5"	49,5cm / 19,5"	45,7cm / 18"	48,3cm / 19"	50,8cm / 20"	47cm / 18,5"	49,5cm / 19,5"	52,1cm / 20,5"	48,3cm / 19"	50,8cm / 20"	53,3cm / 21"	49,5cm / 19,5"	52,1cm / 20,5"	54,6cm / 21,5"		

Height available +/- 0,65cm (1/4") (see stem bolt and fork combination at the end of the order form)

CASTERS												
9	CA01	Composite wheel w/ PU tire (1" width)	3" <input type="checkbox"/>	4" <input type="checkbox"/>	5" <input type="checkbox"/>	6" <input type="checkbox"/>	7" <input type="checkbox"/>	8" <input type="checkbox"/>		NC	<input type="checkbox"/>	
10	CA10	Pneumatic casters 6" x 1 ¼"								£140	<input type="checkbox"/>	
11	CA11	Pneumatic casters 8" x 1 ¼"								£140	<input type="checkbox"/>	
12	CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "								£70	<input type="checkbox"/>
13	CA15		5" x 1 ½ "								£70	<input type="checkbox"/>
14	CA16		6" x 1 ½ "								£70	<input type="checkbox"/>
15	CA17		4" x 1 ½ "								£120	<input type="checkbox"/>
16	CA18	Newton UltraCaster <i>Soft-Roll</i> Aluminum	5" x 1 ½ "								£125	<input type="checkbox"/>
17	CA19		6" x 1 ½ "								£130	<input type="checkbox"/>
CASTER OPTIONS												
18	CA12	Caster Pin Locks (<i>Available with +1" or +2" Stem bolt ONLY</i>) ¹									£250	<input type="checkbox"/>
19	CA13	Carbon fiber Frog Legs suspension forks (<i>Not available with caster pin locks</i>)									£395	<input type="checkbox"/>

1 - Not available with WC-19

C2 - REAR Seat-to-floor Height

WHEEL SIZES	Height available		Specify desired height: _____
20" (451) wheel <input type="checkbox"/>	31,8cm / 12,5"	to 43,8cm / 17,25"	
22" (501) wheel <input type="checkbox"/>	34,3cm / 13,5"	to 46,4cm / 18,25"	
24" (540) wheel <input type="checkbox"/>	36,8cm / 14,5"	to 48,9cm / 19,25"	
25" (559) wheel <input type="checkbox"/>	38,1cm / 15"	to 50,2cm / 19,75"	
26" (590) wheel <input type="checkbox"/>	39,4cm / 15,5"	to 51,4cm / 20,25"	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1cm" of selected seat to floor height"

Footrest

SWING-IN/SWING-OUT FOOTREST										
20	FR03	60° Swing-in/swing-out	34,3 - 43,2cm/13,5-17" <input type="checkbox"/>	36,8 - 49,5cm/14,5-19,5" <input type="checkbox"/>	39,4 - 52cm/15,5-20,5" <input type="checkbox"/>	Set at:*			NC	<input type="checkbox"/>
21	FR04	70° Swing-in/swing-out	34,3 - 43,2cm/13,5-17" <input type="checkbox"/>	36,8 - 49,5cm/14,5-19,5" <input type="checkbox"/>	39,4 - 52cm/15,5-20,5" <input type="checkbox"/>	Set at:*			NC	<input type="checkbox"/>
22	FR15	80° Swing-in/swing-out	34,3 - 40,6cm/13,5-16" <input type="checkbox"/>	36,8 - 49,5cm/14,5-19,5" <input type="checkbox"/>	39,4 - 52cm/15,5-20,5" <input type="checkbox"/>	Set at:*			NC	<input type="checkbox"/>
23	FR06	90° swing-in/swing-out with front mount - Standard (no extension) ¹	19 - 28cm (7.5 - 11")			Set at:*			NC	<input type="checkbox"/>
24	FR07	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	24 - 33cm (9.5 - 13")			Set at:*			NC	<input type="checkbox"/>
25	FR08	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	24 - 36,8cm (9.5 - 14.5")			Set at:*			NC	<input type="checkbox"/>
FOOTREST OPTIONS ³										
26	FR09	High mount (<i>footplate attached to the hanger with clamps</i>)	12,7 - 34,3cm/5-13,5" <input type="checkbox"/>			Set at: *			£66	<input type="checkbox"/>
27	FR10	Elevating legrest with calf pad, 35.5 to 55.88 cm (14"-22") length adjustable ⁴	Left <input type="checkbox"/>		Right <input type="checkbox"/>	Set at:	(Pair ELR £400)		£200	<input type="checkbox"/>
28	FR99	Omit Footrest							NC	<input type="checkbox"/>

1 - Available only with two pieces adjustable footplates 2 - May interfere with front casters 3 - Not available with front mount * +/- 1.27 cm (1/2")

4 - The elevating legrest replace the footrest side selected

FOOTPLATES										
29	FP01	Standard - Composite							NC	<input type="checkbox"/>
30	FP02	Newton adjustable angle - Composite							£180	<input type="checkbox"/>
31	FP03	Newton adjustable angle - Aluminum							£200	<input type="checkbox"/>
32	FP04	Newton adjustable angle - Carbon fiber (<i>with rubber inserts</i>)							£365	<input type="checkbox"/>
33	FP06	One piece adjustable angle flip-up - Aluminum							£166	<input type="checkbox"/>
34	FP07	One piece adjustable angle flip-up - Carbon fiber (<i>with rubber inserts</i>)							£380	<input type="checkbox"/>
FOOTPLATE ACCESSORIES										
35	FP08	Length adjustable heel loop (<i>not available with one piece footplate</i>)							£35	<input type="checkbox"/>
36	FP09	Calf strap							£30	<input type="checkbox"/>
37	FP30	Bodypoint Aeromesh Padded calf strap							£95	<input type="checkbox"/>

SWING AWAY RESIDUAL LIMB SUPPORT										
38	FR27	Residual limb support - flat 15,24x20,32cm (6"x8")	Left <input type="checkbox"/>		Right <input type="checkbox"/>		Pair ELR £710		£355	<input type="checkbox"/>
39	FR28	Residual limb support - flat 20,32x24,4cm (8"x10")	Left <input type="checkbox"/>		Right <input type="checkbox"/>		Pair ELR £710		£355	<input type="checkbox"/>
40	FR29	Residual limb support - flat 25,4x25,4cm (10"x10")	Left <input type="checkbox"/>		Right <input type="checkbox"/>		Pair ELR £710		£355	<input type="checkbox"/>
41	FR30	Residual limb support - flat 25,4x35,56cm (10"x14")	Left <input type="checkbox"/>		Right <input type="checkbox"/>		Pair ELR £710		£355	<input type="checkbox"/>
42	FR32	Residual limb support - contour 15,24x20,32cm (6"x8")	Left <input type="checkbox"/>		Right <input type="checkbox"/>		Pair ELR £710		£355	<input type="checkbox"/>
43	FR33	Residual limb support - contour 20,32x24,4cm (8"x10")	Left <input type="checkbox"/>		Right <input type="checkbox"/>		Pair ELR £710		£355	<input type="checkbox"/>

Backrest

FIXED ANGLE AND HEIGHT ADJUSTABLE BACK WITH PUSH HANDLE OPTION										
44	BC06	Straight Carbon fiber back cane	22,9 - 30,5 (9-12") □	30,5 - 38,1 (12-15") □	38,1 - 45,7 (15-18") □	45,7 - 53,3 (18-21") STD □			NC	□
45	BC08	8° Bend Carbon fiber back cane		34,3 - 38 (13 ½-15") □	38,1 - 45,7 (15-18") □	45,7 - 53,3 (18-21") STD □			£45	□
ANGLE AND HEIGHT ADJUSTABLE BACK WITH PUSH HANDLE OPTION										
46	BC09	Straight Carbon fiber back cane		30,5 - 38,1 (12-15") □	38,1 - 45,7 (15-18") □	45,7 - 53,3 (18-21") STD □			£175	□
47	BC10	8° bend Carbon fiber back cane		34,3 - 38 (13 ½-15") □	38,1 - 45,7 (15-18") □	45,7 - 53,3 (18-21") STD □			£220	□
48		Select back cane angle:	85° ¹ □	90° □	95° STD □	100° □	105° □	110° □		
Set height to: ²		(measured from seat sling to top of back upholstery)								
PUSH HANDLE OPTIONS										
49	PH01	Push handles							STD	□
50	PH05	Fold down push handles ^{3,5}							£120	□
51	PH06	clamp on height adjustable push handle ^{3,4,5}	22,9 - 50,8 cm/9-20" STD □		22,9 - 58,4 cm/ 9-23" □				£260	□
52	PH02	Omit push handles							NC	□

1 - Not available with a 8 degree bend cane 2 - If not indicated, default height setting is at the lowest position 3 - 120 kg (265 pounds) weight limit

4 - Not compatible with swing away armrests, 8° bend back canes and ampute axle plates

5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

STROLLER HANDLE AND STABILIZER BAR ¹ (only available with push handle)										
53	PH04	Stroller handle (not available with the HD option)							£140	□
54	PH03	Newton folding stabilizer bar							£220	□

1 - Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY										
55	BH07	Slip-on Active-Mesh back upholstery							NC	□
56	BH08	Adjustable tension Active-Mesh back upholstery							£145	□
57	BH99	Omit Back upholstery							STD	□

SEAT CUSHION										
58	SC01	5,1 cm (2") Cushion							£145	□
59	SC02	7,6 cm (3") Cushion							£162	□
60	SC99	Omit seat cushion							STD	□

Rear wheels

ONE ARM DRIVE ¹ Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims										
61	HR22	One arm drive ^{2,3}	24"	left □	right □				£1 200	□
62	HR98	Omit handrim on opposite side of the One Arm Drive (OAD)							NC	□

REAR WHEELS											
63	RW01	Mag ⁴	20" □	22" □	24" □				£45	□	
64	RW02	Newton One - spoke wheel ⁵	20" □	22" □	24" □				NC	□	
65	RW03	Newton Gravity - ultralight	20" □	22" □	24" □	25" □	26" □		£275	□	
66	RW04	Spinergy Spox		22" □	24" □	25" □	26" □		£600	□	
67		Spinergy Spox spokes colors ⁶	Blk STD □	Wht □	Blue □	Red □	Yel □		NC	□	
68	RW05	Spinergy LX		22" □	24" □	25" □	26" □		£800	□	
69		Spinergy LX spokes colors ⁶	Blk STD □	Pur □	Wht □	Blue □	Red □	Yel □	Pink □	Grn □	Org □
70	RW07	Spinergy CLX Carbon Blade			24" □	25" □			£1 200	□	
71		Spinergy CLX spokes colors ⁶	Blk STD □	Pur □	Wht □	Blue □	Red □	Yel □	Pink □	Grn □	Org □

TIRES Available										
72	TI01	Soft Urethane 3,17cm (1 1/4")			20", 22" 24"				STD	□
73	TI02	Soft Urethane knobby 3,49cm (1 3/8")			22", 24", 25", 26"				£85	□
74	TI04	Soft Urethane - Shox 2,54cm (1")			22", 24", 25", 26"				£99	□
75	TI05	Pneumatic 3,49cm (1 3/8")			20", 22", 24", 25", 26"				NC	□
76	TI08	Pneumatic with Airless insert 3,49cm (1 3/8") ⁷			22", 24", 26"				£70	□
77	TI06	Pneumatic HP - 2,54cm (1")			20", 22", 24", 25", 26"				£85	□
78	TI07	Pneumatic HP - Schwalbe Marathon Plus 2,54cm (1")			24", 25", 26"				£95	□

1 - Offered with handrims on each side of the wheelchair 2 - Only standard handrim position available 3 - Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position 5 - Not compatible with Shox tires

6 - Colour spokes require an additional 10 days lead-time 7 - Only available with Newton One and Newton Gravity rear wheels

Handrims

HANDRIM											
79	HR01	Aluminum anodized							STD	<input type="checkbox"/>	
80	HR61	Aluminum black hard anodized ^{1,3}						<i>New</i>	£55	<input type="checkbox"/>	
81	HR02	Plastic coated ²							£85	<input type="checkbox"/>	
82	HR03	Q-Grip High friction Neoprene coated ^{2,3,4}							£295	<input type="checkbox"/>	
83	HR07	Newton Air Grip ³	Black <input type="checkbox"/>						£85	<input type="checkbox"/>	
84	HR	Newton Air Grip ³	(08) Sunkissed Orange <input type="checkbox"/>	(09) Burgundy <input type="checkbox"/>	(10) Sapphire Blue <input type="checkbox"/>	(11) Monster Green <input type="checkbox"/>	(12) Charcoal <input type="checkbox"/>	(13) Fuchsia <input type="checkbox"/>	(14) Steel Blue <input type="checkbox"/>	£175	<input type="checkbox"/>
			(15) Ferrari Red <input type="checkbox"/>	(16) Acid Green <input type="checkbox"/>	(17) White <input type="checkbox"/>	(18) Sky Blue <input type="checkbox"/>					<input type="checkbox"/>
86	HR05	Surge LT ^{2,3,4}	(1-1/8" Oval handrim w/ Gripton high friction strip)							£355	<input type="checkbox"/>
87	HR04	Natural Fit ^{2,3,4}							£270	<input type="checkbox"/>	
88	HR99	Omit handrim							NC	<input type="checkbox"/>	

1 - Not available in 26" 2 - Only standard handrim position available 3 - Not available with Mag wheels or OAD 4 - Not available in 20"

Handrims option

Assembly position for aluminum and Newton Air Grip handrims											
89	HR29	Standard Handrim position	Full tab							STD	<input type="checkbox"/>
90	HR30	Narrow Handrim position	Short tab 2 positions							£20	<input type="checkbox"/>
91	HR31	Super Narrow Handrim	Short tab 2 positions STD <input type="checkbox"/>		Short tab 1 position <input type="checkbox"/>				£20	<input type="checkbox"/>	
Thumb piece ⁴ selection for Natural Fit handrims											
92	HR23	Standard grip							£75	<input type="checkbox"/>	
93	HR27	Super grip							£120	<input type="checkbox"/>	
94	HR099	Omit thumb piece							NC	<input type="checkbox"/>	


4 - Compatible with Surge, Surge LT handrims

AXLE ¹										
95	AX01	Quick release axle							STD	<input type="checkbox"/>
96	AX02	Permanent axle							NC	<input type="checkbox"/>
97	AX03	Quad release axle							£108	<input type="checkbox"/>
AXLE PLATE										
98	AX04	Amputee adjustable axle plate ²							£170	<input type="checkbox"/>

1 - Not required with one arm drive 2 - May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER										
99	CM00	0° camber							NC	<input type="checkbox"/>
100	CM03	3° camber							NC	<input type="checkbox"/>
101	CM06	6° camber (not available with amputee axle plate)							£130	<input type="checkbox"/>

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

X =	Most stable										Most responsive		
	2,5 1"	3,2 1/4"	3,8 1/2"	4,4 3/4"	5,1 2"	5,7 1/4"	6,3 1/2"	7 3/4"	7,6 3"	8,3 1/4"	8,9 1/2"	9,5 3/4"	10,1 4"
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

If not indicated, default setting is at 3.8 cm/1.5"

Possible conflict between a seat depth under 40.6 cm (16") and a center of gravity over 6.3 cm (2.5").


Possible conflict between a 66 cm (26") rear wheel and a center of gravity over 6,3 cm (2,5").

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

WHEEL LOCK										
102	WL08	Newton Wheel lock - push to lock							STD	<input type="checkbox"/>
103	WL02	Pull to lock							NC	<input type="checkbox"/>
104	WL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)							£135	<input type="checkbox"/>
105	WL13 WL14	6" removable extension handle							£66	<input type="checkbox"/>
106	WL05	Scissor Lock							£120	<input type="checkbox"/>
107	WL06	Attendant Lock (additional set of wheel lock mounted to the back)							£175	<input type="checkbox"/>
108	WL99	Omit wheel locks							NC	<input type="checkbox"/>

Armrests									
FLIP BACK ARMREST (Convertible to "T" armrest) ^{1,4}			Height Adjustment			Removable function			
109	AR01	10" Desk length arm pad	20 - 30cm (8 - 12") STD <input type="checkbox"/>	25 - 36cm (10 - 14") <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/>	STD	£310	<input type="checkbox"/>
110	AR02	14" Full length arm pad	20 - 30cm (8 - 12") STD <input type="checkbox"/>	25 - 36cm (10 - 14") <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/>	STD	£310	<input type="checkbox"/>
"T" ARMREST ¹			Height Adjustment						
111	AR03	10" Desk length arm pad	20 - 30cm (8 - 12") STD <input type="checkbox"/>	25 - 36cm (10 - 14") <input type="checkbox"/>				£285	<input type="checkbox"/>
112	AR04	14" Full length arm pad	20 - 30cm (8 - 12") STD <input type="checkbox"/>	25 - 36cm (10 - 14") <input type="checkbox"/>				£285	<input type="checkbox"/>
SWING AWAY ARMREST ^{2,4}			Height Adjustment						
113	AR05	Tubular arm pad	20,3 - 27,9cm (8 - 11")					£195	<input type="checkbox"/>
114	AR06	10 "Desk length arm pad	20,3 - 27,9cm (8 - 11")					£195	<input type="checkbox"/>
GEL ARM PAD ³									
115	AR12	36cm (14") Full Length Arm Pad - Gel Ovation						£133	<input type="checkbox"/>
116	AR99	Omit armrest						NC	<input type="checkbox"/>
1 - May have a 1,25 cm adjustment restriction on both limits of the range 2 - Could interfere with some rigid back rests 3 - Not available with swing away armrest; 4 - Not available with attendant lock									
SIDE GUARDS (not available with Flip Back and "T" Armrest)									
117	SG01	Plastic side guards (removable)						£170	<input type="checkbox"/>
118	SG02	Carbon fiber side guards (fixed) ¹						£220	<input type="checkbox"/>
119	SG03	Carbon fiber side guards (removable)						£255	<input type="checkbox"/>
120	SG04	Carbon fiber fender side guards (removable)						£435	<input type="checkbox"/>
121	SG06	Plastic fender side guards (removable) ²						£235	<input type="checkbox"/>
122	SG99	Omit side guard						NC	<input type="checkbox"/>
SIDE GUARDS OPTION									
123	SG07	Side guards clip (pair) ³						NC	<input type="checkbox"/>
1 - Not available with back upholstery 2 - Not compatible with 20" and 22" rear wheels 3 - Not available with fixed carbon fiber side guard									
POSITIONING BELT									
124	SB01	Velcro adjustable	Standard (88,9 cm/35") <input type="checkbox"/>	Longue (152,40 cm/60") <input type="checkbox"/>				£40	<input type="checkbox"/>
125	SB02	Auto buckle	Standard (88,9 cm/35") <input type="checkbox"/>	Longue (152,40 cm/60") <input type="checkbox"/>				£75	<input type="checkbox"/>
126	SB07	Hardware attachment for 2 point seatbelt (1 pair)						NC	<input type="checkbox"/>
127	SB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)						£15	<input type="checkbox"/>
128	SB04	BodyPoint Hip Belt - 2 point not padded	M (142,24 cm/56") <input type="checkbox"/>	L (154,94 cm/61") <input type="checkbox"/>				£70	<input type="checkbox"/>
129	SB05	BodyPoint Hip Belt - 2 point padded	S (129,54 cm/51") <input type="checkbox"/>	M (142,24 cm/56") <input type="checkbox"/>	L (154,94 cm/61") <input type="checkbox"/>			£150	<input type="checkbox"/>
130	SB06	BodyPoint Hip Belt - 4 point padded	S (129,54 cm/51") <input type="checkbox"/>	M (142,24 cm/56") <input type="checkbox"/>	L (154,94 cm/61") <input type="checkbox"/>			£165	<input type="checkbox"/>
131	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	M (73,66 cm/18") <input type="checkbox"/>	L (93,98 cm/24") <input type="checkbox"/>				£250	<input type="checkbox"/>
132	SB99	Omit positioning belt						STD	<input type="checkbox"/>
HEADREST ¹									
133	AC01	Profiled Headrest	S (8,8x15,2 cm/3,5"x6") <input type="checkbox"/>	M (11,4x17,8 cm/4,5"x7") <input type="checkbox"/>	L (14x2,3 cm/5,5"x8") <input type="checkbox"/>			£158	<input type="checkbox"/>
134	AC28	Headrest mounting hardware with removable non folding stabilizer bar						£258	<input type="checkbox"/>
1 - Only available with push handles - Not available with a stroller handle or a folding stabiliser bar									
ACCESSORIES									
135	AC03	Anti tippers						£120	<input type="checkbox"/>
136	AC29	Swing-away anti-tippers ¹	Left <input type="checkbox"/>	Right <input type="checkbox"/>		New		£116	<input type="checkbox"/>
137	AC27	Tip assist ²						£54	<input type="checkbox"/>
138	AC04	Transit Tie Down (Unoccupied)						NC	<input type="checkbox"/>
139	AC05	WC-19 Transit Tie Down w/Q'straint belt ³	Short 50-104 cm/20-41" <input type="checkbox"/>	Long 50-116 cm/20-46" <input type="checkbox"/>				£280	<input type="checkbox"/>
140	AC06	Cane carrier	Left <input type="checkbox"/>	Right <input type="checkbox"/>				£104	<input type="checkbox"/>
141	AC07	Oxygen tank holder - E Type Not compatible with swing away armrest or amputee axle plate	Left					£270	<input type="checkbox"/>
142	AC08	Spoke guards ⁴	Available: 20", 22", 24", 26"					£100	<input type="checkbox"/>
143	AC09	Footrest neoprene impact guard (pair) - Black						£79	<input type="checkbox"/>
144	AC10	Seat pouch						£65	<input type="checkbox"/>
145	AC26	BodyPoint Mobility Bag						£108	<input type="checkbox"/>
146	AC23	Newton Backpack						£100	<input type="checkbox"/>
147	AC12	Velcro Closing Strap						£45	<input type="checkbox"/>
148	AC13	Tool kit						£39	<input type="checkbox"/>
149	AC16	Motion Composites Polo Shirt	W/S <input type="checkbox"/>	W/M <input type="checkbox"/>	W/L <input type="checkbox"/>			£29	<input type="checkbox"/>
			M/S <input type="checkbox"/>	M/M <input type="checkbox"/>	M/L <input type="checkbox"/>	M/XL <input type="checkbox"/>	M/XXL <input type="checkbox"/>	M/XXXL <input type="checkbox"/>	M/XXXXL <input type="checkbox"/>
1 - Must select the pair (left and right) with HD option 2 - Cannot be combined with Anti Tippers 3 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option 4 - Not available with Mag wheels and Spinergy									

FRAME COLOR								
150	FC01	Carbon Fiber Black <input type="checkbox"/>					NC <input type="checkbox"/>	
151	FC	(03) Sunkissed Orange <input type="checkbox"/>	(04) Burgundy <input type="checkbox"/>	(05) Sapphire Blue <input type="checkbox"/>	(06) Monster Green <input type="checkbox"/>		NC <input type="checkbox"/>	
		(07) Charcoal <input type="checkbox"/>	(08) Fuchsia <input type="checkbox"/>	(09) Steel Blue <input type="checkbox"/>	(17) Ferrari Red <input type="checkbox"/>	(18) Acid Green <input type="checkbox"/>		
152	FC99	Custom Color	Specify PPG Paintit color:			Additional 10 days lead-time	£300	<input type="checkbox"/>
TOTAL								



Notes :