

# Carbon Fibre, Ultra Lightweight



Sleek, practical and ultra-lightweight. The all new iGo Lite weighs only 19.8 kg making it our lightest power chair yet! With a built-in USB charger, you'll be able to travel and stay connected with the world around you. It's flight-safe batteries mean the sky is your limit - simply unfold and go, for an adventure of a lifetime.

## Specifications

Weight Capacity	21st 3 lbs (135 kg)
Maximum Speed <sup>1</sup>	Up to 3.7 mph
Overall Weight	19.8 kg (43.6 lbs)
Turning Radius <sup>3</sup>	900 mm (35.43")
Overall Length <sup>3</sup>	960 mm (37.79")
Overall Width <sup>3</sup>	610 mm (24")
Front Wheel	178 mm (7") - solid
Drive Tyres	216 mm (8.5") - solid
Suspension	Limited
Total Weight Without Batteries <sup>2, 4</sup>	18 kg (39.6 lbs)
Heaviest Piece When Disassembled <sup>3</sup>	N/A
Drivetrain	2x 250W motors
Dual Breaking System	Electronic regenerative, disc park brake
Battery Requirements <sup>3</sup>	12 Ah Lithium-Ion battery
Battery Weight <sup>4</sup>	1.8 kg
Battery Charger	2-amp, on board and off-board
Colour	Black

## Features

- Large secure underseat storage compartment
- On-board USB charger
- Adjustable control position
- Stylish Carbon Fibre design
- Simple one touch folding system
- Compact lightweight design
- 12Ah flight safe batteries

Available in:

- Black - Carbon Fibre



<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup> Includes standard seat weight.

<sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

<sup>4</sup> Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.