

Overall Dimension	940x610x960 mm	
Folded Dimension	700x310x610 mm	
Seat Height	500 mm	
Seat width	460 mm	
Backrest Height	480 mm	
N.W. of Wheelchair	19.8kg(including the battery)	
Battery weight	1.8kg	
Maximum Speed	6 km/h	
Braking Distance on Horizontal Road	≤1500 mm	
Minimum Turning Radius	≤900 mm	
Maximum Load Capacity	135kg	
Theoretical Driving Distance	≥15km	
Climbing Ability	≥6°	
Motor Model	24V/250W	
Lithium Battery Specification	24V/12Ah	
Maximum Output Current of Controller	35A	
Maximum Output Current of Charger	2A	
Front Wheel Specification	Solid Tyre	178mm
Rear Wheel Specifications	Solid Tyre	216mm