

Just Fold And Go!



Get ready for the adventure of a lifetime! Feel the anticipation as you wait for your flight in the comfort and style of your i-Go+ power chair. This compact power chair has plenty of storage and folds in just a few steps to make airport check-in effortless. When your plane touches down at your destination, you can just unfold and go!



Specifications

Weight Capacity	17st 12 lbs (113 kg)
Maximum Speed ¹	Up to 3.6 mph
Ground Clearance ³	79 mm (3")
Turning Radius ³	616 mm (24.25")
Overall Length ³	1016 mm (40")
Overall Width ³	597 mm (23.5")
Castors	203 mm (8") - solid
Drive Tyres	305 mm (12") - pneumatic
Suspension	Limited
Total Weight Without Batteries ^{2, 4}	24 kg (54 lbs)
Heaviest Piece When Disassembled ³	N/A
Drivetrain	2x motors
Dual Breaking System	Electronic regenerative, disc park brake
Battery Requirements ³	18.2 Ah Lithium-Ion battery 10.4 Ah Lithium-Ion travel battery
Battery Weight ⁴	Standard battery: 3.4 kg Travel battery: 2.4 kg
Battery Charger	2-amp, off-board
Colour	Silver

Features

- Large secure underseat storage compartment
- Unique comfort seating
- Hi-torque front suspension
- Mesh cup holder
- Simple one touch folding system
- Compact lightweight design
- Optional flight safe batteries available

Available in:

- Silver



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

² Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.