



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

NOVEMBER 2023

Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## Action 3 Junior Large

### Standard Features

- Max user weight 80kg (12.5 stone)
- Crash tested ISO 7176-19
- Self Propulsion or Attendant Propulsion
- Choice of 4 seat widths 12" - 15" (305mm - 380mm)
- Choice of 2 seat depths (depth adjustable upholstery)
- Choice fo 3 seat heights 16" - 18" (410mm - 460mm)
- Black nylon upholstery
- Folding backrest set at 16" (400mm)
- Width and height adjustable backrest upholstery
- 80° Detachable swing away legrest hangers
- Flip up, angle and depth adjustable footplates
- Removable and height adjustable desk armrests
- Adjustable heel straps
- Choice of 7 colours (with 8 special colour options)
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and lap belt (adjustable both sides)
- User tool kit



**Standard Base price: £647**

Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

| Code                 | Description                             |                       |  |
|----------------------|---|-----------------------|--|
| <b>BASE PRODUCT*</b> |   |                       |  |
| CNC0001              | ACTION <sup>3</sup> JUNIOR              | <input type="radio"/> |  |
| CNCSPEC              | SPECIAL                                 | <input type="radio"/> | POA  |
| <b>PROPULSION*</b>   |   |                       |  |
| CNC0100              | SELF PROPULSION                         | <input type="radio"/> | (NC 1055, 1170) £647   |
| CNC0110              | TRANSIT PROPULSION, WITH RH TIPPING AID | <input type="radio"/> | (NC 0280, 1020, 1030, 1045, 1140, 1141, 1180, 1190, 1195, 1275, 1390, 1430, 1440) £647 |

| <b>SEAT WIDTH*</b>          |  |  |                                     |
|-----------------------------|--|--|-------------------------------------|
| CNC0150                     |  | 12" SW   | (NC 0520, 0570) ○                   |
| CNC0151                     |  | 13" SW   | ○                                   |
| CNC0152                     |  | 14" SW   | ◆                                   |
| CNC0153                     |  | 15" SW   | ○                                   |
| <b>SEAT DEPTH*</b>          |  |  |                                     |
| CNC0170                     |  | ADJ SD 11" - 13"   | ○                                   |
| CNC0171                     |  | ADJ SD 13" - 15"   | ◆                                   |
| <b>SEAT HEIGHT*</b>         |  |  |                                     |
| CNC0260                     |  | 18" SH   | ○                                   |
| CNC0270                     |  | 17" SH   | ○                                   |
| CNC0280                     |  | 16" SH   | (NC 1430, 1440) ○                   |
| <b>SEAT UPHOLSTERY*</b>     |  |  |                                     |
| CNC0335                     |  | NYLON WITH HOOK AND LOOP STRIPS SEAT UPHOLSTERY, WIDTH AND DEPTH ADJUSTABLE                                | ◆                                   |
| <b>SEAT CUSHION*</b>        |  |  |                                     |
| CNC33333                    |  | WITHOUT SEAT CUSHION   | ○                                   |
| CNC0362                     |  | BASIC SEAT CUSHION (ESSENTIAL JUNIOR)  | £50                                 |
| <b>BACKREST*</b>            |  |  |                                     |
| CNC0370                     |  | FIXED 0° BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400MM OR 17" / 430MM  | ○                                   |
| CNC0380                     |  | FIXED 7° BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400MM OR 17" / 430MM  | ○                                   |
| CNC0390                     |  | RECLINING (RATCHET) 0/30°<br>5 HEIGHTS 13" - 17"   | (NC 1050, 1190, 1195) £111          |
| CNC0395                     |  | RECLINING (GAS) 0/30° 16"  | (NC 1050, 1190, 1195) £184          |
| CNC0400                     |  | FOLDING BACKREST 3 HEIGHTS 15" / 370MM, 16" / 400 OR 17" / 430MM   | ◆                                   |
| CNC0401                     |  | FOLIDNG BACKREST 20"   | ○                                   |
| CNC0402                     |  | FOLDING BACKREST 2 HEIGHTS 15" / 20"<br>REAR SET PUSH HANDLES  | (NC 0480, 0855, 0902, 1426) £140    |
| CNC0420                     |  | FOLDING BACKREST 16"<br>HEIGHT ADJUSTABLE PUSH HANDLES   | (NC 0855) £76                       |
| CNC0425                     |  | FIXED 7° BACKREST 16"<br>HEIGHT ADJUSTABLE PUSH HANDLES  | (NC 0855) £76                       |
| <b>BACKREST UPHOLSTERY*</b> |  |  |                                     |
| CNC0465                     |  | WITHOUT BACK UPHOLSTERY  | (NC 0855) ○                         |
| CNC0475                     |  | WIDTH AND HEIGHT ADJUSTABLE NYLON SLING TYPE   | ◆                                   |
| CNC0480                     |  | CONTOUR (TENSION) ADJUSTABLE NYLON TYPE  | (NC 0855) £105                      |
| <b>LEGREST HANGER RH*</b>   |  |  |                                     |
| CNC0505                     |  | 90° STRAIGHT RH REMOVABLE HANGER UL STD<br>125 TO 260 MM   | (NC 0600) £29                       |
| CNC0510                     |  | 80° SHORT RH SWING IN/OUT HANGER 0°<br>STANDARD POSITION; UL STD 220 TO 320 MM<br>BLACK COLOUR ONLY        | (NC 0561, 0562, 0563) ◆             |
| CNC0511                     |  | 80° SHORT RH SWING IN/OUT HANGER 30°<br>ABDUCTION POSITION; UL STD 220 TO 320 MM<br>BLACK COLOUR ONLY      | (NC 0560, 0562, 0563, 0611, 0881) ○ |
| CNC0512                     |  | 80° LONG RH SWING IN/OUT HANGER 0°<br>STANDARD POSITION; UL STD 300 TO 450 MM                              | (NC 0560, 0561, 0563) ○             |
| CNC0513                     |  | 80° LONG RH SWING IN/OUT HANGER 30°<br>ABDUCTION POSITION; UL STD 300 TO 450 MM                            | (NC 0560, 0561, 0562, 0611, 0881) ○ |
| CNC0520                     |  | ELEVATING LEGREST (0° TO 80°) RH SWING<br>IN/SWING OUT HANGER (INCLUDING CALFPAD);<br>UL STD 275 TO 375 MM | (NC 0600, 0611, 0880, 0886) £89     |

| LEGREST HANGER LH* |  |   |     |
|--------------------|--|---|-----|
| CNC0555            |  | 90° STRAIGHT LH REMOVABLE HANGER UL STD 125 TO 260 MM (NC 0620)   | £29 |
| CNC0560            |  | 80° SHORT LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY                                   | ◆   |
| CNC0561            |  | 80° SHORT LH SWING IN/OUT HANGER 30° ABDUCTION POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY (NC 0611, 0881)                 | ○   |
| CNC0562            |  | 80° LONG LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM  | ○   |
| CNC0563            |  | 80° LONG LH SWING IN/OUT HANGER 30° ABDUCTION POSITION; UL STD 300 TO 450 MM (NC 0611, 0881)                                    | ○   |
| CNC0570            |  | ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL ST 275 TO 375 MM (NC 0611, 0620, 0880, 0886) | £89 |
| FOOTREST RH*       |  |   |     |
| CNC0600            |  | FLIP UP RH FOOTPLATE (NC 0875, 0886)  | ○   |
| CNC0605            |  | FLIP UP DEPTH AND ANGLE ADJ RH  | ◆   |
| CNC0611            |  | ONE PIECE FOOTBOARD (NC 0620, 0625, 0875)   | ○   |
| FOOTREST LH*       |  |   |     |
| CNC0620            |  | FLIP UP LH FOOTPLATE (NC 0875, 0886)  | ○   |
| CNC0625            |  | FLIP UP DEPTH AND ANGLE ADJ LH  | ◆   |
| CNC66666           |  | NO FOOTPLATE REQUIRED (PLEASE SELECT IN COMBINATION WITH OPTION 0611)   | ○   |
| ARMREST RH*        |  |   |     |
| CNC0640            |  | REMOVABLE RH (N°1) DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED (NC 0745, 0746, 0750, 0755)      | ◆   |
| CNC0645            |  | REMOVABLE RH (N°1) JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED (NC 0740, 0746, 0750, 0755)    | ○   |
| CNC0646            |  | REMOVABLE RH (N°1) NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED (NC 0740, 0745, 0750, 0755)  | ○   |
| CNC0650            |  | MUDGUARD RH HEIGHT ADJUSTABLE (155 TO 220 MM) (NC 0740, 0745, 0746, 0755, 0910, 1201, 1471, 1472)                               | £35 |
| CNC0655            |  | SKIRTGUARD RH JUNIOR (NC 0740, 0745, 0746, 0750, 0910, 1471, 1472)  | ○   |
| ARMREST LH*        |  |   |     |
| CNC0740            |  | REMOVABLE LH (N°1) DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED                                  | ◆   |
| CNC0745            |  | REMOVABLE LH (N°1) JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED                                | ○   |
| CNC0746            |  | REMOVABLE LH (N°1) NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED                              | ○   |
| CNC0750            |  | MUDGUARD LH HEIGHT ADJUSTABLE (155 TO 220 MM) (NC 0910, 1201, 1471, 1472)   | £35 |
| CNC0755            |  | SKIRTGUARD LH JUNIOR (NC 0910, 1471, 1472)  | ○   |

| OTHER SEATING OPTIONS |  |   |      |
|-----------------------|--|---|------|
| CNC0855               |  | BACKREST EXTENSION (NC 0860, 0865)  | £72  |
| CNC0860               |  | HEADREST STANDARD MULTI - ADJUSTABLE (NC 0865, 0902)  | £134 |
| CNC0865               |  | HEADREST ROUND MULTI - ADJUSTABLE (NC 0902)   | £151 |
| CNC0880               |  | ADJUSTABLE CALF STRAP   | £11  |
| CNC0881               |  | FRONT STABILISER KIT + 100 MM   | ○    |
| CNC0875               |  | BODYPOINT ANKLE HUGGERS (PAIR) SIZE S - UNMOUNTED   | £114 |
| CNC0886               |  | FOOTPLATE ADJUSTMENT OPTIMIZER KIT  | £21  |
| CNC0902               |  | WIDTH ADJUSTABLE BACK BRACE   | £41  |
| CNC0910               |  | WITHOUT AUTOLOCK SAFETY DEVICE ON ARMRESTS N°1  | ○    |
| COLOUR*               |  |   |      |
| CNC0952               |  | PEARL WHITE   | ○    |
| CNC0953               |  | MATT BLACK  | ○    |
| CNC0954               |  | MATT RED  | ○    |
| CNC0955               |  | INTENSE YELLOW  | ○    |
| CNC0956               |  | POLAR BLUE  | ○    |
| CNC0957               |  | CANDY PINK  | ○    |
| CNC0959               |  | NAVY BLUE   | ○    |
| CNC0991               |  | SPECIAL COLOUR (ON DEMAND ONLY)   | £169 |
| SWO0001               |  | SLATE GREY - 3 WEEK LEAD TIME   | £169 |
| SWO0002               |  | ELECTRIC GREEN - 3 WEEK LEAD TIME   | £169 |
| SWO0003               |  | HYPNOTIC PURPLE - 3 WEEK LEAD TIME  | £169 |
| SWO0004               |  | ANTHRACITE BLACK - 3 WEEK LEAD TIME   | £169 |
| SWO0005               |  | DARK RED - 3 WEEK LEAD TIME   | £169 |
| SWO0006               |  | PEARL GREY - 3 WEEK LEAD TIME   | £169 |
| SWO0007               |  | VIBRANT ORANGE - 3 WEEK LEAD TIME   | £169 |
| FRONT CASTORS*        |  |   |      |
| CNC0999               |  | 5" X 1" BLACK SOLID TYRE  | ◆    |
| CNC1001               |  | 5" X 1" FLASHING LIGHTS SOLID TYRE  | £46  |
| CNC1007               |  | 5" X 1.3/4" BLACK SOLID TYRE  | ○    |
| REAR WHEEL SIZE*      |  |   |      |
| CNC1010               |  | 12" (NC 1040, 1045, 1130, 1140, 1141, 1180, 1253, 1254, 1255, 1256, 1380, 1390)                                     | ○    |
| CNC1020               |  | 20" (NC 1040, 1055, 1050, 1170)   | ○    |
| CNC1030               |  | 22" (NC 1045, 1055, 1170)   | ○    |
| REAR WHEELS*          |  |   |      |
| CNC1040               |  | STANDARD SELF-PROPEL REAR WHEEL (BLACK HUB, COLORLESS RIM AND ZINC PLATED SPOKES) (NC 1280)                         | ○    |
| CNC1045               |  | REINFORCED STANDARD SELF-PROPEL REAR WHEEL (COLORLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES) (NC 1280) | ○    |
| CNC1050               |  | DRUM BRAKE SELF-PROPEL 22" or TRANSIT 12" (NC 1270, 1275)   | ○    |
| CNC1055               |  | TRANSIT COMPOSITE (BLACK PLASTIC HUB AND RIM) (NC 1280)   | ○    |
| TYRES*                |  |   |      |
| CNC1120               |  | SOLID BLACK   | ○    |
| CNC1130               |  | PNEUMATIC PROFILE BLACK   | ○    |

| <b>HANDRIM - A SELECTION MUST BE MADE FOR SELF PROPULSION</b> |                          |  |                      |
|---|--------------------------|--|----------------------|
| CNC1140   | <input type="checkbox"/> | ALUMINIUM PRE-ANODIZED HANDRIM                         | ○                    |
| CNC1141   | <input type="checkbox"/> | HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC              | £75                  |
| <b>REAR WHEEL AXLE*</b>                                       |                          |  |                      |
| CNC1170   | <input type="checkbox"/> | FIXED AXLE   | ○                    |
| CNC1180   | <input type="checkbox"/> | Q/R AXLE   | ○                    |
| <b>REAR WHEEL POSITION*</b>                                   |                          |  |                      |
| CNC1190   | <input type="checkbox"/> | STANDARD REAR WHEELS POSITION                          | ○                    |
| CNC1195   | <input type="checkbox"/> | MEDIUM REAR WHEEL POSITION                             | ○                    |
| CNC1200   | <input type="checkbox"/> | PASSIVE REAR WHEEL POSITION                            | ○                    |
| CNC1201   | <input type="checkbox"/> | DUAL AMPUTEE REAR WHEEL POSITION (NC 1275, 1280)       | £68                  |
| <b>CHASSIS OPTION</b>   |                          |  |                      |
| CNC1253   | <input type="checkbox"/> | SPOKE GUARDS CLEAR                                     | £81                  |
| CNC1254   | <input type="checkbox"/> | SPOKE GUARDS PIXEL                                     | £100                 |
| CNC1255   | <input type="checkbox"/> | SPOKE GUARDS SNOWFLAKES                                | £100                 |
| CNC1256   | <input type="checkbox"/> | SPOKE GUARDS COLOR ME                                  | £100                 |
| <b>BRAKES*</b>  |                          |  |                      |
| CNC1270   | <input type="checkbox"/> | PUSH TO LOCK BRAKES                                    | ◆                    |
| CNC1275   | <input type="checkbox"/> | PULL TO LOCK BRAKES                                    | ○                    |
| CNC1280   | <input type="checkbox"/> | ATTENDANT HUB BRAKES (NC 1426)                         | £147                 |
| <b>TRANSPORTATION KIT*</b>                                    |                          |  |                      |
| CNC1345   | <input type="checkbox"/> | TRANSPORTATION BRACKETS                                | ◆                    |
| <b>SAFETY OPTIONS*</b>  |                          |  |                      |
| CNC1362   | <input type="checkbox"/> | FLAT-MOUNT WITH SWIVEL BUCKEL<br>MAX USER WEIGHT 50 Kg | ◆                    |
| CNC1364   | <input type="checkbox"/> | FLAT-MOUNT WITH REHAB LATCH                            | ◆                    |
| CNC1370   | <input type="checkbox"/> | REAR ANTI-TIPPERS WITH TIPPING AID (NC 1430, 1440)     | £35                  |
| CNC1380   | <input type="checkbox"/> | REFLECTORS KIT   | £7                   |
| CNC1390   | <input type="checkbox"/> | BRAKE EXTENSION PAIR                                   | £14                  |
| <b>OTHER OPTIONS</b>  |                          |  |                      |
| CNC1426   | <input type="checkbox"/> | STROLLER PUSH HANDLES                                  | £89                  |
| CNC1430   | <input type="checkbox"/> | RH TIPPING AID   | £14                  |
| CNC1440   | <input type="checkbox"/> | LH TIPPING AID   | £14                  |
| CNC1471   | <input type="checkbox"/> | RH SWING AWAY ACRYLIC TRAY                             | £196                 |
| CNC1472   | <input type="checkbox"/> | LH SWING AWAY ACRYLIC TRAY                             | £196                 |
| <b>TOTAL PRICE</b>  |                          |  | <input type="text"/> |

**Non Conforming Options:**

|         |   |
|---------|---|
| CNC0280 | 16" SH  |
| CNC0480 | CONTOUR (TENSION) ADJUSTABLE NYLON TYPE   |
| CNC0520 | ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL STD 275 TO 375 MM      |
| CNC0560 | 80° SHORT LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY             |
| CNC0561 | 80° SHORT LH SWING IN/OUT HANGER 30° ABDUCTION POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY           |
| CNC0562 | 80° LONG LH SWING IN/OUT HANGER 0° STANDARD POSITION; UL STD 300 TO 450 MM                                |
| CNC0563 | 80° LONG LH SWING IN/OUT HANGER 30° ABDUCTION POSITION; UL STD 300 TO 450 MM                              |
| CNC0570 | ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL ST 275 TO 375 MM       |
| CNC0600 | FLIP UP RH FOOTPLATE  |
| CNC0611 | ONE PIECE FOOTBOARD   |
| CNC0620 | FLIP UP LH FOOTPLATE  |
| CNC0625 | FLIP UP DEPTH AND ANGLE ADJ LH  |
| CNC0740 | REMOVABLE LH (N°1) DESK, STD ARMPAD (60 MM) ADJ. 180 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED            |
| CNC0745 | REMOVABLE LH (N°1) JUNIOR, STD ARMPAD (60 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED          |
| CNC0746 | REMOVABLE LH (N°1) NARROW ARMPAD JUNIOR (50 MM) ADJ. 150 TO 250 MM AUTOLOCK SAFETY DEVICE INCLUDED        |
| CNC0750 | MUDGUARD LH HEIGHT ADJUSTABLE (155 TO 220 MM)   |
| CNC0755 | SKIRTGUARD LH JUNIOR  |
| CNC0855 | BACKREST EXTENSION  |
| CNC0860 | HEADREST STANDARD MULTI - ADJUSTABLE  |
| CNC0865 | HEADREST ROUND MULTI - ADJUSTABLE   |
| CNC0880 | ADJUSTABLE CALF STRAP   |
| CNC0881 | FRONT STABILISER KIT + 100 MM   |
| CNC0875 | BODYPOINT ANKLE HUGGERS (PAIR) SIZE S - UNMOUNTED   |
| CNC0886 | FOOTPLATE ADJUSTMENT OPTIMIZER KIT  |
| CNC0902 | WIDTH ADJUSTABLE BACK BRACE   |
| CNC0910 | WITHOUT AUTOLOCK SAFETY DEVICE ON ARMRESTS N°1  |
| CNC1020 | 20"   |
| CNC1030 | 22"   |
| CNC1040 | STANDARD SELF-PROPEL REAR WHEEL (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES)                        |
| CNC1045 | REINFORCED STANDARD SELF-PROPEL REAR WHEEL (COLORLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES) |
| CNC1050 | DRUM BRAKE SELF-PROPEL 22" or TRANSIT 12"   |
| CNC1055 | TRANSIT COMPOSITE (BLACK PLASTIC HUB AND RIM)   |
| CNC1130 | PNEUMATIC PROFILE BLACK   |
| CNC1140 | ALUMINIUM PRE-ANODIZED HANDRIM  |
| CNC1141 | HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC   |
| CNC1170 | FIXED AXLE  |
| CNC1180 | Q/R AXLE  |
| CNC1190 | STANDARD REAR WHEELS POSITION   |
| CNC1195 | MEDIUM REAR WHEEL POSITION  |
| CNC1201 | DUAL AMPUTEE REAR WHEEL POSITION  |
| CNC1253 | SPOKE GUARDS CLEAR  |
| CNC1254 | SPOKE GUARDS PIXEL  |

|         |                            |
|---------|----------------------------|
| CNC1255 | SPOKE GUARDS SNOWFLAKES    |
| CNC1256 | SPOKE GUARDS COLOR ME      |
| CNC1270 | PUSH TO LOCK BRAKES        |
| CNC1275 | PULL TO LOCK BRAKES        |
| CNC1280 | ATTENDANT HUB BRAKES       |
| CNC1380 | REFLECTORS KIT             |
| CNC1390 | BRAKE EXTENSION PAIR       |
| CNC1426 | STROLLER PUSH HANDLES      |
| CNC1430 | RH TIPPING AID             |
| CNC1440 | LH TIPPING AID             |
| CNC1471 | RH SWING AWAY ACRYLIC TRAY |
| CNC1472 | LH SWING AWAY ACRYLIC TRAY |



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):