



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

NOVEMBER 2023
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3NG Comfort

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Self Propulsion or Attendant Propulsion
- Choice of 4 seat widths (16" - 20") 405mm to 505mm
- Choice of 2 seat depths (16" or 18") 400mm or 450mm
- Choice of 3 seat height (16" - 18") 410mm to 460mm
- Black nylon seat upholstery
- Reclining (Ratchet) backrest 0° - 30° with Laguna shell 20° (510mm)
- Comfort Black multi stretch PU coated fabric Passad 2 backrest
- Elevating legrest hangers 0° to 75° including calf pad
- Flip up, angle and depth adjustable footplates
- Removable, height and depth adjustable comfort armrests
- Black multi stretch PU coated fabric trunk supports and headrest
- Choice of 8 colours (with 7 special colour options)
- Push or Pull to lock brakes (configuration dependant)
- Comfort transportation kit and lap belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

Standard Base price: £1,096



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
CEG0005	ACTION 3NG COMFORT (NC 0852)	<input type="radio"/>
CEGSPEC	SPECIAL	POA
PROPULSION*		
CEG0100	SELF PROPULSION (NC 1062, 1074)	£1,096
CEG0110	TRANSIT PROPULSION (NC 0290, 0300, 1068, 1150, 1180, 1260, 1275, 1390)	£1,096

SEAT WIDTH*			
CEG0190		16" / 405 MM SEAT WIDTH	○
CEG0200		17" / 430 MM SEAT WIDTH	◆
CEG0210		18" / 455 MM SEAT WIDTH	○
CEG0230		20" / 505 MM SEAT WIDTH	○
SEAT DEPTH*			
CEG0240		16" / 400 MM SEAT DEPTH (NC 0879)	○
CEG0250		18" / 450 MM SEAT DEPTH (NC 0878)	◆
SEAT HEIGHT*			
CEG0280		18" / 460 MM SEAT HEIGHT	◆
CEG0290		17" / 435 MM SEAT HEIGHT (NC 1010, 1062)	○
CEG0300		16" / 410 MM SEAT HEIGHT (NC 1010, 1062, 1370, 1440, 1441)	○
SEAT UPHOLSTERY*			
CEG0340		NYLON WITH HOOK AND LOOP STRIPS	◆
CEG0345		TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM	○
CEG0372		CUSHION COMFORT MISTRAL 2 LAYER (INCLUDING SEAT UPHOLSTERY 3 VELCROS) (NC 0852, 0869, 0871, 0873, 0874, 0876)	○
BACKREST*			
CEG0441		RECLINING (RATCHET) 0/30° LAGUNA SHELL 20"	◆
BACKREST UPHOLSTERY*			
CEG0490		COMFORT BACKREST LAGUNA BLACK MULTI STRETCH PU COATED 20" / 510 MM HEIGHT	○
CEG0491		COMFORT BACKREST PASSAD 2 BLACK MULTI STRETCH PU COATED 20" / 510 MM HEIGHT	◆
LEGREST HANGER RH*			
CEG0520		ELEVATING LEGREST RH SWING IN/SWING OUT (NC 66666)	◆
CEG0545		AMPUTEE LEGREST RH (NC 0610, 1111)	£133
LEGREST HANGER LH*			
CEG0570		ELEVATING LEGREST LH SWING IN/SWING OUT (NC 55555)	◆
CEG0595		AMPUTEE LEGREST LH (NC 0630, 2222)	£133
FOOTREST RH*			
CEG0610		FLIP UP DEPTH AND ANGLE ADJ RH	◆
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) (NC 0545)	◆
CEE66666		NO FOOTREST	○
FOOTREST LH*			
CEG0630		FLIP UP DEPTH AND ANGLE ADJ LH	◆
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) (NC 0595)	◆
CEE55555		NO FOOTREST	○
ARMREST RH*			
CEG0691		REMOVABLE RH (N°4) COMFORT ARM	◆
CEG0732		REMOVABLE RH (N°4) HEMI ARM (NC 1481)	£59
ARMREST LH*			
CEG0791		REMOVABLE LH (N°4) COMFORT ARM	◆
CEG0832		REMOVABLE LH N°4 HEMI ARM (NC 1481)	£59

SEAT CUSHION			
CEG0869		MATRX FLO-TECH LITE CUSHION - Low Risk	£121
CEG0871		MATRX FLO-TECH CONTOUR VISCO CUSHION - High Risk	£169
CEG0873		MATRX FLO-TECH IMAGE CUSHION - High Risk	£288
CEG0874		MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk	£398
CEG0876		MATRX LIBRA CUSHION - Very High Risk	£410
SEAT CUSHION DEPTH			
CEG0878		FOR 400 MM SEAT DEPTH	○
CEG0879		FOR 450 MM SEAT DEPTH	○
OTHER SEATING OPTIONS			
CEG0852		SEAT RIGIDIZER	£39
CEG0881		TRUNK SUPPORT COMFORT BLACK MULTI STRETCH PU COATED FABRIC (PAIR) ONLY FOR LAGUNA SHELL	◆
CEG0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC	◆
CEG0937		BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED	£114
CEG0938		BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED	£114
COLOUR*			
CEG0961		SLATE GREY (RAL7024 FINE TEXTURE)	○
CEG0964		PEARL WHITE	○
CEG0965		ELECTRIC GREEN	○
CEG0966		MATT BLACK	○
CEG0969		RED MATT	○
CEG0972		WATER BLUE	○
CEG0974		NAVY BLUE	○
CEG0976		DARK PURPLE	○
CEG0980		ANTHRACITE BLACK	○
CEG0991		SPECIAL COLOUR ON DEMAND ONLY	£169
SWO3333		DARK RED - SPECIAL	£169
SWO4444		PEARL GREY - SPECIAL	£169
SWO5555		INTENSE YELLOW	£169
SWO6666		VIBRANT ORANGE	£169
SWO7777		POLAR BLUE	£169
SWO8888		CANDY PINK	£169
CASTOR HOUSING*			
CEG0995		STANDARD	◆
CEG0996		LIGHT WEIGHT	£28
FRONT CASTORS*			
CEG1010		8" X 1.1/4" BLACK SOLID TYRE	◆
CEG1050		6" X 1.1/4" BLACK SOLID TYRE	○
REAR WHEEL SIZE*			
CEG1062		12" (NC 1071, 1072, 1084, 1088, 1140, 1150, 1180, 1260, 1390)	○
CEG1068		24"	◆

REAR WHEELS*			
CEG1071	<input type="checkbox"/>	STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATED SPOKES) (NC 1170, 1280)	◆
CEG1072	<input type="checkbox"/>	SELF-PROPEL REINFORCED REAR WHEEL (NC 1280)	£23
CEG1073	<input type="checkbox"/>	DRUM BRAKE SELF-PROPEL 22 AND 24" or TRANSIT 12" (NC 1201, 1270, 1275)	○
CEG1074	<input type="checkbox"/>	TRANSIT BLACK (NC 1088, 1260, 1280)	○
TYRES*			
CEG1082	<input type="checkbox"/>	SOLID BLACK 1.3/8" / 35MM	◆
CEG1084	<input type="checkbox"/>	PNEUMATIC PROFILE BLACK 1.3/8" / 35MM	£21
CEG1088	<input type="checkbox"/>	PNEUMATIC LOW PROFILE BLACK 1" / 25MM	○
HANDRIM			
CEG1140	<input type="checkbox"/>	ALUMINIUM PRE-ANODIZED HANDRIM	◆
CEG1150	<input type="checkbox"/>	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC	£34
REAR WHEEL AXLE*			
CEG1170	<input type="checkbox"/>	FIXED AXLE	○
CEG1180	<input type="checkbox"/>	Q/R AXLE	◆
REAR WHEEL POSITION*			
CEG1200	<input type="checkbox"/>	PASSIVE REAR WHEELS POSITION	◆
CEG1201	<input type="checkbox"/>	DUAL AMPUTEE REAR WHEEL POSITION (NC 1275, 1280)	£72
PUMP			
CEG1230	<input type="checkbox"/>	PUMP AND SUPPORT	£14
CHASSIS OPTION			
CEG1241	<input type="checkbox"/>	FRONT WHEELS EXTENSION	£22
CEG1260	<input type="checkbox"/>	SPOKE GUARDS PAIR 24" BLK	£78
BRAKES*			
CEG1270	<input type="checkbox"/>	PUSH TO LOCK BRAKES	◆
CEG1275	<input type="checkbox"/>	PULL TO LOCK BRAKES	○
CEG1280	<input type="checkbox"/>	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)	£231
COMFORT TRANSPORTATION KIT*			
CEG1346	<input type="checkbox"/>	COMFORT TRANSPORTATION KIT (NC 1363, 1364)	◆
SAFETY OPTIONS			
CEG1361	<input type="checkbox"/>	BODYPOINT 2-PT POSTURE BELT M (NC 1362, 1363, 1364)	£83
CEG1362	<input type="checkbox"/>	BODYPOINT 2-PT POSTURE BELT L (NC 1363, 1364)	£83
CEG1363	<input type="checkbox"/>	BODYPOINT 4-PT POSTURE BELT MEDIUM (NC 1364)	£98
CEG1364	<input type="checkbox"/>	BODYPOINT 2-PT POSTURE BELT LARGE	£98
CEG1370	<input type="checkbox"/>	REAR ANTI-TIPPERS PAIR (NC 1375, 1440, 1441)	○
CEG1375	<input type="checkbox"/>	REAR ANTI-TIPPERS PAIR WITH TIPPING AID (NC 1440, 1441)	£72
CEG1390	<input type="checkbox"/>	BRAKE EXTENSION PAIR	£14
OTHER OPTIONS			
CEG1430	<input type="checkbox"/>	CANE HOLDER	£23
CEG1440	<input type="checkbox"/>	TIPPING AID RH	£14
CEG1441	<input type="checkbox"/>	TIPPING AID LH	£14
CEG1481	<input type="checkbox"/>	COMFORT ACRYLIC TRAY	£102
SWO10001	<input type="checkbox"/>	ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	£107

TOTAL PRICE	
--------------------	--

Non Conforming Options:

CEG0290	17" / 435 MM SEAT HEIGHT
CEG0300	16" / 410 MM SEAT HEIGHT
CEG0545	AMPUTEE LEGREST RH
CEG0595	AMPUTEE LEGREST LH
CEG0610	FLIP UP DEPTH AND ANGLE ADJ RH
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE66666	NO FOOTREST
CEG0630	FLIP UP DEPTH AND ANGLE ADJ LH
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE55555	NO FOOTREST
CEG0869	MATRX FLO-TECH LITE CUSHION - Low Risk
CEG0871	MATRX FLO-TECH CONTOUR VISCO CUSHION - High Risk
CEG0873	MATRX FLO-TECH IMAGE CUSHION - High Risk
CEG0874	MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk
CEG0876	MATRX LIBRA CUSHION - Very High Risk
CEG0878	FOR 400 MM SEAT DEPTH
CEG0879	FOR 450 MM SEAT DEPTH
CEG0852	SEAT RIGIDIZER
CEG1010	8" X 1.1/4" BLACK SOLID TYRE
CEG1062	12"
CEG1068	24"
CEG1071	STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES)
CEG1072	SELF-PROPEL REINFORCED REAR WHEEL
CEG1074	TRANSIT BLACK
CEG1084	PNEUMATIC PROFILE BLACK 1.3/8" / 35MM
CEG1088	PNEUMATIC LOW PROFILE BLACK 1" / 25MM
CEG1140	ALUMINIUM PRE-ANODIZED HANDRIM
CEG1150	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC
CEG1170	FIXED AXLE
CEG1180	Q/R AXLE
CEG1201	DUAL AMPUTEE REAR WHEEL POSITION
CEG1260	SPOKE GUARDS PAIR 24" BLK
CEG1270	PUSH TO LOCK BRAKES
CEG1275	PULL TO LOCK BRAKES
CEG1280	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)
CEG1362	BODYPOINT 2-PT POSTURE BELT L
CEG1363	BODYPOINT 4-PT POSTURE BELT MEDIUM
CEG1364	BODYPOINT 2-PT POSTURE BELT LARGE
CEG1370	REAR ANTI-TIPPERS PAIR
CEG1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID
CEG1390	BRAKE EXTENSION PAIR
CEG1440	TIPPING AID RH

CEG1441	TIPPING AID LH
CEG1481	COMFORT ACRYLIC TRAY



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):