



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JUNE 2023  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Rea Azalea Transit Tall

### Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Seat Height 500mm
- Depth Adjustable Seat 430mm-500mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +25°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Wide Pad Armrests with Autolock
- Headrest, Sideway Adjustable, Pearl Silver Frame Colour
- Chassis Extender with Solid Black Front Castors
- User Manoeuvred Tyre Brakes, Pelvic Belt, Adjustable both sides
- Anti Tippers, Front Transportation Brackets



**Standard Base price: £2,551**

Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
BNG9004	Rea Azalea Transit Tall	£2,551
BNG8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
<b>PROPULSION*</b>		
BNG0020	Transit 16" wheels with hub brakes	◆
<b>SEAT DIMENSIONS*</b>		
BNG0036	Tall	◆

SEAT WIDTH*				
BNG0040		390 mm/15" (Small chassis)	(NC 0244, 0245, 0246, 0247)	○
BNG0050		440 mm/17" (Medium chassis)	(NC 0241, 0242, 0246, 0247)	○
BNG0060		490 mm/19" (Large chassis)	(NC 0206, 0241, 0242, 0243, 0244)	○
SEAT WIDTH EXTENSION*				
BNG0065		+ 50 mm pre-adjustment on armrests and legrests		○
BNG0075		+ 100 mm pre-adjustment on armrests and legrests		○
BNG0085		Without seat width extension	(NC 0206)	○
SEAT DEPTH*				
BNG0105		480-550 mm / 19" -22" (Azalea Tall - Incl Chassis Extension)		◆
SEAT HEIGHT*				
BNG0125		500 mm / 20" (Azalea Tall)		◆
SEAT ANGLE*				
BNG0127		Negative Seat Angle -1 To +23°	(NC 0360)	○
BNG0000		Standard Seat Angle +1° to +25°		◆
SEAT TYPE*				
BNG0130		Depth Adjustable		◆
SEAT CUSHION*				
BNG0165		Cushion Flosshape seat upholstery	(NC 1111)	£57
BNG0167		Cushion Vicair Adjuster O2 10cm (Black)	(NC 0201, 0202)	£547
BNG0180		Without Seat Cushion, Seat Plate only	(NC 0201, 0202)	○
SEAT UPHOLSTERY*				
BNG0201		Multi-stretch Black PU coated fabric		◆
BNG0202		Black Spacer TR36		○
BNG1111		No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion		○
BACKREST WIDTH EXTENSION*				
BNG0205		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required	(NC 0221, 0223, 0225)	○
BNG0206		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required	(NC 0220, 0221, 0223, 0225)	○
BNG0207		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required		○
BACKREST*				
BNG0210		Flex 3 Backrest, Fixed Plate (Width Adjustable)	(NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0343)	◆
BNG0220		Flex 3 Backrest Tension (Width Adjustable)	(NC 0233, 0234, 0235, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0343)	£102
BNG0221		Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0305, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0920)	£494
BNG0223		Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0305, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0920)	£546
BNG0225		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)	(NC 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0920)	£353
BNG0230		Backrest Tubes Without Backrest	(NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0321, 0322, 0345, 0740, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0920)	○
MATRIX BACKREST HEIGHT				
BNG0233		Matrx Backrest Height 410 mm / 16"		○
BNG0235		Matrx Backrest Height 460 mm / 18"		○
BNG0234		Matrx Backrest Height 510 mm / 20"		○

<b>MATRIX BACKREST WIDTH</b>			
BNG0241		Matrix Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")	○
BNG0242		Matrix Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")	○
BNG0243		Matrix Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")	○
BNG0244		Matrix Backrest Width 430mm / 17" (Chest width 340mm / 13")	○
BNG0245		Matrix Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")	○
BNG0246		Matrix Backrest Width 480mm / 19" (Chest width 390mm / 15")	○
BNG0247		Matrix Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")	○
<b>BACKREST CUSHION*</b>			
BNG0240		Cover for tension adjustable upholstery (NC 2222)	○
BNG0250		Cushion, Laguna (Lateral Support) (NC 2222)	○
BNG0260		Cushion, Mistral 2 (Waist Support) (NC 0322, 2222)	£144
BNG0270		Cushion, Passad 2 (Shoulder Support) (NC 2222)	◆
BNG0300		Without Backrest Cover/ Cushion (NC 0305, 0321, 0322)	○
<b>WIDER BACKREST COVER</b>			
BNG0305		Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected	○
<b>BACKREST UPHOLSTERY*</b>			
BNG0321		Multi-stretch PU coated fabric	◆
BNG0322		Black Spacer TR36	○
BNG2222		No upholstery, select with Matrix, and Without backrest options	○
<b>PUSH HANDLES*</b>			
BNG0335		Push Handles connected	◆
BNG0343		Push Handles connected,with Rea head/neckrest attachment	£72
BNG0345		Push Bar, angle adjustable	£102
<b>SEAT TILT AND BACKREST REGULATION*</b>			
BNG0350		Assistant operated gas piston seat tilt and backrest recline (NC 5555)	◆
BNG0360		Electrical operated seat tilt and backrest recline (NC 1263)	£616
<b>TILT LOCK*</b>			
BNG1263		Tilt Lock (locking of pre-set angle)	◆
BNG5555		No Tilt Lock, Select with 0360 Electric operated seat and backrest	○
<b>CENTRAL LEGREST</b>			
BNG0370		Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444 (NC 0189, 0376, 0378, 0380, 0390, 0410, 0430, 0440, 0460, 0475, 0476, 0477, 0720, 0930, 0940, 0941, 0942, 0952, 0954, 0946, 1290)	£303
<b>LEGREST HANGER ATTACHMENT*</b>			
BNG0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function	£100
BNG0376		High position 50mm higher for leg length compensation	○
BNG0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest (NC 0390, 0410, 0440, 0460)	○

LEGREST HANGER LH*			
BNG0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH	(NC 0475, 0480, 0490, 0942, 0945, 6666) ○
BNG0390		Legrest Hanger fixed 80 deg, swing in/swing out LH	(NC 0928, 0929, 0930) ◆
BNG0410		Legrest Amputee LH	(NC 0475, 0476, 0477, 0480, 0490, 0510, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956) £137
BNG0420		Without Legrest, LH	(NC 0475, 0476, 0477, 0480, 0490, 0510, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956) ○
LEGREST HANGER RH*			
BNG0430		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH	(NC 0475, 0540, 0550, 0942, 0945, 6666) ◆
BNG0440		Legrest Hanger fixed 80 deg, swing in/swing out RH	(NC 0931, 0932, 0940) ○
BNG0460		Legrest Amputee RH	(NC 0475, 0476, 0477, 0540, 0550, 0570, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956) £137
BNG0470		Without Legrest, RH	(NC 0475, 0476, 0477, 0540, 0550, 0570, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956) ○
FOOTREST LH*			
BNG0475		Footrest one piece, depth and angle adjustable short tube	(NC 0477, 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940, 0952, 0954, 0956) £118
BNG0476		Footrest one piece, depth and angle adjustable long tube	(NC 0477, 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940, 0952, 0954, 0956) £118
BNG0477		Footbox, angle adjustable, padded with short tube	(NC 0480, 0490, 0510, 0520, 0540, 0550, 0570, 0580, 0930, 0940, 0941, 0942, 0945, 0952, 0954, 0956) £543
BNG0480		Footrest lower part, fixed with short tube and heel strap, LH	(NC 0930, 0941, 0952, 0954, 0956) ○
BNG0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.	○
BNG0510		Footrest lower part, fixed with long tube and heel strap, LH.	(NC 0930, 0941, 0952, 0954, 0956) ○
BNG0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNG3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○
FOOTREST RH*			
BNG0540		Footrest lower part, fixed with short tube and heel strap, RH	(NC 0940, 0941, 0952, 0954, 0956) ○
BNG0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.	○
BNG0570		Footrest lower part, fixed with long tube and heel strap, RH	(NC 0940, 0941, 0952, 0954, 0956) ○
BNG0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNG4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BNG0600		Armrest Ducktail Pad LH, height adjustable with autolock	(NC 1232) ○
BNG0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only	(NC 1232) ◆
BNG0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock	(NC 1230, 1231, 1232, 1250, 1260) £89
BNG0610		Armrest Padded LH, height adjustable with autolock	(NC 1230, 1231, 1232, 1250, 1260) £64
BNG0615		Armrest Hemiplegic LH, height adjustable with autolock	(NC 0695, 1230, 1231, 1232, 1240, 1245, 1250, 1260) £179
BNG0620		Without Armrest LH Select when choosing BNG1232 Half Tray LH, Narrow Armrest is included	(NC 0695, 1230, 1231, 1240, 1245, 1250, 1260) ○

ARMREST RH*				
BNG0630		Armrest Ducktail Pad RH, height adjustable with autolock	(NC 1233)	○
BNG0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only	(NC 1233)	◆
BNG0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock	(NC 1230, 1231, 1233, 1250, 1260)	£89
BNG0640		Armrest Padded RH, height adjustable with autolock.	(NC 1230, 1231, 1233, 1250, 1260)	£64
BNG0645		Armrest Hemiplegic, RH, height adjustable with autolock	(NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260)	£179
BNG0650		Without Armrest RH Select when choosing BNG1233 Half Tray RH, Narrow Armrest is included	(NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST OPTIONS. Side support cover 0680, must be selected if 0690 Foam for side support is required				
BNG0680		Side Support Cover, removable with tools (wipeable cover for the sideguard)	(NC 0695)	£51
BNG0690		Foam for side support cover (Internal cushioned pad, which is only possible in combination with 0680)	(NC 0695)	£51
BNG0695		Pad for sideguard, removable with hook and loop fastening		£51
SEAT OPTIONS				
BNG0710		Incontinence Cover		£40
BNG0720		Pommel Abduction Cushion		£118
BNG0725		Tilt Scale for seat tilt		£44
HEADREST OPTIONS*				
BNG0730		Headrest, Sideway Adjustable	(NC 0742, 0744, 0746, 0750, 0770, 0782)	◆
BNG0750		Headrest with cheek support, sideway adjustable Fixed hardware included	(NC 0740, 0742, 0744, 0746, 0740, 0770, 0782)	£107
BNG0731		Matrx Elan Standard Pad (254 mm / 10")	(NC 0740, 0770, 0782, 0790, 8888)	£124
BNG0732		Matrx Elan Standard Pad (355 mm / 14")	(NC 0740, 0770, 0782, 0790, 8888)	£124
BNG0733		Matrx Elan Occipital Pad (229 mm / 9")	(NC 0740, 0770, 0782, 0790, 8888)	£124
BNG0734		Matrx Elan Occipital Pad (305 mm / 12")	(NC 0740, 0770, 0782, 0790, 8888)	£124
BNG0735		Matrx Elan 4-Point Pad (229W x 177H mm / 9"W x 7"H)	(NC 0740, 0770, 0782, 0790, 8888)	£232
BNG0736		Matrx Elan 4-Point Pad (280W x 254H mm / 11"W x 10"H)	(NC 0740, 0770, 0782, 0790, 8888)	£232
BNG9991		Without headrest	(NC 0742, 0744, 0746)	○
HEADREST HARDWARE*				
BNG0742		Matrx Elan Hardware	(NC 0740, 0770, 0790)	£184
BNG0744		Matrx Elan Extended Hardware (Extra link)	(NC 0740, 0770, 0790)	£284
BNG0746		Matrx Elan LOXX Hardware	(NC 0740, 0770, 0790)	£166
BNG8888		Without Elan headrest hardware (Only select with option 0730 & 0750 Headrests)		○
OTHER HEADREST OPTIONS				
BNG0740		Cover for headrest		£31
BNG0770		Neckrest		○
BNG9992		Without neckrest		○
BNG0782		Alternative headrest/neckrest mounting bracket		£23
BNG0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)		○
TRUNK SUPPORT LH				
BNG0800		Trunk Support Multi Adjustable with standard arm LH	(NC 0900)	○
BNG0820		Trunk Support Multi Adjustable with extended height arm LH	(NC 0900)	○
BNG0880		Trunk Support Swingaway LH	(NC 0900)	£105
BNG0895		Trunk Support Multi Adjustable Swingaway LH	(NC 0900)	£113

TRUNK SUPPORT RH				
BNG0810		Trunk Support Multi Adjustable with standard arm RH	(NC 0900)	○
BNG0830		Trunk Support Multi Adjustable with extended height arm RH	(NC 0900)	○
BNG0890		Trunk Support Swingaway RH	(NC 0900)	£105
BNG0896		Trunk Support Multi Adjustable Swingaway RH	(NC 0900)	£113
BACKREST OPTIONS				
BNG0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included		◆
BNG0915		Tilt Scale for backrest recline		£44
BNG0920		Wedges (One Pair)		£11
ELEVATING LEGREST OPTIONS*				
BNG0928		Calf pad cover, for LH elevating legrest		○
BNG0929		Knee pad for LH elevating legrest		◆
BNG0931		Calf pad cover, for RH elevating legrest		○
BNG0932		Knee pad for RH elevating legrest		◆
BNG6666		Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)		○
LEGREST OPTIONS				
BNG0930		Footplate Extender LH	(NC 0941)	£34
BNG0940		Footplate Extender RH	(NC 0941)	£34
BNG0941		Cover Footplate, Black Plush Fabric		£54
BNG0942		Calf Strap Padded		£73
BNG0945		Calf Strap (for One piece footrest)		○
BODYPOINT ANKLE HUGGERS				
BNG0952		Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)		£114
BNG0954		Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)		£114
BNG0956		Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)		£114
CHASSIS COLOUR*				
BNG0944		Silver Pearl		◆
BNG0946		Azurite (Dark Blue)		○
BNG0947		Sand (Beige)		○
BNG0948		Electric Green		○
BNG0949		Happy Red		○
CHASSIS SAFETY*				
BNG0950		Chassis extender, castors extend 50mm forward for added stability		◆
CASTOR ATTACHMENT*				
BNG0960		Fixed castor attachment		◆
FRONT CASTORS*				
BNG1015		200x32 mm (8" x 1 1/4") Black Solid Tyre		£21
BNG1025		200x45 mm (8" X 1 3/4") Black Solid tyre		◆
REAR WHEEL SIZE*				
BNG1047		16" / 405 mm Transit rear wheel with fixed attachment		◆
REAR WHEEL TYPE*				
BNG1076		Hub Brake 16" Transit with Black Plastic Hub and Rim		○



<b>REAR WHEEL TYRE*</b>		
BNG1092		Solid Black tyres <span style="float: right;">◆</span>
<b>BRAKE*</b>		
BNG1170		User manoeuvred tyre brakes, includes extended handles <span style="float: right;">◆</span>
<b>PELVIC BELT*</b>		
BNG1200		Pelvic belt <span style="float: right;">◆</span>
BNG1201		Bodypoint 2-Point Centre-pull Push Button Medium (Mounted to chair) <span style="float: right;">£83</span>
BNG1202		Bodypoint 2-Point Centre-pull Push Button Large (Mounted to chair) <span style="float: right;">£83</span>
BNG1203		Bodypoint 4-Point Centre-pull Push Button Medium (Mounted to chair) <span style="float: right;">£98</span>
BNG1204		Bodypoint 4-Point Centre-pull Push Button Large (Mounted to chair) <span style="float: right;">£98</span>
<b>BODYPOINT HARNESS</b>		
BNG1206		Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted) <span style="float: right;">£109</span>
BNG1207		Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted) <span style="float: right;">£109</span>
<b>TRANSPORTATION BRACKETS*</b>		
BNG1225		Transportation Brackets, mounted to the front of the wheelchair <span style="float: right;">◆</span>
<b>OTHER OPTIONS</b>		
BNG1230		Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260) <span style="float: right;">£138</span>
BNG1231		Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260) <span style="float: right;">£138</span>
BNG1232		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260) <span style="float: right;">£103</span>
BNG1233		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260) <span style="float: right;">£103</span>
BNG1240		Swing away table without elbow pads (NC 1237, 1250, 1260) <span style="float: right;">£152</span>
BNG1245		Swing away table with elbow pads (NC 1237, 1260) <span style="float: right;">£152</span>
BNG1235		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray) <span style="float: right;">£51</span>
BNG1237		Half table tray padded cushion <span style="float: right;">£46</span>
BNG1250		Locking device for table tray <span style="float: right;">£15</span>
BNG1260		Mounting kit for table (if ordered separately) <span style="float: right;">£9</span>
BNG1290		Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment) <span style="float: right;">£55</span>
BNG1300		Harness attachment , must be selected in combination with Head/Neckrest <span style="float: right;">£85</span>
BNG1310		Attachment for positioning belt <span style="float: right;">£32</span>

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

BNG0201	Multi-stretch Black PU coated fabric
BNG0202	Black Spacer TR36
BNG1111	No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion
BNG0206	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BNG0220	Flex 3 Backrest Tension (Width Adjustable)
BNG0221	Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)
BNG0223	Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)
BNG0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BNG0233	Matrx Backrest Height 410 mm / 16"
BNG0235	Matrx Backrest Height 460 mm / 18"
BNG0234	Matrx Backrest Height 510 mm / 20"
BNG0241	Matrx Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")
BNG0242	Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")
BNG0243	Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")
BNG0244	Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")
BNG0245	Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")
BNG0246	Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")
BNG0247	Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")
BNG0240	Cover for tension adjustable upholstery
BNG0250	Cushion, Laguna (Lateral Support)
BNG0260	Cushion, Mistral 2 (Waist Support)
BNG0270	Cushion, Passad 2 (Shoulder Support)
BNG0305	Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected
BNG0321	Multi-stretch PU coated fabric
BNG0322	Black Spacer TR36
BNG2222	No upholstery, select with Matrx, and Without backrest options
BNG0335	Push Handles connected
BNG0343	Push Handles connected,with Rea head/neckrest attachment
BNG0345	Push Bar, angle adjustable
BNG0360	Electrical operated seat tilt and backrest recline
BNG1263	Tilt Lock (locking of pre-set angle)
BNG5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
SWO0189	Standard position with 0-30° Legrest Abduction function
BNG0376	High position 50mm higher for leg length compensation
BNG0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BNG0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BNG0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNG0410	Legrest Amputee LH
BNG0430	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH
BNG0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNG0460	Legrest Amputee RH
BNG0475	Footrest one piece, depth and angle adjustable short tube
BNG0476	Footrest one piece, depth and angle adjustable long tube



BNG0477	Footbox, angle adjustable, padded with short tube
BNG0480	Footrest lower part, fixed with short tube and heel strap, LH
BNG0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.
BNG0510	Footrest lower part, fixed with long tube and heel strap, LH.
BNG0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNG0540	Footrest lower part, fixed with short tube and heel strap, RH
BNG0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.
BNG0570	Footrest lower part, fixed with long tube and heel strap, RH
BNG0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNG0695	Pad for sideguard, removable with hook and loop fastening
BNG0720	Pommel Abduction Cushion
BNG0750	Headrest with cheek support, sideway adjustable Fixed hardware included
BNG0742	Matrx Elan Hardware
BNG0744	Matrx Elan Extended Hardware (Extra link)
BNG0746	Matrx Elan LOXX Hardware
BNG8888	Without Elan headrest hardware (Only select with option 0730 & 0750 Headrests)
BNG0740	Cover for headrest
BNG0770	Neckrest
BNG0782	Alternative headrest/neckrest mounting bracket
BNG0790	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BNG0800	Trunk Support Multi Adjustable with standard arm LH
BNG0820	Trunk Support Multi Adjustable with extended height arm LH
BNG0880	Trunk Support Swingaway LH
BNG0895	Trunk Support Multi Adjustable Swingaway LH
BNG0810	Trunk Support Multi Adjustable with standard arm RH
BNG0830	Trunk Support Multi Adjustable with extended height arm RH
BNG0890	Trunk Support Swingaway RH
BNG0896	Trunk Support Multi Adjustable Swingaway RH
BNG0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNG0920	Wedges (One Pair)
BNG0928	Calf pad cover, for LH elevating legrest
BNG0929	Knee pad for LH elevating legrest
BNG0931	Calf pad cover, for RH elevating legrest
BNG0932	Knee pad for RH elevating legrest
BNG6666	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)
BNG0930	Footplate Extender LH
BNG0940	Footplate Extender RH
BNG0941	Cover Footplate, Black Plush Fabric
BNG0942	Calf Strap Padded
BNG0945	Calf Strap (for One piece footrest)
BNG0952	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)
BNG0954	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)
BNG0956	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)

BNG0946	Azurite (Dark Blue)
BNG1230	Fixed Table Tray without elbow pads includes mounting kit
BNG1231	Fixed Table Tray with elbow Pads includes mounting kit
BNG1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNG1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNG1240	Swing away table without elbow pads
BNG1245	Swing away table with elbow pads
BNG1235	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BNG1237	Half table tray padded cushion
BNG1250	Locking device for table tray
BNG1260	Mounting kit for table (if ordered separately)
BNG1290	Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):