

Invoice to:
 Name:
 Street:
 Town: Postal code:
 Tel: Fax:

Deliver to:
 Name:
 Street:
 Town: Postal code:
 Tel: Fax:

Order date:
 Purchase order:

Marked for:
 Sunrise order no:

Jay® ION Features
 Medium risk skin protection cushion
 Contoured moulded foam cushion with viscoelastic top layer
 Water Resistant Dartex inner cover
 Microclimatic outer cover as standard
2 year guarantee, 6 months on cover

Base Price **£180**

E-mail: orders@sunmed.co.uk

Fax: 0845 605 66 89

Maximum user weight 136kg (SW 340 mm - 500 mm)
 226kg (SW 540 mm - 600 mm)

Application:

Client with mild assymmetric posture and mild postural support needs. Clients with moderate pressure sore and low shear risk (medium risk skin protection).

Provides mild-moderate lateral and medial thigh support.

CLIENT MEASUREMENT

SEAT WIDTH

SEAT DEPTH	cm	SEAT WIDTH										
		34	38	40	42	44	46	50	52	54	56	60
JION__38	38	◇	◇	◇	◇	◇	■	■	■	■	■	■
JION__40	40	◇	○	○	◇	◇	◇	■	■	■	■	■
JION__42	42	○	○	○	○	◇	◇	■	■	■	■	■
JION__44	44	◇	○	○	○	○	◇	◇	◇	■	■	■
JION__46	46	◇	◇	○	○	○	○	◇	◇	■	■	■
JION__50	50	■	◇	◇	○	○	○	◇	◇	◇	◇	◇
JION__52	52	■	◇	◇	◇	○	○	◇	◇	◇	◇	◇
JION__54	54	■	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇
JION__56	56	■	■	◇	◇	◇	◇	◇	◇	◇	◇	◇
JION__60	60	■	■	■	◇	◇	◇	◇	◇	◇	◇	◇
JION__62	62	■	■	■	■	◇	◇	◇	◇	◇	◇	◇

Price

£180.00

Product Code

JION ____

- 5 day lead time
- ◇ 15 day lead time
- Not available

Example		Width 42cm	Depth 44cm	JION4244 = Complete Cushion
---------	--	------------	------------	------------------------------------

SPARE CUSHION COVER				Unit Price	Item Code
CJION____EU	<input type="checkbox"/>	Spare Incontinence Resistant Inner Cover		£70.00	CJION____EU
CJION____MEU	<input type="checkbox"/>	Spare Microclimatic Outer Cover		£90.00	CJION____MEU

JAY ION	£180
Option Price	£
VAT	£
Total Price	£



JAY ION CUSHION

Stability and positioning with skin integrity protection

Product Features:

- Medium Skin protection risk
- Contoured foam base with anatomical seat well and femoral loading
- Integral lateral foam inserts for enhanced lateral hip and thigh stability
- Soft memory foam top layer for even pressure distribution, heat dissipation and comfort
- Incontinence resistant inner cover
- Microclimatic outer cover as standard



JAY UNION CUSHION

JION3438	JAY Ion Cushion - 34x38	£180
JION3440	JAY Ion Cushion - 34x40	£180
JION3442	JAY Ion Cushion - 34x42	£180
JION3444	JAY Ion Cushion - 34x44	£180
JION3446	JAY Ion Cushion - 34x46	£180
JION3838	JAY Ion Cushion - 38x38	£180
JION3840	JAY Ion Cushion - 38x40	£180
JION3842	JAY Ion Cushion - 38x42	£180
JION3844	JAY Ion Cushion - 38x44	£180
JION3846	JAY Ion Cushion - 38x46	£180
JION3850	JAY Ion Cushion - 38x50	£180
JION3852	JAY Ion Cushion - 38x52	£180
JION3854	JAY Ion Cushion - 38x54	£180
JION4038	JAY Ion Cushion - 40x38	£180
JION4040	JAY Ion Cushion - 40x40	£180
JION4042	JAY Ion Cushion - 40x42	£180
JION4044	JAY Ion Cushion - 40x44	£180
JION4046	JAY Ion Cushion - 40x46	£180
JION4050	JAY Ion Cushion - 40x50	£180
JION4052	JAY Ion Cushion - 40x52	£180
JION4054	JAY Ion Cushion - 40x54	£180
JION4056	JAY Ion Cushion - 40x56	£180
JION4238	JAY Ion Cushion - 42x38	£180
JION4240	JAY Ion Cushion - 42x40	£180
JION4242	JAY Ion Cushion - 42x42	£180
JION4244	JAY Ion Cushion - 42x44	£180
JION4246	JAY Ion Cushion - 42x46	£180
JION4250	JAY Ion Cushion - 42x50	£180
JION4252	JAY Ion Cushion - 42x52	£180
JION4254	JAY Ion Cushion - 42x54	£180
JION4256	JAY Ion Cushion - 42x56	£180
JION4260	JAY Ion Cushion - 42x60	£180
JION4438	JAY Ion Cushion - 44x38	£180
JION4440	JAY Ion Cushion - 44x40	£180
JION4442	JAY Ion Cushion - 44x42	£180
JION4444	JAY Ion Cushion - 44x44	£180
JION4446	JAY Ion Cushion - 44x46	£180

JION4450	JAY Ion Cushion - 44x50	£180
JION4452	JAY Ion Cushion - 44x52	£180
JION4454	JAY Ion Cushion - 44x54	£180
JION4456	JAY Ion Cushion - 44x56	£180
JION4460	JAY Ion Cushion - 44x60	£180
JION4462	JAY Ion Cushion - 44x62	£180
JION4640	JAY Ion Cushion - 46x40	£180
JION4642	JAY Ion Cushion - 46x42	£180
JION4644	JAY Ion Cushion - 46x44	£180
JION4646	JAY Ion Cushion - 46x46	£180
JION4650	JAY Ion Cushion - 46x50	£180
JION4652	JAY Ion Cushion - 46x52	£180
JION4654	JAY Ion Cushion - 46x54	£180
JION4656	JAY Ion Cushion - 46x56	£180
JION4660	JAY Ion Cushion - 46x60	£180
JION4662	JAY Ion Cushion - 46x62	£180
JION5040	JAY Ion Cushion - 50x40	£180
JION5042	JAY Ion Cushion - 50x42	£180
JION5044	JAY Ion Cushion - 50x44	£180
JION5046	JAY Ion Cushion - 50x46	£180
JION5050	JAY Ion Cushion - 50x50	£180
JION5052	JAY Ion Cushion - 50x52	£180
JION5054	JAY Ion Cushion - 50x54	£180
JION5056	JAY Ion Cushion - 50x56	£180
JION5060	JAY Ion Cushion - 50x60	£180
JION5062	JAY Ion Cushion - 50x62	£180
JION5244	JAY Ion Cushion - 52x44	£180
JION5246	JAY Ion Cushion - 52x46	£180
JION5250	JAY Ion Cushion - 52x50	£180
JION5252	JAY Ion Cushion - 52x52	£180
JION5254	JAY Ion Cushion - 52x54	£180
JION5256	JAY Ion Cushion - 52x56	£180
JION5260	JAY Ion Cushion - 52x60	£180
JION5262	JAY Ion Cushion - 52x62	£180
JION5444	JAY Ion Cushion - 54x44	£180
JION5446	JAY Ion Cushion - 54x46	£180
JION5450	JAY Ion Cushion - 54x50	£180
JION5452	JAY Ion Cushion - 54x52	£180
JION5454	JAY Ion Cushion - 54x54	£180
JION5456	JAY Ion Cushion - 54x56	£180
JION5460	JAY Ion Cushion - 54x60	£180
JION5462	JAY Ion Cushion - 54x62	£180
JION5650	JAY Ion Cushion - 56x50	£180
JION5652	JAY Ion Cushion - 56x52	£180
JION5654	JAY Ion Cushion - 56x54	£180
JION5656	JAY Ion Cushion - 56x56	£180
JION5660	JAY Ion Cushion - 56x60	£180
JION5662	JAY Ion Cushion - 56x62	£180
JION6050	JAY Ion Cushion - 60x50	£180
JION6052	JAY Ion Cushion - 60x52	£180
JION6054	JAY Ion Cushion - 60x54	£180
JION6056	JAY Ion Cushion - 60x56	£180
JION6060	JAY Ion Cushion - 60x60	£180
JION6062	JAY Ion Cushion - 60x62	£180