

Invoice to:
 Name:
 Street:
 Town: Postal code:
 E-Mail:
 Tel: Fax:

Deliver to:
 Name:
 Street:
 Town: Postal code:
 E-Mail:
 Tel: Fax:

Order date: _____ **Marked for:** _____





J3 Backrest Range
 Complete backrest package
 Easily adjustable for assessment and changing needs
 Wide range of positioning and Whitmyer head support options
 No loss of seat depth, when installed in most rear set position
 Backs and options fully transit tested
 Positioning components: 1 large, 1 small lumbar insert included with MT, UT and SH
 1 small lumbar included with LT

Email : orders@sunmed.co.uk
Fax: 0845 605 66 89

Maximum user weight
 Compact Hardware 136Kg
 Heavy Duty Compact Hardware 225Kg

1. Choose contour depth (SC, MC, DC or MDC)
2. Select backrest height range (LT; MT; UT or SH)
3. Select height within range (Short, Medium, Tall)
4. Select backrest size based on client's chest width, cross check with outside width of wheelchair backposts.
5. Add needed postural support options, see last 2 pages



					
		depth 5 cm	depth 8 cm	depth 15/18 cm*	depth 15cm/18 cm*
		Shallow Contour (SC)	Mid Contour (MC)	Deep Contour (DC)	Mid Deep Contour (MDC)
JAY J3 Back Size		31 36 41 46 51	31 36 41 46 51	36 41 46 51	31 36 41 46 51
Backrest Height		£705			£870
Lower Thoracic (LT)	Short (S) - 17	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Med. (M) - 20	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Tall (T) - 24	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Backrest Height		£705			£870
Mid Thoracic (MT)	Short (S) - 30	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Med. (M) - 34	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Tall (T) - 38	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
X Deep MT	Tall (T) X - 38**				<input type="checkbox"/>
Backrest Height		£735		£870	£870
Upper Thoracic (UT)*	Short (S) - 42	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Med. (M) - 46	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Med. (M) - 46**			<input type="checkbox"/>	<input type="checkbox"/>
X Deep UT	Tall (T) - 50	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Backrest Height		£735		£870	£870
Shoulder Height (SH)*	Short (S) - 53	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Med. (M) - 57	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Tall (T) - 61	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

*Laterals, headrest & harnesses available for UT & SH heights

**Extra Deep option = 18 cm contour depth



Invoice to:
 Name:
 Street: Postal code:
 Town:
 E-Mail: Fax:
 Tel:

Deliver to:
 Name:
 Street: Postal code:
 Town:
 E-Mail: Fax:
 Tel:

Order date:
 Purchase order:

Marked for

Product Delivery times 5 days 15 days not available

- Mild lateral support
- Contour depth: 5 cm
- Max. user weight: 136 kg (Compact hardware)
225 kg (HD hardware)
- Laterals, headrest & harnesses available for UT & SH heights



LT and MT backs delivered with Compact Hardware
 UT and SH backs delivered with Heavy Duty Hardware

Jay J3 Back Size	31	36	41	46	51
Chair widths (measured from outside to outside of backcanes)	31-36 cm	36-41 cm	41-46 cm	46-51 cm	51-56 cm
Client max trunk width	28 cm	33 cm	38 cm	43 cm	48 cm

↓ Backrest Height ↓		£705				
Lower Thoracic (LT)	Short (S) - 17	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	Medium (M) - 20	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	Tall (T) - 24	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

↓ Backrest Height ↓		£705				
Mid Thoracic (MT)	Short (S) - 30	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	Medium (M) - 34	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	Tall (T) - 38	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

↓ Backrest Height ↓		£735				
Upper Thoracic (UT)	Short (S) - 42	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	Medium (M) - 46	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	Tall (T) - 50	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

↓ Backrest Height ↓		£735				
Shoulder Height (SH)	Short (S) - 53	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	Medium (M) - 57	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	Tall (T) - 61	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

*Free of charge HD hardware option available for 46cm and 51 cm widths

Invoice to:
 Name:
 Street: Postal code:
 Town:
 E-Mail: Fax:
 Tel:

Deliver to:
 Name:
 Street: Postal code:
 Town:
 E-Mail: Fax:
 Tel:

Order date:
 Purchase order:

Marked for

Product Delivery times 5 days 15 days not available

- Moderate lateral support
- Contour depth: 8 cm
- Max. user weight: 136 kg (Compact hardware)
225 kg (HD hardware)
- Laterals, headrest & harnesses available for UT & SH heights



LT and MT backs delivered with Compact Hardware
 UT and SH backs delivered with Heavy Duty Compact Hardware

Jay J3 Back Size	31	36	41	46	51
Chair widths (measured from outside to outside of backcanes)	31-36 cm	36-41 cm	41-46 cm	46-51 cm	51-56 cm
Client max trunk width	24.5 cm	29.5 cm	34.5 cm	39.5 cm	44.5 cm

Backrest Height		31	36	41	46	51
Lower Thoracic (LT)	Short (S) - 17					
	Medium (M) - 20					
	Tall (T) - 24					

Backrest Height		705				
Mid Thoracic (MT)	Short (S) - 30	<input type="radio"/> J3MCMTS31	<input type="radio"/> J3MCMTS36	<input type="radio"/> J3MCMTS41	<input type="radio"/> J3MCMTS46	<input type="radio"/> J3MCMTS51
					<input type="radio"/> J3MCMTS46HD (*)	<input type="radio"/> J3MCMTS51HD (*)
	Medium (M) - 34	<input type="radio"/> J3MCMTM31	<input type="radio"/> J3MCMTM36	<input type="radio"/> J3MCMTM41	<input type="radio"/> J3MCMTM46	<input type="radio"/> J3MCMTM51
					<input type="radio"/> J3MCMTM46HD (*)	<input type="radio"/> J3MCMTM51HD (*)
		<input type="radio"/> J3MCMTT31	<input type="radio"/> J3MCMTT36	<input type="radio"/> J3MCMTT41	<input type="radio"/> J3MCMTT46	<input type="radio"/> J3MCMTT51
					<input type="radio"/> J3MCMTT46HD (*)	<input type="radio"/> J3MCMTT51HD (*)

Backrest Height		£735				
Upper Thoracic (UT)	Short (S) - 42	<input type="radio"/> J3MCUTS31	<input type="radio"/> J3MCUTS36	<input type="radio"/> J3MCUTS41	<input type="radio"/> J3MCUTS46	<input type="radio"/> J3MCUTS51
	Medium (M) - 46		<input type="radio"/> J3MCUTM36	<input type="radio"/> J3MCUTM41	<input type="radio"/> J3MCUTM46	<input type="radio"/> J3MCUTM51
	Tall (T) - 50		<input type="radio"/> J3MCUTT36	<input type="radio"/> J3MCUTT41	<input type="radio"/> J3MCUTT46	<input type="radio"/> J3MCUTT51

Backrest Height		£735				
Shoulder Height (SH)	Short (S) - 53		<input type="radio"/> J3MCSHS36	<input type="radio"/> J3MCSHS41	<input type="radio"/> J3MCSHS46	<input type="radio"/> J3MCSHS51
	Medium (M) - 57			<input type="radio"/> J3MCSHM41	<input type="radio"/> J3MCSHM46	<input type="radio"/> J3MCSHM51
	Tall (T) - 61				<input type="radio"/> J3MCSHT46	<input type="radio"/> J3MCSHT51

*Free of charge HD hardware option available for 46cm and 51 cm widths

J3 Deep Contour (DC)

Invoice to:
 Name:
 Street: Postal code:
 Town:
 E-Mail: Fax:
 Tel:

Deliver to:
 Name:
 Street: Postal code:
 Town:
 E-Mail: Fax:
 Tel:

Order date:
 Purchase order:


Marked for

Product Delivery times 5 days 15 days not available

- Deep lateral support
- Contour depth: 15 cm or 18cm
- Max. user weight: 136 kg (compact hardware) 225 kg (HD hardware)
- Laterals, headrest & harnesses available for UT & SH heights



All DC backs delivered with Heavy Duty Compact Hardware as standard



Jay J3 Back Size	
Chair widths (measured from outside to outside of backcanes)	
Client max trunk width	

36	41	46	51
36-41 cm	41-46 cm	46-51 cm	51-56 cm
28 cm	33 cm	38 cm	43 cm

↓ Backrest Height ↓	
Lower Thoracic (LT)	Short (S) - 17
	Medium (M) - 20
	Tall (T) - 24

↓ Backrest Height ↓	
Mid Thoracic (MT)	Short (S) - 30
	Medium (M) - 34
	Tall (T) - 38

↓ Backrest Height ↓	
Upper Thoracic Height (UT)	Short (S) - 42
	Medium (M) - 46
	Medium Extra Deep* (T) - 46
	Tall (T) - 50

870			
<input type="radio"/> J3DCUTS36	<input type="radio"/> J3DCUTS41	<input type="radio"/> J3DCUTS46	<input type="radio"/> J3DCUTS51
<input type="radio"/> J3DCUTM36	<input type="radio"/> J3DCUTM41	<input type="radio"/> J3DCUTM46	<input type="radio"/> J3DCUTM51
	<input type="radio"/> J3XDCUTM41	<input type="radio"/> J3XDCUTM46	<input type="radio"/> J3XDCUTM51
<input type="radio"/> J3DCUTT36	<input type="radio"/> J3DCUTT41	<input type="radio"/> J3DCUTT46	<input type="radio"/> J3DCUTT51

↓ Backrest Height ↓	
Shoulder Height (SH)	Short (S) - 53
	Medium (M) - 57
	Tall (T) - 61

870			
<input type="radio"/> J3DCSHS36	<input type="radio"/> J3DCSHS41	<input type="radio"/> J3DCSHS46	<input type="radio"/> J3DCSHS51
	<input type="radio"/> J3DCSHM41	<input type="radio"/> J3DCSHM46	<input type="radio"/> J3DCSHM51
		<input type="radio"/> J3DCSHT46	<input type="radio"/> J3DCSHT51

*Countour depth = 18 cm

Invoice to:		Deliver to:	
Name:		Name:	
Street:	Postal code:	Street:	Postal code:
Town:		Town:	
E-Mail:	Fax:	E-Mail:	Fax:
Tel:		Tel:	


Order date:	Marked for:
Purchase order:	

Product Delivery times

○ 5 days	◇ 15 days
----------	-----------

- Contour depth: 15 cm (Std) 18 cm (Extra Deep)
- Max. user weight: 136kg with Compact Hardware, 225 kg with Heavy Duty Compact Hardware
- Laterals, headrest & harnesses available for UT & SH heights

MT backs delivered with Compact Hardware
 UT and SH backs delivered with Heavy Duty Compact Hardware



JAY J3 Back Size
Chair widths (measured from outside to outside of backcanes)
Client max trunk width

31	36	41	46	51
31-36 cm	36-41 cm	41-46 cm	46-51 cm	51-56 cm
21 cm	28 cm	33 cm	38 cm	43 cm

↓ Backrest Height ↓	
Lower Thoracic (LT)	Short (S) - 17
	Medium (M) - 20
	Tall (T) - 24

870				
	○	○	○	○
	J3MDCLTM36	J3MDCLTM41	J3MDCLTM46	
	○	○	○	○
	J3MDCLTT36	J3MDCLTT41	J3MDCLTT46	J3MDCLTT51
				◇
				J3MDCLTT51HD(*)

↓ Backrest Height ↓	
Mid Thoracic (MT)	Short (S) - 30
	Medium (M) - 34
	Tall (T) - 38
	Tall Extra Deep (T)** - 38

870				
○	○	○	○	○
J3MDCMTS31	J3MDCMTS36	J3MDCMTS41	J3MDCMTS46	J3MDCMTS51
			◇	◇
			J3MDCMTS46HD(*)	J3MDCMTS51HD(*)
○	○	○	○	○
J3MDCMTM31	J3MDCMTM36	J3MDCMTM41	J3MDCMTM46	J3MDCMTM51
			○	○
			J3MDCMTM46HD(*)	J3MDCMTM51HD(*)
○	○	○	○	○
J3MDCMTT31	J3MDCMTT36	J3MDCMTT41	J3MDCMTT46	J3MDCMTT51
			◇	◇
			J3MDCMTT46HD(*)	J3MDCMTT51HD(*)
		○	○	○
		J3XMDCMTT41	J3XMDCMTT46	J3XMDCMTT51
			○	○
			J3XMDCMTT46HD(*)	J3XMDCMTT51HD(*)

↓ Backrest Height ↓	
Upper Thoracic (UT)	Short (S) - 42
	Medium (M) - 46
	Medium Extra Deep** (M) - 46
	Tall (T) - 50

870				
○	○	○	○	○
J3MDCUTS31	J3MDCUTS36	J3MDCUTS41	J3MDCUTS46	J3MDCUTS51
	○	○	○	○
	J3MDCUTM36	J3MDCUTM41	J3MDCUTM46	J3MDCUTM51
		○	○	○
		J3XMDCUTM41	J3XMDCUTM46	J3XMDCUTM51
	○	○	○	○
	J3MDCUTT36	J3MDCUTT41	J3MDCUTT46	J3MDCUTT51

↓ Backrest Height ↓	
Shoulder Height (SH)	Short (S) - 53
	Medium (M) - 57
	Tall (T) - 61

870				
○	○	○	○	○
J3MDCSHS31	J3MDCSHS36	J3MDCSHS41	J3MDCSHS46	J3MDCSHS51
		○	○	○
		J3MDCSHM41	J3MDCSHM46	J3MDCSHM51
			○	○
			J3MDCSHT46	J3MDCSHT51

*Free of charge HD hardware option available for 46 cm and 51 cm widths

**Contour depth of 18 cm

JAY J3 Back - ACCESSORIES



ADDITIONAL / DIFFERENT MOUNTING BRACKETS							Unit Price				
Order both hardware and mounting brackets if changing complete mounting hardware on Backrest											
Compact Hardware											
J3EUCHRDW	<input type="checkbox"/>	Compact hardware (pair) - clamps for back tube					£40				
J3EUCHRDWS	<input type="checkbox"/>	Mounting brackets (pair) - parts which attach directly to the backrest shell					£40				
Heavy Duty Hardware											
J3EUHDCHRDW	<input type="checkbox"/>	New HD hardware (pair) - clamps for back tube					£40				
J3EUHDCHRDWS	<input type="checkbox"/>	New brackets (pair) for HD hardware -parts that go on the backrest shell					£40				
Old style Hardware											
J3EUSHS	<input type="checkbox"/>	OLD Standard hardware (pair) - available only as a spare (short, for LT and MT backs) (until Jan 2021)					£40				
J3EUSHST	<input type="checkbox"/>	OLD Standard hardware (pair) - available only as a spare (tall, for UT and SH backs) (until Jan 2021)					£40				
ACCESSORY BAG											
					cm width	32	36	41	46	51	Unit Price
ABJ3BCFww	Accessory bag				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£50
HEADREST (CRASH TESTED) - choose 1 mount and 1 pad to have a complete headrest with pad							Unit Price				
Mount											
J3EUHDXH	<input type="checkbox"/>	Whitmyer AXYS Mount and Bracket for J3 Back					£120				
J3EUHDFH	<input type="checkbox"/>	Whitmyer Cobra Flip Back Mount and Bracket for J3 Back					£250				
Pad											
P1000L	<input type="checkbox"/>	Whitmyer Plush Pad 254 mm (replaces discontinued J3EUHDLSP)					£60				
P300L	<input type="checkbox"/>	Whitmyer Contoured Cradle Pad 250 mm (replaces discontinued J3EUHDLCP)					£90				
Accessory											
J3EUUHM	<input type="checkbox"/>	Universal Headrest Mount (required to use any other headrest not listed in this order form)					£55				
LATERAL SUPPORT (only available where laterals are indicated on the first page)							Unit Price				
		<- Left Hand Lateral		Right Hand Lateral ->							
	Fixed	10 cm Depth	15 cm Depth	10 cm Depth	15 cm Depth						
10 cm Height	<input type="checkbox"/>	J3LAT1010FXLH	<input type="checkbox"/>	J3LAT1510FXLH	<input type="checkbox"/>	J3LAT1010FXRH	<input type="checkbox"/>	J3LAT1510FXRH	£65		
12 cm Height	<input type="checkbox"/>	J3LAT1012FXLH	<input type="checkbox"/>	J3LAT1512FXLH	<input type="checkbox"/>	J3LAT1012FXRH	<input type="checkbox"/>	J3LAT1512FXRH	£65		
15 cm Height	<input type="checkbox"/>	J3LAT1015FXLH	<input type="checkbox"/>	J3LAT1515FXLH	<input type="checkbox"/>	J3LAT1015FXRH	<input type="checkbox"/>	J3LAT1515FXRH	£65		
	Swing-away	10 cm Depth	15 cm Depth	10 cm Depth	15 cm Depth						
10 cm Height	<input type="checkbox"/>	JLATSA1010LH	<input type="checkbox"/>	JLATSA1510LH	<input type="checkbox"/>	JLATSA1010RH	<input type="checkbox"/>	JLATSA1510RH	£110		
12 cm Height	<input type="checkbox"/>	JLATSA1012LH	<input type="checkbox"/>	JLATSA1512LH	<input type="checkbox"/>	JLATSA1012RH	<input type="checkbox"/>	JLATSA1512RH	£110		
15 cm Height	<input type="checkbox"/>	JLATSA1015LH	<input type="checkbox"/>	JLATSA1515LH	<input type="checkbox"/>	JLATSA1015RH	<input type="checkbox"/>	JLATSA1515RH	£110		
POSITIONING KIT											
					cm width	36	41	46	51	Unit Price	
J3SPNK1ww	Spine Align Mini Fit Kit I		(See back of the order form for content)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£40	
J3SPNK2ww	Spine Align Mini Fit Kit II		(See back of the order form for content)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£60	
J3SPNKEU	<input type="checkbox"/>	JAY J3 Spine Align Assessment Kit			(See last page for content)					£190	
J3HWKEU	<input type="checkbox"/>	JAY J3 Back Hardware Box Kit			(See last page for content)					£600	
CHEST STRAP (only available where indicated on the overview page)							Small	Medium	Large	<larg	Unit Price
J3HRNPASw	Padded Anterior Strap				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		£220	
J3HRNCOHw	Centre open harness				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		£220	
SPARES (separate to product - separate box or delivery)											
					CM width	36	41	46	51	Unit Price	
J3SPNSLww	Small Lumbar		(1 unit included as std with product LT, MT, UT & SH)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£10	
J3SPNLLww	Large Lumbar Support		(1 unit included as std with product MT, UT & SH)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£10	
J3SPNWEww	Wedge				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£10	
							Unit price	Qty	Total Price		
J3SPNCRS	<input type="checkbox"/>	Curved Shim					£10				
J3SPNFLS	<input type="checkbox"/>	Flat Shim					£10				
J3SPNPLP	<input type="checkbox"/>	Posterior/Lateral Pelvic					£15				
J3SPNLLS	<input type="checkbox"/>	Lateral Shim					£5				
J3SPNSBB	<input type="checkbox"/>	Small Bean Bag					£10				
J3SPNLBB	<input type="checkbox"/>	Large Bean Bag					£10				
Jay 3 Back - Basic Price							£				
Option price							£				
VAT.							£				
Total Price							£				

JAY J3 Back - ASSESSMENT KITS & BOXES



Spine Align Mini Fit Kit I			Quantity
J3SPNK1ww	Spine Align Mini Fit Kit I (Sizes available: 36cm; 41cm; 46cm and 51cm)		in the kit
J3SPNWEww	Wedge (Size will be as define for the Spine Align Mini Fit Kit I)		3
J3SPNLLS	Lateral Shim		2
J3SPNPLP	Posterior/Lateral Pelvic		2
Spine Align Mini Fit Kit II			Quantity
J3SPNK2ww	Spine Align Mini Fit Kit II (Sizes available: 36cm; 41cm; 46cm and 51cm)		in the kit
J3SPNWEww	Wedge (Size will be as define for the Spine Align Mini Fit Kit I)		5
J3SPNLLS	Lateral Shim		2
J3SPNPLP	Posterior/Lateral Pelvic		2
J3SPNLLww	Large Lumbar Support (Size will be as define for the Spine Align Mini Fit Kit II)		1
J3SPNFLS	Flat Shim		1
J3SPNCRS	Curved Shim		1
JAY J3 Back Spine Align Assessment Kit			Quantity
J3SPNKEU	JAY J3 Spine Align Assessment Kit		in the kit
J3SPNLL36	Large Lumbar Support (always 36cm, as Jay 3 Back contains lumbar in correct size)		2
J3SPNSL36	Small Lumbar (always 36cm, as Jay 3 Back contains lumbar in correct size)		2
J3SPNWE36	Wedge (always 36cm, as Jay 3 Back contains wedge in correct size)		16
J3SPNCRS	Curved Shim		2
J3SPNFLS	Flat Shim		2
J3SPNPLP	Posterior/Lateral Pelvic		2
J3SPNLLS	Lateral Shim		2
J3SPNSBB	Small Bean Bag		2
J3SPNLBB	Large Bean Bag		2
JAY J3 Back Hardware Box			Quantity
J3HWKEU	JAY J3 Back Hardware Box		in the kit
J3EUAHB	Backrest depth extension		2
J3EUHDXH	Headrest Standard (bracket & mount for J3 Back)		1
J3EUHDFH	Headrest Flip Back (bracket & mount for J3 Back)		1
J3EUHDLSP	Headrest Standard Pad for J3 Back		1
J3EUHDLCP	Headrest Contour Pad for J3 Back		1
J3EULATFB	J3 Fixed Lateral Bracket		2
J3EULATFP	J3 Fixed Lateral Bracket (hardware kit)		2
JEULATSABL	J3 S/A Lateral Bracket Left		1
JEULATSABR	J3 S/A Lateral Bracket Right		1
J3EULATSAP	J3 S/A Lateral Bracket (hardware kit)		2
to creat your complete Lateral use the Fix or the S/A bracket (above) and add the lateral pad of your choice			
	10 x 10 Lateral pad		2
	15 x 12 Lateral pad		2
	15 x 15 Lateral pad		2
J3HRNCOM	Chest strap harnesses M		1
J3HRNPSM	Chest strap padded anterior strap M		1
J3HRNPSX	Chest strap padded anterior strap XL		1
J3EUUHM	Universal Headrest mount		1