

# JAY Zip Back

Clinical seating, designed for kids



The JAY Zip backrest was designed through the use of anthropometric data so it is specially sized to meet a paediatric client's unique needs. It is an ultra lightweight, anti-microbial backrest that features X-static® and Aquaguard™ cover technologies as well as adjustable, one-step release hardware to accommodate a wide range of mounting locations and user requirements.

## Versatile Hardware

- Easy one-step release
- Fits on a wide range of tubing: 1.9cm (0.7") with inserts; 2.2 - 2.3cm (0.8 - 0.9") with inserts or 2.5cm (1") without insert
- One tool adjustments (and the tool is provided!)

## Convenient & Anti-microbial

### Outer Cover

- Machine washable and quick to dry
- Silver impregnated X-static® stretch fabric naturally inhibits microbial growth
- Extra layer of microclimatic spacer fabric for comfort
- Safe non-toxic material

### Inner Cover

- Easy to wipe off
- Aquaguard™ zipper resists moisture
- Anti-wicking seam thread protects foam base



One step release

<b>Backrest width:</b>	<b>20, 25, 30, 35cm (8, 10, 12, 14")</b>	<b>Back angle:</b>	<b>+20° to -20°</b>
<b>Backrest height:</b>	<b>15, 20, 25, 30cm (6, 8, 10, 12")</b>	<b>Width adjustability:</b>	<b>+5cm (2")</b>
<b>Transport weight:</b>	<b>1.1kg</b>	<b>Depth adjustability:</b>	<b>7.5 cm (3")</b>
<b>Max. user weight:</b>	<b>75kg</b>	<b>Accessories: Headrest, lateral support, spine align positioning components, chest strap, vanity flap and spare cover</b>	



### Sunrise Medical Limited

Thorns Road  
 Brierley Hill  
 West Midlands  
 DY5 2LD  
 Tel.: +44 (0) 845 605 66 88  
 Fax: +44 (0) 845 605 66 89  
 www.SunriseMedical.co.uk