

# JAY Zip Cushion

A comfortable, clinically effective cushion designed uniquely for kids.



The JAY Zip cushion was designed through the use of anthropometric data so it is specially sized to meet a paediatric client's unique needs. It is a dual layered contoured foam base that provides low/moderate skin protection and mild positioning in a comfortable, effective and convenient cushion.

## Comfortable

- Dual layered contoured foam base
- 3DX™ fabric spacer in the outer cover promotes air movement
- Soft and stretchable outer cover (black as standard, also available in pink and blue)



## Convenient & Antimicrobial

### Outer Cover

- Machine washable and quick to dry
- Silver impregnated X-static® stretch fabric naturally inhibits microbial growth
- Extra layer of microclimatic spacer fabric for comfort
- Safe non-toxic material



Outer and inner cover

### Inner Cover

- Easy to wipe off
- Aquaguard™ zipper resists moisture
- Anti-wicking seam thread protects foam base



Dual layered contoured foam base




Easy washing and cleaning

Cushion Width		20	22	24	26	28	30	32	34	36	38	40
Cushion Depth	20	■	■	■								
	22	■	■	■	■							
	24	■	■	■	■	■						
	26	■	■	■	■	■	■					
	28	■	■	■	■	■	■	■				
	30	■	■	■	■	■	■	■	■			
	32	■	■	■	■	■	■	■	■	■		
	34	■	■	■	■	■	■	■	■	■	■	
	36	■	■	■	■	■	■	■	■	■	■	■
	38	■	■	■	■	■	■	■	■	■	■	■
	40	■	■	■	■	■	■	■	■	■	■	■
	42	■	■	■	■	■	■	■	■	■	■	■
	44	■	■	■	■	■	■	■	■	■	■	■
46	■	■	■	■	■	■	■	■	■	■	■	

These sizes are available, please see your order form for modifications or special sizes.

Height (in cm)		Maximum Height
Front	Rear	6
5	4.5	

Cushion weight	0.9kg
Max. user weight	75kg



### Sunrise Medical Limited

Thorns Road  
 Brierley Hill  
 West Midlands  
 DY5 2LD  
 Tel.: +44 (0) 845 605 66 88  
 Fax: +44 (0) 845 605 66 89  
 www.SunriseMedical.co.uk