

## JOKER R2



<b>Seat width</b>	33 - 36 - 39 - 42 - 45 - 48 cm
<b>Seat depth</b>	35 - 37,5 - 40 - 42,5 - 45 - 47,5 - 50 cm
<b>Seat height</b>	Front from 44 to 53 cm (adjustable) Rear from 36 to 47 cm (adjustable)
<b>Footplate distance</b>	From 35 to 49 cm (adjustable)
<b>Backrest height</b>	From 24,5 to 47 cm (adjustable)
<b>Front frame angle</b>	Aluminium 90°
<b>Total length</b>	88 cm approx. (std conf.)
<b>Total width</b>	Seat width + 17 cm (with 0° camber)
<b>Backrest angle</b>	From 74° to 94° to the seat (adjustable)
<b>Camber</b>	0° - 2° - 4°
<b>Weight without rear wheels</b>	Starting from 6,0 kg approx.
<b>Wheelchair weight</b>	Starting from 9,0 kg approx.
<b>Maximum load</b>	125 kg