

Sales Account Information	Distributor/Dealer Information
Name:	Order No.
Email:	Order Date:
Phone:	Company Name:

LEON F KAMELEON


Size 1 **Seat width: 250 ~ 300 mm**



Size 2 **Seat width: 350 ~ 400 mm**



anatomic measurements



User's name:	
User's weight (kg):	
User's height (mm) :	
① Top of head (mm) :	
② Top of shoulders (mm) :	
③ Seat pan to elbow (mm) :	
④ Elbow to hand (mm) :	
⑤ Low Leg Length (mm):	
⑥ Depth (mm):	
⑦ Width (mm):	
⑧ Chest width (mm):	
⑨ Chest height (mm) :	
⑩ Shoulder width (mm):	

Standard Base Price: £8,750

Please make selections by checking the corresponding boxes.




◆ Standard feature ○ Alternative options (Free of charge)

+\$ Paid Options (pricing valid for purchase with whole chair only)

Chassis: Leon FWD		
Colour		
<input type="checkbox"/>	Yellow	○
<input type="checkbox"/>	Rose Red	○
<input type="checkbox"/>	Navy Blue	○
<input type="checkbox"/>	Diamond Black	○

Speed			
	<input type="checkbox"/> 6 km/h * <i>Gyro is optional</i>	◆	
	<input type="checkbox"/> 8 km/h with Gyro * <i>Gyro is needed, please select gyro</i>	○	
	<input type="checkbox"/> 10 km/h with Gyro * <i>Gyro is needed, please select gyro</i>	○	
	<input type="checkbox"/> 12 km/h with Gyro * <i>Gyro is needed, please select gyro</i>	○	
Gyro (only available in combination with R-Net)			
	<input type="checkbox"/> No Gyro	◆	
	<input type="checkbox"/> Gyro * <i>Required for 8 km/h 10 km/h and 12 km/h</i>	£450	
Batteries (lead acid)			
	<input type="checkbox"/> 2x 45 Ah MK Batteries	◆	
	<input type="checkbox"/> 2x 73.6 Ah MK Batteries	£300	
Note: Battery box dimension 348.5 x 265 x 225 mm (LxWxH)			
Charger			
	<input type="checkbox"/> Charger 6A (<i>suggested for batteries 50 Ah</i>)	◆	
	<input type="checkbox"/> Charger 8A (<i>suggested for batteries 70 - 80 Ah</i>) For TUV Certificate	£75	
Lights			
	<input type="checkbox"/> No lights	◆	
	<input type="checkbox"/> Lights	£325	
Driving Wheel and Caster			
Rim	<input type="checkbox"/> Black rims	◆	
Colour	<input type="checkbox"/> Black tyre	◆	
Front tyre (Driving wheel)	<input type="checkbox"/> Air filled	◆	
	<input type="checkbox"/> PU tyre	£50	
Rear tyre (Caster)	<input type="checkbox"/> Air filled	◆	
	<input type="checkbox"/> PU tyre	£20	
Chassis Accessories (multiple choice)			
	<input type="checkbox"/> Bumper wheels on rear end	£50	
	<input type="checkbox"/> Dahl docking station interface	£100	
Seat Mounting & Function			
	<input checked="" type="checkbox"/> Powered seat fix * <i>coming soon</i>	◆	
	<input type="checkbox"/> Powered seat tilt, 45° <i>coming soon</i>	TBD	
	<input type="checkbox"/> Powered seat lift & tilt	£2,490	
Seat Base (colour is equal to chassis colour)			
Seat frame size	<input type="checkbox"/> Size 1 seat frame, standard width 250 mm (<i>choose seat configuration of Size 1 seccion</i>)	○	
	<input type="checkbox"/> Size 2 seat frame, standard width 350 mm (<i>choose seat configuration of Size 2 seccion</i>)	○	
Seat	Seat width	Seat depth	
Size 1	<input type="checkbox"/> 250 mm <input type="checkbox"/> 275 mm <input type="checkbox"/> 300 mm	<input type="checkbox"/> 250 mm	○
		<input type="checkbox"/> 270 mm	○
		<input type="checkbox"/> 290 mm	○
		<input type="checkbox"/> 310 mm	○
		<input type="checkbox"/> 330 mm	○
		<input type="checkbox"/> 350 mm	○
		<input type="checkbox"/> 370 mm	○
		<input type="checkbox"/> 390 mm	○
Size 2	<input type="checkbox"/> 350 mm <input type="checkbox"/> 375 mm <input type="checkbox"/> 400 mm	<input type="checkbox"/> 290 mm	○
		<input type="checkbox"/> 310 mm	○
		<input type="checkbox"/> 330 mm	○
		<input type="checkbox"/> 350 mm	○
		<input type="checkbox"/> 370 mm	○
		<input type="checkbox"/> 390 mm	○
		<input type="checkbox"/> 410 mm	○
		<input type="checkbox"/> 430 mm	○

Seat Cushion			
Upholstery			
	<input type="checkbox"/> No cushion	○	
	<input type="checkbox"/> 3D mesh	◆	
<i>* Seat cushion size is based on seat width and seat depth, which you selected above for seat.</i>			
Seat Supports			
Seat Foam Wedge (underneath the seat cushion)			
	<input type="checkbox"/> None	(£75)	
	<input type="checkbox"/> Foam wedges	◆	
Knee Side Support (side supports for knee outside (set) is combined with power legrest, no need to select here)			
	<input type="checkbox"/> None	◆	
	<input type="checkbox"/> Side supports for knee outside (set) PU foam	£175	
	<input type="checkbox"/> Side supports for knee inside (set) PU foam	£175	
Hip Supports			
	<input type="checkbox"/> None	◆	
	<input type="checkbox"/> Seat mounted hip supports (set) PU foam	£250	
Seat cushion	Foam wedges	Knee side supports	Hip supports
Backrest			
Frame size	<input type="checkbox"/> Size 1 Backres frame, standard width 250 mm	○	
	<input type="checkbox"/> Size 2 Backrest frame, standard width 350 mm	○	
Mounting	<input type="checkbox"/> Biomechanical backrest	◆	
Reclining function	<input type="checkbox"/> Manual recline	◆	
	<input type="checkbox"/> Power recline <i>* coming soon</i>	TBD	
Top panel set	<input type="checkbox"/> Standard panel	◆	
	<input type="checkbox"/> Shoulder panel	£50	
Middle panel set	<input type="checkbox"/> Standard panel	◆	
	<input type="checkbox"/> Lateral support panel with lateral supports included	£550	
Bottom panel set	<input type="checkbox"/> Standard panel	◆	
Backrest Cushion			
Upholstery			
	<input type="checkbox"/> No cushion	○	
	<input type="checkbox"/> 3D mesh	◆	
<i>Note: Backrest cushion size is based on the backrest panels, which you selected above for backrest.</i>			
Armrest			
Armrest System			
	<input type="checkbox"/> Fixed armrest, backrest mount	◆	
	<input type="checkbox"/> Biomechanical armrest, backrest mount	£190	
<i>Note: Biomechnaical armrest stays in a comfortable position when the backret is reclined.</i>			
Armrest Pad			
	<input type="checkbox"/> Armpad 215 mm (fits to Size 1)	○	
	<input type="checkbox"/> Armpad 300 mm (fits to Size 2)	○	

Legrest (choose Central System or Separate System)		
Central System		
Adjustment	<input type="checkbox"/> Mechanical adjustable central legrest	◆
	<input type="checkbox"/> Powered adjustable central legrest	£825
Footplate size	<input type="checkbox"/> One-piece footplate 315*216 mm (fits to Size 1)	○
	<input type="checkbox"/> One-piece footplate 385*230 mm (fits to Size 2)	○
Electronics & Controls (Choose Side Control or Tray Table Control for the user)		
Electronics		
	<input type="checkbox"/> VR2 electronics 90 A (<i>max. 2 power functions and speed of 6 km/h</i>)	◆
	<input type="checkbox"/> R-net electronics 120 A (<i>more than 2 power functions, gyro and special controls</i>)	£1,160
Side Control		
Joystick type	<input type="checkbox"/> VR2 LED	◆
	<input type="checkbox"/> LED JSM (<i>R-Net</i>)	£160
	<input type="checkbox"/> CJSM (<i>R-Net</i>)	£325
Mounting position	<input type="checkbox"/> Left side mount	○
	<input type="checkbox"/> Right side mount	◆
Mounting type	<input type="checkbox"/> Fixed mounting	◆
	<input type="checkbox"/> Parallel swing away	£160
Headrest		
Headrest Model		
	<input type="checkbox"/> None	◆
	<input type="checkbox"/> Lectus Standard 3D mesh	£225
	<input type="checkbox"/> Lectus Standard Dartex (water repellent)	£240
	<input type="checkbox"/> Contoured headrest 3D mesh	£240
	<input type="checkbox"/> Multi adjustable headrest 3D mesh	£350
Headrest Accessories		
	<input type="checkbox"/> None	◆
	<input type="checkbox"/> Headrest rail (out of centre positioning)	£50
		
Lectus standard headrest	Contoured headrest	Multi adjustable headrest
Pelvic Positioning Belts		
	<input checked="" type="checkbox"/> Pelvic positioning belt, clip <i>*required item</i>	◆
Positioning Belts		
	<input type="checkbox"/> None	◆
	<input type="checkbox"/> Vest 4-point	£95
Accessories		
	<input type="checkbox"/> None	◆
	<input type="checkbox"/> Rear basket for size 1	£160
Photo		
	Vest 4-point belt	